# 851B SPECTRUM ANALYZER DISPLAY SECTION

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquines@maunitron.co.uk

## OPERATING AND SERVICE MANUAL





TEST EQUIPMENT

### OPERATING AND SERVICE MANUAL

### MODEL 851B

# SPECTRUM ANALYZER DISPLAY SECTION

SERIALS PREFIXED: 526 -

S. No 526-01388

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:-01844-351694 Fax:-01844-352554
Email:- enquiries@mauritron.co.uk

Copyright HEWLETT-PACKARD COMPANY 1965 1501 PAGE MILL ROAD, PALO ALTO, CALIFORNIA, U.S.A.

02148-1


Printed: JUNE 1965

# MANUAL CHANGES

MODEL 851B

#### SPECTRUM ANALYZER - DISPLAY SECTION

Manual Serial Prefixed: 526-Manual Printed: June 65



MAKE ALL CORRECTIONS IN THIS MANUAL ACCORDING TO ERRATA BELOW, THEN CHECK THE FOLLOWING TABLE FOR YOUR INSTRUMENT SERIAL PREFIX (3 DIGITS) OR SERIAL NUMBER (8 DIGITS) AND MAKE ANY LISTED CHANGE(S) IN THE MANUAL.

► NEW ITEM.

SERIAL PREFIX OR NUMBER	MAKE MANUAL CHANGES	SERIAL PREFIX OR NUMBER	MAKE MANUAL CHANGES
526-00101 & above	Errata 1	► 526-00678 & above	Errata 1, Change 1, 2, 3, 4
526-00126 & above	Errata 1, Change 1		
526-00226 & above	Errata 1, Change 1, 2		
►526-00536 & above	Errata 1, Change 1, 2, 3		

ERRATA 1 Page 1-0, Table 1-1:

Change Vertical Display to read:

Vertical Display (7 cm full scale deflection):

Mode

Scale Factor

LINEAR

Relative Voltage/cm

SQUARE LOGARITHMIC Relative Power/cm 10 db/cm calibrated

over 0 to 60 db on

CRT display

For Service Manuals Contact MAURITRON TECHNICAL SERVICES 8 Cherry Tree Rd, Chinnor

Accuracy Oxon OX9 4QY

±3% full scaleTel:-01844-351694 Fax:-01844-352554 ±5% full scale\* Email:-enquiries@mauritron.co.uk

<±0.1 db/db but not more than ±2 db over full calibrated 60 db CRT display range\*

Change Power to read: 115 or 230 volts  $\pm 10\%$ , 50 to 400 cps, less than 55 watts.

#### Add:

External Sweep:

Input: 0 to +15 volt external signal (from 10K ohm source impedance) results in full 10 cm CRT horizontal trace. BNC female connector on rear panel, direct-coupled.

Blanking: -5 volt external blanking signal required to blank retrace. BNC female connector on rear panel.

Output Signals: Vertical and horizontal signals applied to CRT are available for external applications. Rear panel BNC female connectors. IF TEST POINT (20 Mc) also provided; rear panel BNC female connector.

Vertical: 0 to approximately -4 volts, open circuit; 4700 ohms source impedance.

Horizontal: 10 volts p-p  $\pm 0.3$  volt, open circuit; sweep approximately symmetrical about 0 volts. Source impedance 4700 ohms.

Page 3-3, Figure 3-2, changes:

- 7. change "negative 5- to 10-volt pulse," to "negative 4- to 10-volt pulse,"
- 10. change first three lines:

J10: signal to CRT, sampled at output of video detector following 20MC IF Amplifier, and just ahead of Vertical Amplifier;

12. change to read:

J8: sweep voltage, sampled just ahead of Horizontal Amplifier; 10 volts ±0.3V peak-to-peak open circuit, 4700 ohms impedance; BNC female.

Note: VERT and HORIZ outputs will drive high-impedance X-Y recorder to obtain an X-Y plot of spectrum displayed on CRT.

<sup>\*</sup>Except pulse spectra on 1MC IF bandwidth

0

- ERRATA 1 Page 5-6, Table 5-7, Logarithmic, change to read:

  (Cont'd) Logarithmic: < ±0.1 dB dB but not more than ±2 dB over full calibrated 60 dB CRT display range (except pulse spectra on 1MC IF bandwidth).

  Steps e, f, g, and h: delete adjustment of signal level.
  - ▶ Page 5-7, Table 5-7, 2. I. F. BANDWIDTH ACCURACY, Step a. change to read: a. Set VERT Display to LIN. Find 2-Gc BWO signal; see Paragraphs 5-92 through 5-95.
  - ► Page 5-8, Table 5-7, 3. I. F. SENSITIVITY, Step d, add: VERT DISPLAY . . . . . . LIN
  - ▶ Page 5-11, Paragraph 5-19. LOG.

Change specification to: =0.1 dB dB but not more than +2 dB over full calibrated 60 dB CRT display range.

Change Step e to read: Step I. F. GAIN through the rest of its positions without readjusting signal level. Limits are given in Table 5-8.

- ▶ Page 5-11, Paragraph 5-22, change Step a to read:
   a. Connect Attenuator 355D between 851 and 8551, set VERT DISPLAY to LIN, and find...5-100.
- ► Page 5-12, Paragraph 5-26, Step b, add: VERT DISPLAY . . . . . LIN

Page 5-18, VERTICAL DISPLAY: 7th line, change A11R13 to A11R14 9th line, change A11R14 to A11R13

- ▶ Page 5-28, Paragraph 5-94. Substitute following procedure for that given in the Manual:
  - 5-94. 10KC, 3KC, and 1KC ALIGNMENT CHECKS.
  - 5-95. Signals for the three narrower bandwidth filters (10, 3, and 1 kc) pass through two double-tuned crystal filters. The four tuning coils are tapped; change of bandwidth is obtained by using different taps. The same filters are used for all three bandwidths; accurate adjustment of the 10-kc bandwidth should bring the 3-kc and 1-kc bandwidths within specifications. After adjustment of the 10-kc bandwidth, bandwidth is verified at the 3-kc and 1-kc settings.
  - 5-96. IF bandwidth alignment is not a simple technique. While tuning for correct IF bandwidth, remember:
  - a. Ideally, all adjustments should be made simultaneously. Since this is impossible, it will be necessary to repeat the adjustments more than once to obtain the best tuning of the four filters.
  - b. Final adjustment should be the compromise which obtains the best characteristics for all four filters. Do not attempt to adjust the filters unless one or more are out of specifications.

#### 5-97. EQUIPMENT REQUIRED.

Ref No.	Equipment Required	No.			
10*	VHF Attenuator (355D)	1			
D**	Coax Term. w/BNC males (10503A)	2			
G**	GC plastic tuning wand	1			
K**	Screw-holding screwdriver	1			
*Table 5-1 **Table 5-2					

#### 5-98. SIGNAL SOURCE CALIBRATION.

5-99. To check the bandpass characteristics of the narrower IF filters the 851 sweep width must be narrow enough that the IF bandwidth can be determined accurately at the half-power points. This may be done by applying a signal to the 20MC IF which is swept in synchronism with the 851 sweep. Such a signal can be derived from the second harmonic of the 8551's 10MC Reference Oscillator.

remarkation of the confidence of the confidence

### ERRATA 1 (cont'd)

5-99A. When the 8551 is stabilized, the BWO is phase-locked to a 10MC Reference Oscillator. For a BWO frequency of 4 Gc and a spectrum width of 1 Mc/cm the Reference Oscillator would be swept 2.5 kc/cm (BWO locked to the 400th harmonic of 10 Mc; 1 Mc/cm divided by 400 is 2.5 kc/cm). [NOTE: For a 4-Gc BWO frequency, FREQUENCY(GC) must be at .01-2 or 1.8-4.2; (at these settings n = 1).] If the output of the Reference Oscillator is connected to a narrow-band IF filter tuned to 20 Mc, the IF will pass only the second harmonic of the 10MC Reference Oscillator. This is 20 Mc swept at 5 kc/cm. Other values of sweep width may similarly be derived:

Spectrum Width	20MC Sweep Width (851 Display)
1 Mc/cm	5 kc/cm
300 kc/cm	1.5 kc/cm
100 kc/cm	500 cycles/cm
30 kc/cm	150 cycles/cm

#### 5-100. MEASUREMENT SETUP.

5-101. Use the 8551 10MC Reference Oscillator as the signal source for the narrower IF bandwidth alignment procedures. See Figure 5-8A for test setup and Paragraph 5-97 for recommended equipment.

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk




Figure 5-8A. Measurement Setup, 10KC IF Bandwidth Alignment and Checks

#### ERRATA 1 (Cont' d)

#### 5-102. 10 KC PROCEDURE, INITIAL SETUP.

a. Set the 355D to 40 dB.

b. Make the following settings:

8551

LINE. . . . . . . . . . . . . . . . . . STANDBY\*

SPECTRUM WIDTH . . . . . . . . . . . 1 MC/CM

\*Note: LINE remains in STANDBY throughout the procedure.

8551A only

FREQUENCY TUNING . . . . . . . . . STABILIZE\*\*

8551B only

TUNING SELECTOR . . . . STABILIZED NORMAL\*\*
STABILIZATION . . . . . . . . . . . STABILIZED\*\*

\*\*Note: Control setting only; do not perform stabilization procedure.

851

 BASE LINE CLIPPER
 max ccw

 SYNC.
 INT

 I. F. BANDWIDTH
 10 KC

 VERT DISPLAY
 LIN

 SWEEP TIME
 3 MILLISEC/CM

 SWEEP TIME VERNIER
 CAL

 INTENSITY
 about 3 o' clock

 IF GAIN
 .30 + 0

 IF VERNIER
 ccw

c. Check alignment of the base-line trace with the horizontal axis. If necessary, adjust VERT POS and TRACE ALIGN to bring base-line trace exactly parallel with and on the graticule base line.

#### 5-103. 10 KC ALIGNMENT PROCEDURE.

- a. Adjust 8551 TUNE to center the display on the 851. Adjust IF GAIN VERNIER for a maximum vertical deflection of exactly 7.0 cm.
- b. Bandwidth tuning adjustments are inside the RF Circuit Assembly casting (see Figure 5-22); location of adjustments is marked on the cover. Access holes, covered with removable plug-in buttons, are provided in the casting cover. Unless Balance Adj capacitor A2A3C5 or A2A4C8 has been replaced, do not remove the casting cover.

#### Note

It is not likely that capacitor A2A3C5 or A2A4C8 will require replacement. However, if either has to be replaced, before removing it, note degree of mesh between stator and rotor. When installing replacement capacitor, set it to approximately the mesh of original capacitor. After installing and presetting replacement capacitor, fasten cover to casting with five or six of the 26 screws which hold the casting cover in place.

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk

Perform the rest of the 10KC alignment procedure with the cover in place on the casting.


c. Adjust 1-10KC Bandwidth Adj capacitors A2A3C4, A2A3C2, A2A4C5, and A2A4C9 for maximum bandwidth.

#### Note

In tuning capacitor A2A3C4, A2A3C2, or A2A4C5 through its tuning range it will be found there are two points which give vertical deflection peaks. Since there is little difference between the amplitude of the two peaks, it is difficult to distinguish which is the correct tuning region. If correct IF bandwidth tuning cannot be obtained on one peak, try the the other. Correct IF bandwidth tuning can only be obtained

ERRATA 1 (Cont' d) when the adjustment of each capacitor is made in its true tuning region. Maximum bandwidth is usually obtained by tuning off the peak slightly.

- d. Adjust Imped Adj A2A3L3 and Frequency Adj A2A4C7 for maximum vertical deflection.
- e. Center display with TUNE and set maximum vertical deflection to exactly 7.0 cm with IF GAIN. See Figure 5-8B. Display should be 2 cm wide at 5 cm amplitude (half-power points) (sweep width of 851 display is 5 kc/cm). If not within  $\pm 20\%$  of the correct bandwidth (1.6 to 2.4 cm at 5 cm amplitude) repeat steps c through e until the correct bandwidth is obtained.



For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk

Figure 5-8B. Optimum Bandpass Characteristics, IF Bandwidth Adjustments

#### 5-104. 3KC AND 1KC BANDWIDTH CHECK PROCEDURE.

- a. Set IF BANDWIDTH to 3KC and SPECTRUM WIDTH to 300 KC/CM (which gives 851 display sweep width of  $1.5\ kc/cm$ ).
  - b. Adjust IF GAIN and IF VERNIER for maximum vertical deflection of exactly 7.0 cm.
  - c. Width of display at 5.0 cm axis should be between 1.6 and 2.4 cm. See Figure 5-8B.
- d. Set IF BANDWIDTH to 1KC and SPECTRUM WIDTH to 100 KC/CM (which gives 851 display sweep width of 500 cycles/cm).
  - e. Adjust IF GAIN and IF VERNIER for vertical deflection of exactly 7.0 cm.
  - f. Width of display at 5.0 cm axis should be between 1.6 and 2.4 cm. See Figure 5-8B.

#### Note

If 1KC bandwidth appears too wide, recheck tuning of Freq Adj A2A4C7 (Paragraph 5-103, step d).

g. If capacitor A2A3C5 or A2A4C8 was replaced, and casting cover is only partly secured, fasten in place with all 26 screws. For final adjustment (Paragraph 5-105), the cover must be tightly fastened to the casting.

#### Note

A screw-holding screwdriver is recommended for turning the screws.

#### 5-105. FINAL 1 - 10KC BANDWIDTH ADJUSTMENT.

5-106. Set IF BANDWIDTH to 10KC and SPECTRUM WIDTH to 1 MC/CM. Recheck bandwidths (Paragraphs 5-103, 5-104), making adjustment if necessary, until all bandwidths are within specifications.

#### Note

Cover must be fastened down tightly during final adjustment.

#### ERRATA 1 (Cont' d)

Page 5-31, Paragraph 5-119b, change A11R13 to A11R14 Paragraph 5-119c, change A11R14 to A11R13

Page 5-36, Table 5-23A: Change Short Ckt Current for Triplett 630 from "32 ma" to "3.25 ma" for R x 100 range "3.25 ma" to "325  $\mu$ a: for R x 1K range

Page 5-52, Figure 5-19, change PREFIX ALL DESIGNATIONS WITH A1 to PREFIX ALL DESIGNATIONS WITH A12.

Page 5-55, Figure 5-24: In lower middle, add asterisk to A2A4R4.

Page 5-56, Figure 5-25, change R13 designation from 40DB LOG CALIB to 60DB LOG CALIB R14 designation from 60DB LOG CALIB to 40DB LOG CALIB

Page 5-57, Figure 5-27, VERT DISPLAY Switch:
Upper left corner of schematic, change A11C1 from 220 to 300.
Change R13 designation from 40DB to 60DB and R14 designation from 60DB to 40DB.
Add A11CR5 from cable W1, with cathode to ground.
Lower right corner, change A2C2 from 33 to 3.3 pf.

Page 5-59, Figure 5-29:
Upper right corner, add asterisk to 1000-ohm R19.
Add to Notes: \* = Factory selected; average value shown.

Page 5-61, Figure 5-33:
Lower left corner, add asterisk to A6R35, and change value to 39K.
Lower right corner, add asterisk to R59.
Upper right corner, change A6R66 from 33K to 68K.
Lower middle, change A6R43 to 470.
Add to NOTES: \* = factory-selected value; part may be omitted.

Table 6-1, change to read:

A6R43 0683-4715 R:fxd comp 470 ohm 5% 1/4W (Note: A6R43 is 1000 ohms in instruments with serials below 526-00201: however, 470 ohms is the proper replacement for all instruments.)

A6R66 0684-6831 R:fxd 68K 1/4W

A11C1 0140-0225 C:fxd 300 pf 1% 300 vdcw A11CR5 1901-0033 Diode V1 5083-9010 Electron tube: cathode ray, P2 phosphor V1 5083-9011 Electron tube: cathode ray, P7 phosphor V1 5083-9012 Electron tube: cathode ray, P31 phosphor

CHANGE 1 526-00126 & above Page 5-59, Figure 5-29, Vertical Amplifier Schematic: Lower left corner, change A7C2 from 0.47  $\mu f$  to 2.2  $\mu f$ .

Table 6-1, change to read: A7C2 0180-0155 C:fxd 2.2 μf. For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:-01844-351694 Fax:-01844-352554
Email:-enquiries@mauritron.co.uk

CHANGE 2 526-00226 & above Page 5-59, Figure 5-29, Vertical Amplifier Schematic: Lower left corner, add A7R26, 1000 ohms, in series with A7R5.

Table 6-1, add: A7R26 0687-1021 R:fxd 1000 ohm 1/2W. ► CHANGE 3 526-00536 & above

Page 5-63, Figure 5-35, HV Power Supply Schematic: Upper left, place asterisk on A8R1.

Add Note: \* = Factory-selected value; average value shown

► CHANGE 4 526-00678 & above

Page 5-55, Figure 5-24, IF Bandwidth Switching Circuits: Upper left, change A2A2R2 from 1500 to 1200 ohms.

Page 5-57, Figure 5-26, VERT DISPLAY Switch, etc. Schematic: Lower center, change A2A7R4 from 51 to 100 ohms.

Table 6-2, change to read:
A2R2 0683-1225 R:fxd comp 1200 ohm 5% 1/4W
A7R4 0683-1015 R:fxd comp 100 ohm 5% 1/4W

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk

#### TABLE OF CONTENTS

Se	ction	Page	Se	ction		Page
I	GENERAL INFORMATION	1-1	V	MAINTENANO	CE (Cont'd)	Lago
	1-1. Introduction				ntal Amplifier Checks and	
	1-3. Description	1-1			stments	. 5-2
	1-5. Applications	1-2			oration	
	1-7. Cathode-Ray Tube				p Time Calibration	
					EP TIME Vernier Check .	
П	INSTALLATION	2-1			e Sweep and Sweep	
	2-1. Initial Inspection	2-1		0	aplitude Checks	. 5-2
	2-3. Mechanical Check	2-1		5-61. Synchro	onization and Output Checks	5 5 2
	2-4. Performance Check	2-1		5-69. CRT Ch	necks	5 - 5 - 20
	2-6. Claim for Damage		•	5-75. Vertica	l Amplifier Checks and	. 0-20
	2-8. Connections	2-1		Adius	stments	5_94
	2-10. Power Requirements	2-1		5-77. Verti	ical Calibration	. 5-24
	2-15. Establishing Figure-of-Merit Rating	2-1		5-82. 1MC I	F. Bandwidth Alignment	. 5-2-
	2-18. Rack Mounting	2-1			Check	. 5-25
	0			5-87 100KC	I. F. Bandwidth Alignment	. 5-2.
					Check	. 5-27
Ш	OPERATION	3-1		5-91 Final 11	MC and 100KC Bandwidth	. 3-4
	3-1. Introduction	3-1		Ading	stments	. 5-27
	3-7. Description	3-1		5-94 10KC 3	3KC, and 1KC I.F.Band-	. 3-21
	3-9. I. F. Bandwidth	3-1		width	Alignment and Checks	5 20
	3-12. Display	3-1		5-99. 10KC	Procedure	. 5-29
	3-15. Oscilloscope	3 – 1			and 1KC Bandwidth Checks	
	3-20. Operating Instructions				1-10KC Bandwidth Adjustm	
	3-22. Photographic Techniques			5-107 AUTO S	ELECT Check	
	•			5-111 I F Sen	sitivity Check	• 5-3U
				5-116 VERT D	DISPLAY Checks and	. 5-51
IV	PRINCIPLES OF OPERATION	4_1			tments	E 21
	4-1. Introduction	<u>1-1</u> 41		5-117. Prelin	minary Check	. 0-01
	4-5. Operation of Horizontal Sweep	7-1		5-118. LOG I	Display	. 5-31
	Generator	4.4		5-121. SQ Di	splay	. 0-31
	4-7. External				Display Linearity Check	. 0-01
	4-13. Internal	4 5		5-123. Einel I	F. Bandwidth Adjustments.	. 5-32
	4-15. Operation of Vertical Display	4-5		5-124. Final I.	r. Bandwidth Adjustments.	. 5-32
	4-16. Current-Controlled Attenuator.	4-5		D-120. Meast	arement Setup and Initial	- 00
	4-20. Square Mode of Operation	4-5		5-127. Crysta	ocedure	. 5-32
	4-24. Logarithmic Mode of Operation .	4-J		Dole	al Filter (10KC, 3KC, 1KC	
	2 21. 20gartamire mode of Operation .	4-0			ance	
					Bandpass Filter Adjustment	
V	MAINTENANCE	5-1		5-131. Poubles 5-132. Locali	shooting	. 5-33
	5-1. Introduction			5-134. Parts	ization	. 5-33
	5-5. Content	5_1		5 126 Tabletine	Location	. 5-33
	5-6. Cover and Side Panel Removal	5-1			Trouble in Transistor	- 00
	5-7. Performance Checks	5-1		5 140 In Cincu	its	. 5-33
	5-9. Operational Checks	5-1			it Testing of Transistors	<b>5</b> 00
	5-14. Vertical Display Accuracy Check	5-3			l, Q5, Q6	
	5-20. I.F. Bandwidth Accuracy Check	5-11		5-142. Replacer	ment of Cathode-Ray Tube	. 5-36
	5-24. I. F. Input Sensitivity Check	5-12		5 150 Damasist	tor Replacement	. 5-36
	5-27. I. F. GAIN Set Accuracy Check	5 12		5-150. Removin	ng I. F. GAIN Switch Assy A	1 5-39
	5-31. Sweep Rate Accuracy Check	5 12		5-152. Removin	ng Bandpass Filter Assy A1	2 5-41
	5-35. 851B Inspection and Adjustment	0-10		5 156 Description	g Switches	. 5-41
	Test Card	5-14		5-150 Removin	g Assemblies A3 and A5.	. 5-41
	5-37. Checks and Adjustments			5-108. Etched C	Circuits	5-41
	5-39. Preliminary Adjustment Procedure.	5 10		บ-103. Transist	ors	. 5-43
	5-40. LV Power Supply Adjustments	5 10	777	DEDI ACEAST	E Di Des	_
	5-43. HV Power Supply Check	2 20	VI	TEPLACEABLI	E PARTS	6-1
	5-46. Focus Check and Adjustment	5 91		6-1. Introduct	tion	6-1
	5 20. 20005 Oncor and Adjustinent	U-41		<b>6-5.</b> Ordering	Information	. 6-1

#### LIST OF ILLUSTRATIONS

Numb	per Title	Page	Numb	er Title J	Page
1-1.	Model 851B Spectrum Analyzer -	Ü	5-12.	Cathode-Ray Tube, Parts and Connections	3
	Display Section	1-1		Involved in Removal	5_37
2-1.	Installation Connections, Model 851/8551		5-13.	Top View of 851 Display Section	5-38
	Spectrum Analyzer		5-14.	Bottom View of 851 Display Section	5-39
2-2.	Rack Mounting Procedure, 851B			Under Side of Transistor/Transformer	
3-1.	Front Panel Controls, Connectors, and			Deck, Identification of Q3-Q6	
	Indicators	3-2			5-41
3-2.	Switches and Connectors, Rear Panel	3-3	5-16.	Transistor Biasing and Typical Oper-	
3-3.	Initial Operating Procedure for 10-Mc				5-42
	to 10-Gc Inputs	3-4	5-17.	CRT Shield Assembly, Parts	
3-4.	Photographic Procedure	3-5		W 1 1 4 4 4 1 1 1 1 1 1 1 1 1 1 1 1 1 1	5-50
3-5.	Optimum 851 I.F. BANDWIDTH Settings.	3-6	5-18.	I. F. GAIN (DB) Switch Assembly A1,	
4-1.	Model 851B/8551A Spectrum Analyzer,			~	5-52
	Block Diagram	4-2	5-19.	Bandpass Filter Assembly A12,	
4-2.	Model 851B Spectrum Analyzer Display			Component Identification	5-52
	Section, Block Diagram	4-3	5-20.	20MC I. F. Input and Attenuator	
4-3.	Horizontal Sweep Generator, Block			0.1 II 0=4=	5-53
	Diagram	4-4	5-21.	RF Circuit Assembly Boards A2A2,	
4-4.	Typical Dynamic Resistance-vs-Current				5-54
	Curve, Hot Carrier Diodes	4-5	5-22.	RF Circuit Assembly A2, Top Cover	
4-5.	VERT DISPLAY Switch at SQ,				5-54
	Simplified Schematic	4-6	5-23.	RF Circuit Assembly Boards A2A5	
4-6.	VERT DISPLAY Switch at LOG,			and A2A1	5-54
	Simplified Schematic	4-7	5-24.		5-55
5-1.	Measurement Setup, Check and Adjustment	nt	5-25.		5-56
	of LV Power Supply	5-20	5-26.	RF Circuit Assembly Boards A2A6,	
5-2.	CRT Base as Seen from Rear, CRT Pro-			A2A7, and Rear of Casting	5-56
	tective Cover Removed	5-21	5-27.	VERT DISPLAY Switch, Current-	
5-3.	Measurement Setup, Sweep Calibration .	5-22		Controlled Attenuator, and 20MC	
5-4.	Pin-cushioning and Barrelling Defined .	5-24		I. F. Amplifier Schematics, 851B	5-57
5-5.	Measurement Setup, Calibration of 851		5-28.	Vertical Amplifier A7 Board	5-58
	Vertical Amplifier	5-25	5-29.	Vertical Amplifier Schematic, 851B	5-59
5-6.	Setup for 851 Video Bandwidth		5-30.	SWEEP TIME Switch A10S1	5-60
	Measurements	5-26	5-31.	Sweep and Horizontal Amplifier A6 Board	5-60
5-7.	Measurement Setup, 1MC and 100KC		5-32.	SYNC Switch S2	5-60
	I.F.Bandwidth Alignment and Checks.	5-27	5-33.	Sweep and Horizontal Amplifier	
5-8.	Measurement Setup, 10KC I. F. Bandwidth			Schematic, 851B	5-61
	Alignment and Checks	5-28	5-34.	HV Power Supply A8 Board	5-62
5-9.	Measurement Setup, Final I.F. Bandwidth		5-35.	HV Power Supply Schematic, 851B	5-63
	Adjustments	5-32	5-36.		5-64
o-10.	851B Spectrum Analyzer Display Section,		5-37.	LV Power Supply Schematic, 851B !	5-65
	Top View, Top Cover Removed	5-34	5-38.	I. F. BANDWIDTH Switch A4S1	5-66
0-11.	851B Spectrum Analyzer Display Section,		5-39.	851B SWEEP TIME Switch A10S1,	
	Bottom View, Bottom Cover Removed	5-35		Schematic	5_67

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk



Model 851B List of Tables

#### LIST OF TABLES

Numb	er Title	Page	Numbe	r Title Pa	age
1-1.	Specifications	1-0	5-2 <b>3</b> .	Typical Data for Out-of-Circuit	_
1-2.	Accessories Supplied	1-2		Transistor Resistance Measurement . 5	5-36
1-3.	Accessories Available	1-2	5-23A.	Safe Ohmmeter Ranges for Transistor	
1-4.	Options			Resistance Measurements	5-36
			5-23B.	Connection Point, Q3, Q4, Q5, Q6	
5-1.	Test Equipment Required	5-0			5-36
5-2.	Test Accessories Required	5-2	5-23C.	Recommendations, Component Removal 5	5-40
5-3.	Assemblies vs Schematics, 851B			Adjustments Required after Component	
5-4.	Controls, Switches, and Connectors vs			Replacement	5-40
	Schematics	5-3	5-23E.	Etched Circuit Soldering Equipment 5	
5-5.	Chassis Parts Locater			Connections, RF Circuit Assembly A2,	
5-6.	Operational Checks	5-5		· · · · · · · · · · · · · · · · · · ·	5-44
5-7.	Performance Check Test Card, 851B		5-24A.	Connections, RF Circuit Assembly A2,	
5-8.	VERT DISPLAY Accuracy Check	5-11		Boards A2A6, A2A7	5-45
5-9.	I. F. Bandwidth Accuracy Checks (1MC,		5-25.	Waveform Chart Model 851 Spectrum	
	100KC, 10KC)	5-11		Analyzer Display Section 5	5-46
5-10.	I. F. Bandwidth Accuracy Checks, 3KC		5-25A.		5-51
	and 1KC	5-12	5-26.	Connections, Sweep and Horizontal	
5-11.	I. F. Input Sensitivity Check	5-12		Amplifier Assembly A6, 851B 5	5-58
5-12.	Sweep Rate Accuracy Check	5-14	5-27.	Connections, Vertical Amplifier	
	851B Check and Adjustment Test Card .			Assembly A7 5	5-58
	LV Power Supply Measurement Data		5-28.	Connections, SWEEP TIME Switch	
5-15.	Resistances to Ground, LV Power Supply,			A10S1, 851B 5	5-62
	Reference Data	5-20	5-29.	Connections, HV Power Supply A8 Board 5	5-62
5-16.	HV Power Supply Voltages	5-21	5-30.	Connections, SYNC Switch S2, 851B 5	5-64
5-17.	Sweep Time Calibration	5-22	5-31.	Connections, LV Power Supply A9 Board 5	5-64
	Data for Video Bandwidth Check		5 - 32.	Connections, I. F. BANDWIDTH Switch	
	Switch Settings for AUTO SELECT Check			A4S1 5	5-66
			6-1.	Reference Designation Index, General . 6	3-2
	SQ Display Linearity Check Data			Reference Designation Index, Assy A2 . 6	5-14
-22.	LIN Display Linearity Check Data	5-32	6-3.	Replaceable Parts 6	3-20

Section I Table 1-1

Table 1-1. Specifications, 851B Display Section

with the same the same for all is in the health and the same with the same in the same and the same in the same

#### DISPLAY CHARACTERISTICS

Vertical Display:

Linear, square (power), or logarithmic Dynamic Range: Linear, 70:1; square, 70:1; log, 60 db Accuracy: Linear, ±3% of full scale; square, ±5% of full scale\*; log, ±2 db\*.

#### I. F. Bandwidth:

Manual: Bandwidths of 1, 3, 10, 100 Kc, and 1 Mc can be selected.

Auto Select: One of the above bandwidths automatically selected for best resolution of a CW signal with each combination of spectrum width and sweep rate.

Bandwidth Accuracy: Individual bandwidths are calibrated within ±20%; bandwidth repeatability and stability typically better than ±3%.

#### I. F. Input:

Center Frequency: 20 Mc
Input Impedance: 50 ohms, nominal
Input Required for 6-cm Vertical Display:
1-Mc bandwidth, -62 to -53 dbm
100-Kc bandwidth, -75 to -60 dbm
10-Kc bandwidth, -95 to -80 dbm
3-Kc bandwidth, -95 to -80 dbm
1-Kc bandwidth, -86 to -71 dbm

#### Maximum CW Input Signal: -14 dbm

- I. F. GAIN Set: Two-section attenuator provides 0 to 80 db attenuation in 1-db steps. One section provides 0 to 70 db attenuation in 10-db steps; the other, 0 to 10 db in 1-db steps. Vernier provides continuous adjustment between 1-db steps.
- I. F. GAIN Set Accuracy: 70-db section, ±0.5 db; 10-db section, ±0.1 db

Sweep Rate: Six calibrated rates from 3 msec/cm to 1 sec/cm in a 1, 3, 10 sequence. Vernier provides continuous adjustment between calibrated rates and extends slowest rate to at least 3 sec/cm.

Sweep Rate Accuracy: ±3%

#### Sweep Synchronization:

Internal: Sweep free runs

Line: Sweep synchronized with power-line

frequency

External: Sweep synchronized with externally

applied signal +3 to +15 volts peak

Single Sweep: Sweep actuated by panel pushbutton

\* Except pulse spectrums on 1-Mc I. F. bandwidth

#### GENERAL

Output Signals: Vertical and horizontal signals applied to CRT are available for external monitoring. Vertical: 0 to -4 volts; output impedance,  $4700\,\Omega$  Horizontal: 10 volts peak-to-peak,  $\pm 0.3$  volt open circuit; sweep approximately symmetrical about zero; output impedance,  $4700\,\Omega$ 


Cathode-Ray Tube: 7.5-kv post-accelerator tube with P2 medium-persistence phosphor and internal graticule. Light blue filter supplied. Other phosphors optional.

Internal Graticule: Parallax-free 7 x 10 cm, marked in cm squares with 2-mm subdivisions on major vertical and horizontal axes.

RFI: Conducted and radiated leakage limits are below those specified in MIL-I-6181D and MIL-I-16910.

Power: 115 or 230 volts  $\pm 10\%$ , 50 to 1000 cps, approximately 25 watts

#### Dimensions:



Weight: Net 34 lb (15 kg); shipping 45 lb (20, 3 kg)

Accessories Furnished: 7-1/2 foot (2290 mm) power cable; rack mounting kit; joining-bracket kit for bonding Model 851 to Model 8551

#### Options:

- 07. P7 phosphor in lieu of P2 (amber filter supplied), no additional charge.
- 31. P31 phosphor in lieu of P2 (green filter supplied), no additional charge.

# SECTION I GENERAL INFORMATION

#### 1-1. INTRODUCTION.

1-2. The Model 851B is the Display Section of the Hewlett-Packard Spectrum Analyzer; the RF Section is the Model 8551. Together, the 851/8551 provide an Analyzer which can display up to 2 Gc of spectrum. Analyzer input range is from 10 Mc to 42 Gc, its dynamic range is 60 db, its sensitivity is at least -65 dbm, its image separation is 4 Gc, and its functions are calibrated. The RF Section, which is a receiver that electronically scans the input signal, is described in a separate Operating and Service Manual.

#### 1-3. DESCRIPTION.

- 1-4. The 851B Display Section is an oscilloscope with an unusually wide dynamic range. The display is amplitude vs frequency, and range is such that a fundamental and harmonics down as far as 60 db can be viewed simultaneously. In addition to the features offered by other good oscilloscopes -- features such as calibrated sweep times, calibrated gain, choice of synchronizing voltages, vertical and horizontal positioning, focus adjustments, and intensity variation -- the 851B provides additional facilities which widen the scope of the Analyzer as an electronic tool.
- a. Choice of Amplitude Calibrations. The 851B display can be made 1) proportional to voltage (linear), proportional to power (square), or 3) proportional to the log of the input voltage. Use of logarithmic calibration is what makes it possible to view amplitude variations as great as 60 db in the same display.

- b. Choice of Calibrated I. F. Bandwidths. The 851B provides five calibrated I. F. bandwidths: 1 Mc, 100 Kc, 10 Kc, 3 Kc, and 1 Kc. The narrower bandwidths provide greater resolution and the wider bandwidths viewing of a broader band of frequencies.
- c. Automatic Selection of Optimum I. F. Bandwidth. Characteristics of the display are a function of the width of frequency band swept (determined by the setting of the SPECTRUM WIDTH switch in the 8551 RF Section), the I. F. bandwidth, and the sweep speed. For automatic selection of optimum bandwidth for selected sweep speeds and width of band swept, the 851B provides an AUTO SELECT position on the I. F. BANDWIDTH switch.
- d. Facilities for Making Superior Oscillograms. The 851B CRT has an internal graticule, thus providing a parallax-free presentation. In addition, when the internal graticule is illuminated by ultra-violet light, resulting photographs are of unusually fine quality. (The hp 196B Oscilloscope Camera includes a source of ultra-violet light.) Another feature of the Display Section is base-line blanking capability. This feature is useful both when viewing low-level signals or when making an oscillogram since features of interest are clearer when base line glow is blanked. For photographic convenience, the CRT includes a bezel for mounting a camera, and a SINGLE SWEEP lamp which indicates completion of the sweep.




Figure 1-1. Model 851B Spectrum Analyzer - Display Section

#### 1-5. APPLICATIONS.

1-6. Some of the many applications of the Hewlett-Packard Spectrum Analyzer are discussed in Application Notes 63 and 63A. One of these applications is analyzing short RF pulses. Because of features contributed by the 851B -- a choice of I. F. bandwidths, dynamic range of at least 60 db, and calibrated sweep speeds, the Analyzer is a valuable tool in pulse work. Short RF pulses (tens of nanoseconds) have previously been difficult to analyze in the frequency domain because of limitations in dynamic range and the I. F. bandwidth of available analyzers. The 1-Mc I. F. bandwidth of the 851 gives 11 db of additional dynamic range when measuring short pulses, additional by

comparison to a hypothetical system having 80-Kc I. F. bandwidth and equal CW dynamic range. With the 851B calibrated sweep times, pulse repetition rate can be determined directly from the display, obviating the need for measuring repetition rate externally.

#### 1-7. CATHODE-RAY TUBE WARRANTY.

1-8. The cathode-ray tube (CRT) supplied with the 851B is guaranteed against electrical failure by Hewlett-Packard for one year from the date of sale. Warranty claim and adjustment procedures for the CRT are given on the warranty at the rear of this manual. Use this form and follow claim instructions exactly when returning a CRT for warranty adjustment.

Table 1-2. Accessories Supplied

distinct of the second second

Part Number	Name	Description
8120-0078	Power cable	Standard 3-conductor 7-1/2 foot NEMA power cable
5060-0216	Joining Bracket Kit	Plates and hardware for bonding 851 to 8551
5060-0076	Rack Mounting Kit	Parts and hardware for mounting 851 in 19-inch rack

#### Table 1-3. Accessories Available

Model Number	Name	Description
8442	1 KC Bandwidth Crystal Filter	For use ahead of 851B I. F. input; extremely good skirt selectivity pass band less than 10 Kc 60 db down

#### Table 1-4. Options

Number	Description
07	In lieu of P2 phosphor, P7 long-persistence phosphor and amber filter supplied; no additional charge.
31	In lieu of P2 phosphor, P31 medium-persistence phosphor and green filter supplied; no additional charge.

# SECTION II

#### 2-1. INITIAL INSPECTION.

#### 2-2. MECHANICAL CHECK.

2-3. If damage to the shipping carton is evident, ask that carrier's agent be present when instrument is unpacked. Inspect instrument for mechanical damage such as scratches, dents, or broken knobs. Also check the cushioning material for signs of severe stress.

#### 2-4. PERFORMANCE CHECK.

2-5. The electrical performance of the 851B should be verified as soon as possible after receipt. Performance checks suitable for incoming inspection are given in Paragraphs 5-7 through 5-34.

#### 2-6. CLAIM FOR DAMAGE.

2-7. If the 851B is mechanically damaged or fails to meet specifications on receipt, notify the carrier and the nearest Hewlett-Packard office immediately. (A list of sales and service offices is at the back of this manual.) Retain the shipping carton and the padding material for the carrier's inspection. The field office will arrange for the repair or replacement of your instrument without waiting for the claim against the carrier to be settled.

#### 2-8. CONNECTIONS.

- 2-9. Connect the two Sections of the Analyzer:
- a. Place the Model 851 Display Section on the Model 8551 RF Section.
- b. A power cable is supplied with the 851, and five cables are supplied with the 8551. Connect cables as shown in Figure 2-1.
- c. To obtain the best common ground for the two Sections, strap the Sections together with the plates provided in the Joining Kit (supplied). Bonding instructions are supplied with the Kit.

#### 2-10. POWER REQUIREMENTS.

2-11. The 851B is designed to operate from either a 115- or 230-volt 50- to 1000-cycle source, and requires approximately 25 watts. However, when used as the Display Section of the Analyzer the line input for the 851B is in the 8551 RF Section, and power is extended to the 851B by external cable. The two Sections of the Analyzer require approximately 300 watts and a nominally 115- or 230-volt 50- to 60-cycle source.

- 2-12. Both Sections are equipped with input transformers. Primary windings on each input transformer can be connected in series or in parallel; changing from one type of connection to the other is by operation of a slide switch (115/230) located on the rear panel (see Figure 3-2). Always check the setting of the slide switches in both Sections before plugging the Analyzer into a power source; the setting of the 115/230 switch must agree with the voltage of the power source. Refer to Figure 2-1 for sequence of the plug-in procedure. (Sequence for turn-on is given in Figure 3-3.)
- 2-13. The fuse installed at the factory is for 115-volt operation. When operating from 230 volts, use a fuse of the value shown adjacent to the 230-volt position of the slide switch.
- 2-14. To protect operating personnel, the National Electrical Manufacturers' Association (NEMA) recommends that instrument panel and cabinet be grounded. The Analyzer is equipped with a three-conductor power cable; the third conductor is the ground conductor, and when the cable is plugged into an appropriate receptacle, the instrument is grounded. The offset pin on the power cable three-prong connector is the ground connection. To preserve the protection feature when operating the instrument from a two-contact outlet, use a three-prong to two-prong adapter, and connect the green lead on the adapter to ground.


# 2-15. ESTABLISHING FIGURE-OF-MERIT RATING.

- 2-16. Immediately following initial inspection, it is good practice to establish a figure-of-merit rating for your 851B Display Section. Throughout the life of the components for which the checks establish a rating, the figure of merit can be used for comparison purposes to determine whether the circuits are performing at the level of excellence they had when the instrument was shipped from the factory.
- 2-17. A figure of merit for the 851B can be established by performing a sensitivity check. Procedure is given in Paragraph 5-24. In the table below, space is provided for recording the figures obtained.

#### 2-18. RACK MOUNTING.

2-19. Procedure for rack-mounting the 851B is indicated in Figure 2-2.

<u></u>	put		ettings	Gi - C H I	Cable Hand	Power Input for 6-cm
Freq	Point	I. F. BW	I. F. GAIN (DB)	Sig Gen Used	Cable Used	Vertical Deflection
20 Mc	I. F. INPUT	1 MC 100 KC 10 KC 3 KC 1 KC	70 + 10 I. F. VERNIER, max cw	606A	MAURITRON TEC 8 Cherry Tre Oxon ( Tel:- 01844-35169	lanuals Contact HNICAL SERVICES e Rd, Chinnor DX9 4QY Fax: 01844-352554



more and the same of the selection of th

- 1. Model 908A Coaxial Termination: install in LO OUTPUT - TERMINATE in  $50 \Omega$ .
- 2.  $\triangle$  Connect 2GC OUTPUT to 2GC I. F. INPUT.
- 3. Connect 200MC INPUT to OUTPUT.
- Connect SWEEP OUTPUT to SWEEP INPUT.
- Connect I. F. OUTPUT to I. F. INPUT.
- Connect 851 CONTROL to 8551 CONTROL.

- $^{10}, \ \triangle$  Connect SWITCHED LINE OUTPUT to 11.  $\bigtriangleup$  LINE INPUT.
- 12. SET LINE to OFF.
- 13.8551 line voltage switch: set for nominal voltage of power source (set with blade of screwdriver); check that fuse is value marked adjacent to selected setting.
- 14.851 line voltage switch: set to same setting as set at 8551 line voltage switch; check that fuse is proper value for voltage set.
- 15. LINE INPUT: connect to 115/230V 50/60 cps 300-watt source.





Figure 2-1. Installation Connections, Model 851/8551 Spectrum Analyzer (sheet 2 of 2)



Completing on the Control of the State of the Control of the Contr

the way to the work which will and the state of the state

Figure 2-2. Rack-Mounting Procedure, 851B

# SECTION III OPERATION

#### 3-1. INTRODUCTION.

- 3-2. The Model 851/8551 Spectrum Analyzer is a triple-conversion superheterodyne scanning receiver with a visual amplitude-vs-frequency output. Information obtained by the 8551 RF Section is displayed on the 851 CRT. Analyzer controls are calibrated and thus considerable information can be read directly from the display; calibration accuracies are given in Table 1-1.
- 3-3. Basic step-by-step procedures for putting the Analyzer into operation are given in this Section of the Manual. Information on spectrum analysis and applications of the 851/8551 Spectrum Analyzer are provided in hp Application Note 63.
- 3-4. Operating the Analyzer requires both the Display and RF sections. Instructions for the 851B will make more sense if instructions for the 8551 RF Section are included also, and so the Operating Plates, Figures 3-3 and 3-4, give instructions that include both instruments. However, always remember that instructions given in this Manual are incomplete in that they do not discuss limitations on input signal level. Before turning on the Analyzer, therefore, refer also to the operating instructions (Section III) in the Manual for the 8551.
- 3-5. Front panel controls are identified and briefly described in Figure 3-1, and rear panel connectors and switches are identified in Figure 3-2; initial turnon instructions are given in Figure 3-3, and photographic procedures in Figure 3-4.
- 3-6. Optimum I. F. BANDWIDTH setting for selected SPECTRUM WIDTH and SWEEP TIME settings is given in Figure 3-5. As used here, optimum is defined as the narrowest bandwidth which does not attenuate the signal because of limitations in the rise time of the 20-Mc I. F. Amplifier. (The SPECTRUM WIDTH switch is on the 8551, and determines width of band swept by the 8551Local Oscillator.)

#### Note

With I. F. BANDWIDTH at AUTO SELECT, optimum bandwidth is automatically selected.

#### 3-7. DESCRIPTION.

3-8. The 851B Spectrum Analyzer Display Section includes a 20-Mc I. F. amplifier with five calibrated bandwidths, shaping circuits which provide a choice of amplitude calibration, and a cathode-ray tube and associated circuits, one of which is a calibrated SWEEP TIME switch.

#### 3-9. I. F. BANDWIDTH.

3-10. SELECTABLE. Bandwidth of the 20-Mc I. F. Amplifier is 1 Mc. However, by means of selectable

precision filters, bandwidth of the 20-Mc I. F. Amplifier can be narrowed to 100 Kc, 10 Kc, 3 Kc, or 1 Kc; selection is made with the I. F. BANDWIDTH switch.


3-11. DETERMINES RESOLUTION. Display resolution is determined by the setting of I. F. BANDWIDTH. The signal shown on the CRT can be considered as a presentation of the spectrum as seen through a moving window. How much of spectrum can be seen at any one instant is the ratio of 20 MC I. F. Amplifier bandwidth to spectrum width displayed. For example, if the I. F. BANDWIDTH switch is set at 1 Mc and the SPECTRUM WIDTH switch is set at 100 MC/CM (spectrum being examined, 1 Gc), one-thousandth of the spectrum can be seen at any instant as the horizontal sweep voltage moves the "window" across the CRT. The rate at which the "window" moves across the graticule is set by the SWEEP TIME switch; the shape of the "window" is the passband characteristic of the 20 MC I. F. Amplifier.

#### 3-12. DISPLAY.

- 3-13. SHAPING. In addition to the usual oscilloscope controls and circuits, the 851B Display Section provides circuits for shaping the signal derived from the 20-Mc I. F. signal. By means of the shaping circuits, in addition to the conventional linear display (proportional to signal voltage), the display can be made proportional to signal power (square) or to the logarithm (level in db) of the signal. Signal shaping (together with the wide dynamic range of the Analyzer) makes it possible to show signals of widely-varying amplitudes on the same display; for example, signals at -30 dbm and -90 dbm can be viewed simultaneously. Choice of display ratio is made with the VERT DIS-PLAY switch.
- 3-14. SCALES. There are two calibrations on the 851 CRT graticule. One is numerical (0.2/cm), the other is logarithmic (10 db/cm). With VERT DIS-PLAY at LIN or SQ, use the numerical calibration; at LOG, use the DB scale.

#### 3-15. OSCILLOSCOPE.

- 3-16. TUBE. The CRT used in the 851 is a 5-inch tube with an internal (parallax-free) graticule. Unless otherwise ordered, the tube is supplied with a medium persistence (P-2) phosphor and light blue filter.
- 3-17. TIME BASE. Time-base range is from 3 milliseconds/centimeter to 1 second/centimeter, and is selected with the SWEEP TIME switch. The time base can be synchronized with an internal or externally-supplied signal, or line frequency. For photographic use, a single-sweep mode is provided; selection of mode is made with the SYNC switch.
- 3-18. ALIGN AND BASE LINE CLIPPER CONTROLS. In addition to the usual oscilloscope controls, the




and the second s

8 Cherry Tree Rd, Chinnor Oxon OX9 4QY Tel:- 01844-351694 Fax:- 01844-352554 Email:- enquiries@mauritron.co.uk

- 1. CRT:
  - a. Scale used with VERT DISPLAY at SQ or LIN.
    b. Markers used in signal identification.
  - b. Markers used in signal-identification technique. (See Manual for 8551 RF Section.)
  - c. Scale used with VERT DISPLAY at LOG.
- 2. To blank base line, turn clockwise.
- 3. Adjusts brightness of trace.
- 4. Adjusts focus of trace.
- 5. To obtain one nonrecurring sweep, set SYNC at SINGLE SWEEP, and depress pushbutton.
- 6. Lights when single sweep starts, goes out when single sweep ends.
- 7. SYNC:
  - a. SINGLE SWEEP: sets up internal connections for single-sweep operation.
  - b. LINE, INT, EXT: sets up internal connections for type of sync voltage selected. For EXT operation, input (SYNC INPUT) is on rear panel.
- 8. From six sweep rates, selects time base for horizontal sweep. At EXT, sets up internal conditions required when using sweep voltage supplied from external source. VERNIER provides continuous adjustment between calibrated steps. Note: At


EXT, apply sweep voltage to SWEEP input on rear panel, and compatible blanking pulse to BLANKING input, also on rear panel.

- Selects vertical calibration:
   LIN: amplitude proportional to voltage
   SQ: amplitude proportional to power
   LOG: amplitude proportional to logarithm of
   input signal; level indicated in db.
- 10. I. F. BANDWIDTH switch: 1 KC to 1 MC: manual selection of I. F. bandwidth AUTO SELECT: automatic selection, for CW signals, of optimum I. F. bandwidth for chosen SPECTRUM WIDTH (on 8551) and SWEEP TIME settings.
- 11. Vernier for I. F. GAIN (DB): as vernier is turned cw, up to 1 db of additional gain is inserted.
- 12. Controls I. F. input attenuator, in 10-db and 1-db steps.
  Highest gain setting: 70 + 10 outer control at 70, inner at 10
  Max atten setting: both controls at 0
- 13. Position adjustments: HORIZ POS shifts trace to right or left; VERT POS shifts trace up or down.
- 14. Adjusts trace angularly, enabling operator to align trace with graticule horizontal axes.



- 1.  $\triangle$  J4: power cable connects here. Cable supplied with 851B.
- 2. S1: Line Voltage slide switch: controls powersupply input connections. Always check that switch is set for nominal voltage of external power source.
- 3. Fuseholder: rating of fuse is marked at the Line Voltage switch setting which corresponds to voltage of power source.
- 4. O J1: cable for carrying 20-Mc I. F. signal from RF Section connects here. Cable supplied with RF Section.
- J7: cable carrying sweep voltage from 851B to RF Section connects here. Cable supplied with RF Section.
- 6. J6: input for external sync signal; requirements: positive-going pulse of between 3 volts peak and 15 volts peak.
- 7. J3: input for externally-supplied blanking voltage; requirements: negative 5- to 10-volt pulse, width of which is compatible with retrace time of external sweep voltage used.

- 8. J2: input for externally-supplied sweep voltage; requirements: 0 to approximately +15V saw-tooth from 10,000-ohm source.
- CRT protective cover; may be removed for servicing and/or tube replacement.
- 10. J10: signal to CRT, sampled at output of video detector in 20 MC I. F. Amplifier just ahead of Vertical Amplifier; 0 to -4 volts open circuit, 4700 ohms impedance; BNC female. With high-impedance earphones, output can be used to monitor modulated signals tuned in on Analyzer.
- 11.J5: for sampling 20-Mc I. F. signal just ahead of Video Detector; BNC female.
- 12. J8: sweep voltage, sampled just ahead of Horizontal Amplifier; 10 volts ±3V peak-to-peak open circuit, 4700 ohms impedance; BNC female. Note: With appropriate amplifier, VERT and HORIZ outputs will drive an X-Y recorder to obtain an X-Y plot of spectrum displayed on CRT.
- 13. J9: 14-conductor cable connects here; carries ±15 vdc to RF Section, and SWEEP TIME/SPECTRUM WIDTH connections required for I. F. bandwidth AUTO SELECT operation. Cable supplied with RF Section.




The first state which will be seen a second state of the second state of the second state of the second state of

त्तरः व्यवस्थान्त्रकानः स्थानकान्तर्यः स्थाप्तः । द्वार्षः । वस्य विकासन्तरः स्थाप्तरः स्थाप्तः । वस्य विकासन्

- LINE: set to STANDBY; light should glow, fan turn. After about 10 seconds, base line should appear on CRT. <u>CAUTION</u>: Before setting LINE to ON (step 12), be sure fan is turning; BWO can be damaged if fan is not operating.
- 2. Set: SIGNAL IDENTIFIER . . . . . . OFF FREQUENCY TUNING . . . . COARSE SPECTRUM WIDTH VERNIER . . . CAL
- 3. ATTENUATOR (DB): Set to 60. Be sure switch seats in detent.
- 4. Set range switch, FREQUENCY (GC).
- 5. Set I. F.
  - a. 200 MC if FREQUENCY (GC) is at 1.8 4.2b. 2 GC for all other inputs.
- Connect signal under investigation; range: 10 Mc to 10 Gc. <u>CAUTION</u>: Input must not exceed 1 watt average with <u>ATTEN</u> (DB) set at 60; see Manual for 8551 RF Section.

- 8. INTENSITY: Set at about 3 o'clock.
- 9. If input frequency is unknown, set TUNE to about 3 on LOCAL OSC FREQ scale; if known, set TUNE for Frequency Scale reading near that of input signal.
- 10. Set SPECTRUM WIDTH for 200MC/CM.
- 11. Set I. F. GAIN (DB) at 70 or 80, I. F. VERNIER max ccw.
- 12. Set LINE to ON; indicator light will glow. When TUNE is at 2 Gc, large signal will appear; this is the BWO signal used for self-check, and can be ignored.
- 13. ATTENUATOR (DB): Adjust to bring signal on CRT. <u>CAUTION</u>: Attenuator must be set for sufficient loss to reduce input signal to 1 mw or less.
- 14. Set I. F. GAIN for convenient signal-to-noise ratio, and adjust TUNE to center signal on CRT.
- 15. Adjust SPECTRUM WIDTH for best detail in region of interest. If display shifts position, check COARSE POS BALANCE (see Paragraph 5-125, 851A/8551A Manual). (COARSE POS controls range-to-range position.)



Note: Before taking pictures, read Paragraph 3-22.

- 1. Perform steps 1 through 14 of initial operating procedure, Figure 3-3; refer to information on input signal levels in 8551 Manual.
- If 8551 SPECTRUM WIDTH is set for 1 MC/CM or less, stabilize the Analyzer; see Figure 3-5, 8551 Manual.
- Align trace, and adjust SWEEP TIME, SPEC-TRUM WIDTH, I. F. GAIN for detail of interest.
- 4. Install and align Oscilloscope Camera. Load film pack; Polaroid ASA 3000 is recommended.
- 6. Adjust INTENSITY so fast transients in waveform almost disappear.

- 7. Adjust FOCUS for finest trace.
- 8. Depress ultra-violet light pushbutton on camera; determine exposure required when graticule is illuminated by ultra-violet light and no sweep is present. For Polaroid's ASA 3000, about 1/5 second at F11 is recommended.
- 9. After setting camera, depress ultra-violet light pushbutton, and hold it on until phosphor glows. Photograph the graticule.
- 10. Release ultra-violet light.
- 11. Set camera for photographing the trace. Exposure required is a function of sweep speed and intensity.
- 12. Open camera shutter, depress SINGLE SWEEP pushbutton, and watch SINGLE SWEEP light. When light goes out, close shutter.

***			

	ı ( ]		1.	
			,	
	1 1)	,	,	1μ C
		1.	tr c	1    1
	1.1.1	,	16-1	l je l
	• · · ·	<u> </u>	161	164
	CM		101	186
	4.841	0		3 <b>K</b> G
	1 1 8/8	of .		4 M E
Ì	sorta	M 1/4		រង្សម
	1.7.41.2	. A		iúgs.
	,	l q		112 <b>14181</b>
ı		r i	'	n 10

The state of the s

The second secon

The second secon

#### 1 JU OPERATORS INSTRUCTIONS

r the region of a real pulling to a met a definite process of the design of the second of the second

#### 11011

do timita i proposali, continuata i proposali, continu

#### E > > | MICHOGRAPHIC TECHNIQUES

to the first terms grations to he with the control of the possibilities of a possibilities of the possibilities of

t t

a contra i to the contra contra contra contra contra

And the state of t



For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-35254
Email:- enquiries@mauritron.co.uk

Section IV Paragraphs 4-1 to 4-4

# SECTION IV PRINCIPLES OF OPERATION

#### 4-1. INTRODUCT!ON.

- 4-2. The Model 851/8551 Spectrum Analyzer receives and scans signals in the 10-Mc to 42-Gc range, and displays the amplitude of signal components as a function of frequency. In the same display the Analyzer can present signals which vary up to 2 Gc in frequency and up to 60 db in amplitude. The Analyzer is shown in block diagram form in Figure 4-1.
- RF SECTION. The scanning receiver is the 8551 RF Section; it is a triple-conversion, narrowband superheterodyne receiver. Tuning and width of frequency range displayed are controlled by circuits in the RF Section. Amplitude-vs-frequency information over a selected portion of the frequency spectrum is obtained by sweeping the first local oscillator between 2 and 4 Gc. Sweep voltage to drive the first Local Oscillator (a backward-wave oscillator) is supplied by the 851 Display Section, and is the same voltage used to drive the 851 CRT horizontal plates; this arrangement maintains frequency calibration of the CRT display. The 8551 heterodynes the input signal and its components to an I.F. of 20 Mc. The 20-Mc I.F. is carried by external jumper cable to the 851 Display Section.
- 4-4. DISPLAY SECTION. Figure 4-2 is a block diagram which relates the Display Section to its schematic diagrams.
- a. I. F. Input Attenuator. In the Display Section, the 20-Mc I. F. is first applied to an attenuator which is calibrated as a gain switch [I. F. GAIN (DB)]. The I. F. attenuator includes a two-section switch (see Figure 5-16) which inserts attenuation from 0 to 80 db in 1-db steps. I. F. VERNIER, which provides up to 1-db adjustment between steps, and which is associated with the I. F. GAIN switch, is not in the input circuit; it is in the final 20-Mc I. F. Amplifier, and attenuates the signal just ahead of the Video Detector (see Figure 5-27).
- b. Bandwidth Switching Circuits. Following the I. F. attenuator is the first 20-Mc I. F. Amplifier, bandwidth of which is selected with the I. F. BAND-WIDTH switch. Amplifier bandwidth determines the resolution of the display; the narrower the bandwidth, the more detailed the presentation of the frequency distribution of signal and components. Variation in 20-Mc I. F. Amplifier bandwidth is provided by a set of bandpass filters which are connected into the Amplifier circuit by relays (see Figure 5-24). The relays are controlled by the I. F. BANDWIDTH switch; with no relays energized, 20-Mc I. F. Amplifier bandwidth is 10 Kc. In addition to the 10-Kc bandwidth, 1-Kc, 3-Kc, 100-Kc, and 1-Mc bandwidths can be selected. Or an automatic bandwidth selection system can be used which selects that I. F. bandwidth which will provide optimum display for whatever SPECTRUM WIDTH and SWEEP TIME settings are selected. (As used

here, optimum is the narrowest bandwidth which does not attenuate the signal because of limitations in the rise time of the 20-Mc I. F. Amplifier.) With I. F. BANDWIDTH at AUTO SELECT, current to operate the bandwidth-switching relays is brought through contacts on the SWEEP TIME switch and on the 8551 SPECTRUM WIDTH switch as well as through contacts on I. F. BANDWIDTH; connections for this mode of operation are carried in the inter-unit CONTROL cable.

#### Note


The narrower the I.F. bandwidth, the more the incoming signal is attenuated if the tuning is too fast.


c. Display Calibration. The Analyzer Display Section provides a choice of display calibration, thus making it possible to present in the same display responses with a difference of up to 60 db. The switch from linear calibration to logarithmic (or square) calibration is made with the VERT DISPLAY switch. The signal is converted to the proper ratio by the Current-Controlled Attenuator (see Figure 5-27); this circuit includes diodes which shunt the signal path. The diodes act as variable resistors whose resistance depends on the bias current supplied. For a linear display, bias current is fixed. For log or square displays, bias current is changed by shaping the output voltage of the video amplifier and applying it to the current-controlled attenuator. The shaping circuits are discussed in more detail in Paragraph 4-15.

#### d. Video Signal.

- (1) From the Current-Controlled Attenuator, the 20-Mc I. F. signal passes to the final 20-Mc I. F. Amplifier, the Video Detector (A2A7T1 and A2A7CR1-A2A7CR4), and the Video Amplifier. These circuits are shown in Figure 5-27.
- (2) Two outputs from the 20-Mc-to-Video signal path are provided: 1) I. F. TEST POINT, which samples the 20-Mc I. F. just ahead of the Video Detector, and 2) VERT OUTPUT, which samples the Video signal at the output of the Video Amplifier, just ahead of the Vertical Amplifier.
- (3) In the Vertical Amplifier (Figure 5-29), the video signal is applied to the base of A7Q8 in differential amplifier A7Q8-A7Q9; the signal on the base of A7Q9 is determined by the setting of VERT POS adjust R8. Vert Gain adjust A7R15 in the emitter circuit of differential amplifier A7Q8-A7Q9 is adjusted when calibrating the Vertical Amplifier (see Paragraph 5-80). A variable RC network in the emitter circuit of the final stage of the Vertical Amplifier, A7Q6-A7Q7, provides further filtering of the video signal before it is applied to the CRT vertical deflection plates. Selection of capacitor for the RC filter is controlled by I.F.BANDWIDTH.

The state of the s





ecters the transfer to a man defendance when the state of the state of

Figure 4-2. Model 851B Spectrum Analyzer Display Section, Block Diagram

marie term traffice and arrangements to profess more in any strategy the research

- e. Sweep and Horizontal Amplifier. This circuit (see Figure 5-33) generates the sweep voltage and blanking signal.
  - (1) Rate of sweep generated is determined by RC networks connected into the sweep-generating circuit by SWEEP TIME switch A10S1. Connections for the type of SYNC voltage selected are set up by SYNC switch S2. Operation of the Sweep Generator is discussed in Paragraph 4-5.
  - (2) Sweep voltage for driving the first Local Oscillator (BWO) in the 8551 RF Section is taken at the output of the Sweep Generator, just ahead of the Horizontal Amplifier. This voltage appears at SWEEP OUTPUT on the rear panel, and is carried by inter-unit SWEEP cable to the 8551. Sweep voltage, sampled at the same point in the circuit, appears also at the HORIZ output connector on the rear panel.
  - (3) Horizontal Amplifier A6Q16-A6Q18 is a differential amplifier. Sweep voltage is applied to the base of A6Q16, and the signal on the base of A6Q18 is determined by the setting of HORIZ POS adjust R9. Horiz Gain adjust A6R54 in the collector circuit of the Amplifier is adjusted

- when calibrating the Horizontal Amplifier (see Paragraph 5-50). Amplified sweep voltage is applied to the CRT horizontal deflection plates.
- (4) The blanking signal, taken from the emitter of A6Q6, is amplified by A7Q3 on the Vertical Amplifier Board (see Figure 5-29) before it is applied to the CRT.
- (5) Via contacts at the EXT position of SWEEP TIME, sweep voltage from a suitable external source, such as one of the hp 690 Sweep Oscillators, can be applied to the Horizontal Amplifier to drive CRT horizontal plates. The CRT requires a sawtooth voltage of from 0 to +15 volts. Inputs for sweep voltage and compatible blanking signal are on the 851B rear panel.

## 4-5. OPERATION OF HORIZONTAL SWEEP GENERATOR.

4-6. The Horizontal Sweep Generator is shown in block-diagram form in Figure 4-3; the schematic is Figure 5-33.

#### 4-7. EXTERNAL.

Marie Britant and the automated the Landing Comment of the land

4-8. With SYNC at EXT, a positive input signal will cause Schmitt Trigger A6Q1, A6Q2 to generate a




Figure 4-3. Horizontal Sweep Generator, Block Diagram

negative trigger pulse at the collector of A6Q1. The negative pulse is coupled through A6CR2 to the base of A6Q7, one half of a bistable multivibrator. The negative pulse turns A6Q7 off and A6Q8 on.

- 4-9. When A6Q7 turns off, it turns off A6Q6. This cuts off the blanking signal, allowing the horizontal sweep to be seen on the CRT if the base line clipper has not biased off the CRT.
- 4-10. When A6Q8 turns on, Pulse Amplifier A6Q9 conducts, turning off -5V Clamp A6Q10. When A6Q10 turns off, its emitter goes positive, biasing off Reset Switch A6Q15. When the Reset Switch cuts off, the Sweep Capacitor (A10C7-A10C12) in the base circuit of A6Q11 starts to charge, applying a negative-going signal to the base of A6Q11.
- 4-11. Transistors A6Q11, A6Q12, A6Q13, and A6Q14 form a Miller Integrator. The output of A6Q14 is a positive-going ramp. The positive ramp voltage is fed back to the Sweep Capacitor (A10C7-A10C12). As the Sweep Capacitor charges negatively on its bottom plate, the top of the capacitor is going positive. The result is that the voltage drop across R35 and A6R29-A6R34 is almost constant as the Sweep Capacitor charges. If the voltage drop is constant, the current through the resistors is constant. This same current is the charging current for the Sweep Capacitor. If the charging current is constant then the capacitor is charging linearly and a linear ramp voltage out of the Miller Integrator is the result.
- 4-12. The positive ramp fed back to the Sweep Capacitor also goes to the base of A6Q4, one half of the Sweep Reset Multivibrator. The signal is amplified and appears as a positive-going voltage at the collector of A6Q3. The signal is coupled through A6CR1 and delivered to the base of A6Q7. When the sweep voltage from the Miller Integrator circuit reaches a predetermined level, A6Q7 starts to conduct. The conduction of A6Q7 cuts off A6Q8. This causes Pulse Amplifier A6Q9 to cut off, turning on the -5V Clamp, A6Q10. When A6Q10 conducts it turns on Reset Switch A6Q15 which discharges the Sweep Capacitor, ending the sweep. At the time A6Q7 turns on it turns on A6Q6, blanking the CRT during retrace.

#### 4-13. INTERNAL.

- 4-14. With SYNC at INT, operation of the sweep circuit is essentially the same except that no external trigger is needed to turn off A6Q7.
- a. With SYNC at INT, the Reset Capacitor (A10C1-A10C6) is connected through A6R13 to the -15V supply. As the Reset Capacitor charges negatively, the voltage is coupled through A6R12 and A6CR2 to the base of A6Q7. This triggers A6Q7 and starts the sweep. Sweep termination is the same as when operating from an external trigger.
- b. The one other difference in operation is that the conduction of A6Q8 also turns on A6Q5, discharging the Reset Capacitor until the end of sweep. At the end of sweep when the Reset Multivibrator flips back, A6Q5 is cut off, allowing Reset Capacitor (A10C1-A10C6) to charge negatively again and restart sweep.

#### 4-15. OPERATION OF VERTICAL DISPLAY.

#### 4-16. CURRENT-CONTROLLED ATTENUATOR.

- 4-17. Between the first and second 20-Mc I. F. Amplifiers, the 20-Mc I. F. is passed through the Current-controlled Attenuator (see Figure 5-27). The attenuating element is a network of hot carrier diodes which shunt the signal path.
- 4-18. Hot carrier diodes are used because they have very low shunt capacity and a very predictable dynamic resistance-vs-current characteristic. This predictable characteristic makes it possible to design shaping circuits which will give the desired attenuation characteristics in the LOG and SQuare modes of operation.
- 4-19. Figure 4-4 shows a dynamic resistance-vs-current curve for a hot carrier diode. As current through the diodes increases, dynamic resistance decreases. In the Current-Controlled Attenuator, lower diode resistance causes more signal shunting, i.e., more attenuation of the signal.




Figure 4-4. Typical Dynamic Resistance-vs-Current Curve, Hot Carrier Diodes

#### 4-20. SQUARE MODE OF OPERATION.

- 4-21. The purpose of the SQuare mode of operation is to convert the voltage (linear) indication on the CRT to a voltage-squared display representing power. If two signals are present with a voltage ratio of 2:1 (VERT DISPLAY at LIN), at SQ they will appear on the CRT as signals with a 4:1 ratio.
- 4-22. To achieve this change in the display, the amount of current to the Current-Controlled Attenuator must decrease with an increase in signal level.
- 4-23. The video signal into the VERT DISPLAY Switch Assembly is negative-going. A negative signal on the base of A11Q2 (see Figure 4-5) will increase its conduction. This will decrease current through A11CR1, AllCR2, and AllQl; this current flows through the hot carrier diodes in the Current-Controlled Attenuator. Larger signals will cause a much greater decrease in current through the Attenuator diodes than small signals. Shaping circuit characteristics are such that any increase in signal level will cause the square of the increase to appear on the CRT. For example, as the signal goes from 1 to 2 in voltage, the decrease in shaping circuit current is such that four times as much signal gets through the Attenuator. In general, any increase in signal level will cause the square of the increase to appear on the CRT.




Figure 4-5. VERT DISPLAY Switch at SQ, Simplified Schematic

#### 4-24. LOGARITHMIC MODE OF OPERATION.

4-25. The purpose of the LOG mode of operation is to convert the incoming voltage to bias current of such value that the resulting display is proportional to the log of the input voltage. To achieve such a display, larger signals must cause much greater attenuation than small signals; that is, as signal level increases, a much greater amount of current must flow through the diodes in the Attenuator. With VERT DISPLAY at LOG, current out of the shaping circuit is such that gain through the Current-Controlled Attenuator is logarithmic; that is, for each 10 db of change

in signal level there is a 1-centimeter change in signal display.

4-26. Refer to the simplified schematic of VERT DISPLAY at LOG, Figure 4-6. The video signal fed back to A11Q1 is negative-going. As A11Q1 conducts more, diodes A11CR3 and A11CR4 are biased on. When they conduct, they decrease the emitter resistance of A11Q1, increasing the gain. This causes proportionately more current to flow through the Attenuator hot carrier diodes on large signals than on small. The shaping circuit in the emitter of A11Q1 is designed to provide a logarithmic gain through the Current-Controlled Attenuator.

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:-01844-351694 Fax:-01844-35254
Email:-enquiries@mauritron.co.uk




Figure 4-6. VERT DISPLAY Switch at LOG, Simplified Schematic

Table 5-1. Test Equipment Required

Ref No.	Instrument Type	Critical Specifications	Recommended Equipment
1	Oscilloscope Low-Frequency	Sensitivity: 0.1 mv/cm	hp 130C 200 μv/cm Oscilloscope
2	DC Voltmeter	Accuracy: 0.05% Input Impedance: 10.2 megohms Automatic Range Selection Range: to 150V	hp 3440A Digital Voltmeter with hp 3442A Automatic Range Selector Plug-In
3	Transformer for varying input voltage	Range: 103 to 253 VAC at approx 1/2 amp Voltmeter Range: 103 to 127 volts Voltmeter Accuracy: ±1 volt	General Radio Type W10MT3A Superior Electric UC1M
4	Clip-On DC Milliammeter	Accuracy: ±0.1 ma ±3% of FS	hp 428B
5	DC Voltmeter	Accuracy: ±2% of FS Input Resistance: 100 megohms Can accommodate voltage-divider probe	hp 410C Electronic Voltmeter
6	DC Voltage Divider	Accuracy: ±5% Division Ratio: 100:1 Input Resistance: 10,000 megohms Max Volts: 6000	hp 11045A DC Voltage Divider
7	Electronic Counter	Frequency: 200 Mc Accuracy: 5 parts in 10 <sup>8</sup> ±1 count Multiple period averaging feature	hp 5245L Electronic Counter and hp 5253B Frequency Converte
8	Low-Frequency Oscillator	Frequency Range: 1 cps to 350 cps, continuously variable Output: 5 volts peak Distortion: less than 0.5% above 5 cps	hp 202C
9	HF Signal Generator	Output Frequency: 50 Kc to 20 Mc Frequency Accuracy: ±1% Output: at least 3 volts into 50 ohms Modulating capability with external modu- lating voltage input Meter which monitors generator output level	hp 606A
10	VHF Attenuator	To 60 db, in 10-db steps, at 2 Gc	hp 355D
11	UHF Signal Generator	Frequency: 2 Gc	hp 8614A
12	Precision 10-db/ Step Attenuator	Accuracy at 20 Mc: 0-10 db, ±0.02 db 10-20 db, ±0.03 db 20-70 db, ±0.03 db + 0.03 db/20 db	hp H25-355D VHF Attenuator
13	Precision 1-db/ Step Attenuator	Accuracy at 20 Mc: ±0.02 db	hp H25-355C VHF Attenuator

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk

### SECTION V MAINTENANCE

#### 5-1. INTRODUCTION.

- 5-2. Information required to maintain the 851 Spectrum Analyzer Display Section in working condition is provided in this section. Type of information covered is summarized briefly in Paragraph 5-5.
- 5-3. "Right" and "left", "backward" and "forward" as used in this section are with respect to the instrument as seen by the operator when he is facing the front panel and the instrument is in its normal (upright) position.
- 5-4. Unless otherwise specified, test procedures assume the 851 is connected to a 115 or 230-volt, 50- to 1000-cycle, approximately 25-watt source.

#### 5-5. CONTENT.

- a. Performance Checks. Tables 5-6 and 5-7; Paragraphs 5-9 through 5-34.
- b. Checks and Adjustments. Procedures are given in brief form in Table 5-13 and, in more detail, in Paragraphs 5-37 through 5-130.
- c. Required Test Equipment. Instruments required for tests and adjustments are listed and briefly specified in Table 5-1; each instrument is given a number reference. Accessories required for the procedures are listed in Table 5-2; each accessory is given a letter reference. Ahead of each section of the procedure, equipment required is listed by its reference number or letter.
- d. Schematics. Inserted at the end of Section V, ahead of the Replaceable Parts Lists (Section VI).
  - e. Waveforms. Table 5-25.
- f. Assembly and Component Identification. screening on the instrument identifies Assemblies and parts plainly and thoroughly. In addition, the following aids are provided.
  - (1) Location of Assemblies is called out in Figures 5-10 and 5-11.
  - (2) Table 5-3 lists Assemblies numerically and gives the schematic or schematics on which each is shown.
  - (3) Table 5-4 lists controls, switches, and connectors alphabetically and gives the schematic on which each is shown.
  - (4) Table 5-5 lists each chassis-mounted component numerically by reference designation, and references the schematic on which each is shown, and either tells where the part is located or references an illustration which calls out the
  - (5) Each Board-mounted component is called out on a picture of the Board. In the main, these illustrations face the schematic in which the

Assembly appears. For help in locating the illustration for a given Board, refer to the List of Illustrations in the front of the Manual. Paragraph 5-134 includes suggestions on how to proceed when it is necessary to locate a part.

g. Troubleshooting Information, Disassembly Instructions. Paragraphs 5-131 through 5-157.

### 5-6. COVER AND SIDE PANEL REMOVAL.

- a. Equipment Required: Phillips driver No. 2.
- b. Top Cover Removal.
- (1) Remove four phillips head screws (6-32 x7/16").
- (2) Slide cover to rear, and off instrument.
- c. Side Panel Removal. After removing the top cover, from each panel remove the four phillips head screws (6-32 x 3/16"). The side panels lift off.
  - d. Bottom Plate Removal.
    - (1) Remove the four phillips head screws (6-32  $\times$  7/16").
    - (2) Push plate to rear, and off instrument.

#### 5-7. PERFORMANCE CHECKS.

Operational checks for incoming or routine inspection are given in Table 5-6, and procedures for verifying that the 851 meets specifications are given in Table 5-7. Both sets of procedures are spelled out in greater detail in Paragraphs 5-9 through 5-34. Both Tables reference the more detailed procedures as an aid in case brevity has obscured clarity. Table 5-7 is in test-card form, briefly describes test sequences, and provides space for recording measurement results.

#### 5-9. OPERATIONAL CHECKS.

- 5-10. INTENSITY CONTROL.
  - a. Set 8551 LINE to ON.
  - b. Set INTENSITY maximum cw.
- c. Turn SWEEP TIME through its range, and watch display for retrace.

There should be no retrace at any setting of of SWEEP TIME.

d. Set INTENSITY maximum ccw. No trace should be visible.

#### 5-11. BASE LINE CLIPPER.

- a. Perform steps 1 through 14 of initial turn-on procedure (see Figure 3-3), using any signal from 10 Mc to 10 Gc.
  - b. Set I. F. GAIN (DB) at 70.

Table 5-2. Test Accessories Required

Ref No.	Instrument Type	Critical Specifications	Recommended Equipment
A	Cable Assembly (2 each)	Shielded 50-ohm cable terminated with dual banana plugs	hp 11000A
В	Cable Assembly	Shielded 50-ohm cable, dual banana plug to alligator clips	hp 11037A
C	Cable Assembly (3 each)	RG-58C/U, BNC male to dual banana plug	hp 11001A
D	Cable Assembly	RG-58C/U, BNC male to BNC male	hp 10503A
Е	Adapter	BNC female to dual banana plug adapter	hp 10111A
F	BNC Tee	BNC male to 2 BNC females	UG-274A/U, hp 1250-0072
G	Plastic tuning wand	Approx 7" long x 3/8" diam plastic	Modified* General Cement #GC8721
Н	Cable Assembly	Shielded coax, type N male to type N male, 3 feet long	Special hp 11500A
J	Adapter	BNC male to male	UG-491A/U
K	Screwholding Screwdriver		Quick Wedge 1734-XM or 736-50

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk

Table 5-3. Assemblies vs Schematics, 851B

Assy No.	Designation	Schematic	Figure No.	
A1	I. F. GAIN Switch Assy	I. F. Input and Attenuator	5-20	
A2	RF Circuit Assembly	I. F. BANDWIDTH Switching Circuits VERT DISPLAY Switch, Current- Controlled Attenuator, and 200MC I. F. Amplifier Assemblies	5 -24 5 -27	
A3	100-Kc Bandpass Filter Assy	I. F. BANDWIDTH Switching Circuits	5-24	
A4	I. F. BANDWIDTH Switch Assy	I. F. BANDWIDTH Switching Circuits I. F. BANDWIDTH Switch A4S1	5-24 5-38	
A5	100-Kc Bandpass Filter Assy I. F. BANDWIDTH Switching Circuits		5-24	
A6	Sweep and Horizontal Amplifier Assy		5-33	
A7	Vertical Amplifier Assy Vertical Amplifier		5-29	
A8	High-Voltage Board High-Voltage Power Supply		5-35	
A9	Low-Voltage Board	Low-Voltage Power Supply	5-37	
A10	SWEEP TIME Switch Assy	Sweep and Horizontal Amplifier SWEEP TIME Switch A10S1	5-33 5-30	
A11	VERT DISPLAY Switch Assy	VERT DISPLAY Switch, Current- Controlled Attenuator, and 20MC I. F. Amplifier Assemblies	5 - 27	
A12	Bandpass Filter Assy	I. F. Input and Attenuator	5-20	

le atrols, Switches, and (

Name	Ret Desig	
BASE LINE CLIFT	R7	Vertica
BLANKING Input	J3	Sweep &
CRT	<b>V</b> 1	HV Pow
FOCUS	R4	HV 1k (w)
HORIZ Output	18	Sweep &
HORIZ POS	R9	Sweep 8
I. F. BANDWID116	A4S1	Switch I
L. F. GAIN (DB - 86	A1S1 A1S2	L.F. Input K
I. F. INPUT	J1	I. F. Input /
I. F. TEST POINT	J5	VERT DISE Attenuato
I.F. VERNIFR	R10	VERT DISE Attenuato
INTENSITY	R1	HV Power
LINE INPUT	J4	LV Power
Line Switch (11)	S1	LV Hower
SINGLE SWEEP SATE	S3	Sweep & H
SWEEP OUT PET	J7	Sweep & H
SWEEP TIME Switch	A10S1	Sweep & H
SWEEP TIME VERNIER	A10R1	Sweep & He
SYNC Input	J6	Sweep & H
SYNC Switch	S2	Sweep & H
TRACE ALIGN	R5	HV Power
VERT DISPLAY Switch	A11S1	VERT DISE Attenuato
VERT Output	19	Vertical Ar
VERT POS	R8	Vertical Amplifier

- c. Adjust level at Signal Generator for a 7-cm display. Center display with TUNE.
  - d. Turn BASE LINE CLIPPER maximum cw. At least the bottom 2 cm of the display should be blanked.

#### 5-12. FOCUS CONTROL.

5-13. The FOCUS control is within specifications if focus is obtained somewhere within -90° and +90° of FOCUS travel where  $0^{\circ}$  is defined as the position the control has with its white arrow vertical.

#### 5-14. VERTICAL DISPLAY ACCURACY CHECK.

5-15. EQUIPMENT REQUIRED

Ref No.	Egri
11*	UHF Signa
	8551 RF
H**	<b>Cable</b> terms of <b>(1150</b> 0A)
*T:	able 5-1

5-16. PRELIMINAR (Fig. 1994) turned on (see Figura any other frequency set

SPECTRUM WITH SWEEP TIME
I. F. BANDWIDTH

Table 5-5. Chassis Parts Locater

Circuit				Location (Fig. 1	40 )
Desig	Name		Photogra	ph   Tig. 1	Schematic
A1	I. F. GAIN (DB) Switch Assy		5-11		5-20
A2	RF Circuit Assy		5-11		5-24,5-27
A3	100KC Bandpass Filter			under A9 board	5-24
A3L1	100KC Bandwidth Adj		5-10		5-24
A4S1	I. F. BANDWIDTH Switch		5-10		5 <b>-3</b> 8
A5	100KC Bandpass Filter			under A9 board	5-24
A5L1	100KC Bandwidth Adj		5-10		5-24
A6	Sweep & Horizontal Amplifi	er Assy	5-10		5 -33
A7	Vertical Amplifier Assy		5-11		5 <b>-2</b> 9
A8	HV Power Supply Assy		5-10		5 -35
<b>A</b> 9	LV Power Supply Assy		5-10		5-37
A10R1	VERNIER (Sweep Time Swit	ch)	5-11		5 -33
A10S1	SWEEP TIME Switch Assy		5-11		5-33
A11	VERT DISPLAY Switch Ass	у	5-10		5-27
A12	Bandpass Filter Assy		5-10		5-20
C1	p/o Line-input filter			on inside of rear panel	5-37
C3	+100V supply filter		5-10		5-37
C4	-24V supply filter		5-10		5-37
C5	-24V supply filter		5-10		5-37
C6	Hold-Sweep-On Capacitor; or	n S3	5-10		5-33
C7	+15V supply filter			on inside of rear panel	5-37
C8	-15V supply filter	For Service Manua	als Contact	on inside of rear panel	5-37
L1	p/o Line-input filter	MAURITRON TECHNIC 8 Cherry Tree Rd	Chinnor	on inside of rear panel	5-37
L2	p/o Line-input filter	Oxon OX9 Tel:- 01844-351694 Fax:	IOV	on inside of rear panel	5-37
L3	Couples TRACE ALIGN to C	R7Email:- enquiries@ma	uritron.co.uk	on inside of CRT shield see Figure 5-17	; 5-35
L4	p/o 24V filter			on side of A2 casting; see Figure 5-26C	5 -37
Q1	Drives Step-up Xfmr A8T1	ì	5-10, 5-	11	5-35
Q2	Drives Step-up Xfmr A8T1		5-10, 5-1	1	5-35
Q3	LV Series Regulator		5-10		5-37
Q4	LV Series Regulator		5-10		5-37
Q5	LV Series Regulator		5-10		5-37
Q6	LV Series Regulator		5-10		5-37
R1	INTENSITY control		5-11		5-35
R2	Int Level Adj	ĺ	5-10		5-35
R3	Astig Adj	ł	5-10, 5-1	1	5 -35
R4	FOCUS control		5-10		5-35
R5	TRACE ALIGN		5-10		5 -35
R6	Isolation for Vert Ampl A7	ļ		on inside of rear panel	5-29
R7	BASE LINE CLIPPER		5-11		5 <b>-29</b>
R8	VERT POS		5-11		5 <b>-29</b>

Table 5-5. Chassis Parts Locater (cont'd)

Circuit		Location (Fig. N	
Desig	Name	Photograph	Schematic
R9	HORIZ POS	5-11	5 -33
R10	I. F. VERNIER	5-11	5-37
R11	p/o +15V filter	on inside of rear pane	1 5-37
R12	p/o -15V filter	on inside of rear pane	1 5-37
S1	115/230 slide switch	under left cover plate	5-37
S2	SYNC switch	5-10	5-33
S3	SINGLE SWEEP switch	5-10	5 -33
T1	Line Transformer	5-10	5 - 37
V1	CRT	5-10	5 <b>-3</b> 5
W1	20-Mc I. F. input	5-10	5-20
w2	From Horizontal Amplifier to CRT (D1, D2)	5-10	5-33
W3	From Horizontal Amplifier to SWEEP OUTPUT	5-11	5-33
W4	CRT Post Accelerator lead; from HV Supply A8 to CRT V1	5-10	5 - 37

Table 5-6. Operational Checks

	Note: Operatio	nal Checks are made with 851 connecte	ed to the 8551.
Par. Ref	Control Under Check	Procedure	Proper Performance
5-10	INTENSITY	LINE ON INTENSITY max cw Turn SWEEP TIME through range watching for retrace.	No retrace at any SWEEP TIME setting
		INTENSITY max ccw	No trace visible
5-11	BASE LINE CLIPPER	Perform initial turn-on (Fig. 3-3), 1 Gc input I.F.GAIN 70 Adjust signal level for 7.0 cm display Set BASE LINE CLIPPER . max cw	At least bottom 2 cm of display should blank
5-12	FOCUS	Set FOCUS with white arrow vertical; this is 0°. Set FOCUS to -90°, then to +90°	Focus should be obtained between -90° and +90°.

Table 5-7. Performance Check Test Card, 851B

Ref	Procedure		Min	Act.	Мах
5-14	1. VERTICAL DISPLAY ACCURACY:			<del></del>	
	a. Equipment Required: Stable Sig Gen (8614A) 8551 RF Section				
	b. SPECTRUM WIDTH 1 Mc/cm SWEEP TIME 3 ms/cm I.F.BANDWIDTH AUTO SELECT				
	c. Perform initial turn-on (Fig. 3-3), 1 Gc input	For Service	e Manuals Conta	act	
5-17	Linear: ±3% of full scale	8 Cherry	TECHNICAL SERVI Tree Rd, Chinno	CES or	
	a. VERT DISPLAY LIN Inner I. F. GAIN 10 Outer I. F. GAIN for low-noise base- line trace (about 50)	Tel:- 01844-35	n OX9 4QY 694 Fax:- 01844-35 iries@mauritron.co.	32554 uk	
	b. Align trace base exactly w/graticule base line	e.			
	<ul> <li>c. Adjust ATTENUATOR (DB) and output of Sig C for 7.0 cm display.</li> </ul>	Gen			
	d. Set Inner I.F.GAIN to .4	cm	3.3		3.7
5-18	Square: ±5% of full scale				0
	a. VERT DISPLAY SQ Inner I. F. GAIN 10 Outer I. F. GAIN for low-noise base- line trace				
	b. Align trace base exactly w/graticule base line	٠.			
	<ul> <li>c. Adjust ATTENUATOR (DB) and output of Sig G for 7.0 cm display.</li> </ul>	en			
	d. Set Inner I. F. GAIN to 7 -3 db	cm	3.15		3.85
5-19	Logarithmic: ±2 db				
	a. VERT DISPLAY LOG  Inner I. F. GAIN 0  Outer I. F. GAIN 70				
	<ul><li>b. Adjust Sig Gen and ATTENUATOR (DB) for 7.0-cm display.</li></ul>				
	c. I. F. GAIN 60	cm	5.8		6.2
	d. Adjust input sig level for 6.0 cm I.F.GAIN 50	cm	4.8		5.2
	e. Input sig level for 5.0 cm I.F.GAIN 40	cm	3.8		4.2
į	f. Input sig level for 4.0 cm I.F.GAIN 30	cm	2.8		3.2
	g. Input sig level for 3.0 cm I.F.GAIN 20	cm	1.8		2.2
	h. Input sig level for 2.0 cm I.F.GAIN 10	cm	0.8		1.2
	<ul> <li>i. I. F. GAIN outer control. 10</li> <li>I. F. GAIN inner control. 10</li> <li>Sig Gen for 2.0-cm display</li> <li>j. I. F. GAIN outer control. 0</li> </ul>	cm	0.8		1.2

Table 5-7. Performance Check Test Card, 851B (cont'd)

Ref	Table 5-7. Performance Check Test Card, 851	7		
		Min	Act.	Max
	2. I.F.BANDWIDTH ACCURACY:  Individual bandwidths Bandwidth repeatabilit better than ±3%.  Equipment Required:  VHF Attenuator (355D) Cable term. w/BNC males (10503A) 8551 RF Section  I.F. CABLE SUPPLIED W/8551	y and stabi	ted within ±20 lity typically	0%.
			851	B - A - 7
5-22	1MC, 100KC, and 10KC Bandwidths			
	a. Find 2-Gc BWO signal: see Paragraphs 5-92 thru 5-95.			
	b. SWEEP TIME 3 ms/cm c. I.F.BANDWIDTH 1 MC			
	SPECTRUM WIDTH 1 Mc/cm			
	d. Adjust VHF Atten and I. F. GAIN for 7.0-cm display.			
	e. Read display at 5.0 cm. cm	0.8		1.2
	f. I.F.BANDWIDTH 100KC SPECTRUM WIDTH 100 Kc/cm		<del> </del>	
	g. Adjust for 7.0-cm display, read at 5.0 cm. cm	0.8		1.2
	h. I.F.BANDWIDTH 10KC SPECTRUM WIDTH 10 Kc/cm			
	1. Adjust for 7.0-cm display, read at 5.0 cm. cm	0.8		1.2
	<ol> <li>Return to each setting of I. F. BANDWIDTH, and note bandwidth at 5.0 cm.</li> </ol>			
	k. Each should be within $\pm 3\%$ of recorded bandwidth.			
	Maximums cm	0.77	<del> </del>	1.23
-23	3KC and 1KC Bandwidths			
	a. Calibrate SPECTRUM WIDTH: see Paragraph 5-98.			
	<ul> <li>b. Set SPECTRUM WIDTH to 10 Kc/cm</li> <li>(by calibration, scale actually is 1 Kc/cm)</li> </ul>			
	<ul> <li>c. 3KC</li> <li>(1) I. F. BANDWIDTH . 3KC</li> <li>SWEEP TIME 3 ms/cm</li> <li>(2) Adjust I. F. GAIN and Ext Atten for 7.0-cm display, read at 5.0 cm.</li> </ul>			
	(3) Switch I. F. BANDWIDTH to any other setting, then back to 3KC; bandwidth should be within ±3% of that noted in step (2).	2.4		3.6
	Maximums cm d. 1KC (1) I. F. BANDWIDTH 1KC SWEEP TIME 10 ms/cm	2.31		3.69
	(2) Adjust Atten for 7.0-cm display, read at 5.0 cm. cm	0.8		1.2
	(3) Switch I. F. BANDWIDTH to any other setting, then back to 1KC; bandwidth should be within ±3% of that noted in step (2).			
	Maximums cm	0.77		1.23

Table 5-7. Performance Check Test Card, 851B (cont'd)

was to the war and he was a second of the state of the st

S. I.F.INPUT SENSITIVITY:	Ref	Procedu	re		Min	Act.	Max
Bandwidth   Limits (dbm)		3. I.F.INPUT SENSITIVITY:	Input required*	for 6-cm ve	ertical disp	olay	
100 Kc		Bandwidth			•	•	
### A Comparing Required:    Sig Gen W/calibrated power output (606A)   Cable term. w/BNC males (10503A)		100 Kc 10 Kc 3 Kc	-75 to -60 -95 to -80 -95 to -80				
Sig Gen w/BNC males (10503A)							
c. Set Sig Gen for 20 Mc; connect to 851 I. F. INPUT.  d. I. F. GAIN 70 + 10		Sig Gen w/calibrated I					
d. I. F. GAIN 70 + 10 I. F. VERNIER fully cow I. F. BANDWIDTH 1 IMC e. Set Sig Gen output for 6.0 cm 951 display. f. Read output at Sig Gen, taking into consideration loss through input cable  g. Perform steps e and f for other I. F. BANDWIDTH settings.  100KC dbm  10KC  3KC  -95  -80  3KC  -95  -80  1KC  4. I.F.GAIN SET ACCURACY: 70-db section: ±0.5 db 10-db section: ±0.1 db  a. Equipment Required: Precision 10-db/step Attenuator (hp H25-355D) Precision 10-db/step Attenuator (hp H25-355C) Signal Generator (506A) Adapter, BNC m-to-m (UG-491A/U) 2 - coax cables term. w/BNC males (10503A)  SIGNAL GENERATOR  H25-355C  MAUNITRON TECHNICAL SERVICES 8 Chernical Maunitron Columbia SERVICES 8 Chernical Maunitron Columbia SERVICES 8 Chernical Maunitron Columbia SERVICES 8 Chernical Services 9 Chernical Services 8 Chernical Services 8 Chernical Services 8 Chernical Services 8 Chernical Services 9 Chernical Services 8 Chernical Services 8 Chernical Services 8 Chernical Services 9 Chernical Services 8 Chernical Services 9 Chernical Services 8 Chernical Services 8 Chernical Services 8 Chernical Services 9 Chernical Services 8 Chernical Services 8 Chernical Services 8 Chernical Services 9 Chernical Services 8 Chernical Services 8 Chernical Services 9 Chernical Services 8 Chernical Services 9 Chernical	5-26	b. Connect 851 to 115/230	7 60/1000 cps.				
I. F. VERNIER		c. Set Sig Gen for 20 Mc;	connect to 851 I.F. INF	TU.			
f. Read output at Sig Gen, taking into consideration loss through input cable  g. Perform steps e and f for other I. F. BANDWIDTH settings.  100KC dbm -75 -60  10KC -95 -80  3KC -95 -80  1KC -95 -80  1KC -86 -71  4. I.F.GAIN SET ACCURACY: 70-db section: ±0.5 db 10-db section: ±0.1 db  a. Equipment Required: Precision 10-db/step Attenuator (hp H25-355D) Precision 1-db/step Attenuator (hp H25-355D) Signal Generator (606A) Adapter, BNC m-to-m (UG-491A/U) 2 - coax cables term. w/BNC males (10503A)  SIGNAL GENERATOR H25-355D		I. F. VERNIER	fully ccw				
10   10   10   10   10   10   10   10		e. Set Sig Gen output for 6					
1. F. BANDWIDTH settings.  100KC dbm -75 -60  10KC -95 -80  3KC -95 -80  1KC -86 -71  4. I.F.GAIN SET ACCURACY: 70-db section: ±0.5 db 10-db section: ±0.1 db  a. Equipment Required: Precision 10-db/step Attenuator (hp H25-355C) Signal Generator (606A) Adapter, BNC m-to-m (UG-491A/U) 2 - coax cables term. w/BNC males (10503A)  SIGNAL GENERATOR H25-355C 10503A 851  For Service Manuals Contact MAURITROLIFICAL SERVICES 8 Cheny Tree NG. Chinnor Oxon OX9 4QY Tel-0184-31634 Fax: 01844-32554 Email: enquiries@mauritron.co.uk					-62		-53
10KC  3KC  1KC  -95  -80  1KC  4. I.F.GAIN SET ACCURACY: 70-db section: ±0.5 db 10-db section: ±0.1 db  a. Equipment Required: Precision 10-db/step Attenuator (hp H25-355D) Precision 1-db/step Attenuator (hp H25-355C) Signal Generator (606A) Adapter, BNC m-to-m (UG-491A/U) 2 - coax cables term. w/BNC males (10503A)  SIGNAL GENERATOR  H25-355D  ADAPTER  UG-491A/U  10503A  For Service Manuals Contact MURITHON TECHNICAL SERVICES 8 Cherry Tree Rd, Chinnor Oxon OX9 40Y  Tel: 0184-351894 Fac: 01844-352554 Email: enquiries@mauritron.co.uk		g. Perform steps e and f for I. F. BANDWIDTH settin					
3KC  1KC  -95  -86  -71  4. I.F.GAIN SET ACCURACY: 70-db section: ±0.5 db 10-db section: ±0.1 db  a. Equipment Required:  Precision 10-db/step Attenuator (hp H25-355D)  Precision 1-db/step Attenuator (hp H25-355C)  Signal Generator (606A)  Adapter, BNC m-to-m (UG-491A/U)  2 - coax cables term. w/BNC males (10503A)  SIGNAL GENERATOR  H25-355C  H25-355D  ADAPTER  J0503A  ADAPTER  J0503A  ADAPTER  J0503A  ADAPTER  J0503A			100KC	dbm	-75		-60
4. I.F.GAIN SET ACCURACY: 70-db section: ±0.5 db 10-db section: ±0.1 db  a. Equipment Required: Precision 10-db/step Attenuator (hp H25-355D) Precision 1-db/step Attenuator (hp H25-355C) Signal Generator (606A) Adapter, BNC m-to-m (UG-491A/U) 2 - coax cables term. w/BNC males (10503A)  SIGNAL GENERATOR H25-355C H25-355D H25-355D ONO OX9 4QY Tel: 01844-351634 Fax: 01844-35254 Email: enquiries@mauritron.co.uk			10KC		-95		-80
4. I.F.GAIN SET ACCURACY: 70-db section: ±0.5 db 10-db section: ±0.1 db  a. Equipment Required: Precision 10-db/step Attenuator (hp H25-355D) Precision 1-db/step Attenuator (hp H25-355C) Signal Generator (606A) Adapter, BNC m-to-m (UG-491A/U) 2 - coax cables term. w/BNC males (10503A)  SIGNAL GENERATOR  H25-355D  ADAPTER UG-491A/U  10503A  ADAPTER UG-491A/U  10503A  For Service Manuals Contact MAUNITRON TECHNICAL SERVICES 8 Cherry Tree Rd, Chinnor Oxon OX9 40Y Tel: 01844-351894 fax: 01844-35254 Email: enquiries@mauritron.co.uk			3KC		-95	<del></del>	-80
10-db section: ±0.1 db  a. Equipment Required: Precision 10-db/step Attenuator (hp H25-355D) Precision 1-db/step Attenuator (hp H25-355C) Signal Generator (606A) Adapter, BNC m-to-m (UG-491A/U) 2 - coax cables term. w/BNC males (10503A)  SIGNAL GENERATOR  H25-355D  H25-355D  ADAPTER  OUTPUT  AD			1KC		-86		-71
Precision 10-db/step Attenuator (hp H25-355D) Precision 1-db/step Attenuator (hp H25-355C) Signal Generator (606A) Adapter, BNC m-to-m (UG-491A/U) 2 - coax cables term. w/BNC males (10503A)  SIGNAL GENERATOR  H25-355C  Oxon OX9 4QY Tel: 01844-351894 Fax: 01844-352554 Email: enquiries@mauritron.co.uk		4. I.F.GAIN SET ACCURACY	o ab acciton. Io				
5-29 b. Set Sig Gen for 20 Mc.		Precision 10-db/step A Precision 1-db/step A Signal Generator (606A Adapter, BNC m-to-m 2 - coax cables term.  SIGNAL GENERATOR	MAURITR 8 Che ( Tel:- 01844	ON TECHNICAL SE Prry Tree Rd, Chii Oxon OX9 4QY -351694 Fax:- 0184	RVICES nnor 1-352554		
,	5-29	b. Set Sig Gen for 20 Mc.					

Table 5-7. Performance Check Test Card, 851B (cont'd)

Ref		Proced	lure		Min	Act.	Max
	4. I. F. GAIN S	ET ACCURACY	(cont'd):				
	I. I. SY SW VE	F. GAIN inner co F. GAIN outer co F. BANDWIDTH NC VEEP TIME ERT DISPLAY	ntrol. 0 100K LINE 3 MII SQ				
	1	tternal Attens for					
	i	Sig Gen for 6.0		3.77			
5-30a	1	ternal Attenuato	~	•			
		AIN (DB) to 10.	- 101 -0 45 100	c:	m 5.2		6.8
	h. If nece 851 di	essary, adjust Si splay.	g Gen for 6.0-	-cm			
	setting	I. F. GAIN positions is changed, if n for 6.0-cm 85.	necessary adj	n set of ust			
		Ext Atten	I. F. GAIN	7			
		20	20	cr	n 5.2		6.8
		30	30	•	5.2		6.8
		40	40		5.2		6.8
		50	50		5.2		6.8
		60	60		5.2		6.8
		70	70		5.2		6.8
	j. Set I. F and Ex	C. GAIN inner and ternal Attenuator	outer control	s to 0,			
	k. Adjust	Sig Gen for 6.0-	cm 851 displa	у.			
5-30b	m. Set Ext inner o	ernal Atten for I control to 1.	-db loss, I.F				
	n. If nece	ssary, adjust Sig	Gen for 6 0-	cm 851 dienlay	5.8		6.2
	p. Other l	.F.GAIN positio					
	oor urs	Ext Atten	I. F. GAIN	]			
		2	2	cm	5.8		6.2
		3	3		5.8		6.2
		4	4		5.8		6.2
		5	5		5.8		6.2
		6	6		5.8		6.2
		7	7		5.8		6.2
		8	8		5.8		6.2
		9	9		5.8	<del></del>	6.2
		10	10		5.8		6.2

Table 5-7. Performance Check Test Card, 851B (cont'd)

a. Equation in the second seco	EEP TIME VE BANDWIDTH NC	ired: nter (5245L) (202C) or (606A) BNC male, br BNC males (1  5245L  A  202C OUTPUT  MODULATED  ERNIER	anana plug (11001 10503A)  606A MOD RET 20MC  CAL 1MC EXT about 3 volts 10 PERIOD AVER 20 Mc -20 dbm EXT DC fully cw	851A-c-6			
b. SW I. F SYI c. LF Con Sig O M d. SW LF	EEP TIME VE Con SELECT.  OD SELECT.  OD AMPLITU  EEP TIME .  OSC	THER (5245L) (202C) OF (606A) PENC male, by PENC males (1)  S245L  A 202C OUTPUT  MODULATED  ERNIER	CAL 1MC EXT about 3 volts 10 PERIOD AVER 20 Mc -20 dbm EXT DC fully cw 3 ms/cm 333 cps (reading	851A-c-6			
I. F SYI C. LF Coo Sig O M M d. SW.	EEP TIME VET BANDWIDTH NC	A 202C SYNC OUTPUT  MODULATED  ERNIER	CAL 1MC EXT about 3 volts 10 PERIOD AVER 20 Mc -20 dbm EXT DC fully cw 3 ms/cm 333 cps (reading	RAGE			
I. F SYI C. LF Coo Sig O M M d. SW.	C.BANDWIDTH NC Oscillator . Inter Gen utput OD SELECT . OD AMPLITU	f	1MC EXT about 3 volts 10 PERIOD AVEI 20 Mc -20 dbm EXT DC fully cw 3 ms/cm 333 cps (reading	į			
Cot Sig O M M d. SW LF	Inter	DE	10 PERIOD AVER 20 Mc -20 dbm EXT DC fully cw 3 ms/cm 333 cps (reading	į			
d. SW LF	EEP TIME . Osc		3 ms/cm 333 cps (reading	of 30 ms			
e. Not	o 'minus la sur sur sur s		/	ŀ			
•		les displayed. times:		cycles	9.7		10.3
	LF Osc set for cps	*Counter Reading	SWEEP TIME Setting				
	100	100 ms	10 ms/cm	cycles	9.7	T-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1-1	10.3
	33.3	300 ms	30 ms/cm		9.7		10.3
	10	1 sec	0.1 sec/cm		9.7	<del></del>	10.3
	3.33	3 sec	0.3 sec/cm				10.3
-		L	<u> </u>	1	9.1	<del></del>	10.3
Ĺ				۱			
		MA Tel:	NRITRON TECHNICAL SER 8 Cherry Tree Rd, Chinr Oxon OX9 4QY - 01844-351694 Fax:- 01844-3	VICES nor 352554			
	f. For	LF Osc set for cps  100  33.3  10  3.33	for cps Reading  100 100 ms  33.3 300 ms  10 1 sec  3.33 3 sec  1 10 sec  * Set for 10 Period Av  Tel:	LF Osc set   *Counter   SWEEP TIME   for cps   Reading   Setting	LF Osc set	LF Osc set for cps	LF Osc set for cps

#### 5-17. LINEAR.

Accuracy Specification: ±3% of full scale

- a. Set VERT DISPLAY to LIN, inner I.F.GAIN control to 10.
- b. Adjust I. F. GAIN outer control for low-noise base-line trace (set to about 50).
- c. Adjust VERT POS and TRACE ALIGN to align trace base exactly with graticule base line.
- d. Adjust 8551 ATTENUATOR (DB) and Signal Generator output for 7.0-cm 851 display.
  - e. Set inner I. F. GAIN control to 4 (attenuate 6 db). Display should be no higher than 3.7 cm and no lower than 3.3 cm.

#### 5-18. SQUARE.

Accuracy Specification: ±5% of full scale

- a. Set VERT DISPLAY to SQ, inner I. F. GAIN control to 10.
- b. Adjust I. F. GAIN outer control for low-noise base-line trace, bottom of which should coincide with first horizontal axis.
- c. Adjust ATTENUATOR (DB) and Signal Generator output for 7.0-cm 851 display.
  - d. Set inner I. F. GAIN control to 7 (attenuate 3 db). Display should be no higher than 3.85 cm and no lower than 3.15 cm.

## 5-19. LOG.

Accuracy Specification: ±2 db

- a. Set VERT DISPLAY to LOG, I.F. GAIN to 70+0.
- b. Adjust 8551 ATTENUATOR (DB) and level at Signal Generator for 7.0-cm 851 display.
  - c. Set I.F.GAIN outer control to 60. Display should be within 5.8 to 6.2 cm.
- d. If display does not coincide exactly with 6.0-cm graticule line, at Signal Generator readjust signal level for coincidence.
- e. Proceed in same manner for other I. F. GAIN positions; see Table 5-8. If necessary, readjust signal level for trace coincidence with graticule line.

Table 5-8. VERT DISPLAY Accuracy Check

Set Display	I. F. GAIN	Setting	VERT DISPLAY
to (cm)	From	То	Limits (cm)
7.0	70	60	5.8 to 6.2
6.0	60	50	4.8 to 5.2
5.0	50	40	3.8 to 4.2
4.0	40	30	2.8 to 3.2
3.0	30	20	1.8 to 2.2
2.0	20	10	0.8 to 1.2

f. With I. F. GAIN outer control at 10, set inner control to 10, and adjust level at Signal Generator for 851 2.0-cm display. Set outer control to 0. Display should be within 0.8 to 1.2.

#### 5-20. I.F.BANDWIDTH ACCURACY CHECK.

Specification: Individual bandwidths are calibrated within ±20%. Bandwidth repeatability and stability typically better than ±3%.

### 5-21. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
10*	VHF Attenuator	1
D**	Shielded coax cable term. with BNC males (10503A)	1
	8551 RF Section	
*Tal	ole 5-1 **Table 5-2	

- 5-22. 1MC, 100KC, AND 10KC BANDWIDTHS.
- a. Connect Attenuator 355D between 851 and 8551, and find 2-Gc BWO signal; see Paragraphs 5-96 through 5-100.
- b. Check bandwidths as indicated in Table 5-9, in each case recording actual bandwidth in cm.
- c. Switch to any other setting of I. F. BANDWIDTH, then back to setting for bandwidth under test. Bandwidth should be within ±3% of recorded bandwidth (no < 0.77, no > 1.23 cm).

Table 5-9. I. F. Bandwidth Accuracy Checks (1MC, 100KC, 10KC)

	Settings					Record
I.F. BW	SPECT WIDTH	SWEEP TIME	Adjust Atten* for Display of	Read Display at	Spec Limits, BW (cm)	Actual BW (cm)
1MC	1 Mc/cm	3 ms/cm	7.0 cm	5.0 ст	0.8 - 1.2	
100KC	100 Kc/cm	3 ms/cm	7.0 cm	5.0 cm	0.8 - 1.2	
10KC	10 Kc/cm	3 ms/cm	7.0 cm	5.0 cm	0.8 - 1.2	

\*I.F.GAIN (DB) and external Attenuator (355D)

# MAURITRON TECHNICAL SERVICES 8 Cherry Tree Rd, Chinnor

Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk

Section V Paragraphs 5-23 to 5-29

Model 851B

Table 5-10. I.F. Bandwidth Accuracy Checks, 3KC and 1KC

	Settings					T
I.F.BW	SPECT WIDTH	SWEEP TIME	Adjust Atten** for Display of	Read Display at	Spec Limits, BW (cm)	Record Actual BW (cm)
3KC	10* Kc/cm	3 ms/cm	7.0 cm	5.0 cm	2.4 - 3.6	
1KC	10* Kc/cm	10 ms/cm	7.0 cm	5.0 cm	0.8 - 1.2	

<sup>\*</sup>Calibrated to 1 Kc/cm

#### 5-23. 3KC AND 1KC BANDWIDTHS.

- a. Calibrate SPECTRUM WIDTH to obtain increased resolution; see Paragraph 5-103.
  - b. Check as indicated in Table 5-10.
- c. Switch to any other setting of I.F.BANDWIDTH, then back to setting for bandwidth under test.
  - At 3KC, bandwidth should be between 2.31 and 3.69 At 1KC, bandwidth should be between 0.77 and 1.23

## 5-24. I.F.INPUT SENSITIVITY CHECK.

Specification - 6-cm Ve	Input Required* for ertical Display
1-Mc bandwidth	-62 to -53 dbm
100-Kc bandwidth	-75 to -60 dbm
10-Kc bandwidth	-95 to -80 dbm
3-Kc bandwidth	-95 to -80 dbm
1-Kc bandwidth	-86 to -71 dbm
	and I. F. VERNIER full

### 5-25. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
9*	Signal Generator set for 20 Mc; Generator must have calibrated power output (606A)	1
D**	Shielded coax term. with BNC males (10503A)	1
*Tabl	.e 5-1 **Table 5-2	•

#### 5-26. PROCEDURE.

a. Connect Signal Generator to I.F. INPUT on 851 rear panel. Check Line Switch 115/230 is set for voltage of power source, and connect 851 directly to 115/230V power source.

## b. Set 851:

DCC UUI,						
I. F. GAIN (DB)						. 70 + 10
I. F. VERNIER						. fully ccw
I. F. BANDWIDT	Η					1MC

c. At Signal Generator, adjust for  $6.0\ cm$  display on 851.

- d. At Signal Generator, read output signal level, and take into consideration loss through input cable. Limits are given in Table 5-11.
- e. Perform steps c and d at other settings of I. F.  ${\tt BANDWIDTH.}$

Table 5-11. I. F. Input Sensitivity Check

I. F. BANDWIDTH Setting	Input Signal Level Limits* (dbm)
1 MC	-62 to -53
100KC	-75 to -60
10KC	-95 to -80
3KC	-95 to <b>-80</b>
1KC	-86 to -71

<sup>\*</sup>For 6-cm deflection with I. F. GAIN at 80 db and I. F. VERNIER full counterclockwise

## 5-27. I.F. GAIN SET ACCURACY CHECK.

i					
	Specification	- 70-db	section:	±0.5 db	
		10-db	section:	±0.1 db	

#### 5-28. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
12*	Precision 10-db/step Attenuator (H25-355D)	1
13*	Precision 1-db/step Attenuator (H25-355C)	1
9*	Signal Generator (606A)	1
D**	Coax cable term. w/BNC males (10503A)	2
J**	Adapter, BNC male-to-male (UG-491A/U)	1
*T	able 5-1 **Table 5-2	

5-29. SETUP. Accuracy of I. F. input attenuator I. F. GAIN (DB) is excellent and therefore cannot be checked

<sup>\*\*</sup>I. F. GAIN (DB) and External Attenuator (355D)

without special equipment of exceptional accuracy. The H25-355C and H25-355D are calibrated at 20 Mc to give required accuracy.

- a. Set Signal Generator for 20 Mc.
- b. Check that  $851\ 115/230V$  Line Switch is set for voltage of power source, and connect 851 directly to  $115/230V\ 50-1000$  cps source.
  - c. Set 851: I. F. GAI

I. F. GA	N inner	c	or	ıtr	ol									0
I. F. GAI	N outer	: c	or	ıtr	ol									0
I.F.BA	MDWID!	ГΗ	Ι.									10	<b>10</b> C	ζC
SYNC .												]	LII	VΕ
SWEEP	TIME						3	N	III	٠L	IS]	EΟ	:/c	M
VERT D	ISPLA	Y											. 8	SQ

- d. Connect H25-355D to H25-355C via Adapter UG-491A/U. Set both External Attenuators for 0.
- e. Connect one Attenuator to Signal Generator RF Output and other to I. F. INPUT on 851 rear panel.

#### 5-30. PROCEDURE.

### a. I. F. GAIN Outer Control.

- (1) With External Attenuator and I. F. GAIN outer control both set for 0, adjust level at Signal Generator for 6.0-cm trace on 851 CRT.
- (2) Set External Attenuator for 10-db loss; set I. F. GAIN outer control to 10. Trace should be within 5.2 and 6.8 cm. Note: With signal reference at 6.0 cm, ±0.8 cm is approximately ±0.5 db.
- (3) If necessary, adjust level at Signal Generator to return reference trace to 6.0 cm.
- (4) Check other I. F. GAIN positions in same manner, turning both External Attenuator and I. F. GAIN in 10-db steps to 70 db.

At each 10-db change, trace should be within 5.2 and 6.8 cm.

Note: If necessary, readjust signal level at 10-db change to maintain reference at 6.0 cm.

### b. I. F. GAIN Inner Control.

(1) Set I. F. GAIN inner control							0
I. F. GAIN outer control							0
External Attenuator .							0
Signal Generator	for	6	. 0	cr	n	tra	ace
•			on	8	51	С	RT

(2) Set External Attenuator for 1-db loss; set I. F. GAIN inner control to 1.

Trace should be within 5.8 and 6.2. Note: With signal reference at 6.0 cm,  $\pm 0.2$  cm is approximately  $\pm 0.1$  db.

- (3) If necessary, adjust level at Signal Generator to return reference trace to 6.0 cm.
- (4) Check other I. F. GAIN inner control positions in same manner, maintaining reference trace at 6.0 cm.

At each 1-db change, trace should be within 5.8 and 6.2 cm.

#### 5-31. SWEEP RATE ACCURACY CHECK.

Specification - Sweep Rate Accuracy: ±3%

#### 5-32. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
7*	Electronic Counter (5245L)	1
8*	Low-Frequency Oscillator (202C)	1
9*	Signal Generator (606A)	1
C**	Shielded cable, BNC male to banana plug (11001A)	1
D**	Shielded cable, term. with BNC males (10503A)	1
*Ta	ble 5-1 **Table 5-2	

5-33. SETUP. Connect as indicated in Figure 5-3, except that 851 can be connected directly to 115/230V, 50/1000-cycle source.

5-34. PROCEDURE. Instructions assume use of equipments shown in Figure 5-3.

a. Set 851:

SWEEP TIME VERNIER.				CAL
I. F. BANDWIDTH				1MC
SYNC				EXT

- b. Set Low-Frequency Oscillator for output of about 3 volts (AMPLITUDE control at about 90).
- c. Set Counter FUNCTION selector for 10 PERIOD AVERAGE.
  - d. Set Signal Generator controls:
    RANGE, FREQUENCY. . . for 20 Mc output
    ATTENUATOR, VERNIER for -20 dbm output
    MODULATION SELECTOR . . . . EXT DC
    MODULATION AMPLITUDE . . . . fully cw
- e. With equipments connected as shown in Figure 5-3, output of Oscillator 202C, monitored by Counter, is modulating the 20-Mc output of the Signal Generator. Output of Signal Generator is displayed on 851 CRT.
  - f. To check sweep-rate accuracy specification:
  - (1) Set SWEEP TIME to 3 MILLISEC/CM.
  - (2) Adjust Oscillator 202C for output of precisely 333 cps (reading of 30 ms with Counter set for 10 period average).

Sweep rate is within specifications if 9.7 - 10.3 cycles appear on display.

Note: Period of 333-cycle signal is 0.003 second.

g. Check other SWEEP TIME positions, using procedure indicated in Table 5-12.

With settings as specified in Table 5-12, sweep rate is within specifications if 9.7 - 10.3 cycles appear on display.


# TEST CARD

1			
	1		
	1000		
	300 mis		
	1 sec	9.1	
	3 sec	73	
	10 sec	!	


	100	1	

		Tabib	9		1 *	
Ref	Seq				· ·	
INSTRU	MENT OFF				•	
	1	<b>Mec</b> has	;			
5-39	2	Presta	n de la companya de			
		N. Baser Madelle.	17.5		For Service Me	Krieter   Sterich
			6 7 F ( ) (1417)		8 Cherry Tree	NICAL DERMANS S
			The Maria		Oxon O:	(9.4 <b>0Y</b>
		). ).	(ME CLAS (ME VER)	w.	Tat- 01844-351694   Email:- enquinees	mau <b>ntron.co.uk</b>
	ŀ	<b>7</b>	. VIII - V I D	X.		
		· 1	PUATE I			
		;	9.54 (C.) 41			
INSTRU	MENT ON					
	3 1				1150	2:
5-40		LV POWEE	71 H ·			
	4	1	$\chi_{(0)} = \chi_{(0)} \circ \gamma$	vite	-14 3	
			10.3 × 10.0 × 10	. vdc		and difference on a
			$P_{A}: \mathbb{R}^{n}$	P-1/		strationing as a
	5			vdc	-14 )	
			1 ( 	vdc.		
				12 F		
	6	* 1		vdc	-92.5	p. 100 m.
1			100 p	) <b></b>		
			i de la companya de l	mv p-p		11.00
		HV POWER	of S	• •		
		NA NOW F	MAN STATE			
		v i	to the second			
			411 (1)			
		. '	$\frac{N_{\rm bol} \sim N_{\rm bol}}{100000000000000000000000000000000000$			
	7		b.	vdc	2340	

				Record			
Ref	Seq	Operation		Min	Act.	Max	
		HV POWER SUPPLY (cont'd)					
	8	INTENSITY & FOCUS cw cathode V change	vdc			60	
	9	post accelerator V	vdc	4750		5450	
- 10		FOCUS, Astigmatism, & Pattern					
5-46	10	Shape adjusts produce sharp spot FOCUS produces sharp spot at all INTENSITY levels					
		HORIZONTAL AMPLIFIER					
5-49	11 12	Adj TRACE ALIGN Calibration: cw rotation of HORIZ POS moves trace right - center trace Adj Horiz Gain A6R54 for 10-cm deflection HORIZ POS trace movement - each direction	0.5		1.5		
		TIME BASE					
5-51	13	Sweep Calibration: SWEEP TIME VERNIER - CAL					
		MOD FREQ 10 PERIOD (cps) MEAS SWEEP TIME AI	DJ				
		100 100 ms 10 ms A6 33.3 300 ms 30 ms A6	3R29 3R30 3R31				
		3.33 3 sec .3 sec A6	SR32 SR33 SR34				
5-54e	14	Sweep Linearity: SWEEP TIME: 3 ms/cm Distance between successive positive mod peaks	cm	0.8		1.2	
5-55	15	SWEEP TIME VERNIER:  Mod frequency: 100 cps SWEEP TIME: 3 ms/cm SWEEP TIME VERNIER: full ccw 1-cycle waveform width	cm			1.0	
5-58	16	SINGLE SWEEP & sweep amplitude	V p-p	9.7	<del></del>	10.3	
5-61	17 18 19 20 21	SYNC & OUTPUT CHECKS  EXT SYNC: 1 cps - 5 Kc  VERT OUTPUT check  HORIZ OUTPUT check  LINE sync check  BASE LINE CLIPPER blanks at least  lower 2 cm of trace					
5-69	22	CRT CHECKS	N	<b>AURIT</b> RON TE	Manuals Contact CHNICAL SERVICES ee Rd, Chinnor	<b>3</b>	
	23	Check TRACE ALIGN		Oxon OX9 4QY Tel: 01844-351694 Fax:- 01844-352554 Email:- enquiries@mauritron.co.uk			
		Check 1 Kc pattern for uniform focus at normal intensity					
	24	Blanking: no retrace all sweep speeds					



				Record			
Ref	Seq	Operation	Min	Act.	Max		
5-76	25	VERTICAL AMPLIFIER  Calibration:  a. cw rotation of VERT POS moves trace upward  b. Align trace with base line - no input c. Set 20-Mc Sig Gen for 4.0 ±0.1 vdc detected input to Vert Ampl d. Adjust A7R15 for 6.0 cm defl e. Repeat steps b & d until both conditions are met.					
5-81	26	Video Bandwidth:  I. F. BANDWIDTH: 1MC  Set 50 Kc Vertical Amplifier input for 7.0 cm defl Increase freq to 1.2 Mc, same input V Vert defl cm	5.0				
		Other I. F. BANDWIDTH positions 1 Kc ref 100KC I. F. BANDWIDTH Kc 10KC 3KC 1KC	160 32 9.6 3.2		240 48 14.4 4.8		
		I.F.BANDWIDTH VERT DISPLAY - LIN		<del></del> -			
5-82	27	1 Mc Alignment: Adj Detect. Tune A2A7T1 Tune Imped Adj A2A6L11 Preset A12C1, A12C3, A12L2					
5-86	28	1MC BANDWIDTH check Mc	1.4		2.2		
5-89	30	100 Kc Alignment: Tune 100 Kc BW adj in A3 & A5 Tune Imped Adj A2A2L1 100KC BANDWIDTH check: Kc	80		120		
-91, -j	31	Final 1MC and 100KC Bandwidth Adjusts  Connect 851 to 8551  Perform 1-14, Fig. 3-3  SWEEP TIME 3 ms/cm  I.F.BANDWIDTH . 1MC  Input signal - anywhere between  10 Mc and 5 Gc  Set SPECTRUM WIDTH 1 Mc/cm  VERT DISPLAY LOG  Check symmetry and if necessary  readjust A12C1, A12C3, A12L1.  Set SPECTRUM WIDTH 100 Kc/cm  I.F.BANDWIDTH 100KC  Check symmetry and if necessary  readjust A12C1, C3, L2.  Switch back and forth between LIN and LOG  and between 1MC and 100KC I.F.BANDWIDTH			120		
-94 -101	32	readjusting as required for best compromise on amplitude and symmetry while keeping bandwidths within specifications.  1-3-10 Kc Alignment: SPECTRUM WIDTH . 10 Kc I. F. BANDWIDTH 10 Kc					



<del></del>				Record		
Ref	Seq	Operation	Min	Act.	Max	
		I.F.BANDWIDTH (cont'd)				
5-101		1-3-10 Kc Alignment (cont'd):  If desired, remove RF Ckt Assy A2 cover.  Preset Bal Adj A2A3C5 & A2A4C8 at 1/3 mes	n.			
		Tune 1-10 Kc BW Adj A2A3C4, A2A3C2, A2A4C5, & A2A4C9 for max BW				
		Adj Imped Adj A2A3L3 and Freq Adj A2A4C7 for max defl Set I. F. GAIN for 7.0 cm display. Check bandwidth at 5.0 cm cm Readjust until bw at 5.0 cm vert defl is 1 cm	1		1	
5-102	33	3KC & 1KC BANDWIDTH checks: See Calib, Paragraph 5-103.	-			
5-104		BW at 5.0 cm with 7.0 cm max defl 3KC cm 1KC Reinstall Assy A2 cover	2.4 0.8		3.6 1.2	
5-105	34	Recheck 10 Kc BW adj: SPECTRUM WIDTH - 10 Kc I. F. BANDWIDTH - 10 Kc				
		BW at 5.0 cm with 7.0 cm max defl cm	0.8		1.2	
5-107	35	AUTO SELECT CHECK  Optimum BW is automatically selected For check see Paragraph 5-110 and Table 5-19				
5-111		I. F. SENSITIVITY				
5-114	36	I. F. GAIN - 80 db I. F. GAIN VERNIER - full cw Power input for 6-cm vert defl				
		I. F. BW				
		1 Mc dbm	-53		-62	
		100 Kc	-60		-75	
		10 Kc	-80		-95	
		3 Kc	-80		-95	
		1 Kc	-71		-86	
5-115	37	Noise Level				
		I. F. BW For Service Manuals Con			0.45	
		1 Mc MAURITRON TECHNICAL SERV 8 Cherry Tree Rd, Chinn				
		100 Kc Oxon OX9 4QY  Tel:- 01844-351694 Fax:- 01844-3  Email:- enquiries@mauritron.or			0.45	
		10 Kc			0.45	
		3 Kc			0.45	
		1 Kc			0.45	

Table 5-13. 851B Check and Adjustment Test Card (cont'd)

Dof	See 1	<b>^</b>		Record			
Ref	Seq	Operation	Min	Act.	Max		
5-116		VERTICAL DISPLAY		-			
		I. F. GAIN - 20 db I. F. BW - 1 Mc					
5-117c	38	Check trace alignment					
5-118	39	LOG Display Adj: Adj input for 1 cm defl					
		Increase I. F. GAIN to 40 db Adj A11R13 for 3 cm defl					
		Increase I. F. GAIN to 60 db Adj A11R14 for 5 cm defl					
		Increase I. F. GAIN to 70 + 10 db Adj A11R20 for 7 cm defl					
		I. F. GAIN VERNIER - full cw Gain decrease: db Reset VERNIER full ccw	1				
5 -120	40	LOG Display Linearity:  Decrease I. F. GAIN in 10-db steps  Trace should decrease 1 cm/step  Error at each cm div cm  Repeat 40 for other I. F. bandwidths			±0.2		
5-122	41	SQ Display Adj: Adj input for 7.0 cm defl with 20 + 10 db I. F. GAIN Decrease I. F. GAIN 6 db Adj A11R2 for 1.75 cm defl					
-123	42	SQ Display Linearity: I. F. GAIN		Vert Defl - cr	<u>n</u>		
		-3 db -6 db	3.15 1.40		3.85 2.10		
		Repeat 42 for other I.F.BWs					
	43	LIN Display Linearity: I. F. GAIN		Vert Defl - cn	<u>n</u>		
į		-6 db -12 db	3.29		3.71		
ļ			1.54		1.96		
-124		Repeat 43 for other I. F. BWs FINAL I. F. BANDWIDTH ADJUSTS					
		(with 8551)					
-127	44	Crystal Filter Balance: VERT DISPLAY - LOG					
		60-Mc input -'7-cm defl Tune A2A3C5 & A2A4C8 for best sym- metrical display & best compromise:					
		8551 851 SPECTRUM I.F. WIDTH BANDWIDTH					
		300 Kc/cm 10 Kc 100 Kc/cm 3 Kc 30 Kc/cm 1 Kc					

Table 5-13. 851B Check and Adjustment Test Card (cont'd)

				Record	
Ref	Seq	Operation	Min	Act.	Max
		FINAL I.F.BANDWIDTH ADJUSTS (with 8551) (cont'd)			
5-128	45	1 Mc Bandpass Filter Adjustments: VERT DISPLAY - LIN SPECTRUM WIDTH - 1 Mc/cm I. F. BW - 100 Kc Center display on CRT I. F. BW - 1 Mc Adj 851 A12C1, A12C2, & 8551 A9A2L2 for symmetrical display and 1-Mc BW			
5-130		I. F. BW - 100 Kc Max ampl on CRT is at same freq as 1 Mc bandwidth			

## 5-37. CHECKS AND ADJUSTMENTS.

- 5-38. Procedures for checking and adjusting the 851B are provided in Paragraphs 5-39 through 5-130.
- a. Most of the procedures call for the use of other equipment; only those instructions pertinent to the procedure are given -- for full operating instructions use the Manual supplied with the instrument.
- b. Unless specified otherwise, the 851B is not connected to the 8551, but is powered separately. Procedures assume a 115-volt line.
- c. Removal of cover plates is simple; instructions are given in Paragraph 5-6 and will not be referenced again.
- d. When making a thorough check of the instrument, it is recommended that procedures be performed in the order presented.

## 5-39. PRELIMINARY ADJUSTMENT PROCEDURE.

a. Equipment Required.

Ref No.	Equipment	No.
3	Variable Transformer	1

- b. Remove 851 top cover plate.
- c. Set 115V/230V slide switch on rear panel to 115V.
- d. Set front panel controls:

  BASE LINE CLIPPER ... max ccw
  INTENSITY ... 1 SEC/CM
  SWEEP TIME ... 1 SEC/CM
  SWEEP TIME VERNIER ... max ccw
  SYNC ... INT
  VERT DISPLAY ... LIN
  I. F. BANDWIDTH ... 100 KC
  I. F. GAIN ... 30 DB
- e. Set Int Level R2 max cw. (R2 is located to right of cathode-ray tube toward rear of instrument; see Figure 5-1.)

f. Set Variable Transformer to minimum. Connect 851 to power source through Variable Transformer, and increase Transformer voltage slowly to 115 volts.

## 5-40. LV POWER SUPPLY ADJUSTMENTS.

## 5-41. EQUIPMENT REQUIRED.


Ref No.	Equipment	No.
3*	Variable Transformer, set for 115V input	1
1*	Oscilloscope (130C)	1
2*	Digital DC Voltmeter (3440A)	1
A**	Shielded cable, banana plug to banana plug (11000A)	1
B**	Shielded cable, banana plug to alligator clips (11037A)	1
*7	Table 5-1 **Table 5-2	

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel: 01844-351694 Fax: 01844-352554
Email: enquiries@mauritron.co.uk

## 5-42. PROCEDURE.

- a. Connect Digital Voltmeter and vertical input of Oscilloscope in parallel using shielded cables.
- b. Locate Low-Voltage Power Supply Board A9 (see Figure 5-1).
  - (1) Make ground connection at A9C6 (point A); Figure 5-36).
  - (2) Measure and adjust the low-voltage power supplies in accordance with Table 5-14. (Normal resistances to ground are given for reference in Table 5-15.)

in the secondary



m mei la distributa (metilikusis ist singa <del>kandandana</del> metiributa daki sikuran metiributa mendilikan kecamatan daki sikuran singa mendilikan kecamatan daki sikuran singa mendilikan kecamatan daki sikuran daki sik

Figure 5-1. Measurement Setup, Check and Adjustment of LV Power Supply

Table 5-14. LV Power Supply Measurement Data

	Supply	Meas Point	Ref Fig. 5-36	Adjust	115V Line (vdc)	Reg (max	to 126.5V ine Max Ripple (mv p/p)
	+ 15vdc	+A9C6	<b>B</b>	A9R14	+15±0.1	±0.03	1.5
	- 15vdc	-A9C10	<b>Q</b>	A9R29	-15±0.1	±0.03	6.0
-	+100vdc	+A9C1	<b>P</b>		100±7.5	±3.0	75

Table 5-15. Resistances to Ground, LV Power Supply, Reference Data

Supply	Meas Point	Ref (Fig.5-36)	Normal Resistance* (ohms)
+15 vdc	+A9C6	B	> 300
-15 vdc	-A9C10	Û	> 27
+100 vdc	+A9C1	<b>1</b>	> 9000

\*As measured with electronic volt-ohmmeter such as hp 410B, 410C, or 412A

## 5-43. HV POWER SUPPLY CHECK.

## 5-44. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
3*	Variable Transformer, set for 115V input	1
4*	Clip-On DC Milliammeter (428B)	1
5*	DC Voltmeter (410C)	1
6*	DC Voltage Divider (11045A)	1
*T	able 5-1	

## 5-45. PROCEDURE.

a. <u>Instrument Condition</u>. Top cover plate, left side plate, and <u>CRT</u> protective cover (on rear of instrument) are removed.

## b. Cathode Current.

- (1) Wiring to CRT base is shown in Figure 5-2.

  Note gray wire designated (A); clip Milliam-meter probe around gray conductor (A).
- (2) Turn INTENSITY max cw, and adjust Int Level R2 (to right of CRT, see Figure 5-1) for 0.5 ma.
- (3) Turn INTENSITY fully ccw, and check that beam is extinguished (no cathode current).

STRAPS FILAMENT TO CATHODE CURRENT FOR CRT CATHODE FROM FILAMENT HV PWR SUP A8 GRA GRA GRA TO GRID, FROM HV PWR SUP A8, VIA INTENSITY (RI) FROM HV PWR SUP A8, VIA FOCUS (R4) GROUND FROM +100VDC VIA ASTIG (R3)

Figure 5-2. CRT Base as Seen from Rear, CRT Protective Cover Removed

c. Cathode and Post-Accelerator Voltages. Using a dc Voltmeter and Voltage-divider Probe, check as shown in Table 5-16. Location of A8 board is shown in Figure 5-1.

Table 5-16. HV Power Supply Voltages

Voltage Checked	Settings	Meas Point	Ref (Fig. 5-34)	Test Limits (vdc)
Cathode	INTENSITY, full ccw	-A8C7	A	2500 ±190
	INTENSITY, full cw FOCUS, full cw			
Post Acceler	same	Junction A8R3, A8C4	11	

\*The change in cathode voltage should not exceed 60 vdc.

#### 5-46. FOCUS CHECK AND ADJUSTMENT.

## 5-47. CHECK.

- a. With the 851B connected to a source of power (either the Model 8551 or through a Variable Transformer set for 115V), adjust FOCUS for a sharp spot on the CRT.
- b. Turn INTENSITY through its range and adjust FOCUS to maintain a sharp spot at all INTENSITY levels.
- 5-48. ADJUSTMENT. If a sharp spot is not obtained, adjust FOCUS, Astig adjust R3, and Pattern Shape adjust A7R22 for a sharp spot. Then check that a sharp spot can be obtained as INTENSITY is turned

through its range. See Figure 5-1 for location of Astig adjust R3 and Pattern Shape adjust A7R22.

# 5-49. HORIZONTAL AMPLIFIER CHECKS AND ADJUSTMENTS.

## 5-50. CALIBRATION.

- a. Adjust TRACE ALIGN so that horizontal trace is parallel to horizontal axis.
- b. Rotate HORIZ POS (R9); trace should move to right. Adjust HORIZ POS to center trace on graticule.
- c. Adjust Horiz Gain adjust A6R54 (see Figure 5-31) for 10 cm of horizontal deflection.
- d. Turn front panel HORIZ POS adjust full cw and note how far trace moves, then turn HORIZ POS full ccw and note trace movement; trace should move at least 1.0  $\pm 0.5$  cm each direction.
  - e. Center trace.

#### 5-51. SWEEP TIME CALIBRATION.

### 5-52. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
3*	Variable Transformer, set for 115V	1
7*	Electronic Counter (5245L)	1
8*	Low-Frequency Oscillator (202C)	1
9*	Signal Generator (606A)	1
C**	Shielded cable, BNC male to dual banana plug (11001A)	3
D**	Shielded cable, BNC male to BNC male (10503A)	1
*Ta	able 5-1 **Table 5-2	

- 5-53. MEASUREMENT SETUP. Connect as indicated in Figure 5-3.
- 5-54. PROCEDURE. The following instructions assume use of equipment shown in Figure 5-3.
  - a. Set 851 controls:

    SWEEP TIME VERNIER. . . . . . CAL

    I, F, BANDWIDTH . . . . . . . 1 MC

    SYNC . . . . . . . . . . . . EXT
- b. Set Low-Frequency Oscillator for output of about 3 volts (AMPLITUDE control at about 90).
- c. Set Counter FUNCTION selector for 10 PERIOD AVERAGE.
  - d. Set Signal Generator controls:
    RANGE, FREQUENCY. . . . for 20 Mc output
    ATTENUATOR, VERNIER . for -20 dbm output
    MODULATION SELECTOR . . . . . EXT DC
    MODULATION AMPLITUDE . . . . full cw
- e. With equipments connected as shown in Figure 5-3, the output of Oscillator 202C, monitored by Counter, is modulating the 20-Mc output of the Signal




Figure 5-3. Measurement Setup, Sweep Calibration

Generator. Output of the Signal Generator is displayed on the 851 CRT. To check the linearity of the sweep generator output and to calibrate the 3 MILLISECOND/CM position of the 851 SWEEP TIME switch:

- (1) Set SWEEP TIME to 3 MILLISEC/CM.
- (2) Adjust Oscillator 202C for an output of precisely 333 cps (reading of 30 ms with counter set for 10 period average).
- (3) Adjust A6R29 (see Figure 5-31) so modulation peaks are precisely aligned with graticule vertical lines. The distance between each successive modulation peak on the display should not exceed 1 ±0.2 cm. This checks sweep linearity; if linearity is good in one position of SWEEP TIME it will be good in all others.
- f. Set SWEEP TIME to 10 MILLISEC/CM, and Oscillator output for precisely 100 cps. Adjust A6R30 to align the first and last modulation peak with the first and tenth vertical lines on graticule.
- g. Follow the same procedure at other positions of SWEEP TIME, using data given in Table 5-17.

## 5-55. SWEEP TIME VERNIER CHECK.

5-56. SETUP. To check that the SWEEP TIME VERNIER has the proper range, use the setup used for SWEEP TIME calibration (see Paragraph 5-51 and Figure 5-3).

Table 5-17. Sweep Time Calibration

		The Land Cult		
Mod Freq (cps)	Counter Reading*	SWEEP TIME Setting	851 Adjust	Fig. Ref
333	30 ms	3 MILLISEC/CM	A6R29	5-31
100	100 ms	10 MILLISEC/CM	A6R30	
33.3	300 ms	30 MILLISEC/CM	A6R31	
10	1 sec	.1 SEC/CM	A6R32	
3.33	3 sec	.3 SEC/CM	A6R33	
1	10 sec	1 SEC/CM	A6R34	↓
*Fo	r 10-period	average		

## 5-57. PROCEDURE.

- a. Set SWEEP TIME. . . . . 3 MILLISEC/CM SWEEP TIME VERNIER . . . . . . CAL
- b. Set Low-Frequency Oscillator for precisely 100 cps (Counter reading of 100 ms when set for 10 PERIOD AVERAGE).
- c. Rotate SWEEP TIME VERNIER full ccw, and note period of one cycle as displayed on 851 CRT; width of cycle waveform should be less than 1.0 cm.

# 5-58. SINGLE SWEEP AND SWEEP AMPLITUDE CHECKS.

#### 5-59. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
3*	Variable Transformer, set for 115V	1
2*	Digital Voltmeter (3440A)	1
C**	Shielded cable, dual banana plug to BNC male (11001A)	1
*7	Table 5-1 **Table 5-2	

#### 5-60. PROCEDURE.

- a. On the 851, set

  SWEEP TIME . . . . . . . . . 1 SEC/CM

  SWEEP TIME VERNIER. . . . . . . . cew

  SYNC . . . . . . . . . . SINGLE SWEEP
- b. Connect SWEEP OUTPUT (on rear of 851) to Digital Voltmeter. Note reading obtained.
- c. Depress SINGLE SWEEP button on 851 front panel.
- d. Note that single sweep is obtained, and note maximum positive voltage indicated by Voltmeter. Sweep amplitude should be  $10.0 \pm 0.3$  volts.

### 5-61. SYNCHRONIZATION & OUTPUT CHECKS.

## 5-62. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
3*	Variable Transformer set for 115V	1
1*	Oscilloscope (130C)	1
8*	Low-Frequency Oscillator (202C)	1
9*	Signal Generator (606A)	1
A**	Shielded cable, dual banana plug to dual banana plug (11000A)	1
C**	Shielded cable, BNC male to dual banana plug (11001A)	2
D**	Shielded cable, BNC male to BNC male (10503A)	1
*T	able 5-1 **Table 5-2	

5-63. MEASUREMENT SETUP. Similar to that shown in Figure 5-3 except that Oscilloscope replaces Counter, and is connected as noted in Paragraphs 5-65, 5-66, and 5-67.

#### 5-64. EXTERNAL SYNC CHECK.

- a. On 851, set SYNC to EXT.
- b. Set Low-Frequency Oscillator for output of 6 volts peak-to-peak.
  - c. Set Signal Generator output attenuator to -20 dbm.

- d. Vary Oscillator output frequency from 1 cps to 5 Kc and, changing 851 sweep time as required, observe signal displayed on 851. Signal displayed should be stable from 1 cps to 5 Kc.
- 5-65. VERTICAL OUTPUT CHECK.
- a. Connect Oscilloscope to VERT OUTPUT on 851 rear panel.
  - b. Check that signal is displayed on Oscilloscope.
- 5-66. HORIZONTAL OUTPUT CHECK.
- a. Connect Oscilloscope to HORIZ OUTPUT on 851 rear panel.
  - b. Check that signal is displayed on Oscilloscope.
- 5-67. LINE SYNC CHECK.
  - a. On 851, set SYNC to LINE.
- b. Set Oscilloscope input for dc coupling, sync on line.
- c. Connect 851 SWEEP OUTPUT to Oscilloscope. Display of sweep signal should remain in synchronization.
- 5-68. BASE LINE CLIPPER CHECK. Rotate BASE LINE CLIPPER full cw; trace on at least lower 2 cm of 851 CRT should blank.

#### Note

At high INTENSITY levels, it is normal for trace to defocus slightly when BASE LINE CLIPPER is set cw.

## 5-69. CRT CHECKS.

## 5-70. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
8*	Low-Frequency Oscillator (202C)	1
9*	Signal Generator (606A)	1
C**	Shielded cable, banana plug to BNC male (11001A)	2
D**	Shielded cable, BNC male to BNC male (10503A)	1
*	Table 5-1 **Table 5-2	


- 5-71. MEASUREMENT SETUP. Make connections between Low-Frequency Oscillator, Signal Generator, and 851 as indicated in Figure 5-3.
- 5-72. ALIGNMENT. Before starting to check the CRT, make sure horizontal trace is parallel to horizontal axis of graticule; if not, readjust TRACE ALIGN.

### 5-73. PATTERN DISTORTION AND RESOLUTION.

- a. 100% modulate Signal Generator at 20 Kc using Low-Frequency Oscillator as modulating voltage source.
- b. Set  $851\ I.F.GAIN$  for 6 cm of vertical deflection on  $851\ CRT.$

Section V Paragraphs 5-74 to 5-81





BARRELLIN

Figure 5-4. Pin-cushioning and Barrelling Defined

- c. Check pattern for excessive barrelling or pincushioning (see Figure 5-4); if present, adjust A7R22, Pattern Shape Adj on Vert Ampl Bd A7 (see Figure 5-1), for best compromise (minimum average distortion of horizontal and vertical edges of pattern).
- d. Decrease Low-Frequency Oscillator output frequency to 1 Kc; at normal intensity, focus should be uniform throughout the 6 x 10 cm screen area.

### 5-74. BLANKING.

- a. Set INTENSITY full cw.
- b. Observe trace on all sweep speeds. No retrace should be seen.

# 5-75. VERTICAL AMPLIFIER CHECKS AND ADJUSTMENTS.

## 5-76. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
3*	Variable Transformer, set for 115V	1
2*	Voltmeter with automatic range finder (3440A & 3442A)	1
9*	Signal Generator (606A)	1
8*	Low-Frequency Oscillator (202C)	1
B**	Shielded cable, dual banana plug to alligator clips (11037A)	1
D**	Shielded cable, BNC male to BNC male (10503A)	1
E**	Adapter, BNC female to dual banana plug (10111A)	1
A**	Shielded cable term. w/dual banana plugs (11000A)	
*7	Table 5-1 **Table 5-2	

## 5-77. VERTICAL CALIBRATION.

## 5-78. VERTICAL POSITION.

- a. Rotate VERT POS adjust cw; trace should move upward.
- b. With no input, align trace with base line of graticule.

## 5-79. CALIBRATION SETUP.

a. Turn instrument so it is resting on top plate. Remove bottom plate.

b. Measurement setup is shown in Figure 5-5. Connect Voltmeter (digital with automatic range finding capability) at feed-thru terminal (Video Out) at output of RF Circuit Assembly A2. This terminal projects through the casting that encloses Assembly A2, and is identified in Figure 5-6.

## 5-80. CALIBRATION PROCEDURE.

- a. Set Signal Generator: RANGE, FREQUENCY... for 20-Mc output ATTENUATOR, VERNIER. for 4.0 ±0.1 vdc detected input to Vertical Amplifier (as read on Digital Voltmeter)
- b. Adjust Vert Gain Adj A7R15 (see Figure 5-5) for 6.0 cm vertical deflection on CRT.
- c. Disconnect Signal Generator; trace should return to graticule base line.
  - If trace does not return to base line, again adjust VERT POS to align trace with base line.
  - (2) Then again perform calibration procedure (steps a, b, and c).
  - (3) Continue until requirements of both steps b and c are met.

#### Note

Since VERT POS and Vert Gain interact, it may be necessary to repeat adjustments several times.

#### 5-81. VIDEO BANDWIDTH.

- a. Disconnect coaxial lead at feed-thru terminal coming out of casting which houses RF Circuit Assembly A2. Connect Signal Generator to this coaxial lead (this lead is Vertical Amplifier Input Cable A2W6 and is identified in Figure 5-6).
  - b. Set 851 I.F.BANDWIDTH to 1 Mc.
  - c. Set Signal Generator:
    RANGE, FREQUENCY. . . for 50 Kc output
    ATTENUATOR . . . . . . . . 3 VOLT range
- d. Adjust Signal Generator output (use VERNIER on hp 606A) for 7.0 cm vertical deflection on 851 CRT, and note reading of Signal Generator output meter.
- e. Increase frequency of Signal Generator output to 1.2 Mc, and adjust level of output to obtain same meter reading as was noted in step d.
- f. Observe vertical deflection on 851 CRT. Vertical deflection should exceed 5.0 cm.
- g. Replace Signal Generator with Low-Frequency Oscillator such as the hp 202C. Monitor amplitude of Oscillator output with an Oscilloscope.
  - h. Set 851 I. F. BANDWIDTH at 100KC.
  - (1) Set Oscillator for 1 Kc output, and adjust output level (use AMPLITUDE on 202C) for 7.0 cm of vertical deflection on 851 CRT. On monitoring Oscilloscope, note amplitude at which 7.0 cm deflection was obtained.
  - (2) Increase Oscillator output frequency until 851 CRT vertical deflection decreases to 5.0 cm.




Figure 5-5. Measurement Setup, Calibration of 851 Vertical Amplifier

While increases: Oscillator frequency, adjust output level is essessary to maintain signal amplitude same as noted on Oscilloscope in step (1).

- (3) The vertical deflection decrease to 5.0 cm should take place at 200 ±40 Kc.
- i. Using a 1 Kc signal while obtaining a 7.0 cm vertical deflection, check other positions of I. F. BAND-WIDTH using procedure given in step h. Frequency at which deflection should decrease to 5.0 cm is given in Frequency column of Table 5-18.

Table 5-18. Data for Video Bandwidth Check

I. F. BANDWII Setting	отн	Frequency
100 KC		200 ±40 Kc
10 KC		40 ± 8 Kc
3 KC	:	12 ±2.4 Kc
1 KC		4 ±0.8 Kc

j. Disconnect Oscillator from Vertical Amplifier Input Cable A2W6, and reconnect Cable to RF Circuit Assembly A2 feed-thru terminal.

# 5-82. 1MC 1.F.BANDWIDTH ALIGNMENT AND CHECK.

82.1

44

## 5-83. FULLPMENT REQUIRED

Ref No.	Equipment		No.
7*	er (5245L)	·	1
9*	Signal Generator (606A)		1
D**	Shielded cable term w BN males (10503A)	(C	1
F**	BNC tee, male and 2 femal (UG-274A U)	les	1
G**	Plastic tuning wand		1
	Adapter, male type N to fe BNC (UG-201A U) hp 125	male 0-0067	1
	*Table 5-1 **Table 5	5 - 2	

#### 5-84. SETUP AND INITIAL SETTINGS

a. Setup. See Figure 5-7 Casting while RF Circuit Assembly A2 is on bottom of 8 that which houses Bandpass Filter Assemin top of 851, on right side near front panto both is required in this procedure, local Assemblies and adjustments are dalled ou

365

6 - 5

đ




Figure 5-6. Setup for 851 Video Bandwidth Measurements

5-7 and 5-8. With top and bottom covers removed, rest 851 on its right side. Access holes in cover plates of Assembly castings are provided for adjustments called for in alignment procedure; holes are covered by removable plug-in buttons.

b	Settings	s:											
	SYNC.	<i>.</i> .										. 1	NT
	SWEEP	TIME					3	ΜI	LI	IS	Œ	C/	СМ
	I. F. BA	NDWID	TF	ł								1	MC
	VERT I												

c. <u>Initial Procedure</u>. Set Signal Generator for 20 Mc at -10 dbm.

## 5-85. 1 MC ALIGNMENT.

a. With modified GC plastic Tuning Wand, adjust 851 Detector Tune T1 (A2A7T1) for maximum deflection on CRT.

## Note

Two peaks are present; adjust for maximum deflection of the highest.

b. Tune L11 Imped Adj (A2A6L11) and A12 adjustments C1, C3, L2 for maximum vertical deflection.

### Note

Adjustment of A12C1, A12C3, and A12L2 presets them; final adjustment of A12C1, A12C3, and A12L2 is made with 851 connected to 8551 (see Paragraph 5-93b).

## 5-86. 1MC BANDWIDTH CHECK.

a. Set Signal Generator output level for 7.0 cm 851 display.

- b. While watching 851 display, decrease frequency at Signal Generator until 851 display amplitude is 5.0 cm. Note Counter reading.
- c. Still watching 851 display, increase frequency and observe display go through maximum and return to 5.0 cm. Note Counter reading.

Frequency difference between the two readings should be within 1.4 and 2.2 Mc.

# 5-87. 100KC I.F.BANDWIDTH ALIGNMENT AND CHECK.

#### 5-88. SETUP AND INITIAL SETTINGS.

- a. Use setup indicated in Figure 5-7. 100KC Bandpass Filter Assemblies A3 and A5 are located toward rear of 851 on right side, beneath Low-Voltage Power Supply A9; access to adjustments A3L1 and A5L1 is through holes in the A9 Board; see Figure 5-36.
- b. Use same initial procedure as given in Paragraph 5-84, subparagraph c, changing control settings as follows:

SPECTRUM WIDTH				1	00	KC/CM
I.F.BANDWIDTH.						100KC

5-89. 100KC ALIGNMENT. Tune L1 Imped Adj (A2A2L1 in RF Circuit Assembly A2) and 100KC Bandpass Filter adjustments A3L1 and A5L1 for maximum vertical deflection.

## 5-90. 100KC BANDWIDTH CHECK.

- a. Set Signal Generator output level for 7.0 cm 851 display.
- b. While watching 851 display, decrease frequency at Signal Generator until 851 display amplitude is 5.0 cm. Note Counter reading.
- c. Still watching 851 display, increase frequency and observe display go through maximum and return to  $5.0\ cm$ . Note Counter reading.

Frequency difference between the two readings should be within 80 and 120 Kc.

# 5-91. FINAL 1MC AND 100KC BANDWIDTH ADJUSTMENTS.

5-92. SETUP. Connect the 851 to the 8551; see Figures 2-1 and 5-8.

Set 8551	:														
LINE .												SI	'A	ND	BY
SIGNAL	IDE	NTI	FIE	ER										0	FF
SPECTR	NUS	WII	TI	ΙV	E.	RN	ΙI	EF	S					_	AL
FREQUE	ENC	Y (G	C)			•								. 0	1 - 2
I.F															GC.
FREQUE	ENC	Y T	UN]	IN(	3							. (	CO	ΑF	RSE
ATTENU	JAT(	OR (	(DE	3)					•						60
Set 851:															
SYNC.														. 1	INT
SWEEP.		F	•		•	•		•	3	ΜI	$\mathbf{L}$	LIS	SE	c/	CM
I. F. BAN															MC
VERT D														. 1	LIN




Figure 5-7. Measurement Setup, 1MC and 100KC I. F. Bandwidth Alignment and Checks

Harman harmalitetter, miler is before tertable a




Figure 5-8. Measurement Setup, 10KC I. F. Bandwidth Alignment and Checks

- b. Set Signal Generator for some frequency above 10 Mc, such as 50 Mc, and at some level less than 1 watt.
- c. Perform steps 1 through 14 of Initial Operating Procedure, Figure 3-3, using settings given in stepa.

## 5-93. PROCEDURE.

- a. Set SPECTRUM WIDTH . . . . . 1 MC/CM VERT DISPLAY . . . . . . . . LOG
- b. Check symmetry, and if necessary readjust A12C1, A12C3, and A12L2 for best symmetry.
  - c. Set SPECTRUM WIDTH to 100KC/CM.
- d. Check symmetry, and if necessary readjust A12C1, A12C3, and A12L2 in Bandpass Filter Assembly A12 for best symmetry.
- e. Switch VERT DISPLAY back to LIN, and if necessary readjust A2A2L1, A3L1, and A3L5 for best compromise on symmetry and amplitude while keeping bandwidth within  $100 \pm 20$  Kc.
- f. Recheck at LOG, then LIN, readjusting as required.
- g. Return VERT DISPLAY to LIN, set I.F.BAND-WIDTH to 1 MC, and SPECTRUM WIDTH to 1 MC/CM.

- h. Check display for maximum amplitude, correct bandwidth, and symmetry, readjusting Impedance Adj A2A6L11 if required.
- i. Set VERT DISPLAY at LOG and, if required, readjust A12C1, A12C3, and A12L2 for best symmetry and amplitude.
- j. Set I. F. BANDWIDTH to 100KC and SPECTRUM WIDTH to 100 KC/CM, and again check the 100 Kc filter bandpass characteristics.
- k. Continue readjusting as required to obtain the best compromise on amplitude and symmetry with VERT DISPLAY at LIN and LOG for both 1 Mc and 100 Kc filters while keeping respective bandwidths within specifications.

# 5-94. 10KC, 3KC, AND 1KC I.F.BANDWIDTH ALIGNMENT AND CHECKS.

5-95. Signals for the three narrower bandwidths (1, 3, and 10 Kc) pass through four tuned filters. The coil (A2A3L1, A2A3L2, A2A4L1, A2A4L2) in each filter is tapped; change in bandwidth is obtained by using different taps for each bandwidth. The same four filters, however, are used for the three bandwidths, and therefore accurate adjustment of the 10-Kc bandwidth should bring the 3-Kc and 1-Kc bandwidths within specifications. After bandwidth is set with I. F. BANDWIDTH at 10KC, bandwidth is checked

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554

Tel:- 01844-351694 Fax:- 01844-35255 Email:- enquiries@mauritron.co.uk Section V Paragraphs 5-96 to 5-101

at the 3KC and 1KC settings. I. F. bandwidth alignment is not a simple technique. While tuning for correct I. F. bandwidth, remember:

- a. Ideally, all adjustments should be made simultaneously. Since this is impossible, it will be necessary to repeat adjustments more than once to obtain best tuning of the four filters.
- b. Final adjustment should be the compromise which obtains best characteristics for all four filters.

### 5-96. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.
10*	VHF Attenuator (355D)	1
D**	Shielded coax cable term. w/ BNC males (10503A)	1
G**	GC plastic tuning wand	1
K**	Screwholding screwdriver	1
*	Table 5-1 **Table 5-2	

#### 5-97. MEASUREMENT SETUP.

- 5-98. Use the 8551 RF Section as the signal source for the narrower I.F. bandwidth alignment procedures. See Figure 5-8.
- a. Connect the 8551 to 851 as shown in Figure 2-1, except for the 200MC cable.
- b. Insert a VHF Attenuator, such as the 355D, between 200MC OUTPUT and INPUT on the rear of the 8551. The cable which straps 200MC OUTPUT to INPUT can be used for the connection to one 355D port; use a coaxial cable such as the 10503A to make connection between the other 355D port and 200MC INPUT.

## 5-99. 10KC PROCEDURE.

#### 5-100. INITIAL.

8551

- a. Set the 355D to 30 db.
- b. Perform steps 1 through 6 of Turn-On Procedure, Figure 3-3, using the following settings:

<u> </u>															
LINE .								•	•	•		ST.	ľΑ.	ND.	ΒY
SIGNAL	<b>DEN</b>	TIE	ľΕ	R										.0	FF
SPECTR	UM V	VID'	TH									1	M	C/C	CM
SPECTR	UM V	VID'	ΤН	V	EI	RN	Œ	R						C.	AL
FREQUE	NCY	(G	C)											. 01	1-2
TUNE.															ale
FREQUE										•					NE
I TELL OF CE		- 0	- 1 -	- ' -	•	•	•	٠	٠	٠	•	•	•		
851															
BASE LI	NE C	LIE	PPI	ΞR								• 1	ma	хс	cw
SYNC .														. I	NT
I. F. BAN	DWI	DTH	Ι.										1	100	KC
VERT D															
SWEEP'	TIME	١.							3	M	L	LI	SE	C/C	СМ
VERNIE	R.													,	
INTENSI															
I. F. GAI															
I.F.VEF	CNIE	₹.	•	•	n	na:	X C	C	V I	m	ını	mı	ım	ga	un)

- c. Set LINE to ON. After a warmup of about five minutes, find and center the 2-Gc feed-through signal. Bringing the signal onto the display probably will require some readjustment of TUNE and insertion of more attenuation; use the 355D -- a setting of 70 db is typical.
- d. After finding and centering the signal, stabilize the Analyzer (see Figure 3-5, 851A/8551A Manual).
- e. Check alignment of the base-line trace with the horizontal axis. If necessary, adjust VERT POS and TRACE ALIGN to bring base-line trace exactly parallel with and on the graticule base line.

#### 5-101. 10KC ALIGNMENT.

- a. Set
  SPECTRUM WIDTH . . . . . . 10KC/CM
  I.F.BANDWIDTH . . . . . . . . 10KC
  SWEEP TIME . . . . . . . 3 MILLISEC/CM
- b. Adjust TUNE to center signal on 851 display. Adjust 851 I. F. GAIN and I. F. VERNIER for a vertical deflection of exactly 7.0 cm.
- c. Bandwidth tuning adjustments are inside the RF Circuit Assembly casting, and location of adjustments is marked on the cover. Access holes, covered with removable plug-in buttons, are provided in the casting cover. However, though it can be done, presetting (step d) two of the capacitors is a little difficult with the casting cover on. If you prefer to remove the casting cover, do so at this point in the procedure. Casting cover is held on by 26 screws with integral washers; use of a screwholding screwdriver is helpful. Figure 5-22 shows the boards in RF Circuit Assembly A2.
- d. Preset BALANCE ADJ capacitors A2A3C5 and A2A4C8 to approximately 1/3 mesh.
- e. Adjust 1-10KC BANDWIDTH ADJ capacitors A2A3C4, A2A3C2, A2A4C5, and A2A4C9 for maximum bandwidth.

#### Note

In tuning capacitor A2A3C4, A2A3C2, or A2A4C5 through its tuning range it will be found there are two points which give vertical deflection peaks. Since there is little difference between the amplitude of the two peaks, it is difficult to distinguish which is the spurious tuning region. However, correct I. F. bandwidth tuning can be obtained only when adjustment of each capacitor is made in the true tuning region. Also, maximum bandwidth usually is obtained by tuning off the peak slightly.

- f. Adjust IMPED ADJ A2A3L3 and FREQUENCY ADJ A2A4C7 for maximum vertical deflection.
- g. Center display with TUNE, and set vertical deflection to 7.0 cm with I. F. GAIN. If display is not 1 cm wide at the 5-cm (half-power) points, again perform steps e through g until a 1-cm bandwidth at the 5-cm points is obtained.

02148-1 5-29

Model 851B

## 5-102. 3KC AND 1KC BANDWIDTH CHECKS.

5-103. CALIBRATION. To read bandwidth at 3KC and 1KC, it is necessary to increase resolution to 1 Kc/cm by calibrating the VERNIER; procedure follows:

and the second of the second of the second second

- b. At the 355D, adjust for vertical display of 4 to 5 cm.
- c. With 0-to-10 I. F. GAIN control, adjust for vertical display of between 6 and 7 cm.
- d. Use I.F. VERNIER to bring signal amplitude to exactly 7.0 cm.
  - e. Center signal with TUNE.
- f. With SPECTRUM WIDTH VERNIER, expand signal until it is 10 cm at the 5.0-cm axis. Since signal width at the 5.0-cm axis was initially adjusted to 1 cm, SPECTRUM WIDTH was set for 10 KC/CM, and the 10-Kc display has been expanded to 10 cm, display scale at this setting of SPECTRUM WIDTH VERNIER is 1 Kc/cm.

### 5-104. BANDWIDTH CHECKS.

- a. Set I. F. BANDWIDTH to 3KC.
- b. Adjust I. F. GAIN and I. F. VERNIER for vertical deflection of exactly  $7.0\ \text{cm}$ .
- c. Width of display at 5.0-cm axis should be between 2.4 and 3.6 cm.
  - d. Set I. F. BANDWIDTH to 1KC.
- e. Adjust I. F. GAIN and I. F. VERNIER for vertical deflection of exactly 7.0 cm.
- f. Width of display at 5.0-cm axis should be between 0.8 and 1.2 cm.

#### Note

If 1-Kc bandwidth appears to be too wide, recheck tuning of FREQ ADJ A2A4C7 (see Paragraph 5-101f).

## 5-105. FINAL 1-10KC BANDWIDTH ADJUSTMENT.

5-106. Set I.F.BANDWIDTH at 10KC, and recheck bandwidth making adjustments if necessary. When maximum vertical deflection is 7.0 cm, at 5.0 cm bandwidth should be within 0.8 and 1.2 cm.

## 5-107. AUTO SELECT CHECK.

5-108. With I.F.BANDWIDTH at AUTO SELECT, the Analyzer automatically selects the I.F. bandwidth which provides optimum operation for whatever combination of 8551 SPECTRUM WIDTH and 851 SWEEP TIME settings is selected.

5-109. Connections to the filters which determine I.F. bandwidth are made through relays. With I.F.BAND-WIDTH at 1KC, 3KC, 100KC, or 1 MC, dc to operate

the relays is applied via contacts on the I. F. BAND-WIDTH switch. With I. F. BANDWIDTH at AUTO SELECT, however, dc to operate the relays is applied via contacts on the 8551 SPECTRUM WIDTH switch and the 851 SWEEP TIME switch. Inter-unit connections required for automatic selection of I. F. bandwidth are carried in the CONTROL cable.

- 5-110. To check that the AUTO SELECT feature is functioning, I.F. BANDWIDTH, SPECTRUM WIDTH, and SWEEP TIME are given the settings known to result in optimum operation, the display is noted, then I.F. BANDWIDTH is set to AUTO SELECT, and the resulting display is compared to the preceding display. To perform this check:
- a. Connect (see Figure 2-1) the 851 to an 8551 known to be in adjustment.
- b. Perform the initial operating procedure, Figure 3-3, using an input signal of less than a watt, 10 Mc or higher in frequency.
  - c. Set
    SPECTRUM WIDTH . . . . . . 10 KC/CM
    SWEEP TIME . . . . . 10 MILLISEC/CM
    I. F. BANDWIDTH . . . . . . . . . . . . 1 KC
  - d. Note display.
  - e. Switch I. F. BANDWIDTH to AUTO SELECT.
- f. Note display; it should be same as display noted in step d.
- g. Follow same procedure for all settings shown in Table 5-19, switching to AUTO SELECT after each change of switch settings.

Table 5-19. Switch Settings for AUTO SELECT Check

AUTO SELECT Check								
I.F. BANDWIDTH	SPECTRUM WIDTH	SWEEP TIME						
1 KC	10 KC/CM 30 KC/CM 100 KC/CM 300 KC/CM 1 MC/CM	10 MILLISEC/CM 30 MILLISEC/CM .1 SEC/CM .3 SEC/CM 1 SEC/CM						
3 KC	10 MC/CM 3 MC/CM 1 MC/CM 300 KC/CM 100 KC/CM 30 KC/CM	1 SEC/CM .3 SEC/CM .1 SEC/CM 30 MILLISEC/CM 10 MILLISEC/CM 3 MILLISEC/CM						
10 KC	300 KC/CM 1 MC/CM 3 MC/CM 10 MC/CM 30 MC/CM 100 MC/CM	3 MILLISEC/CM 10 MILLISEC/CM 30 MILLISEC/CM .1 SEC/CM .3 SEC/CM 1 SEC/CM						
100 KC	200 MC/CM 100 MC/CM 30 MC/CM 10 MC/CM 3 MC/CM	.3 SEC/CM .1 SEC/CM 30 MILLISEC/CM 10 MILLISEC/CM 3 MILLISEC/CM						

### 5-111. I.F.SENSITIVITY CHECK.

## 5-112. EQUIPMENT REQUIRED.

Ref No.	Equipment	No.						
3*	Variable Transformer, set for 115V 1							
9*	Signal Generator (606A)							
D**								
*	Table 5-1 **Table 5-2							

5-113. MEASUREMENT SETUP. Connect Signal Generator to I. F. INPUT (on 851 rear) and line voltage through a variable transformer set for 115V. Set 851:

#### 5-114. SIGNAL LEVEL CHECK.

- a. Set Signal Generator for 20 Mc output.
- b. Set I. F. BANDWIDTH to 1 Mc, and adjust Signal Generator output level to obtain 6 cm of vertical deflection on 851 CRT. The 6-cm deflection should be obtained with Signal Generator output level at -57.5 dbm ±4.5 dbm.
- c. Adjust Signal Generator output level to obtain 6 cm of vertical deflection on 851 CRT at each setting of I.F.BANDWIDTH. Level at which 6 cm deflection should be obtained at each I.F.BANDWIDTH setting is given in Table 5-20.
  - d. Disconnect Signal Generator.

Table 5-20. Data for I.F. Sensitivity Check

nal Level Limits*
to -53 dbm to -60 dbm to -80 dbm to -80 dbm to -71 dbm
;

\* For 6 cm deflection with I. F. GAIN at 80 db and VERNIER full counterclockwise

5-115. NOISE LEVEL CHECK. With no signal connected to I.F. INPUT, switch I.F. BANDWIDTH through all positions. The noise displayed on the CRT should not exceed 0.45 cm at any setting of I.F. BANDWIDTH.

# 5-116. VERT DISPLAY CHECKS AND ADJUSTMENTS.

### 5-117. PRELIMINARY CHECK.

a. Connect to line voltage through Variable Transformer, set for 115V.

b.	Set 851	VERT DISPLAY					LOG
		I.F.GAIN (DB)					20
		I. F. VERNIER					.full cw
		I F BANDWIDT	Н				. 1 MC

c. With no signal input, check that trace aligns with graticule base line.

## 5-118. LOG DISPLAY.

#### 5-119. ADJUSTMENT.

- a. Connect Signal Generator to I. F. INPUT; set for 20 Mc output at level which obtains 1 cm of vertical deflection on 851 CRT.
- b. On 851, increase I.F.GAIN to 40 db, and adjust A11R13 (on VERT DISPLAY switch, see Figure 5-25) for 3 cm of vertical deflection.
- c. Increase I.F. GAIN to 60 db, and adjust A11R14 (Figure 5-25) for 5 cm of vertical deflection.
- d. Increase I. F. GAIN to 80 db (outer control at 70, inner at 10), and adjust potentiometer A11R20 (Figure 5-25) for 7 cm of vertical deflection.
- e. Decrease I. F. GAIN 1 db, and note deflection level. Reset I. F. GAIN to 80 db (70 + 10), and rotate I. F. GAIN VERNIER fully counterclockwise.

Deflection level decrease should exceed 1 db.

f. Reset I.F. GAIN VERNIER fully clockwise.

## 5-120. LINEARITY CHECK.

a. Decrease I.F.GAIN in steps of 10 db, and observe trace.

Each step should lower trace 1.0 cm on CRT and, at each 10-db step, alignment between trace and horizontal line on graticule should be within  $\pm 0.2$  cm.

b. Reset I. F. GAIN to 70 - 10, and repeat step a for all other I. F. bandwidths.

#### 5-121. SQ DISPLAY.

5-122. ADJUSTMENT.

a.	On 851 set:					
	VERT DISPLAY					SQ
	I. F. GAIN (DB)				20	+ 10
	(outer control at 20,					

- b. Set Signal Generator signal level to give 7.0 cm of vertical deflection on 851 CRT.
- c. Decrease I.F.GAIN 6 db, and adjust .710 SQ CALIB A11R2 (Figure 5-25) for 1.75 cm of vertical deflection.
- d. Increase signal level for 7.0 cm vertical deflection; decrease I.F. GAIN 6 db in 3-db steps. See Table 5-21 for vertical deflection limits.

Table 5-21. SQ Display Linearity Check Data

I. F. GAIN (DB) Settings	Step	Vertical Deflection (cm)
30	Ref	7.0
27	-3 db	3.15 - 3.85
24	-6 db	1.40 - 2.10

e. Perform step d at all other I.F. bandwidths.

## 5-123. LIN DISPLAY LINEARITY CHECK.

a.	Set 851								
	VERT DISPLAY							1	LIN
ï	I. F. GAIN (DB)					30	(20	+	10)

and the control of th

- b. Increase signal level for 7.0 cm of vertical deflection.
- c. Decrease I.F.GAIN 12 db in 6-db steps; what the vertical deflection should be at each step is shown in Table 5-22.

Table 5-22. LIN Display Linearity Check Data

I. F. GAIN (DB) Settings	Step	Vertical Deflection (cm)
30	Ref	7.0
24	-6 db	3.29 - 3.71
18	-12 db	1.54 - 1.96

d. Repeat steps b and c at all other I. F. bandwidths.

## 5-124. FINAL I.F.BANDWIDTH ADJUSTMENTS.

## 5-125. EQUIPMENT REQUIRED.

Ref No.	Equipment No.							
11*	UHF Signal Generator (8614A) 8551 (adjusted)	1 1						
H**	* 3-ft shielded coax cable term. w/ 1 type N males (11500A)							
G**	G** GC plastic tuning wand 1							
*Table 5-1								

# 5-126. MEASUREMENT SETUP AND INITIAL PROCEDURE.

a. Setup. Remove top and bottom covers from 851 and top cover from 8551; connect as indicated in Figure 5-9.

# b. <u>Settings:</u> (1) 8551:

\-/																
	LINE.												.5	ЗT	ANI	BY
	SPECT	rRU	M V	VIC	<b>T</b>	H	V	ΞR	N.	ŒF	3				. (	AL
	SPECT	rru	ΜV	VΙΣ	T:	Η							30	0	KC/	'CM
	FREQ	UEN	ICY	ΤĮ	JN	IN	ľG			C	O.A	ιR	SE	: c	or F	INE
	FREQ	UEN	1CY	(G	C)										.0	1-2
	I.F	•													2	GC
	TUNE														1.8	GC
	FREQU	JEN	ICY	ID	Εľ	ľΤ	IF	Έ	R						. (	FF
(2)	851:															
\-/	SWEET	P T	IME							30	M	Ш	<sub>L</sub>	IS	EC/	СМ
	SWEEL	? T	ME	VI	ΞR	N.	Œ	R							´c	AL
	I. F. BA	IND	WII	TI	I										10	KC
	VERT	DIS	PLA	Y											. I	OG
	I. F. GA	ΙN	(DB	)												80
	I. F. GA	III	VE	RN:	Œ.	R									full	cw

## c. Initial Procedure:

(3) Signal Generator:

(1) Follow steps 1-14 of initial operating procedure, Figure 3-3, using above settings.

Output level. . . . . . . . . . . 0 dbm

(2) Switch FREQUENCY TUNING to STABILIZE, and stabilize Analyzer (see Figure 3-5, 8551A Manual).

## 5-127. CRYSTAL FILTER (10KC, 3KC, 1KC) BALANCE.

- a. Adjust Signal Generator output level at 60 Mc for 7 cm of vertical deflection on  $851\ \text{CRT}.$
- b. Tune BALANCE ADJ A2A3C5 and A2A4C8 for a symmetrical display. Capacitor locations are marked on RF Circuit Assembly A2 cover plate, and capacitors can be tuned through access holes in cover plate.
- c. Set I. F. BANDWIDTH and SPECTRUM WIDTH as shown in the following table, and readjust A2A3C5




Figure 5-9. Measurement Setup, Final I.F. Bandwidth Adjustments

and A2A4C8 for best symmetry compromise for the three bandwidths.

851	8551
I. F. BANDWIDTH	SPECTRUM WIDTH
3KC	100 KC/CM
1KC	30 KC/CM

## 5-128. 1MC BANDPASS FILTER ADJUSTMENTS.

5-129. ADJUSTMENT LOCATIONS. Location of Bandpass Filter Assembly A12 in top of 851 is shown in Figure 5-8, and that of Converter Assembly A9 in top of 8551 is indicated in Figure 5-9. Both assemblies are housed in castings, and each provides access to adjustments through holes in top cover plate.

#### 5-130. PROCEDURE.

- a. Set:
  - SPECTRUM WIDTH . . . . . . 1 MC/CM VERT DISPLAY . . . . . . . . . LIN
- b. Set I.F.BANDWIDTH to 100 KC/CM, and center display with TUNE.
- c. Adjust A12C1 and A12C2 in the 851 and A9A2L2 in the 8551 for best symmetry and a 1-Mc bandwidth.

#### Note

A12C1 and A12C2 control bandwidth, and A9A2L2 mainly the frequency at which maximum amplitude occurs.

d. Switch I. F. BANDWIDTH to 100 KC/CM, and check that frequency at which maximum amplitude occurs does not shift. Readjust A9A2L2, as required, to assure maximum amplitude occurs at same frequency with I. F. BANDWIDTH in 100KC and 1MC positions.

## 5-131. TROUBLESHOOTING.

## 5-132. LOCALIZATION.

- 5-133. First use the in-the-cabinet performance checks, Tables 5-6 and 5-7, to localize trouble.
- a. If these checks localize trouble to a particular part of the instrument, first make a visual check for broken leads, overheated resistors, or cold solder joints before making an electrical check. If this inspection yields no information, first check the LV Power Supply (see Paragraph 5-40), and then the part of the circuit that appears to be in trouble (see the appropriate part of Table 5-13).
- b. If the performance checks fail to localize the trouble, check the instrument using procedures given in Table 5-13.

## 5-134. PARTS LOCATION.

- 5-135. The key to parts locations is in the part designation.
- a. If a component is mounted on an Assembly board, the designation is prefixed with the Assembly number, e.g., A1R5. Location of each Assembly is called out in Figures 5-10 and 5-11. In addition, a picture of

each Assembly board is provided, and all components on the Board are identified. For the most part, Board pictures face the schematic in which the Assembly appears. All Board pictures are listed in the List of Illustrations.

b. If a component is mounted on the chassis, the designation has no prefix, e.g., R5. Since these parts are harder to locate, a Locater list, Table 5-5, has been prepared which gives information on how to locate the part.

# 5-136. ISOLATING TROUBLE IN TRANSISTOR CIRCUITS.

5-137. For general data on transistors, see Paragraph 5-163 and Figure 5-16.

## 5-138. IN-CIRCUIT TESTING.

- a. When checking a transistor stage, first determine if the emitter-base junction is forward-biased. Do not place an electronic voltmeter directly across the junction to measure the voltage difference; there could be sufficient loop current between the voltmeter leads to damage the transistor. Instead, measure each voltage separately with respect to a common point (e.g., chassis). If junction is not forward-biased, and power supply voltages are known to be correct, the base-emitter junction may be open (see Paragraph 5-139).
- b. If the emitter-base junction is forward-biased, check for amplifier action by short-circuiting base to emitter while observing collector voltage. The short eliminates base-emitter bias and should cause the transistor to stop conducting. Collector voltage should then shift to near the supply voltage. Any difference is due to leakage current through the transistor and, in general, the smaller this current, the better the transistor. If collector voltage does not change, the transistor either has an emitter-collector short circuit or emitter-base open circuit.
- 5-139. OUT-OF-CIRCUIT TESTING WITH OHM-METER. If a short or open circuit is suspected, remove the transistor from the circuit (see Paragraph 5-145) and use an ohmmeter to measure internal resistance. See Table 5-23 for typical measurement data.

## CAUTION

Most ohmmeters can supply enough current or voltage to damage a transistor. Before using the ohmmeter, check ohmmeter open-circuit voltage and short-circuit current output ON THE RANGE TO BE USED. Open-circuit voltage must not exceed 1.5 volts and short-circuit current must be less than 3 ma. See Table 5-23A for safe resistance ranges for some common ohmmeters.

# 5-140. IN-CIRCUIT TESTING OF TRANSISTORS $\overline{Q3}$ , $\overline{Q4}$ , $\overline{Q5}$ , $\overline{Q6}$ .

5-141. To check base-emitter junction of transistors Q3, Q4, Q5, or Q6, connect Voltmeter as noted in Table 5-23B. Any sensitive high-impedance voltmeter,

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:-01844-351694 Fax:-01844-352554
Email:-enquiries@mauritron.co.uk




Figure 5-10. 851B Spectrum Analyzer Display Section, Top View, Top Cover Removed




Figure 5-11. 851B Spectrum Analyzer Display Section, Bottom View, Bottom Cover Removed

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tet: 01844-351694 Fax:- 01844-352554

Model 851B

Table 5-23. Typical Data for Out-of-Circuit

Transistor Resistance Measurement

		Connect Ohmmeter		Measure	
1	sistor pe	Pos lead to	Neg lead to	Resistance (ohms)	
	Small Signal	emitter	base*	200-500	
	PNP   Signal   emitter   collector	10K-100K			
manium	Power	emitter	base*	30-50	
		emitter	collector	several hundred	
	Small	base emitter	1K-3K		
NPN Silicon	Signal	∞llector	emitter	very high (might read open)	
	Power	base	emitter	200-1000	
		collector	emitter	high, often greater than 1M	

<sup>\*</sup>To check collector, short collector to base; resistance should decrease.

Table 5-23A. Safe Ohmmeter Ranges for Transistor Resistance Measurements

Ohmmeter	Safe Range(s)	Open Ckt Voltage	Short Ckt Current	Color	Polarity
hp 412A	R x 1K R x 10K R x 100K R x 1M R x 10M	1.0V 1.0V 1.0V 1.0V 1.0V	1 ma 100 μa 10 μa 1 μa 0.1 μa	Red Black	+ -
hp 410C	R x 1K R x 10K R x 100K R x 1M R x 10M	1.3V 1.3V 1.3V	0.57 ma 57 μa 5.7 μa 0.5 μa 0.05 μa	Red Black	+
hp 410B	R x 100 R x 1K R x 10K R x 100K R x 1M	1.1V 1.1V 1.1V 1.1V 1.1V	1.1 ma 110 μa 11 μa 1.1 μa 0.11 μa	Black Red	+
Simpson 260	R x 100	1.5V	1 ma	Red Black	+ -
Simpson 269	R x 1K	1.5V	). 82 ma	Black Red	+
•	R x 100 R x 1K	1.5V 1.5V	32 ma 3.25 ma	Varies with Serial Number	
	R x 10 R x 100	1.5V 1.5V	750 μ <b>a</b> 75 μ <b>a</b>		

Table 5-23B. Connection Point, Q3, Q4, Q5, Q6 Base-Emitter Forward Bias Check

	,	Connect VM Between Chassis and		
Xstr	Measurement	Component	Point on Fig. 5-36	
Q3	Base to chassis	A9R1	13	
Q4	Base to chassis	A9R3	5	
	Emitter to chassis	A9R4	4	
Q5	Base to chassis	A9C7	3	
	Emitter to chassis	A9R16	20	
Q6	Base to chassis	A9R16	20	
	Emitter to chassis	A9R21	11	

such as the hp 3440A Digital Voltmeter or hp 412A Precision V-O-A is suitable.

## 5-142. REPLACEMENT OF CATHODE-RAY TUBE.

5-143. REMOVAL. It is recommended that a face mask or goggles and gloves be used when it is necessary to handle the CRT. Perform removal procedures with 851 in normal position. To reach the CRT, remove 851 top cover. Parts mentioned in the following procedure are identified in Figure 5-12 by numbered callouts. To remove:

- a. Disconnect post-accelerator lead (1).
  - (1) The post-accelerator lead connects to the tube by means of a spring-clip arrangement (2), and the connection is protected by a rubber cap (3).
  - (2) Lift edge of cap with screwdriver and, using a pair of long-nose pliers, compress spring contacts as indicated in Figure 5-12. This will free lead-and-spring assembly from recess.
- b. Disconnect the six leads (4) at the neck of the CRT. The lead pins pull straight out; be careful not to bend the pins.
- c. Remove the four screws which hold bezel to front panel; a No. 2 phillips driver is required.
  - d. Loosen clamp (5) at socket of CRT.
- e. The socket (6) is a tight fit; with a screwdriver carefully pry socket loose, and remove it.
- f. While keeping one hand on front face of CRT, carefully slide CRT forward and out of instrument.
- 5-144. INSTALLATION. Reverse removal procedure. Color-coding of leads to CRT is stamped on CRT shield. After installing new tube, perform CRT checks specified in Table 5-23D.

## 5-145. TRANSISTOR REPLACEMENT.

5-146. REMOVAL. For the most part, transistors are to be removed from the front of the circuit board. This can be done safely by a skilled technician; a pointed soldering iron is recommended. Refer to

The state of the s




Figure 5-12. Cathode-Ray Tube, Parts and Connections Involved in Removal

Table 5-23C for recommendations regarding component removal, to Paragraph 5-158 for general information about working on etched circuits, to Table 5-23D for recommendations regarding checks to be made after replacing transistors, and to Table 5-23E for recommended soldering equipment.

#### Note

Do not change an operating voltage or calibration adjustment unless it is either definitely outside specified tolerance, or calibration of a dependent function is unsatisfactory. Improving a marginal adjustment can adversely affect calibration.

5-147. CHASSIS-MOUNTED TRANSISTORS. Transistors Q1 and Q2, which drive the Step-up Transformer in HV Power Supply A8, and Q3, Q4, Q5, Q6,

Series Regulators for LV Power Supply A9, are highcurrent types which require good thermal contact with mounting surfaces for adequate heat dissipation. To assure good thermal contact for a replacement transistor, coat both sides of the black insulator with Dow Corning #5 silicone compound or equivalent before fastening the transistor to the chassis. Dow-Corning #5 compound is available in 8-oz tubes from Hewlett-Packard; order hp stock No. 8500-0059.

5-148. TRANSISTORS Q1, Q2. Location of Q1 and Q2 on the left side of the chassis is called out on Figures 5-10 and 5-11. To test these transistors it is necessary to remove the left (and top) cover plates. To replace Q1 or Q2, it is necessary to remove both the bottom and left-side cover plates. Base, emitter, and collector terminals are identified on the inner side of the deck that Q1 and Q2 are mounted on.

Section V Paragraph 5-149 For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tet: 01844-351694 Fax: 01844-352554
Email:- enquiries@mauritron.co.uk

Model 851B




Figure 5-13. Top View of 851 Display Section

- 5-149. TRANSISTORS Q3, Q4, Q5, Q6. To replace Q3, Q4, Q5, or Q6, it is necessary to gain access to the under side of the deck (see Figure 5-10) on which they are mounted. Procedure follows:
- a. Rest 851 on left side. Remove top, bottom, and right-side covers.
- b. Refer to Figure 5-13. Remove the three screws designated (1); these are 6-32 x 3/8 BH machine screws with integral lockwasher, and are accessible from the rear plate.
- c. Remove screws (2), (3), (4); these are 6-32 x 3/8 BH machine screws, each fitted with a split lockwasher.
- d. Remove the two screws designated (5); these are 8-32 x 3/8 FH machine screws fitted with integral

- lockwashers. Screw heads are accessible from right side of instrument.
- e. Remove three screws  $\bigcirc$  on the deck; these are 6-32 x 3/8 BH machine screws fitted with toothed lockwashers.
- f. Remove the  $6-32 \times 5/16$  stainless steel hex nut designated (1); a 5/16 socket wrench (Spintite) is recommended. This nut secures a  $6-32 \times 1/2$  spade lug connected to a cable clamp. After removing nut, push down on screw so it will drop out of deck hole and hole in A6 Board.
- g. On rear of A6 Board, find screw designated 8 on Figure 5-14. This screw is also a 6-32 x 1/2 spade lug connected to a cable clamp, and is secured to the A6 Board by a nut which is located under Assembly A12.

Longnose pliers can be put on the nut while loosening the screw by going in through an opening on the left side of the instrument.

h. The deck on which the transistors are mounted is now free of its fastenings and can be shifted so the under side can be exposed. One method is to pull it gently out from under the A6 Board, and then turn the deck over. Transistors and their terminals are identified in Figure 5-15.

# 5-150. REMOVING I.F.GAIN SWITCH ASSEMBLY A1.

- 5-151. To check or replace components on the I.F. GAIN switch Assembly, it is necessary to remove the switch and its shield. Proceed as follows:
- a. Rest 851 on right side, and remove bottom and left-side covers.
- b. Remove knobs; each secures to the shaft with an  $8-32 \times 3/16$  setscrew that can be loosened with a No. 8 allen wrench. Loosen locknut under knob with a 1/2" wrench; locknut is a  $3/8-32 \times 1/2$  hex nut.

c. Disconnect cable from A2J1 (on RF Circuit Assembly casting, see Figure 5-14).

and the second s

- d. Disconnect cable from A12J2 (see Figure 5-13). This is a right-angle connector and is a tight fit; if there is difficulty disconnecting it, disconnect the cable from A12J1, remove screws designated (9) which secure A12 to the side casting, and lift Assembly A12 up far enough to disconnect the cable from A12J2.
- e. Remove screw (10) (Figure 5-14) which holds Assembly A1 bracket to left-side casting; this is an  $8-32 \times 1/2$  FH machine screw with integral lockwasher.
- f. Assembly A1 is now free of its fastenings, but clearance is small. Carefully slide A1 shaft out of front panel being ready to slant A1 to the left as soon as panel-clearance permits. Maneuver A1 free of the instrument.
- g. To unfasten the shield, remove five screws (11); these are 6-32 x 1/2 BH machine screws with integral lockwasher.




Figure 5-14. Bottom View of 851 Display Section

Table 5-23C. Recommendations, Component Removal

Component p/o (Assy No.)	Access for Unsoldering	Additional Information
A1	Inside metal shield	
A2	1	Par. 5-150, Fig. 5-18
A3	Front of Boards, inside the casting	Remove 851 bottom plate, A2 casting cover; Figs. 5-11, -21, -23, -2
<b>A</b> 5		Par. 5-156, Fig. 5-26C
A6		Par. 5-156, Fig. 5-26C
AO	From front or rear of Board; remove 851 top and bottom covers	Figs. 5-10, -11, -31
A7	From front of Board; remove 851 bottom cover	Figs. 5-11, -28
A8	From front of Board; remove 851 top and left-side covers	Figs. 5-10, -34
A9	From front of Board; remove 851 top cover	Time 5 40
A10	From switch Assembly; remove 851 bottom cover	Figs. 5-10, -36 Figs. 5-11, -30; Par. 5-15
A11	From switch Assembly; remove 851 top cover	Figs 5 10 95 D
A12	Inside metal shield	Figs. 5-10, -25; Par. 5-15
Q1, 2		Par. 5-152, Fig. 5-19
Q3, 4, 5, 6	1	Par. 5-147, -148; Fig. 5-1
		Par. 5-147, -149; Fig. 5-1

Table 5-23D. Adjustments Required After Component Replacement

Component	Type/Part No.	Function	Adjustment, Par. No.
A2A3Y1	1410-0091	Xtal in 1-10KC BP Filter	1-10KC I. F. Bandwidth Align., Pars 5-82 thru 5-92, 5-128
A2A6CR1- A2A6CR6	1901-0162	Shunt diodes in Current-Controlled Atten	VERT DISPLAY Checks and Adjusts, Pars 5-116 thru 5-12;
A3, A5	00851-6028	100KC BP Filter	100KC I.F. Bandwidth Align., Pars 5-87 thru 5-92, 5-128
A7Q8 A7Q9	2N708	p/o Vertical Ampl	Vertical Calibration, Pars 5-77 thru 5-80
A11CR1 A11CR2	1901-0047	p/o SQuare shaping network	VERT DISPLAY Checks and Adjusts., Pars 5-117, 5-122
A11CR3 A11CR4	1901-0047	p/o LOG shaping network	VERT DISPLAY Checks and Adjusts., Pars 5-117, 5-119
V1	5083-0654	CRT For Service Manuals Contact MAURITRON TECHNICAL SERVICES 8 Cherry Tree Rd, Chinnor Oxon OX9 4QY Tel:- 01844-351694 Fax:- 01844-352554 Email:- enquiries@mauritron.co.uk	CRT current, voltage checks, Pars 5-45 thru 5-48 Horizontal Calib, Linearity Checks, Pars 5-50 thru 5-54e (3 CRT Checks, Pars 5-69 thru 5-74 Vertical Amplifier Checks and Adjusts., Pars 5-75 thru 5-80




Figure 5-15. Under Side of Transistor/ Transformer Deck, Identification of Q3-Q6 Terminals

# 5-152. REMOVING BANDPASS FILTER ASSEMBLY A12.

5-153. A12 components are mounted on a Board inside the casting; see Figure 5-19. To free A12, proceed as follows:

- a. Rest 851 on left side, and remove top and right-side covers.
  - b. Disconnect cable from A12J1 (Figure 5-13).
- c. Remove two screws (9) which hold A12 to the side casting; these are 6-32 x 3/8 FH machine screws with integral lockwashers.
- d. Lift A12 free of the mounting recess and turn it over so top cover can be removed.
- e. Remove four screws which hold top cover on; these are  $4-40 \times 1/4$  RH machine screws with integral lockwashers.

### 5-154. REMOVING SWITCHES.

5-155. Larger knobs secure to the shaft with an 8-32 x 3/16 setscrew which is loosened with a No. 8 allen wrench. The red verniers secure to the shaft with a 6-32 x 1/8 setscrew; loosen with a No. 6 allen wrench. Each shaft is secured to the panel with a  $3/8-32 \times 1/2$  hex nut which takes a 1/2" wrench.

# 5-156. REMOVING ASSEMBLIES A3 AND A5.

5-157. To reach 100KC Bandpass Filter Assembly A3 or A5:

- a. Remove 851 top cover.
- b. A3 and A5 are beneath the LV Power Supply A9 Board (see Figure 5-10), and some of the cabling is beneath the Transformer/Transistor Deck. The A9 Board and the Deck lift as one piece; free them as described in Paragraph 5-149.
- c. Disconnect the two cables connected to the Filter Assembly of interest, and remove the four screws which attach the Assembly to the bottom of the A2

casting. Assembly A3 is fastened with three 6-32 x 3/8 BH machine screws with integral lockwasher and one 6-32 x 1/2 BH machine screw with lockwasher (this screw also secures a 3-terminal tie point to the casting). Assembly A5 is fastened with four 6-32 x 3/8 BH machine screws with integral lockwasher.

#### 5-158. ETCHED CIRCUITS.

- 5-159. The etched circuit boards in the 851 Spectrum Analyzer Display Section are of the plated-through type consisting of metallic conductors bonded to both sides of insulating material. The metallic conductors are extended through the component mounting holes by a plating process. Soldering can be done from either side of the board with equally good results. Table 5-23E lists recommended tools and materials. Following are recommendations and precautions pertinent to etched circuit repair work.
- a. Avoid unnecessary component substitution: it can result in damage to the circuit board and/or adjacent components.
- b. Do not use a high-power soldering iron on etched circuit boards. Excessive heat may lift a conductor or damage the board.
- c. Use a suction device (Table 23E) or wooden toothpick to remove solder from component mounting holes. DO NOT USE A SHARP METAL OBJECT SUCH AS AN AWL OR TWIST DRILL FOR THIS PURPOSE. SHARP OBJECTS MAY DAMAGE THE PLATED-THROUGH CONDUCTOR.
- d. After soldering, remove excess flux from the soldered areas and apply a protective coating to prevent contamination and corrosion. See Table 5-23E for recommendations.

#### 5-160. TRANSISTOR REPLACEMENT.

- a. Do not apply excessive heat; see Table 5-23E for recommended soldering tools.
- b. Use long-nose pliers between transistor and hot soldering iron as a heat sink. The instant solder is melted, use pliers to pull lead free of Board.
- c. When installing replacement transistor, ensure sufficient lead length to dissipate soldering heat by using about the same length of exposed lead as used for original transistor.

## 5-161. COMPONENT REPLACEMENT.

a. Remove defective component from Board.

#### Note

Axial lead components, such as resistors and tubular capacitors, can be replaced without unsoldering. Clip leads near body of defective component, remove component and straighten leads left in board. Wrap leads of replacement component one turn around original leads. Solder wrapped connection, and clip off excess lead.

b. If component was unsoldered, remove solder from mounting holes with a suction desoldering aid (Table 5-23E) or wooden toothpick.

तः १९ १८ (१६) वस्ते स्वाप्त्रक्षेत्रे व्यवस्थातः । १९ विश्वपूर्णकानक्ष्यमुग्ने स्वाप्तिः । सः । अञ्चलकार्यः ।

a single transfer out year

A. TRANSISTOR BIASING					
DEVICE	SYMBOL	CUT OFF	CONDUCTING		
VACUUM TUBE	GRID	+200V -15V	+200v		
N P N TRANSISTOR	BASE EMITTER	OV- (OR-)	+20V +.3V CURRENT CONTROL CURRENT		
PNP TRANSISTOR	COLLECTOR  BASE EMITTER	-20V (OR+)	-20V3V CONTROL CURRENT		

B. AM	1PLIFIER CHAP	RACTERISTICS	
CHARACTERISTIC	COMMON BASE	COMMON Emitter	COMMON COLLECTOR
INPUT Z	30-50 Ω	500-1500 Ω	20-500ΚΩ
OUTPUT Z	300-500ΚΩ	30-50κ Ω	50-1000 Ω
VOLTAGE GAIN	500-1500	300-1000	< 1
CURRENT GAIN	< 1	25-50	25-50
POWER GAIN	20-30 db	25-40 db	10-20 db
For Service Manuals Contact MAURITRON TECHNICAL SERVICES 8 Cherry Tree Rd, Chinnor Oxon OX9 4QY Tel:-01844-351694 Fax:-01844-352554 Email:-enquiries@mauritron.co.uk	INPUT OUTPUT	OUTPUT	-15V INPUT OUTP

Figure 5-16. Transistor Biasing and Typical Amplifier Characteristics

Table 5-23E. Etched Circuit Soldering Equipment

14010 0 2021					
Item	Use	Specification	Item Recommended		
Soldering tool	Soldering Unsoldering	Wattage rating: 47-1/2 - 56-1/2 Tip Temp: 850 - 900°	Ungar #776 Handle with *Ungar #4037 Heating Unit		
Soldering *Tip	Soldering Unsoldering	*Shape: pointed	*Ungar #PL111		
De-soldering aid	To remove molten solder from connection	Suction device	Soldapullt by Edsyn Co. Arleta, California		
Resin (flux) solvent	Remove excess flux from soldered area before application of protective coating	Must not dissolve etched circuit base board material or conductor bonding agent	Freon Acetone Lacquer Thinner Isopropyl Alcohol (100% dry)		
Solder	Component replace- ment Circuit board repair Wiring	Resin (flux) core, high tin content (60/40 tin/lead), 18 gauge (SWG) preferred	//1009		
Protective coating	Contamination, corrosion protection after soldering	Good electrical insulation, corrosion-prevention properties	Humiseal Protective Coating, Type 1B12 by Columbia Technical Corp. Woodside 77, New York		
	<u> </u>	1 Hage #1937 Heat	ing Unit		

<sup>\*</sup>For working on 851 Boards: for general purpose work, use Ungar #1237 Heating Unit (37.5W, tip temp of 750-800°) and Ungar #PL113 1/8" chisel tip.

- c. Shape leads of replacement component to match mounting hole spacing.
- d. Insert component leads into mounting holes, and position component as original was positioned. DO NOT FORCE LEADS INTO MOUNTING HOLES; sharp lead ends may damage plated-through conductor.
- 5-162. ETCHED CONDUCTOR REPAIR. A broken or burned section of conductor can be repaired by bridging the damaged section with a length of tinned copper wire. Allow adequate overlap and remove any varnish from etched conductor before soldering wire into place.

#### 5-163. TRANSISTORS.

- 5-164. The following general information is provided for those who may not have had extensive experience with transistors.
- 5-165. In transistor testing the most important consideration is the base-emitter junction; like the control grid of a vacuum tube, this is the operational

- control point in the transistor. This junction is essentially a solid-state diode, and for the transistor to conduct this diode must be forward-biased.
- 5-166. The transistor symbol (see Figure 5-16) can be used to determine the polarity required to forward-bias the base-emitter junction. Remember that the base material is the middle letter of the transistor type (NPN or PNP). Referring to part A of Figure 5-16, notice that the emitter arrow points toward the N-type material. Thus when the arrow points away from the base (NPN), the base must be positive with respect to the emitter to forward-bias the junction, and when the arrow points toward the base (PNP), the base must be negative with respect to the emitter to forward-bias the junction.
- 5-167. Bias polarity for cutoff and conduction for vacuum tubes as well as transistors is also shown in part A of Figure 5-16. Part B shows simplified versions of the three basic transistor circuits, and gives the amplifier characteristics of each.

<sup>\*\*</sup>Krylon, Inc., Norristown, Pennsylvania

Table 5-24. Connections, RF Circuit Assembly A2, Boards A2A1 through A2A5

Ref No.	Color Code	Connection	Fi Re
A2A	1 Board (00851	-6025) Input Switching Circuit	+-
1 2	wht-blk-orn	20MC input; from A2J1, which accepts cable A1W2 from I. F. GAIN Assy A1 Output, 1MC I. F. bandwidth path; to Output Switching Circuit A2A5, point 30	5-2 5-2
3 4	wht-red-yel wht-red-grn	-24VDC input at A to A2A1K1, from point 28, A2A5 board -24VDC input at B to A2A1K2; from point 26, A2A5	
5	5 Coax Output, 100KC I. F. bandwidth path; to 100KC Bandpass Filter A3 (on rear of casting)		
6	wht	Output, 1-10KC I. F. bandwidth paths; to point 11, 1st 1-10KC BP Flt and Ampl Assy A2A3	
<u>A2A2</u>	2 Board (00851-	6022) 20MC Amplifier	
7	coax	Input 100KC I.F. bandwidth path from 100KC Filter A3, via cable A2W2	5-2
8	coax	Output, 100KC I.F. bandwidth path; to 100KC Filter A5, via cable A2W3	5-2
9	vio	-15VF supply, input; from LV Pwr Supply A9 via A2C1	5-3
10	wht-vio	-15VF supply, output; to Current-Controlled Attenuator A2A6, point 4	5-2
A2A3	Board (00851-6	6023) 1st 1-10KC BP Flt and Ampl Assy	-
11	wht	Input, 1-10KC I.F. bandwidth paths; from point 6, A2A1 Assembly	
12	wht-brn-vio 24VDC input at C to A2A3K1; from point 20, 2nd 1-10KC BP Flt and Amp		5-2
13	wht-blk-vio	-24VDC input at D to A2A3K2; from point 19, A2A4 board	
14	vio	-15VF supply, input; via point 17, A2A4 board	
15	wht	Output, 1-10KC I.F. bandwidth paths; to point 16, A2A4 board	
<u>A2A4</u>	Board (00851-6	024) 2nd 1-10KC BP Flt and Ampl Assy	
16	wht	Input, 1-10KC I.F. bandwidth paths; from point 15, A2A3 board	5-24
17	vio	-15VF supply output to point 14, A2A3 board	5-24
18	wht-orn-blu	-24VDC supply for relay A2A4K1; incoming at E via network A2Z2 and I.F.BANDWIDTH switch	5-24 5-38
19	wht-blk-vio	-24VDC supply outgoing to relay A2A3K2	5-24
20	wht-brn-vio	-24VDC supply outgoing to relay A2A3K1	5-24
21	wht	Output, 1-10KC I.F. bandwidth paths; to point 24, A2A5 Output Switching Circuit board	5-24
22	wht-orn-yel	-24VDC supply for relay A2A4K2; incoming at F via network A2Z3 and I.F.BANDWIDTH switch	5-24 5-38
23	vio	-15VF supply, input; via A2Z1 network and LC filter on rear of A2 casting; from LV Pwr Supply A9	5-37
A2A5 I		26) Output Switching Circuit Assy	
24	wht	Input, 1-10KC I.F. bandwidth paths; from point 21, A2A4 board	5-24
25	wht-orn-vio	-24VDC supply for relay A2A5K1; incoming at G via network A2Z4 and I. F. BANDWIDTH switch	5-24 5-38
26	wht-red-grn	-24VDC supply outgoing to relay A2A1K2	5-24
27	coax	Input, 100KC I.F. bandwidth path; from 100KC BP Flt A5, via cable A2W4	5-24

Table 5-24. Connections, RF Circuit Assembly A2, Boards A2A1 through A2A5 (cont'd)

Ref No.	Color Code	Connection	Fig. Ref
A2A5	Board (00851-6	0026) Output Switching Circuit Assy (cont'd)	
28	wht-red-yel	-24VDC supply outgoing to relay A2A1K1	5-24
29	wht-orn-grn	-24VDC supply for relay A2A5K2; incoming at H via network A2Z5 and I.F.BANDWIDTH switch	5-24 5-38
30	wht-blk-orn	Input, 1MC I.F. bandwidth path; incoming from point 2 on A2A1 board	5-24
31	wht-red-blu	Output, I.F. bandwidth switching circuits; to A2A6 Current-Controlled Attenuator input	5-27

Table 5-24A. Connections, RF Circuit Assembly A2, Boards A2A6, A2A7

Ref No.	Color Code	Connection	Fig. Ref		
A2A6 Board (00851-6021) Current-Controlled Attenuator					
1	wht-red-blu	red-blu 20MC input; from output of I.F. bandwidth switching circuits, point 31, A2A5 board			
2	coax	Control-current input; from VERT DISPLAY switch via cable A2W5	5-27		
3	wht-blk-blu	20MC output, to 20MC I.F. Amplifier (point 5 on A2A7 board)	5-27		
4	wht-vio	-15VF supply; from point 10, A2A2 board	5-24		
A2A7	Board (00851-6	6020) 20MC I.F. Amplifier Assy			
5	wht-blk-blu	20MC input; from point 3, A2A6 board	5-2		
6	coax	Connection to I.F.VERNIER, through feed-through capacitor A2C6 via cable A2W7			
7		-15VF supply; incoming via feed-through capacitor A2C3			
8		Connection to connector J5, I.F. TEST POINT, on rear panel; capacitor A2C2 is in the line to J5.			
9		+15VDC supply; incoming via feed-through capacitor A2C4 and resistor A2R1	1		
10	coax	Video output; to Vertical Amplifier A7 via cable A2W6	5-2		

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk

Table 5-25. Waveform Chart, Model 851 Spectrum Analyzer Display Section (cont'd)

		Torin Chart, Model 85	1 Spectrum Analyzer Dis	play Section (cont'd)
	Analys C. W.		Test Oscilloscope Sensitivity and	
$\vdash$	Analyzer Condition	Test Point	Sweep Speed	Waveform
	I. F. BANDWIDTH3KC SWEEP TIME3MS/CM Input signalCW	ł	CH AND 20MC IF AMPL	
1	VERT DISPLAYSQ Display signal amplitude1.4	Base, AllQ2	2 v/cm 5 ms/cm Ext sync, from 851 HORIZ OUTPUT	
2	VERT DISPLAYSQ Displayed signal amplitude1.4	Emitter, A11Q2	0.2 v/cm 5 ms/cm Ext sync, from 851 HORIZ OUTPUT	
3	VERT DISPLAYLOG Displayed signal amplitude60DB	I. F. Test Point	50 mv/cm Sweep from 851 HORIZ OUTPUT	
4	VERT DISPLAYLOG Displayed signal amplitude60DB	Base, A2A7Q4	2 v/cm Sweep from 851 HORIZ OUTPUT	
		, 851 VERTICA	I L AMPLIFIER	
	SYNCINT I. F. BANDWIDTH1KC SPECTRUM WIDTH10KC, SWEEP TIME3MS/CM Input signalCW	/CM		
5		Input to A7R6 (blanking voltage)	10 v/cm 10 ms/cm	
6		Base, A7Q5 (blanking voltage)	50 v/cm 10 ms/cm	

Table 5-25. Waveform Chart, Model 851 Spectrum Analyzer Display Section (cont'd)

	Table 5-25. Waveform	Chart, Model 851 Spe	ectrum Analyzer Displa	Section (cont'd)
			Test Oscilloscope	
	Analyzer Condition	Test Point	Sensitivity and Sweep Speed	Waveform
<del>                                     </del>			MPLIFIER (cont'd)	Wavelorin
	SYNC INT I. F. BANDWIDTH1KC SPEC WIDTH10KC/CI SWEEP TIME 3MS/CM Input signal CW		WPDIFER (cont a)	
7		Base, A7Q8 (video)	1 v/cm 10 ms/cm	
8		Collector, A7Q8	50 mv/cm 10 ms/cm	
9		Collector, A7Q7	10 v/cm 10 ms/cm	
10	For Service Manuals Contact MAURITRON TECHNICAL SERVICES 8 Cherry Tree Rd, Chinnor Oxon OX9 4QY Tel:- 01844-351694 Fax:- 01844-352554 Email:- enquiries@mauritron.co.uk	Collector, A7Q9	0.1 v/cm Sweep from 851A HORIZ OUTPUT	
11		Collector, A7Q6	10 v/cm Sweep from 851A HORIZ OUTPUT	
	Unless otherwise specified:	SWEEP & HORIZ	AMPLIFIER	
	SYNC LINE SWEEP TIME3MS/CM			
12		Base, A6Q1	10 v/cm 5 ms/cm	
13		Collector, A6Q1	10 v/cm 5 ms/cm	-w-w-

Table 5-25. Waveform Chart, Model 851 Spectrum Analyzer Display Section (cont'd)

		Test Oscilloscope	
Analyzer Condition	Test Point	Sensitivity and Sweep Speed	Waveform
GYDYG	SWEEP & HOR	IZ AMPLIFIER (cont'd	
SYNC LINE SWEEP TIME3MS/CM			
14	Base, A6Q7	2 v/cm 10 ms/cm	
15	Base, A6Q9	2 v/cm 10 ms/cm	
16	Base, A6Q15	5 v/cm 10 ms/cm	
17	Collector, A6Q12	0.5 v/cm 10 ms/cm	
18	Emitter, A6Q14	5 v/cm 10 ms/cm	
19	Connector, A6Q3	5 v/cm 10 ms/cm	
20	Collector, A6Q16	20 v/cm 10 ms/cm	

Table 5-25. Waveform Chart, Model 851 Spectrum Analyzer Display Section (cont'd)

	Table 5-25. Wavelolli	<u> </u>	Test Oscilloscope	
ļ			Sensitivity and	Waveform
<u> </u>	Analyzer Condition	Test Point	Sweep Speed	wavelorm
21	SYNC LINE SWEEP TIME3MS/CM	SWEEP & HORIZ  Collector, A6Q17	20 v/cm 10 ms/cm	
22	SYNCINT SWEEP TIME3MS/CM	Base, A6Q7	5 v/cm 10 ms/cm	
23	SYNCINT SWEEP TIME3MS/CM	Collector, A6Q5	5 v/cm 10 ms/cm	
24		Collector, A6Q13	5 v/cm 10 ms/cm	
		HIGH-VOLTAGE SU	JPPLY AND CRT	•
25		Base, Q1 Q2	5 v/cm 10 μs/cm	WIN
26	For Service Manuals Con MAURITRON TECHNICAL SERV 8 Cherry Tree Rd, Chinn Oxon OX9 4QY Tel:- 01844-351694 Fax:- 01844-3 Email:- enquiries@mauritron.cc	/ICES orCollector, Q1 Q2	10 v/cm 50 μs/cm	MMM
		LOW-VOLTAGE P	OWER SUPPLY	
27		Junction, A9R1 A9R2	2 v/cm 5 ms/cm	مندم
28		Collector, A9Q1	5 v/cm 5 ms/cm	






Figure 5-17. CRT Shield Assembly, Parts Identification


For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:-01844-351694 Fax:-01844-352554
Email:-enquiries@mauritron.co.uk

Table 5-25A. Symbols Used on Schematic Diagrams



For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk








Figure 5-20. 20MC I.F. Input and Attenuator Schematics, 851B






Figure 5-21. RF Circuit Assembly Boards A2A2, A2A3, A2A4



851A-A-12


bly Boards A2A2, A2A3, A2A4



For Service Manuals Contact

MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY

Tel:-01844-351694 Fax:-01844-352554
Email:- enquiries@mauritron.co.uk



SEE TABLE 5-24

Figure 5-21

Figure 5-23. RF Circuit As




Figure 5-22. RF Circuit Assembly A2, Top Cover Removed




Figure 5-21


Figure 5-23. RF Circuit Assembly Boards A2A5 and A2A1

NOTES:

- I. RESISTANCE IN OHMS, CAPACITANCE IN PICOFARADS, INDUCTANCE IN MICROHENRIES, UNLESS OTHERWISE NOTED.
- 2. RELAYS SHOWN DE-ENERGIZED ( I.F. BANDWIDTH AT IOKC)
- 3. VF = FILTERED VOLTAGE (SEE LV PWR SUPPLY DIAG)
- \* A4SI = I.F. BANDWIDTH SWITCH, SEE SWITCH DETAIL, FIG. 5-38
- VOLTAGES MEASURED WITH -hp-410C ELECTRONIC VOLTMETER 100 MEGOHMS INPUT RESISTANCE

OPTIMUM VALUE SELECTED AT FACTORY, AVERAGE VALUE SHOWN O HEX NUT ADJ

I.F. BANDWIOTH	RELAY			
POSITION	ENERGIZED			
IMC	AZAIKI 8 AZA5KZ			
100 KC	AZAIKZ & AZA5KI			
FORC	NONE			
3 KC	A2A3K2 8 A2A4KI			
IKC	A2A3KI & A2A4K2			




20MC

FROM

COFARADS,
FHERWISE NOTED.

OWIDTH AT LOKC)
PPLY DIAG)
VITCH DETAIL, FIG. 5-38
ECTRONIC

I.F. BANDWIDTH POSITION	RELAY ENERGIZED		
I M C	A2AIKI & A2A5K2		
100 K C	A2AIK2 & A2A5KI		
10 K C	NONE		
3 K C	A2A3K2 & A2A4KI		
1 K C	A2A3KI & A2A4K2		



For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tet:-01844-351694 Fax:-01844-352554
Email:-enquiries@mauritron.co.uk

REFERENCE DESIGNATIONS WITHIN ASSEMBLIES ARE ABBREVIATED. ADD ASSEMBLY DESIGNATION AS PREFIX TO FORM COMPLETE DESIGNATION. e.g., RI OF ASSEMBLY AI IS AIRI, AND IS LISTED AIRI IN THE TABLE OF REPLACEABLE PARTS. DESIGNATIONS OF COMPONENTS NOT WITHIN ASSEMBLIES ARE COMPLETE AS SHOWN.





Figure 5-24. I. F. Bandwidth Switching Circuits, 851B








Figure 5-25. VERT DISPLAY Switch A11





ERT DISPLAY Switch A11




For the Armonia of all Services 8 Charn Chinnor
Oxon Ox9 4QY
Tel: 01844-351894 Fax - 01844-352554
Effail: enquinee@mauritron to uk



ह राज्यसम् ५ १४४

3.36

12.









Figure 5-25


Figure 5-26. RF Circuit Assembly Boards A2A6, A2A7, and Rear of Casting







RIO SOO



COPYRIGHT 1965 BY HEWLETT-PACKARD CO

modern and a management of the management of the second and the second of the second o




Figure 5-27. VERT DISPLAY Switch, Current-Controlled Attenuator, and 20MC I.F. Amplifier Schematics, 851B

Table 5-26. Connections, Sweep and Horizontal Amplifier Assy A6

Ref No.	Color Code	Connection	Fig. Ref
1	red	+100 vdc input, from LV Power Supply A9	5-37
2	wht	To CRT (V1) horizontal deflection plate, terminal D2	5-33, 5-35
3	grn	To CRT (V1) horizontal deflection plate, terminal D1	5-33, 5-35
4	wht-grn-blu	To R9, HORIZ POS adjust	5-33
5	wht-grn-gra	From adjustable contact on R9, HORIZ POS adjust	5-33
6	wht-grn-vio	To R9, HORIZ POS adjust	5-33
7	blk	Chassis ground	
8	yel	Blanking signal to Vert Amplifier A7	5-29
9	vio	-15 vdc, from LV Power Supply A9	5-37
10	wht-red	+15 vdc, from LV Power Supply A9	5-37
11	wht-red	+15 vdc, from LV Power Supply A9	5-37
12	vio	-15 vdc, from LV Power Supply A9	5-37
13	coax cable W3	To J7, SWEEP OUTPUT, on rear panel	5-33
14	wht-orn-yel	To J8, HORIZ OUTPUT, on rear panel	5-33
15	wht-red-blu	To fixed contact, rear of wafer 2, SWEEP TIME switch, A10S1	5-30, 5-39
16	wht-red-vio	To contactor, front of wafer 1, SWEEP TIME switch, A10S1	5-30, 5-39
17	wht-brn-gra	To contactor, front of wafer 3, SWEEP TIME switch, A10S1	5-30, 5-39
18	wht	To adjustable contact on SWEEP TIME VERNIER, A10R1 To fixed contact on front of wafer 3, SWEEP TIME switch, A10S1	5-30, 5-33
19	wht-brn-yel		5-30, 5-33
20	wht-brn-blu	To fixed contact on front of wafer 3, SWEEP TIME switch, A10S1	5-30, 5-33
21	wht-brn-red	To fixed contact on front of wafer 3, SWEEP TIME switch, A10S1	5-30, 5-33
22	wht-grn	To junction A10R2, VERNIER A10R1	5-33
23	wht-brn-orn	To fixed contact on front of wafer 3, SWEEP TIME switch, A10S1	5-30, 5-33
24	wht-brn-vio	To fixed contact on front of wafer 3, SWEEP TIME switch, A10S1	5-30, 5-33
25	wht-brn-grn	To fixed contact on front of wafer 3, SWEEP TIME switch, A10S1	5-30, 5-33
26	vio	To junction A10R2, VERNIER A10R1	5-33
27	wht-red-grn	To contactor on front of wafer 2, SWEEP TIME switch, A10S1	5-30, 5-33
28	wht	From contactor on rear of wafer 1, SYNC switch S2	5-32
29	vio	-15 vdc, from LV Power Supply A9	5-37
30	blk	Chassis ground	
31	vio	-15 vdc, from LV Power Supply A9	5-37
32	wht-red	+15 vdc, from LV Power Supply A9	5-37
33	wht-grn	To fixed contact on SINGLE SWEEP switch S3	5-33
34	blk	Chassis ground	
35	blk	Chassis ground	5-33
36	coax	To contactor on front of wafer 1, SYNC switch S2	5-33
37	blu	To SINGLE SWEEP lamp DS1	
38	grn	To wafer 1F on SYNC switch S2	5-32 5-32
39	grn	To wafer 1F on SYNC switch S2	
40	wht-yel-grn	To wafer 6R on SWEEP TIME switch A10S1	5-30, 5-39
41	wht-red-grn	To SWEEP INPUT J2 on rear panel	5-33
42	wht-orn-grn	To J3, BLANKING INPUT	5-33
43	wht-blu-gra	To wafer 6F on SWEEP TIME switch A10S1	5-30
44	grn	To wafer 6R on SWEEP TIME switch A10S1	5-30
45	wht-yel	To wafer 6R on SWEEP TIME switch A10S1	5-30

#### Horizontal Amplifier Assy A6

Fig. Ref

inection

Supply A9	5-37
on plate, terminal D2	5-33, 5-35
on plate, terminal D1	5-33, 5-35
•	5-33
HORIZ POS adjust	5-33
	5-33
er A7	5-29
<b>A</b> 9	5-37
A9	5-37
A9	5-37
A9	5-37
ir panel	5-33
r panel	5-33
2, SWEEP TIME switch, A10S1	5-30, 5-39
SWEEP TIME switch, A10S1	5-30, 5-39
SWEEP TIME switch, A10S1	5-30, 5-39
TIME VERNIER, A10R1	5-30, 5-33
er 3, SWEEP TIME switch, A10S1	5-30, 5-33
er 3, SWEEP TIME switch, A10S1	5-30, 5-33
er 3, SWEEP TIME switch, A10S1	5-30, 5-33
.10R1	5-33
	5-30, 5-33
er 3, SWEEP TIME switch, A10S1	1 ' 1
er 3, SWEEP TIME switch, A10S1	5-30, 5-33
er 3, SWEEP TIME switch, A10S1	5-30, 5-33
.10R1	5-33
;, SWEEP TIME switch, A10S1	5-30, 5-33
r 1, SYNC switch S2	5-32
7 <b>A</b> 9	5-37
	E 27
, A9	5-37 5-37
7 A9	5-33
EEP switch S3	7-33
., SYNC switch S2	5-33
	5-33
	5-32
	5-32
witch A10S1	5-30, 5-39
anel	5-33
<del></del>	5-33
witch A10S1	5-30
witch A10S1	5-30
	5-30
witch A10S1	3-30

Section V Tables 5-26, 5-27 and Figure 5-28




Figure 5-28. Vertic

Table 5-27 Connections

Table 5-26		Table 5-27. Connections,				
Ref No.	Color Code	Connection	Fig. Ref			
1	red	+100VDC regulated from LV Pwr Sup A9	5-37			
2	wht-blk-red	To Int Level Adj R2	5-35			
3	wht	Tc 3rd Anode, CRT V1	5-35			
4	vio	-15VDC reg from LVPS A9	5-37			
5	coax	(D) video signal from 20MC I. F. Ampl Assy A2A7	5-27			
6	wht	To CRT vert deflection plate D3	5-35			
7	grn	To CRT vert deflection plate D4	5-35			
8	blk	Chassis ground	5-37			
9	wht-red	+15VDC regulated from LV Pwr Sup A9	5-37			
10	wht-red-vio	To VERT POS adj R8	5-29			
11	wht-red-gra	To VERT POS adj R8	5-29			






Figure 5-28. Vertical Amplifier A7 Board

Table 5-26

Table 5-27. Connections, Vertical Amplifier A7 Board

					<del></del>		
Ref No.	Color Code	Connection	Fig. Ref	Ref No.	Color Code	Connection	Fig. Ref
1	red	+100VDC regulated from	5-37	12	wht-grn-gra	To BASE LINE CLIP. R7	5-29
		LV Pwr Sup A9		z	wht-orn-yel	-24VDC via I. F. BAND-	5-38
2	wht-blk-red	To Int Level Adj R2	5-35			WIDTH switch	
3	wht	To 3rd Anode, CRT V1	5-35	13	grn	From movable contact,	5-29
4	vio	-15VDC reg from LVPS A9	5-37	ŀ	:	BASE LINE CLIP. R7	
5	coax	(D) video signal from 20MC I. F. Ampl Assy A2A7	5-27	14	yel	Blanking signal from Sweep & Horiz Ampl A6 (emitter of A6Q6)	5-33
6	wht	To CRT vert deflection	5-35	Y	wht-orn-grn	-24VDC via I. F. BW	5-38
		plate D3		w	wht-orn-blu	-24VDC via I. F. BW	5-38
7	grn	To CRT vert deflection	5-35	x	wht-orn-vio	-24VDC via I. F. BW	5-38
8	blk	plate D4 Chassis ground	5-37	15	wht-yel	Blank. sig to HV Pwr Sup A8 (applied to top of IN-	5-35
9	wht-red	+15VDC regulated from	5-37		*	TENSITY divider)	
		LV Pwr Sup A9		16	wht-blk-grn	To BASE LINE CLIP. R7	5-29
10	wht-red-vio	To VERT POS adj R8	5-29	17	wht-grn	To movable contact	5-29
11	wht-red-gra	To VERT POS adj R8	5-29	1		VERT POS adjust R8	









Figure 5-29. Vertical Amplifier Schematic, 851B



Component Identification (see Table 5-28)




Figure 5-31. Sweep and Horizonta

For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:-01844-351694 Fax:-01844-352554
Email:-enquiries@mauritron.co.uk




Figure 5-32. SYNC.




Figure 5-31. Sweep and Horizontal Amplifier A6 Board (see Table 5-26)






Figure 5-32. SYNC Switch S2 (see Table 5-30)








Figure 5-33. Sweep and Horizontal Amplifier Schematic, 851B

Table 5-28. Connections, SWEEP TIME Switch A10S1, 851B

Ref No.	Color Code	Wafer	Connections	Fig. Ref			
1	wht-red-vio	1 F	Contactor; connects to Sweep Capacitors A10C7-A10C12 From junction of A6R35, base of A6Q11, p/o Miller Integrator in Sweep Generator	5-33			
2	blk	1R	Ground for Reset Capacitors A10C1-A10C6; strap				
3	wht-red-grn	2 F	Contactor; connects to Reset Capacitors A10C1-A10C6 From A6R13: when SYNC is at INT, -15VDC comes in over this lead.				
4	wht-red-blu	2R	Ramp voltage from Miller Integrator A6Q11-A6Q14. To tie point for Sweep Capacitors A10C7-A10C12	5-33			
5	wht-brn-gra	3 F	Contactor; connects to leads outgoing to Sweep Time Adjusts From A6R35	5-33			
6	wht-brn-yel		To A6R34, 1 Sec Adj	5-33			
7	wht-brn-orn		To A6R33, .3 Sec Adj  For Service Manuals Contact MAURITRON TECHNICAL SERVICES	5-33			
8	wht-brn-blu		To A6R32, .1 Sec Adj  8 Cherry Tree Rd, Chinnor Oxon OX9 4QY	5-33			
9	wht-brn-vio		To A6R31, 30 Msec Adj <b>Tel:-01844-351694 Fax:-01844-352554</b>	5-33			
10	wht-brn-red		To A6R30, 10 Msec Adj	5-33			
11	wht-brn-grn		To A6R29, 3 Msec Adj	5-33			
12	wht-blk-yel	4F	Contactor; to I. F. BANDWIDTH, wafer 1R	5-38			
13	wht-brn wht-red wht-orn wht-grn		To pins on CONTROL connector J9 pin 1 pin 2 pin 3 pin 4	5-39			
14	wht-blk-blu wht-blu wht-vio	4R	Contactor; to I. F. BANDWIDTH, wafer 1F To pins on CONTROL connector J9 pin 5 pin 6	5-38 5-39			
15	wht-blk-vio wht-gra wht wht-blk	5 F	Contactor; to I. F. BANDWIDTH, wafer 2R To pins on CONTROL connector J9 pin 8 pin 9 pin 9 pin 10	5-38 5-39			
16	wht-blk-grn wht-gra wht wht-blk	5R   	Contactor; to I. F. BANDWIDTH, wafer 2F To pins on CONTROL connector J9 pin 8 pin 9 pin 10	5-38 5-39			
17	yel	6F	From A6Q6 emitter (blanking signal)	5-33			
18	wht-blu-gra		From A6Q20 collector (amplified external blanking signal)				
19	yel	*	Contactor; to Blanking Amplifier A7Q4-A7Q3				
20	grn	6R	Contactor; to base of A6Q16, Horizontal Amplifier				
21	wht-yel-grn		From emitter of A6Q21 in external sweep input circuit				
22	wht-yel		From emitter of A6Q14, Sweep Generator output				
23	wht-grn		From SWEEP TIME VERNIER A10R1, via A6R38, to base of A6Q12 in Sweep Generator Miller Integrator				
24	wht		From Sweep Time Adjusts A6R29-A6R34 to adjustable contact on SWEEP TIME VERNIER A10R1				
25	vio	₩	Returns A10R1 to -15VDC	7			

## FEP TIME Switch AloS1, 851B

Connections	Fig. Ref
nects to Sweep Capacitors A10C7-A10C12 of A6R35, base of A6Q11, p/o Miller Sweep Generator	5-33
set Capacitors A10C1-A10C6; strap	
nects to Reset Capacitors A10C1-A10C6 when SYNC is at INT, -15VDC comes in i.	5-33
from Miller Integrator A6Q11-A6Q14. r Sweep Capacitors A10C7-A10C12	5-33
nects to leads outgoing to Sweep Time Adjusts	5-33
ec Adj	5-33
Sec Adj	5-33
Sec Adj	5-33
Msec Adj	5-33
Msec Adj	5-33
Isec Adj	5-33
.F.BANDWIDTH, wafer 1R	5-38
NTROL connector J9	5-39
.F.BANDWIDTH, wafer 1F NTROL connector J9	5-38 5-39
.F.BANDWIDTH, wafer 2R NTROL connector J9	5-38 5-39
strap, from 5R	
.F.BANDWIDTH, wafer 2F NTROL connector J9	5-38 5-39
itter (blanking signal)	5-33
ollector (amplified external blanking signal)	
lanking Amplifier A7Q4-A7Q3	
ase of A6Q16, Horizontal Amplifier	
of A6Q21 in external sweep input circuit	
of A6Q14, Sweep Generator output	
TIME VERNIER A10R1, via A6R38, to base weep Generator Miller Integrator	
me Adjusts A6R29-A6R34 to adjustable EEP TIME VERNIER A10R1	
to -15VDC	1




Figure 5-34.

•	Table 5	-28	Table 5-29. Conn∈							
	Ref* No.	Color Code	Connection							
	1	yel	To terminal on INTENSITY control R1							
	2	wht-yel	From Vert Ampl A7 (blankir signal)							
	3	wht	From adjustable contact on INTENSITY control R1							
	4 grn 5 red		To grid, pin 3, CRT (V1)							
			To terminal on INTENSITY control R1							
	6	gra	To cathode, pin 2, CRT (V1)							






Figure 5-34. HV Power Supply A8 Board

Table 5-28

Table 5-29. Connections, HV Power Supply A8 Board

Ref* No.	Color Code	Connection	Fig. Ref	Ref* No.	Color Code	Connection	Fig. Ref
1	yel	To terminal on INTENSITY control R1	5-35	7	wht-brn-orn	To terminal on Astig adjust (R3)	5-35 I
2	wht-yel	From Vert Ampl A7 (blanking signal)	5-29	8	orn	To terminal on FOCUS control (R4)	
3	wht	From adjustable contact on INTENSITY control R1	5-35	9	brn	To terminal on FOCUS control (R3)	*
4	grn	To grid, pin 3, CRT (V1)		10	vio	From -15 vdc supply	5-37
5	red	To terminal on INTENSITY control R1		11	red	To CRT post-accelerator anode	5-35
6	gra	To cathode, pin 2, CRT (V1)	*	12	grn	To Int Level adjust (R2)	5-35
			* Figu:	re 5-34			








Table 5-30. Connections, SYNC Switch S2, 851B

tef √o.	Color Code	Connection	Fi Re	
1	coax	To contactor on wafer 1F; from base of A6Q1, p/o Schmitt Trigger in input to Sweep Generator	5-	33
2	wht-grn-blu	Also connects to contactor on wafer 1F: from SINGLE SWEEP switch S3		
3	coax	Cable from SYNC INPUT J6 on rear panel connects here; wafer 1F; EXT position.		
4	blk	Chassis ground; connects to Sweep and Horiz Ampl board, point 30 (see Figure 5-31); wafer 1F, INT position; wafer 1R, SINGLE SWEEP, LINE, and EXT positions.		
5	grn	From A6Q7, p/o Reset Multivibrator, via A6R66; wafer 1F; SINGLE SWEEP position.		
6	grn	To A6Q22, amplifier in SINGLE SWEEP indicator lamp circuit; wafer 1F; SINGLE SWEEP position.		1
7	gra	Conductor from Line Transformer T1 connects here; wafer 1F; LINE position.	5-	37
8	wht	To contactor on wafer 1R; at INT, connects -15V to Reset Capacitors A10C1-A10C6 on SWEEP TIME, via A6R13.	5-	33

## 1 S2, 851B

a	Fig. Ref
3Q1, p/o Schmitt Trigger	5-33
om SINGLE SWEEP	
connects here;	
riz Ampl board, point 30 n; wafer 1R, SINGLE SWEEP,	
A6R66; wafer 1F;	
ndicator lamp circuit;	<b>Y</b>
nects here; wafer 1F;	5-37
s -15V to Reset Capacitors	5-33




Figure 5-36. LV Power Supp

Table 5-31. Connections, LV Power

Table	5-30	Table 5-31. Connec	ctions, I	7 N LOM
Ref* No.	Color Code	Connection	Fig. Ref	Ref* No.
1	red	+100 vdc regulated: to 851 circuits	5-37	15
2	blk	Chassis ground		16 17
3	wht-red-grn	To base of Emitter Foll Q5		18
4	wht-orn-yel	To emitter of Series Reg Q4		19
5	wht-orn-grn	To base of Series Reg Q4		20
6	wht-orn-blu	To coll of Series Reg Q4		
7	wht-red-yel	From T1-15		21
8	wht-red-yel	From T1-13		
9	wht-orn-blu	From T1-14		
10	wht-vio	To coll of Series Reg Q6 and Emitter Foll Q5		22
11	wht-brn-yel	To emitter of Series Reg Q6		
12	wht-blk-blu	To coll of Series Reg Q3		
13	wht-blk-grn	To base of Series Reg Q3		23
14	wht-blk-red	From T1-12	1	
			* Figu	re 5-36






Figure 5-36. LV Power Supply A9 Board

Table 5-30

Table 5-31. Connections, LV Power Supply A9 Board

Ref*	G.1 C-d-	Connection	Fig.	Ref*	Color Code	Connection	Fig Ref
No.	Color Code	Connection	Ref	+		<u> </u>	
1	red	+100 vdc regulated: to	5-37	15	wht-blk-orn	From C3	5-37
		851 circuits		16	wht-blk-vio	From T1-9	
2	blk	Chassis ground		17	wht-blk-vio	From T1-10	
3	wht-red-grn	To base of Emitter Foll Q5		18	blk	Chassis ground	
4	wht-orn-yel	To emitter of Series Reg Q4		19	wht-brn-vio	To junction of C4 and L4	
5	wht-orn-grn	To base of Series Reg Q4		20	wht-brn-grn	To junction of emitter of Q5	
6	wht-orn-blu	To coll of Series Reg Q4				& base of Series Reg Q6	1
7	wht-red-yel	From T1-15		21	vio	-15 vdc regulated:	
8	wht-red-yel	From T1-13				to filter in RF Ckt Assy A2	5-24 5-37
9	wht-orn-blu	From T1-14				to 851 circuits via	5-39
10	wht-vio	To coll of Series Reg Q6				J9-12	
		and Emitter Foll Q5		22	wht-red	+15 vdc regulated:	5-37
11	wht-brn-yel	To emitter of Series Reg Q6			,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	to 851 circuits	
12	wht-blk-blu	To coll of Series Reg Q3				to 8551 circuits via J9-13	5-39
13	wht-blk-grn	To base of Series Reg Q3		23	blk		
14	wht-blk-red	From T1-12	*	23	DIK	Chassis ground to 8551 via J9-14	





REFERENCE DESIGNATIONS WITHIN ASSEMBLIES ARE ABBREVIATED. ADD ASSEMELY DESIGNATION AS PREFIX TO FORM COMPLETE DESIGNATION. e.g., RI OF ASSEMBLY AI IS AIRI, AND IS LISTED AIRI IN THE TABLE OF REPLACEABLE PARTS. DESIGNATIONS OF COMPONENTS NOT WITHIN ASSEMBLIES ARE COMPLETE AS SHOWN.




Figure 5-37. LV Power Supply Schematic, 851B

Table 5-32. Connections, I.F.BANDWIDTH Switch A4S1, 851B


Ref	Color Code	Connection	Fig. Ref
1	wht-orn-yel	-24V to tie point at A2Z3 on rear of A2 RF Circuit casting; energizes I. F. Bandwidth Switching relays A2A4K2, A2A3K1, and filter-switching relay A7K2.	5-24 5-26 5-29
2	wht-blk-yel	To SWEEP TIME switch, wafer 4F, on which connection is made to leads outgoing to CONTROL connector J9*.	5-30
3	wht-orn-blu	-24V to tie point at A2Z2 on A2 casting rear; energizes relays A2A4K1, A2A3K2, A7K1, A7K4.	5-24 5-26 5-29
4	wht-blk-blu	To SWEEP TIME switch, wafer 4R, on which connection is made to leads outgoing to CONTROL connector J9*.	5-30
5	wht-orn-vio	-24V to tie point at A2Z4 on A2 casting rear; energizes relays A2A5K1, A2A1K2, A7K1, A7K3, A7K4.	5-24 5-26 5-29
6	wht-blk-vio	To SWEEP TIME switch, wafer 5F, on which connection is made to leads outgoing to CONTROL connector J9*.	5-30
7	wht-orn-grn	-24V to tie point at A2Z5 on A2 casting rear; energizes relays A2A5K2, A2A1K1, A7K1.	5-24 5-26 5-29
8	wht-blk-grn	To SWEEP TIME switch, wafer 5R, on which connection is made to leads outgoing to CONTROL connector J9*.	5-30
9	wht-vio	-24VDC from LV Power Supply A9	5-37

<sup>\*</sup> For AUTO SELECT operation; via inter-unit CONTROL cable, connection is made to SWEEP TIME switch in 8551 RF Section.

## BANDWIDTH Switch A4S1, 851B

on.

Connection	Fig. Ref
3 on rear of A2 RF Circuit casting; ith Switching relays A2A4K2, A2A3K1, lay A7K2.	5-24 5-26 5-29
, wafer 4F, on which connection is made )NTROL connector J9*.	5-30
2 on A2 casting rear: energizes relays <1, A7K4.	5-24 5-26 5-29
, wafer 4R, on which connection is made NTROL connector J9*.	5-30
4 on A2 casting rear; energizes relays <1, A7K3, A7K4	5-24 5-26 5-29
, wafer 5F, on which connection is made ONTROL connector J9*.	5-30
5 on A2 casting rear; energizes relays %1	5-24 5-26 5-29
n, wafer 5R, on which connection is made DNTROL connector J9*.	5-30
· Supply A9	5-37
unit CONTROL cable, connection is made	



For Service Manuals Contact
MAURITRON TECHNICAL SERVICES
8 Cherry Tree Rd, Chinnor
Oxon OX9 4QY
Tel:- 01844-351694 Fax:- 01844-352554
Email:- enquiries@mauritron.co.uk





Figure 5-38.



SEE NOTE 1

SEE NOTE 2

SEE NOTE 2

SEE NOTE 2

SEE NOTE 4

OF REAR 100 5

OF REA

NOTES

I. SWITCH VIEWED FROM KNOB END, IN MAX CCW POSITION (IKC)

CODE	A4SI POS
1	IKC
2	зкс
3	IOKC
4	IOOKC
5	IMC
6	AUTO SELECT

2. SEE TABLE 5-32 :






Figure 5-38. I. F. BANDWIDTH Switch A4S1





MALE: 6

Tet: 01844 (5) 693 (Email: enquire-sur

YH! SRA

þ

b.\*




Figure 5-39. 851B SWEEP TIME Switch A10S1, Schematic

terminal board

test point

# SECTION VI REPLACEABLE PARTS

#### 6-1. INTRODUCTION.

assembly

motor

- 6-2. This section contains information for ordering replacement parts.
- 6-3. Tables 6-1 and 6-2 list parts in the alpha-nu-merical order of their reference designations and give the Hewlett-Packard stock number and description for each part, together with any applicable notes. Miscellaneous parts not assigned a reference designation are listed at the end of the Table which covers the Assembly with which the part is associated. Reference Designation Index Tables cover the following Assemblies:
  - Table 6-1. Assemblies A1, A3 through A12, and parts mounted on the chassis
  - Table 6-2. RF Circuit Assembly A2
- 6-4. Table 6-3 lists parts in the alpha-numerical order of their hp Stock Numbers, and provides the following information on each part: 1) description of part (see list of abbreviations below), 2) typical manufacturer of part in five-digit code (see code list of manufacturers in Appendix), 3) manufacturer's stock number, 4) total quantity used in instrument (TQ col).

#### 6-5. ORDERING INFORMATION.

- 6-6. To order a replacement part, address order or inquiry to your nearest Hewlett-Packard sales and service office. Addresses of sales and service offices around the world are given at the rear of this manual.
- 6-7. Specify the following for each part: 1) model and complete serial number of instrument, 2) Hewlett-Packard stock number, 3) reference designation, and 4) description.
- 6-8. When ordering from Hewlett-Packard always furnish the hp stock number. The part you receive may not be made by the manufacturer listed but will be electrically and mechanically interchangeable, and performance will be equal. Manufacturer's part number is listed for your convenience should you want to order directly.
- 6-9. To order a part not listed, give complete description of the part and include its function and location.

mechanical part

plug

## REFERENCE DESIGNATORS

MP

misc electronic part

fuse

5	_	1110001	F		_	luse	P	#	piug	TP	=	test point
C	#	capacitor	FL		32	filter	Q	-		v	=	vacuum tube, neon
CP	¥	coupling	J		=	jack	R	×	resistor			bulb, photocell, etc.
CR	=	diode	K		=	relay	RT	=	thermistor	W	=	cable
DL	=	delay line	L		3	inductor	S	=		x	=	socket
DS	=	device signaling (lamp)	M		=	meter	Ť	=	transformer	Ÿ		crystal
		, , , , , , , , , , , , , , , , , , , ,					•		cransior mer	•	_	Cijstai
ABBREVIATIONS												
Α	=	amperes	GE		=	germanium	N/C	=	normally closed	RMO	=	rack mount only
A.F.C	=	automatic frequency control	GL		=	glass	NE		neon	RMS		root-mean-square
AMPL	=	amplifier	GRD	)	=	ground(ed)	NI PL		nickel plate	111.10		root mean square
		•				<b>3</b> (,	N/O		normally open	S-B	=	slow-blow
B.F.O.	=	beat frequency oscillator	H		=	henries	NPO		negative positive zero	SCR	=	screw
		beryllium copper	HEX	•	=	hexagonal	0		(zero temperature	SE		selenium
вн		binder head	HG		=				coefficient)	SECT	=	
BP	=	bandpass	HR		_		NRFR	=	·			i = semiconductor
BRS	=	brass	1114			11041 (b)	MILE IL	-		SI		silicon
BWO	=		IF		_	intermediate freq	NSR	_	field replacement	SIL	-	silver
D.11 O		buckward wave oscillator	IMP	C.		impregnated	Non	-	not separately	SL	_	slide
CCW	_	counter-clockwise	INCI			incandescent			replaceable	SPL		
CER	=		INCI			include(s)	022				==	special
CMO		cabinet mount only	INS			insulation(ed)	OBD		order by description	SST	#	stainless steel
COEF		coefficient	INT			internal	OH		oval head	SR	=	split ring
COM			11/11		=	internal	OX	=	oxide	$\mathtt{STL}$	=	steel
COMP		common				1.17 4000	_					
CONID	=	oomposition.	K		=	kilo = 1000	P	=	F	TA	æ	tantalum
CP	=	connector.					PC		printed circuit	TD	=	time delay
		cadmium plate	LIN			linear taper	PF	=	picofarads =	TGL	=	toggle
CRT		cathode-ray tube				lock washer			10-12 farads	TI	=	titanium
CW	=	clockwise	LOG			logarithmic taper	PH BRZ	, =	phosphor bronze	TOL	Ξ	tolerance
D 20 2			LPF		#	low pass filter	PHL	==	Phillips	TRIM	=	trimmer
DEPC		deposited carbon				•	PIV	=	peak inverse voltage	TWT	=	traveling wave tube
DR	=	drive	M		=	milli = 10 <sup>-3</sup>	P/O	=	part of			-
			MEG	ł	=	$meg = 10^6$	POLY	=	polystyrene	U	=	$micro = 10^{-6}$
	7	electrolytic	MET	FL	M =	metal film	PORC		porcelain			
	=	encapsulated	MFR	Ł	=	manufacturer	POS	=	·	VAR	=	variable
EXT	=	external	MIN	AT	=	miniature	POT	=	potentiometer			dc working volts
			MOM		=	momentary	PP	-	peak-to-peak	. 2011		ac norming rotes
F	=	farads	MTG		=	mounting	PT	_	point	W/	=	with
FH	=	flat head	MY	•	=	•	RECT		rectifier	w		watts
		fillister head	747 7		_	,	RF			ww	_	watts wirewound
FXD		fixed	N		_	nano (10 <sup>-9</sup> )	RH		radio frequency	ww W/O	=	
		11nou	14		-	namo (10 -)	нn	=	round head	w/O	=	without

Section VI Table 6-1

Table 6-1. Reference Designation Index, General

Reference Designation	⊕ Stock No.	Description #	Note
A 1	00851-6002	SWITCH ASSY: IF GAIN (DB)	
A1R1	0727-0036	R*FXD UEPC 71.16 OHM 1% 1/2W	ŀ
A1R2	0727-0042	R1FXD DEPC 96.25 OHM 1/2% 1/2W	1
A1R3	0727-0042	R * F X D U E P C 96 • 25 OHM 1/2% 1/2W	i
AIR4	0727-0062	R*FXD DEPC 247.5 OHM 1/2% 1/2W	
A1R5	0727-0033	R*FXD DEPC 61-11 OHM 1% 1/2W	
A1R6	0727-0033	R*FXD UEPC 61-11 OHM 1% 1/2W	
AIR7	0727-0122	RIFXD DEPC 2.51K OHM 1% 1/2W	İ
AIRS	0727-0025	R:FXD DEPC 51.5 OHM 1% 1/2W	- [
AIR9	0727-0025	R*FXD DEPC 51.5 OHM 1% 1/2W	
AIR10	0727-0005	RIFXU DEPC 5.77 OHM 1/2% 1/2%	
A1R11	0727-0004		
	0727-0094 0727-0094	RIFXD DEPC 670 OHM 1/2% 1/2W	
AIR12		RIFXU DEPC 870 OHM 1/2% 1/2W	1
A1R13	0727-0008	RIFXU DEPC 11.61 OHM 1/2% 1/2W	
A1R14	0727-0074	R*FXD DEPC 436 OHM 1/2% 1/2W	
A1R15	0727-0074	RIFXD DEPC 436 OHM 1/2% 1/2W	
AIR16	0727-0010	RIFXD DEPC 17.61 OHM 1/2% 1/2W	
A1R17	0727-0063	RIFXD DEPC 292.5 OHM 1% 1/2W	
AIR18	0727-0063	RIFXD DEPC 292.5 OHM 1% 1/2W	
A1R19	0727-0017	RIFXD DEPC 37.35 OHM 1/2% 1/2W	
A1R20	0727-0048	RIFXD DEPC 150 OHM 1% 1/2W	
A1R21	0727-0048	R:FXD DEPC 150 OHM 1% 1/2W	
4151	3100-0812	ROTARY SWITCH	
	00851-0014	IF GAIN SWITCH   COVER:SWITCH IF GAIN	
	00851-0015	PLATE COVER IF GAIN	
	00851-0016		
	00851-2027	BRACKET: IF GAIN (DB) SWITCH KNOB: IF GAIN 0-70 DB	
		For Service Manuals Contact	
A152	3100-0812	ROTARY SWITCH MAUHITRON TECHNICAL SERVICES	
	AA.E. 6	IF GAIN SWITCH 8 Cherry Tree Rd, Chinnor	
	00851-0014	COVER I SWITCH IF GAIN Oxon OX9 4QY	
İ	00851-0015	PLATE . COVER IF GAIN Tel: 01844-351694 Fax: 01844-352554	
	00851-2028	KNOB I IF GAIN 0-10 DB Email:- enquiries@mauritron.co.uk	
Alwi	00851-6014	CABLE ASSY:ATTEN INPUT	
_		13-INCH COAX W/BNC MALE PLUG P1	1
A1#2	00851-6013	CABLE ASSYLATTEN OUTPUT	
42	00851-6003	12-INCH COAX W/BNC MALE PLUG P2 RF CIRCUIT ASSEMBLY	
-		SEE TABLE 6-2 FOR LISTING OF COMPONENTS	
	00851-2022	CAVITY FILTER P/O A3	
43	00851-6028	FILTER ASSY. : 100KC BAND-PASS	1
	08551-2083	BUSHING	
A3C1	0160-0822	CIFXO TI 2.2PF 5% 500VDCW	
\3J1	1250-0228	CONNECTORING JACK CHASSIS	
3332	1250-0228	CONNECTORIRE JACK CHASSIS	
3L1	00851-8008	COIL:RF VARIABLE	
451	00851-6007	SWITCH ASSY .: I .F. BANDWIDTH	
	0370-0112	KNOB: I F BANDWIDTH	
			1

<sup>#</sup> See list of abbreviations in introduction to this section

Table 6-1. Reference Designation Index, General (cont'd)

Designation	⊕ Stock No.	Description #	Note
	00851-2022	CAVITY: FILTER, P/O A5	
A5	00851-6028	FILTER ASSY. 1100KC BAND-PASS	•
	08551-2083	BUSHING	
A5C1	0160-0822	C#FXD TI 2.2PF 5% 500VDCW	
A5J1	1250-0228	CONNECTORING JACK CHASSIS	1
A5J2	1250-0228	CONNECTORIRE JACK CHASSIS	
A5L1	00851-8008	COIL:RF VARIABLE	
A6	00851-6038	BOARD ASSY: SWEEP & HORIZ AMPL	
A6C1	0140-0207	CIFXU MICA 330PF 5% 500VDCW	
A6C2	0140-0180	C*FXD MICA 2000 PF 2% 300 VDCW	
A6C3	0160-0174	C:FXD CER 0.470F 80% 25VDCW	
A6C4	0140-0207	CIFXD MICA 330PF 5% 500VDCW	
A6C5	0140-0207	CSFXD MICA 330PF 5% 500VDCW	
A6C6	0150-0121	CIFXD CER 0.1UF 50 VDCW	
A6C7	0170-0064	C:FXD MY 0.47UF 10% 100VDCW	
A6C8	0170-0079	CIFXD MY 0.047UF 20% 50VDCW	
A6CR1	1901-0096	SEMICON DEVICE DIODE SILICON	
A6CR2	1901-0096	SEMICON DEVICE: DIODE SILICON	
A6CR3	1901-0096	SEMICON DEVICE: DIODE SILICON	
A6CR4	1902-0050	SEMICON DEVICE DIONE SI JUNC 8.664 5%	
46Q1	1851-0017	TRANSISTOR (2N1304	
A6Q2	1851-0017	TRANSISTOR: 2N1304	i
A6Q3	1850-0062	TRANSISTOR GERMANIUM ALLOY JUNCTION	
4604	1850-0062	TRANSISTOR GERMANIUM ALLOY JUNCTION	
A6Q5	1850-0062	TRANSISTOR GERMANIUM ALLOY JUNCTION	
A6Q6	1850-0062	TRANSISTOR GERMANIUM ALLOY JUNCTION	
A6Q7	1851-0017	TRANSISTOR: 2N1304	
4698	1851-0017	TRANSISTOR: 2N1304	1
4609	1850-0062	TRANSISTOR: GERMANIUM ALLOY JUNCTION	İ
46010	1851-0017	TRANSISTOR : 2N1304	
46011	1854-0003	TRANSISTURINPN SILICON	
46012	1854-0003	TRANSISTOR INPN SILICON	
46013	1850-0065	TRANSISTOR GERMANIUM 2N1370	
46014	1850-0065	TRANSISTOR:GERMANIUM 2N1370	1
46015	1850-0065	TRANSISTOR GERMANIUM 2N1370	
46016	1854-0022	TRANSISTOR: NPN SILICON	
6017	1854-0022	TRANSISTOR NPN SILICON	
46018	1851-0017	TRANSISTOR 12N1304	1
76019	1854-0003	TRANSISTORINPN SILICON	
16020	1854-0003	TRANSISTORINPN SILICON	
16021	1854-0033	TRANSISTOR SILICON NPN	
16022	1854-0022	TRANSISTOR INPN SILICON	
16R1	0684-1031	RIFXU COMP 10K OHM 10% 1/4W	
16R2	0683-4725	R:FXD COMP 4700 OHM 5% 1/4W	
6R3	0727-0124	RIFXD DEPC 3000 OHM 1% 1/2W	
i			l l

Reference Designation Index, General (cont'd) Table 6-1

Reference Designation		Description #		Note
Debigiant				
	0403 0115	DIEVO COMP DIO OHM EG 1700		
46R4	0683-9115	RIFXD COMP 910 OHM 5% 1/4W		
46R5	0727-0163	R#FXU DEPC 11.88K OHM 1% 1/2W		
46R6	0727-0173	RIFXD DEPC 20K OHM 1% 1/2W		
46R7	0684-1011	R:FXD COMP 100 OHM 10% 1/4#		
		R:FXD DEPC - 10 -1K OHM 1% 1/2#		
16R8	0727-0158	ROPAD DEFCO 1001R OFFE 12 1/2#	İ	
46R9	0727-0136	RIFXD DEPC 5.03K OHM 1% 1/2W		
A6R10	0683-6825	RIFXD COMP 6800 OHM 5% 1/4W		
46R11	0683-4735	RIFXD COMP 47K OHM 5% 1/4W		
	0683-4725	RIFXD COMP 4700 OHM 5% 1/4W		
46R12 46R13	0683-2735	RIFXD COMP 27K OHM 5% 1/4W		
A6R14	0683-1835	RIFXD COMP 18K OHM 5% 1/4W		
46R15	0683-2725	RIFXD COMP 2700 OHM 5% 1/4W		
46R16	0684-4711	RIFXD COMP 470 OHM 10% 1/4W		
A6R17	0683-3025	RIFXD COMP 3000 OHM 5% 1/4W		
46R18	0683-1635	RIFXD COMP 16K OHM 5% 1/4W		
44310	0403 2035	DIEVO COMP SAME OUM EN TAND	}	
A6R19	0683-2435	RIFXD COMP 24K OHM 5% 1/4W RIFXD COMP 24K OHM 5% 1/4W		
A6R20	0683-2435			
46R21	0683-1635	RIFXU COMP 16K OHM 5% 1/4W		
46R22	0683-1025	RIFXD COMP 1000 OHM 5% 1/4#	į	
46R23	0683-2025	R:FXD COMP 2000 OHM 5% 1/4W		
46R24	0683-1825	R:FXD COMP 1800 OHM 5% 1/4%		
		RIFXU COMP 12K OHM 5% 1/4W	1	
46R25	0683-1235			
46R26	0683-1535	RIFXD COMP 15K OHM 5% 1/4W		
46R27	0727-0158	R:FXU DEPC. 10.1K OHM 1% 1/2W		
46R28	0727-0136	R#FXD DEPC 5.03K OHM 1% 1/2W		
44330	3100-0910	DIVAD COME STREET OUM I THE SOR I VIII		
46R29	2100-0910	RIVAR COMP 2X35K OHM LIN 20% 1/4W		
46R30	2100-0910	RIVAR COMP 2X35K OHM LIN 20% 1/4W		
46R31	2100-0910	RIVAR COMP 2X35K OHM LIN 20% 1/4W		
46R32	2100-0910	RIVAR COMP 2X35K OHM LIN 20% 1/4W		
46R33	2100-0910	RIVAR COMP 2X35K OHM LIN 20% 1/4W		
A & D 3/1	2100-0910	RIVAR COMP 2X35K OHM LIN 20% 1/4W		
46R34				
46R35	0758-0051	RIFXD MET FLM 43K OHM 5% 1/2W		
46R36	0683-2225	RIFXD COMP 2.2K OHM 5% 1/4W		
46R37	0683-6825	RIFXD COMP 6800 OHM 5% 1/4W		
A6R38	0758-0002	RIFXU MET FLM 560 OHM 5% 1/2W		
14070	0759 007	DIEVO MET SIM OLOG ONN SA A COM		
16R39	0758-0038	RIFXD MET FLM 9100 OHM 5% 1/2W	į	
46R40	0683-1205	R:FXD COMP 12 OHM 5% 1/4W	1	
46R41	0683-2225	RIFXD COMP 2.2K OHM 5% 1/4W		
46R42	0683-1535	RIFXD COMP 15K OHM 5% 1/4W		
A6R43	0683-1025	R&FXD COMP 1000 OHM 5% 1/4W	For Senion Manual Co.	
	0.03.055	DAGNO CAND VIEW OVER THE ANDRESS	For Service Manuals Contact MAURITRON TECHNICAL SERVICES	
16R44	0683-2735	RIFXD COMP 27K OHM 5% 1/4W	8 Cherry Tree Rd, Chinnor	
16R45	0683-4725	RIFXD COMP 4700 OHM 5% 1/4W	Oxon OX9 4QY	
16R46	0683-1215	RIFXU COMP 120 OHM 5% 1/4W	Tel: 01844-351604 Carris 04044 Carris	
A6R47	0758-0005	RIFXD MET UX 4700 OHM 5% 1/2#	Tel:- 01844-351694 Fax:- 01844-352554	ļ
16R48	0758-0057	RIFXD MET FLM 5600 OHM 5% 1/2W	Email:- enquiries@mauritron.co.uk	
	0727 0102	BIEVO AERO IL DE OUM LO LACE		
16R49	0727-0189	RIFXD DEPC 41.7K OHM 1% 1/2W		
46R50	0687-1041	R:FXO COMP 100K OHM 10% 1/2W		
16R51	0758-0034	R:FXD MET FLM 2400 OHM 5% 1/2W	Ĭ	
46R52	0758-0004	RIFXD MET FLM 2700 OHM 5% 1/2W		
16R53	0758-0043	R#FXU MET FLM 1800 OHM 5% 1/2W		
46R54	2100-0144	R:VAR COMP 250K UHM 30% LIN 2/5%		

Table 6-1. Reference Designation Index, General (cont'd)

Reference Designation	₩ Stock No.	Description #	Note
A6R55	0758-0034	R*FXD MET FLM 2400 OHM 5% 1/2W	
A6R56	0727-0189	RIFXD DEPC 41.7K OHM 1% 1/2W	
A6R57	0758-0022	RIFXD MET FLM 82K OHM 5% 1/2W	
A6R58	0758-0044	R*FXD MET FLM 2200 OHM 5% 1/2W	
A6R59	0758-0012	R*FXD MET FLM 12K OHM 5% 1/2m	
	0.30 0012	FACTORY SELECTED PARTATURES AND AND AND AND AND AND AND AND AND AND	
		FACTORY SELECTED PARTITYPICAL VALUE GIVEN	
A6R60	0758-0044	RIFXU MET FLM 2200 OHM 5% 1/2W	
A6R61	0683-3925	RIFXD COMP 3900 OHM 5% 1/4w	
A6R62	0683-2735	RIFXD COMP 27K OHM 5% 1/4W	
A6R63	0683-5125	RIFXD COMP 5180 OHM 5% 1/4W	
A6R64	0683-1235	RIFXD COMP 12K OHM 5% 1/4W	1
A6R65	0683-1525	RIEVO COMP LEGO DUM EN A MIL	1
A6R66	1 <u> </u>	R:FXD COMP 1500 OHM 5% 1/4W	1
A6R67	0683-3335	RIFXD COMP 33K OHM 5% 1/4W	
	0683-3335	RIFXD COMP 33K OHM 5% 1/4W	İ
A6R68	0683-3935	RIFXD COMP 39K OHM 5% 1/4W	
A6R69	2100-0092	RIVAR COMP 10K OHM 20% LIN 1/5W	
A6R70	0683-1345	RIFXU COMP 130K OHM 5%A1/4W	
A6R71	0683-2025	RIFXD COMP 2000 OHM 5% 1/4W	
A6R72	0683-5125	RIFXD COMP 5100 OHM 5% 1/4W	
A6RT1	0852-0021	RITEMPERATURE COMPENSATING 100 OHM 5% 250	
A6TB1	00851-2005	BLANK PC BUARD. SWEEP & HORIZ AMPL	
A6W1	00851-6032	CABLE ASSY.	
A6#2	00851-6036	8-IN COAX:SYNC SWITCH-SWEEP AMPL A6 CABLE ASSY:HORIZ. OUTPUT TO CRT 12-INCH CABLE TERM. W/FEMALE PIN	
A7	00851-6019	BOARD ASSY VERT. AMPL. ASSY.	
A7C1	0150-0096	CAEND CER O OFFICE ADDITION	
A7C2	0160-0174	C:FXD CER 0.05UF 100VDCW	
A7C3	i	CIFXD CER 0.47UF 80% 25VDCW	
A7C4	0170-0086	CIFXD MY U-22UF 20% 50VDCW	
A7C5	0170-0083 0140-0160	CIFXD MY 0.022UF 20% 50VDCW	
A. 65	0140-0100	CIFXD MICA 3400 PF 5% 500 VDCW	
A7C6	0170-0084	C : FXD MY 0 + C68UF 20% 50VDCW	
A7C7	0170-0018	CIFXD MY 1UF 5% 200VDCW	
A7C8	0150-0121	CIFXD CER 0.18F 50 VDCW	ļ
A7C9	0160-0174	C:FXD CER 0.47UF 80% 25VDCW	1
A7CR1	1002 0005		
A7CR2	1902-0025	SEMICON DEVICE: DIODE SILICON	1
	1901-0096	SEMICON DEVICE DIODE SILICON	
47CR3	1902-0017	SEMICON DEVICE DIODE SI	1
17CR4	1901-0025	SEMICON DEVICE DIODE JUNCTION	İ
A7CR5	1901-0025	SEMICON DEVICE DIODE JUNCTION	1
17CR6	1901-0025	SEMICON DEVICE:DIODE JUNCTION	
A7CR7	1901-0025	SEMICON DEVICE DIOUS HINDS	
17CR8	1901-0025	SEMICON DEVICE:DIODE JUNCTION	
A7CR9	1901-0025	SEMICON DEVICE:DIODE JUNCTION	•
7CR10	1901-0033	SEMICON DEVICE:DIODE 1N629	1
7CR11	1901-0033	SEMICON DEVICE: DIODE SILICON 1N485 B	
17K1	0490-0125	SEMICON DEVICE:DIODE SILICON 1N485 B	
	0470-0123	RELAYIDPDT 1/4A 24VDCICOIL 24VDC	
	L		

Table 6-1. Reference Designation Index, General (cont'd)				
Reference Designation	♠ Stock No.	Description #	Note	
A7K2	0490-0125	RELATIUPDT 1/4A 24VDCICOIL 24VDC		
A7K3	0490-0125	RELAYIUPDT 1/4A 24VDCICOIL 24VDC		
A7K4	0490-0125	RELAYEDPDT 1/4A 24VDCICOIL 24VDC		
A7L1	9140-0137			
A7L2	9140-0137	COIL FXD RF 1 MH		
A7L3	9140-0137	COIL: FXD RF 1 MH COIL: FXD NF 1 MH		
A7Q1	1854-0005	TOANGEEWONEANTOA		
A702	1854-0005	TRANSISTORIZNOOB NPN SILICON		
A7Q3	1854-0022	TRANSISTOR 2N708 NPN SILICON		
A704	1854-0005	TRANSISTORINPN SILICON TRANSISTORIZNOS NPN SILICON		
A7Q5	1854-0022	TRANSISTOR INPN SILICON		
A7Q6	1854-0022			
A707	1854-0022	TRANSISTOR:NPN SILICON TRANSISTOR:NPN SILICON		
A798	1854-0005	TRANSISTORIAN SILICON TRANSISTORIZNOS NPN SILICON		
A709	1854-0005	TRANSISTOR: 2N708 NPN SILICON	1	
A7R1	0686-6225		j	
A7R2	0686-5115	RIFXD COMP 6200 OHM 5% 1/2W RIFXD COMP 510 OHM 5% 1/2W	İ	
A7R3	0686-2025	R:FXD COMP 2000 OHM 5% 1/2W		
A7R4	0686-5625	R*FXD COMP 5600 OHM 5% 1/2W		
A7R5	0686-1825	R:FXÚ COMP 1800 OHM 5% 1/2W	[	
A7R6	0683-5625	R*FXD COMP 5600 OHM 5% 1/4W	•	
A7R7	0686-7525	R:FXD COMP 7500 OHM 5% 1/2W		
A7R8	0687-3931	R:FXD COMP 39K OHM 10% 1/2W		
A7R9	0690-1231	RIFXD COMP 12K OHM 10% 1W		
A7R10	0686-2725	RIFXD COMP 2700 OHM 5% 1/2W		
A7R11	0761-0074	RIFXO MET UX FLM 15K OHM 5% 1W	ļ	
A7R12	0758-0024	RIFXD MET FLM 100 OHM 5% 1/2W		
A7R13	0758-0005	RIFXD MET OX 4700 OHM 5% 1/2W		
A7R14	0758-0017	RIFXD MET FLM 1500 OHM 5% 1/2W		
A7R15	2100-0154	RIVAR COMP 1000 OHM 30% LIN 0.15W		
A7R16	0758-0005	R*FXD MET OX 4700 OHM 5% 1/2W		
A7R17	0758-0024	RIFXD MET FLM 100 OHM 5% 1/2W		
A7R18	0761-0074	RIFXD MET OX FLM 15K OHM 5% 1W		
A7R19	0758-0003	RIFXO MET FLM 1000 OHM 5% 1/2#		
A7R20	0758-0038	RIFXD MET FLM 911K OHM 5% 1/2W		
A7R21	0690-2721	R#FXD COMP 2700 OHM 10% 1W		
A7R22	2100-0095	RIVAR COMP 100K OHM 30% LIN 1/5W		
A7R23	0683-2225	R:FXD COMP 2.2K OHM 5% 1/4w	1	
A7R24	0683-1015	RIFXD COMP 100 OHM 5% 1/4W		
A7R25	0684-1001	R:FXD COMP 10 OHM 10% 1/4W		
A7TB1	00851-2013	BLANK PC BUARD: VERT. AMP.		
A7#1	00851-6037	CABLE ASSEMBLY: VERTICAL OUTPUT TO CRT	For Service Manuals Contact MAURITRON TECHNICAL SERVICES	
	5 - 55,	4-INCH CABLE TERM. W/FEMALE PIN	8 Cherry Tree Rd, Chinnor	
A8	00851-6001	HV POWER SUFPLY ASSY	Oxon OX9 4QY Tel:- 01844-351694 Fax:- 01844-352554	
4861			Email:- enquiries@mauritron.co.uk	
A8C1 A8C2	0170-0064	C:FXD MY 0.47UF 10% 100VDCW		
A8C3	0180-0104 0150-0036	CIFXD ELECT 200UF 15VDCW		
A8C3	5040-0400	CIFXD CER 470 PF 20% 6KV SUPPORTICAPACITOR		

<sup>#</sup> See list of abbreviations in introduction to this section

Table 6-1. Reference Designation Index, General (cont'd)

		e 6-1. Reference Designation Index, General (cont'd)	
Reference	⊕ Stock No.		
Designation	Stock No.	Description #	Note
A8C4	0150-0036	C:FXD CER 470 PF 20% 6KV	
,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	5040-0400	SUPPORT: CAPACITOR	
A8C5	0160-0151	CIFXD CER 4700PF +80%-20% 4000VDCW	
	5040-0401	SUPPORTICATACITOR	
A8C6	0160-2054	C:FXD MY 0.015 UF 10% 3000VDCW	
		OWNER OF TON SOUGHOU	
A8C7	0160-2054	C 1 FXU MY 0.015 UF 10% 3000 VDCW	
A8C8	0180-0089	CIFXD ELECT 10UF-108+100% 150VUCW	
ABC9	0160-0151	CIFXD CER 4700PF +80%-20% 4000VDCW	
	5040-0401	SUPPORTICAPACITOR	
A8C10	0150-0023	CIFXD CER 2000PF 20% 1000VDC%	
40011	21.2.22		
A8C11	0160-2054	C:FXD MY 0.015 UF 10% 3000VDCW	
A8CR1	1001-0143	DECTIFICALLY TOOL	
ABCR2	1901-0142	RECTIFIERISILICON	
ABCR3	1901-0142	RECTIFIER:SILICON	
A8CR4	1901-0142 1901-0142	RECTIFIERISILICON	
AUCKT	1901-0142	RECTIFIER: SILICON	
A8L1	9140-0051	COILIFXD 400 UH	
_ <del>-</del>		1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	}
ASR1	0687-6801	R#FXD COMP 68 OHM 10% 1/2W	
A8R2	0687-1521	RIFXD COMP 1500 OHM 10% 1/2W	
A8R3	0687-2741	R:FXD COMP 270K OHM 10% 1/2#	
A8R4	0687-1231	RIFXO COMP 12K OHM 10% 1/2W	
ABR5	0690-3951	RIFXD COMP 3.9 MEGOHM 10% 1W	]
1004	0400 7054		Ì
A8R6 A8R7	0690-3951	RIFXD COMP 3.9 MEGOHM 10% 1W	
ABR8	0690-3951	RIFXD COMP 3.9 MEGOHM 10% 1W	
ABR9	0690 <b>-</b> 3951 0690 <b>-</b> 3951	RIFXU COMP 3.9 MEGOHM 10% 1W	
ABRIO	0687-1031	RIFXU COMP 3.9 MEGOHM 10% 1W	
A01.10	0087-1031	RIFXD COMP 10K OHM 10% 1/2W	
A8R11	0687-1051	RIFXU COMP 1 MEGOHM 10% 1/2W	
A8R12	0687-4731	R*FXD COMP 47K OHM 10% 1/2%	}
A8R13	0690-1851	RIFXU COMP 1.8 MEGOHM 10% 1W	
A8R14	0690-1851	RIFXO COMP 1.8 MEGOHM 10% 1W	
A8R15	0690-1851	RIFXD COMP 1.8 MEGUHM 10% 1W	
		100 11	
A8R16	0690-1851	RIFXD COMP 1.8 MEGOHM 10% 1W	
A8R17	0690-8241	RIFXD COMP 620K OHM 10% 1W	
		1	
A8T1	9120-0092	TRANSFORMEK:AUDIO	
		STEP-UP	
AOTDI	0000. 0004		
A8TB1	00851-2006	BLANK PC BOARD HV POWER SUPPLY	
49	00851-6017	LOW VOLTAGE BOACH CHOOLY ACCOUNT	
	20071-0017	LOW VOLTAGE POWER SUPPLY ASSY.	
4901	0180-0089	C*FXD ELECT 100F-10%+100% 150VDCW	
4902	0180-0138	C:FXU ELECT 100 UF -10+100% 40VDCW	
1903	0180-0049	CIFXD ELECT 20UF SOVDCW	
1904	0170-0042	C:FXD MY 0.33UF 5% 100VDCW	
N9C5	0170-0064	CIFXD MY 0.47UF 10% 100VDCW	
			-
1906	0180-0097	CIFXD ELECT 47 UF 10% 35VDCW	
1907	0170-0064	C:FXD MY 0.47UF 10% 100VDCW	-
\9C8	0160-0163	CIFXD MY 3300PF 10%	
1969	0180-0119	C:FXD ELECT 1UF -10+100% 25VUCW	1
V9C10	0180-0097	CIFXD ELECT 47UF 10% 35VDCW	
1			
			ı
			1

Table 6-1. Reference Designation Index, General (cont'd)

			No
A9CR1	1901-0029	SEMICON DEVICE DIODE SI 600V	
49CR2	1902-0241	SEMICON DEVICE: DIODE SILICON 100V 5%	i
A9CR3	1901-0045	SENTION DEVICE DIODE SILICON 100V 5%	i
49CR4		SEMICON DEVICE DIODE SILICON	
	1901-0025	SEMICON DEVICE: DIODE JUNCTION	
19CR5	1901-0045	SEMICON DEVICE DIODE SILICON	
19CR6	1901-0025	SEMICON DEVICE: DIODE JUNCTION	
19CR7	1902-0025	SEMICON DEVICE DIODE SILICON	
A9CR8	1902-0017	SENTION DEVICE DIODE SILICON	Ì
		SEMICON DEVICE DIODE SI	į.
A9CR9	1901-0025	SEMICON DEVICE: DIODE JUNCTION	i
19CR10	1901-0025	SEMICON DEVICE DIODE JUNCTION	
9CR11	1902-0017	SEMICON DEVICE:DIODE SI	
9CR12	1901-0049	SENTENN DEVICE DIODE SI	1
9CR13		SEMICON DEVICE: DIODE SILICON	
	1901-0049	SEMICON DEVICE: DIOUE SILICON	1
9CR14	1901-0049	SEMICON DEVICE: DIODE SILICON	Į.
9CR15	1901-0049	SEMICON DEVICE: DIODE SILICON	
9CR16	1902-0017	SEMICON DEVICE: DIODE SI	
9CR17	1901-0025	SEMICON DEVICE TO LODGE TO THE TENT	1
9CR16	1901-0025	SEMICON DEVICE DIODE JUNCTION	}
9CR19		SEMICON DEVICE : DIODE JUNCTION	J
	1901-0025	SEMICON DEVICE: DIODE JUNCTION	
9CR20	1901-0025	SEMICON DEVICE DIODE JUNCTION	
901	1850-0040	TRANSISTOR:GERMANIUM 2N383 PNP	
902	1850-0065	TRANSISTUR : GERMANIUM 201370	
903	1854-0003	TRANSISTOR: NPN SILICON	
904	1850-0040	THE NETSTORIE TRANSPORT	
905	1854-0003	TRANSISTUR:GERMANIUM 2N383 PNP TRANSISTOR:NPN SILICON	
996			
<b>,40</b>	1854-0003	TRANSISTOR: NPM SILICON	
PR1	0758-0012	RIFXU MET FLM 12K OHM 5% 1/2W	
9R2	0757-0817	RIFXU MET FLM 750 UHM 1% 1/2W	
PR3	0687-1231	RIFXU COMP 12K OHM 10% 1/2W	
9R4	0699-0005	R:FXU COMP 2.7 OHM 10% 1W	
R5	0761-0016	RIFXD MET FLM 7500 OHM 5% 1W	
104			
PR6	0687-2211	RIFXU COMP 220 OHM 10% 1/2%	
PR7	0687-1011	R#FXU COMP 100 OHM 10% 1/2%	
R8	0687-5611	RIFXD COMP 560 OHM 10% 1/2W Fors	Service Manuals Contact
R9	0687-4721	R:FXD COMP 4700 OHM 10% 1/2W FORS	POVICE MAINUAIS CONTROL
R10	0687-3311	RIFXD COMP 330 OHM 10% 1/2% MAURI	TRON TECHNICAL SERVICES herry Tree Rd, Chinnor
R11	0607330+	1	Oven OX9 4OY
	0687-3321	RIFXD COMP 3300 OHM 10% 1/2W Tel-018	344-351694 Fax:- 01844-352554
R12	0687-1021	KIFXU COMP 1000 OHM 10% 1/2W Fmai	I:- enquiries@mauritron.co.uk
R13	0812-0027	R&FXD WW 3100 OHM 5% 3W	
R14	2100-0154	RIVAR COMP 1000 OHM 30% LIN U.15#	
R15	0812-0027	R:FXD WW 3100 OHM 5% 3W	
R16	0687-3321		
R17		R#FXO COMP 3300 OHM 10% 1/2W	
	0687-6801	R*FXD COMP 68 OHM 10% 1/2W	
R18	0687-6801	RIFXD COMP 68 OHM 10% 1/2W	
R19	0687-1521	RIFXD COMP 1500 OHM 10% 1/2W	
R20	0687-3921	R*FXD COMP 3980 OHM 10% 1/2W	
R21	0811-0040	RIFXD WW 1 OHM 1% 5W	
R22	0687-5611	RIFYO COMP SAO ONE LOS ACO.	
R23	0686-1225	R*FXO COMP 560 OHM 10% 1/2W	[
R24	_	R*FXD COMP 1200 OHM 5% 1/2W	1
14.7	0686-4725	R:FXD COMP 4700 OHM 5% 1/2W	

<sup>#</sup> See list of abbreviations in introduction to this section

Table 6-1. Reference Designation Index, General (cont'd)

Designation	Stock No.	Description #		Note
A9R25	0687-8221	R#FXD COMP 6200 OHM 10% 1/2W		•
A9R26	0687-5621	R:FXU COMP 5600 OHM 10% 1/2W		
A9R27		RALAD COMP DOOD OHM TOW 1/5#		
	0687-8221	R # F X D C OMP 6200 OHM 10% 1/2W		
A9R28	0758-0047	RIFXD MET FLM 7500 OHM 5% 1/2W		
A9R29	2100-0154	RIVAR COMP 1000 OHM 30% LIN 0-15W		
A9R30	0758-0047	R:FXD MET FLM 7500 OHM 5% 1/2#		
A9TB1	00851-2004	BLANK PC BOARD LOW VOLT POWER SUPPLY		
A10	00851-6039 5040-0218	SWITCH ASSY: SWEEP TIME COUPLING: MECHANICAL P/O SWEEP TIME SWITCH		
A10C1	0170-0038	C:FXD MY 0.22UF 10% 200VDCW		
A10C2	0170-0051	CIFXD MY U.635UF 5% 100VDC#		
A10C3	0180-0101	C:FXD ELECT TA 1.8UF 10% 35VDC#	i	
A10C4	0180-0116	CIFAD ELECT TA 6.8UF 10% 35VDCW	1	
A10C5	0180-0233	C:FXU TA ELECT. 20UF +20-15% 60VUC#		
A10C6	0180-0235	C:FXD TA ELECT. 56UF 20% 75VUC#		
A10C7	0170-0042	C:FXU MY 0.33UF 5% 100VDCW		
A10C8	0180-0230	CIFXD TA ELECT. 1UF 20% 50VDCW		
A1009	0180-0231	C:FXU TA ELECT. 3.5UF +20-154 75VDCW		
A10C10	0180-0232	CIFXU TA ELECT. 10UF 20% 100VDCW		
A10C11	0180-0234	CIFXD TA ELECT. 33UF 20% 75VDCW		
A10C12	0180-0113	C:FXD ELEC) TA 100UF +20-15% 30VUCW		
A10R1	2100-0107	R:VAR COMP 50K OHM 30% 1/3W, VERNIER		
		KNOB : VERNIER FOR SWEEP TIME SWITCH	1	
A10R2	0370-0114 0757-0831	RIFXD MET FLM 4.32K OHM 1% 1/2W		
A1051	3100-1500	SWITCH ROTARY		
	0370-0113	KNOB.		
		SWEEP TIME		
	3130-0041	SWITCH SHIELD		
A11	00851-6006	SWITCH ASSY.: VERT. DISPLAY		
A11C1	0160-0134	CIFXD MICA 220PF 5% 300VDCW		
A11C2	0160-0178	CIFXD MICA 27PF 5% 300VDCW	l	
A11C3	0150-0093	C*FX0 CER U.01UF +80-20 100VDCW		
A11C4	0150-0093	C:FXD CER 0.01UF +80-20 100VDCW		
A11C5	0150-0093	C:FXD CER 0.01 UF +80-20% 100V0CW	For Service Manuals Conta	ıct
Alicri	1901-0047	SEMICON DEVICE:DIODE JUNCTION	MAURITRON TECHNICAL SERVING 8 Cherry Tree Rd, Chinno	PEO
AllCR2	1901-0047	SEMICON DEVICE DIODE JUNCTION	Oxon Oxo any	
Alicr3	1901-0047	SEMICON DEVICE DIODE JUNCTION	Tel:- 01844-351694 Fax:- 01844-35	DEE4
Alicr4	1901-0047	SEMICON DEVICE DIODE JUNCTION	Email:- enquiries@mauritron.co.u	:054 k
AllL1	9140-0118	COIL:FXD 500 UH 5%		
A11Q1	1853-0003	TRANSISTORIPHP SILICON F 50MC MIN		
A1102	1853-0003	TRANSISTORIPHP SILICON F 50MC MIN		
Aliri	0683-1025	RIFXD COMP 1000 OHM 5% 1/4W		
A11R2	2100-0958	RIVAR COMP 10K OHM 20% 0.5W	1	
A11R3	0727-0405	RIFXD DEPE 57.46K OHM 1/2% 1/2W		

Table 6-1. Reference Designar and referal reference

		Table 6-1. Reference Designar and referal re-
Hotelesse		
- Harden '	·	·
		Manufacture and an analysis of the second analysis of the second analysis of the second and an analysis of the second and an analysis of the second and an analysis of the second and an analysis of the second and an analysis of the second and an analysis of the second and an analysis of the second and an analysis of the second and an analysis of the second and an analysis of the second and an analysis of the second and an a
AllR4	1	NOT ASSIGNED
Alir5	07,	RIFXU JEFC IBK OMM 12
AliR6 AliR7	0757 - 557	MITAU MET FEM 32.4K OHM 18 1.77
ALIRE	0757	PITAU TEI FUM 191K OHM IN 12/10
******	0757-0885	RIFXU MET FLM 15.4K CHM IN 1748
AllR9	0757-0889	
ALIRIO	0683-2025	
ALIRII	0683-3625	RIPAU COMP 2000 OHM SE 1244
AllR12	0727-0123	MIPAU COMP 3600 OHM SE 1244
ALIRIJ	2100-0957	HIF XD DEPC 2900 OHM LK 1/2m
	1100-0757	HIVAR COMP SK OHM 20% LIN 1/20#
A11R14	2100-0957	
411R15	0727-0403	
411R16	0727-0126	THE OF THE PROPERTY OF THE PRO
411R17	0727-0178	
A11R18	0727-0101	
		14 175#
A11R19	0727-0398	RIFXU UEPC 3.79K OHM 1/2% 1/2%
AllR20	2100-0956	RIVAR COMP 500 OHM 20% LIN 1/20%
; A1101		200 OUT 504 CIM 1/508
A1151	3100-0815	RCTARY SWITCHTZ-SECTION 3-POSITION
	0370-0112	KNOBIVERTICAL DISPLAY
A11#1	00.5	
	00921-0033	
		17-IN COAREVERT IN (C) TO ALISE
	00851-0017	
A12	00851-6035	COVERTINPUT OF FILTERIALS
	00031-0033	PC BUARD ASSYLLAPUT BANDPASS FILTER
A12C1	0130-0017	
A12C2	0140-0194	CIVAR CEN 8-50 PF N750
A12C3	0130-0017	CIFYU MICA 110 PF 5% 300VDCW
A12C4	0160-0178	CIVAR CEH 8-50 PF N755
A12C5	0140-0197	CIFXD MICA 180 PF 5%
		011 VO WICK 100 PF 2%
A12C6	0140-0204	CIFXU MICA 47 PF 54 NPO 500VUC
		FACTORY SELECTED PARTITYPICAL VALUE SIVEN
A12.0		
A12J1 A12J2	1250-0212	CONNECTORIJACK CHASSIS BNC
~	1250-0149	CONNECTORING JACK CHASSIS RIGHT ANGLE
A12L1	00051-00-1	
A12L2	00851-8001	COILIRE FXD O. JUN
	00031-0002	COILIRE VAR D. JUH MAX
A12TB1	00851-2016	M() ( 0 (3 * 1 to 0 to 1 * a
	2010	BOARDIINPUT B.P. FILTER
C1	0150-0119	CIEXO CEN 1110 OL OL OL
C2		CIFXD CEM 2X(0.01 UF) 20% 250VUC# NOT ASSIGNED
C3	0180-0042	CIFXD ELECT 1200F 350VDCW
	1520-0001	PLATE I MOUNTING ELECTROLATIC CAPACITOR
C4	0180-0047	CIFXD ELECT DOUGH TOVOCH
	1520-0001	PLATE I MOUNTING ELECTHOLYTIC CAPACITOR
C5	0180-0047	しょとえい たにたにし ちひかいた チャッカの 裏
C6	1520-0001	PLH ETMOUNTANG ELECTROLYTTO CAPACITON
C7	0150-0121	しょうさい したさ シャネーンと キャンカー20g もりしりたき
	0180-0059	CIFAD ELEC 10 OF +1.04108 25VUCH
CB	0180-0098	
l	UU-UU78	CIFXO ELECT 100 OF 20% 20VUC#
1		
ĺ		
1		

<sup>#</sup> See list of abbreviations in introduction to this section

Table 6-1. Reference Designation Index, General (cont'd)

Reference Designation	⊕ Stock No.	Description #	Note
051	2140-0018	LAMP&GLOW 1/10W	
	5040-0234	LAMPHOLDER	
	5040-0235	BASE: LAMPHOLDER	1
F1	2110-0016	FUSE CARTRIDGE 0.6 AMP SLO-BLO	
	2440 00##	115V OPERATION	
F1	2110-0044	FUSE:CARTRILGE 0.3 AMP SLO-BLO 230V OPERATION	
J1	-	1.F. INPUT, PART OF CABLE W1	
J2	1250-0171	CUNNECTORIENC JACK SWEEP INPUT	
J3	1250-0171	CONNECTORIBNE JACK	
<b>پ</b> ر	1251-0148	BLANKING INPUT CONNECTOR: FOWER 3 PIN MALE	
		LINE INPUT	
J5	1250-0083	CONNECTORIBAC	
	1200-0081	1.F. TEST POINT BUSHING:INSULATOR NYLON	
	1250-0053	CAPICONNECTOR BNC WITH CHAIN	
J6		SYNC INPUT: PART OF CABLE W2	
J <b>7</b>		SWEEP OUTPUT: PART OF CABLE W3	
J8_	1250-0171	CONNECTORIBNE JACK, HORIZ OUTPUT	
J <del>9</del>	1251-0143	CONNECTOR: FEMALE 14 CONTACTS, CONTROL	
19	1250-0171	CONNECTOR & BNC JACK, VERT OUTPOT	
_1	9140-0082	COILIFAD RF 15 UH	
_2	9140-0082	COIL:FXD RF 15 UH	
L3	5060-0409	COIL:ALIGNMENT   IN CRT TRACE ALIGN CIRCUIT	ļ
_4	9110-0042	CHUKE:FILTER 70 MH 1.0 AMP 1.5 OHM	
MP1	175A-83A	RETAINER: CHT. SHIELD	
MP2	00851-0006	BRACKET:PU#ER SUPPLY	
MP3	00851-0007	SHIELD : HIGH VOLTAGE	
MP4 MP5	00851-0008	COVER:SOCKET BRACKET:SWEEP AND HORIZ AMPL PCBD A6	
MF3	00851-0009	BUNCKET AND HOUSE AND TO THE	
MP6	00851-6008	CRT. SHIELD ASSEMBLY	
MP7	0905-0050	GASKET*FELT BLACK 5/16 INCH THICK EXTRUSION*RUBBER For Service Manuals Contact	
4P8 4P9	4320-0007 5000-0408	COIL BRACKET MAURITRON TECHNICAL SERVICE	<b>S</b>
		8 Cherry Tree Rd, Chinnor	
<b>91</b>	1850-0090	Tel:-01844-351694 Fax:-01844-3525	554
	1200-0092	BUSHING: TRANSISTOR INSULATOR: TRANSISTOR Email: enquiries@mauritron.co.uk	
22	1850-0090	TRANSISTOR GERMANIUM 2N11838 PNP	
<del></del>	1200-0092	BUSHING:TRANSISTOR	
	1200-0076	INSULATOR: TRANSISTOR	
43	1850-0098 1200-0043	TRANSISTOR:GERMANIUM PNP SELECTED INSULATOR:TRANSISTOR ANODIZEU ALUMINUM	
24	1850-0064	TRANSISTOR GERMANIUM 2N1183 PNP	
<b>=</b> 7	1200-0076	INSULATOR : TRANSISTOR	
	1200-0087	CLAMPITRANSISTOR	
	1200-0081	BUSHING: INSULATOR NYLON	

Designation	₩ Stock No	Description #	Note
0.5			11000
Q5	1850-0064	TRANSISTUR GERMANIUM 201183 PNP	
	1200-0076	INSULATOR: TRANSISTOR	
_	1200-0087	CLAMPITRANSISTOR	
96	1850-0098	TRANSISTOR ( C BUANTAN	1
	1200-0045	TRANSISTOR GERMANIUM PNP SELECTED	1
		INSULATORITHANSISTOR ANODIZED ALUMINUM	
R1	2100-0893		İ
	0370-0127	RIVAR 2K(FRONT)750K(REAR)0HM LIN 1/2W	l
	5040-0421	1 WHOD IN LINGTH	
R2	2100-0027	INSULATORIPOTENTIOMETER	
	-100-0027	RIVAR COMP 10K OHM 10% LIN 2W	
R3	2100-0189	I ANT LEVEL ALU	1
	2100-0189	RIVAR COMP 1 MEGOHM 30% LIN 1/4W	
		ASTIG. ADJ.	
R4	3100 00		
•••	2100-0218	RIVAR COMP 1.2 MEGOHM 20% LIN 24	
	07-0	FOCUS	
	0370-0026	KNOD	
k5	5040-0418	INSULATORIPOTENTIOMETER	
N.O	2100-015U	RIVAR 2-SECT TOK OHM 20% LIN 1/4W	
0.4		TRACE ALIGN ADJ.	
R6	0758-0005	RIFXD DEPC. 4700 OHM 5% 1/2W4	
o=		700 UMM 5% 1/2WW	
R7	2100-0893	RIVAR ON FEMANTA TO THE TO	1
	0370-0120	RIVAR 2K(FRONT) 750K(REAR) OHM LIN 2W KNOB: BASE LINE CLIPPER	
R8	2100-0030	BIVAS COMP 100	
		RIVAR COMP 1000 OHM LIN	
२9	2100-0067	VERTICAL PUSITION	
		RIVAR COMP 2500 OHM 10% LIN 1/2#	Į.
		l e e e e e e e e e e e e e e e e e e e	
R10	2100-0019	HORIZONTAL POSITION	
	0370-0020	RIVAR COMP 500 OHM 10% LIN 1/2%	1
211	0727-0004	I NIVO BLACK C. 750 DIA LE VERMES	
(12	0727-0004		
	0727-0004	RIFXD DEPC 5 OHM 1% 1/2W	
.1	3101-0035		
_	2101-0033	SWITCHISLIDE DPDT	
.2	0005	115V/230V	
-	00851-6040	SMITCH ASSY: SYNC	
1	0370-0112	KNOE:	
١	_	SYNC	1
,	3101-0052	SWITCH PUSHBUTTON COST NOOM	
		SWITCH PUSHBUTTON SPST, NORMALLY OPEN SINGLE SWELF	
1	9100-0274	TRANSFURMER:POWER	
		LINE	
	5083-0624	FLECTRON TIESTACE	
	· - <b>- ·</b>	ELECTRON TUEE: CATHODE-RAY P-2 PHOSPHOR	
	00851-2026	THE POSCU WITCH OPTION OF TO 31 SECONDER	
		' * E ' E N * C N   E   a   B)	
	120A-20	USED WITH F-2 PHOSPHOR	[
	5083-0634	BEZEL:CRT.	1 1
	3003 0034	ELECTRON TUBE: CATHODE-RAY P-7 PHOSPHOR	
	1201-634	1 0: 11014 07	
	120A-83A	LIGHT FILTERIAMBER	For Service Manuals Contact
		FCK CRT VI	MAURITRON TECHNICAL SERVICES
	5/107 A -		8 Cherry Tree Rd, Chinnor
1 :	5083-0654	ELECTRON TUBE: CATHODE-RAY P-31 PHOSPHOR	Oxon OXI9 4OV
1		1 0. 12014 31	Tel:- 01844-351694 Fax:- 01844-35255
1 .	120A-83G	LIGHT FILTER: GREEN	Email:- enquiries@mauritron.co.uk
ŀ		FOR CRT VI	- Aminocontaction.ux
1	\ <del>-</del> -		
(	00851-6027	CABLE ASSY. IF INPUT	1
		13-IN COAX W/BNC FEMALE JI BNC MALE PI	1
		MONO FEMALE DI BNC MALE PI	] ]
1			
]	i		
1	1		1 1
ı	I		
	1		

 $<sup>\</sup>ensuremath{\,\#\,}$  See list of abbreviations in introduction to this section

Table 6-1. Reference Designation Index, General (cont'd)

Table 6-1. Reference Designation Index, General (cont'd)				
Reference Designation	֍ Stock No.	Description #	Note	
w2 w3	00851 <b>-</b> 6016	CABLE ASSY. ISYNC INPUT 23-INCH COAX W/BNC FEMALEJACK J6		
W4		CABLE ASSY. SWEEP OUTPUT 24-IN COAX W/BNC FEMALE J7		
xF1	00140-61606 1400-0084	CABLE HIGH VOLFAGE, INCLUDES 270K RES		
XVI	1200-0084 1200-0037 1200-0085 1200-0050	FUSEHOLDER EXTRACTOR POST TYPE  SOCKETICRT TUBE COVER PLATEICRT SOCKET 14-CONTACT PINICRT SOCKET		
L				

Table 6-2. Reference Designation Index, Assembly A2

Reference	Tab.	Reference Designation Index, Assembly A2	
Designation	⊕ Stock No.	Description #	Note
	00851-6003	RF CIRCUIT ASSEMBLY	
		PREFIX DESIGNATIONS IN THIS TABLE WITH A2	
A1	00851-6025 00851-0013	INPUT SWITCHING CIRCUIT ASSY COVERIRE CKT ASSY A2 HOUSING	
A1C1	0150-0095	C:FXD CEN U-0110F +80+20 1000000	
A1C2	0150-0093	C:FXJ CER 0.01UF +80-20 100VUCW	
AIKI AIK2	0490-0125	RELAYIUPUT 1/4A 24VDCICOIL 24VDC	
AILI	9140-0146	RELAYIUPUT 1/4A 24VDC:COIL 24VUC	
A1L2	9140-0146	COIL:FXD RF 10.0 UH	
A1781	00851-2014	BLANK PC BUARD INPUT SHITCHING CIRCUIT	
A2	00851-6042	20MC AMPLIFIER ASSY	
A2C1	0150-0050	C:FXU CER 1000PF 600 VUCH	
A2C2	0150-0093	( CIFXD CER U.01UF +AC-20 IOCVDC-	
A2C3 A2C4	0150-0042	- 1 CIFXU II 4.7 PF 59 500 Vicu	
A2C5	0140-0225 0150-0093	C:FXD MICA 300PF 1% 300VDCW	ł
	0130-0093	C : FXU CER 0.01UF +80-20 100VUC#	
A2C6	0150-0093	C : FXD CER U. 01 UF +80-20 100 VDCW	
A2C7	0150-0050	CIFXD CER 1000PF 600 VDCW	
AZLI	9140-0235	COILIRF TAPPED 0.95-1.8UH	
A2L2	9140-0232	COILIRF TAPPED 0.254UH-0.50UH	
A2L3	9140-0146	COIL FXD KF 10.0 UH	
A2Q1	1850-0153	TRANSISTUR#PNP SM1642	
A2R1	0683-2725	RIFXU COMP 2700 OHM 5% 1/4W	ļ
A2R2	0683-1525	R*FXU COMP 1500 OHM 59 1/44	
A2R3	0683-6825	FACTORY SELECTED PART: TYPICAL VALUE GIVEN R:FXD COMP 6800 OHM 5% 1/4#	- 1
A2R4	0683-2215	RIFXU COMF 220 OHM 5% 1/4W	1
A2R5	0663-1615	R:FXD COMP 180 OHM 5% 1/4W	
A2T81	00851-2007	BLANK PC BUARDIAMPLIFIER (20MC)	
A3	00851-6023	FIRST 1-10KC BANDPASS FILTER ASSY.	
A3C1	0140-0175	i l	
A3C2	0121-0037	CIFXD MICA 39 PF 2% 300 VDCW CIVAR CER 7+25PF	İ
A3C3	0140-0175	CIFXD MICA 39 PF 2% 300 VDCW	ŀ
A3C4	0121-0037	CIVAR CER 7-25PF	ĺ
A3C5	0121-0033	CIVAR AIR 1.4-7.3PF	
A3C6	0160-0179	C:FXD MICA 33PF 5% 300VDCW For Service Manuals Contac	
A3C7	0150-0050	CIFXD MICA 33PF 5% 300VDCW For Service Manuals Contac CIFXD CER 1000PF 600 VDCW MAURITRON TECHNICAL SERVICE	<u>t</u>
A3CB	0150-0093	CIFXD CER U-01UF +80-20 100V0C4 8 Cherry Tree Rd Chinnel	ES
4309	0150-0093	CIFXD CER U-UIUF +30-20 106V0CW Oxon OX94OV	
43010	0150-0050	C * FXD CER 1000PF 600 VDC * Tel:-01844-351894 Fax:-01844-3528 Email:- enquiries@mauritron.co.uk	554
A3C11 A3C12	0150-0093	C * FXU CER 0.01UF +80-20 100VLC	
43012	0150-0093	C:FXD CER 0.01UF +80-20 100VUCW	
3K1	0490-0125	RELAYIDPUT 1/4A 24VDCICOIL 24VDC	
			j
	ļ		
			1

<sup>#</sup> See list of abbreviations in introduction to this section

Table 6-2. Reference Designation Index, Assembly A2 (cont'd)

Reference Designation	₱ Stock No.	Description #	Not
43K2	0490-0125	PELAY-DOOT LAND DUCCHOOT LAND	
		RELAY:DPDT 1/4A 24VDC:COIL 24VDC	
A3L1	00851-8005	CUILIRF	
A3L2	00851-8004	CCILIRF	
A3L3	9140-0235	CUILIRE TAFFED U.95-1.0UH	
A3L4	9140-0150	COIL FAU HE 2.7 UH	
A3L5	9140-0146	CUILIFAD RF 10.0 UH	İ
A3L6	9140-0146	COIL FAD RF 10.0 UH	
A3Q1	1850-0153	TRANSISTURIFNP SM1642	
A3R1	0683-2725	RIFXU COMP 2700 UHM 5% 1/4%	į.
A3R2	0683-6825	R:FAU COMP 6800 UHM 5% 1/4%	
A3R3	0683-1525	RIFXU COMP 1500 0HM 5% 1/4%	[
A3R4	0683-4305	RIFXÙ COMF 43 OHM 5% .25W	
G#11.7	3005 4505	EALTING COME TO OME DA 1208	
A3R5	0683-1015	FACTORY SELECTED PART: TYPICAL VALUE GIVEN RIFXD COMP 100 CHM 5% 1/4%	
A 3 TB 1	00851-2008	BLANK PC BUFRD: FIRST, 1-10KC, BP FILTER	
A3Y1	0410-0091	CRYSTAL: WUARTZ ZUMC	
		ABYL AND AHYL (MATCHED PAIR)	
Α4	00851-6024	SECOND 1-10 KC BANDPASS FILTER ASSY.	
A4C1	(1157, 000)	CARNA CELLA LA CARNA CAR	1
	0150-0093	C:FXD CER 0.010F +80-20 100VDCW	
A4C2	0150-0093	CIFAU CER C.01UF +80-20 100VUCH	
44C3	0150-0093	C : FXD CER 0.010F +80-20 100VDCW	
4404	0150-0050	CIFAU CER ICOUPF 600 VUC#	
4405	0121-0037	CIVAR CER 7-25FF	
A4C6	0140-0175	C:FAU MICA 39 PF 2% 300 VDCW	
A4C7	0130-0017	C:VAR CER 8-50 PF N750	
A4C8	0121-0033	CIVAR AIR 1.4-7.3PF	
A4C9	0121-0037	CIVAR CER 7-25PF	1
A4C10	0140-0175	CIFXW MICA 39 PF 2% 300 VCCW	
	0140-0175	CT AD THEA DY PP 2% SOO VICE	
A4C11	0150-0093	C:FXD CER G.01UF +80-20 100VUC#	
A4C12	0150-0093	C:FXD CER 0.010F +60-20 100VDC#	
A4K1	0490-0125	RELAYIDPUT 174A 24VUCICUIL 24VUC	
44K2	0490-0125	RELAY: OPUT 1/4A 24VDC; COIL 24VDC	
A4L1	00851-8006	(Gu) +85	
		COLLERF	
A4L2	00851-8004	COILIRF	
A4L3	9140-0146	CCILIFXD RF 10.0 UH	}
44L4 44L5	9140-0146	COLLEFXO RF 10.0 UH	
	9140-0146	COLLETAD RF 10.0 UH	
1431	1850-0153	TRANSISTURIENP SM1642	
A4R1	0663-2725	R#FXD COMP 2700 OHM 5% 1/4%	
44R2	0683-1525	RIFXU COMP 1500 OHM 5% 1/4%	
44R3	0663-6525	RIFXU COMP 6800 UHM 5% 1/4%	1
\4R4	0683 <b>~33</b> 15	RIFXU COMP 330 OHM 5% 1/4W	1
4R5	0683-1015	RIFXU COMP 100 OHM 5% 1/4W	
44781	00851-2009	BLANK PC BUARDISECOND 1-10KC BP FILTER	
14Y1	0410-0091	CRYSTAL QUARTZ 26MC	

Table 6-2. Reference Designation Index, Assembly A2 (cont'd)

Reference		ole 6-2. Reference Designation Index, Assembly A2 (cont'd)	
Designation		Description #	Not
			1100
		A3Y1 AND A4Y1 (MATCHED PAIR)	
A5	00851-6026	OUTPUT SWITCHING CIRCUIT ASSY.	
A5C1	1		
A5C2	0150-0093 0150-0093	C:FXD CER 0.01UF +80-20 100VUCW	1
AEV.		0.1 XD CER 0.010F +80-20 100VDCW	
A5K1 A5K2	0490-0125 0490-0125	RELAYIDPUT 1/4A 24VDCICOIL 24VUC	}
- · · · -	049040125	RELAY: UPDT 1/4A 24VDC: COIL 24VDC	
A5L1	9140-0146	COIL: FXD RF 10.0 UH	
A5L2	9140-0146	COIL FXD RF 10.0 UH	
ASTB1	00851-2015	BLANK PC BUARD-OUTPUT SWITCHING CIRCUIT	
A6	00851-6021		
	00831-6021	ASSY. CURRENT-CONTROLLED ATTEN.	
A6C1 A6C2	0160-0179	C:FXD MICA 33PF 5% 300VDCW	
46C3	0150-0093	CIFXU CER U_01UF _80_20 100000	
A6C4	0160-0179 0150-0095	I CTI AD MICH DOPP BY REALINATED	
A6C5	0160-0179	C:FXD CER U.01UF +80-20 100VUCW C:FXD MICA 33PF 5% 300VDCW	
4606			
4607	0140-0192	CIFXÚ MICA 68PF 5% 300VDCW	
1608	0140-0192 0160-0179	LIFAU MICA 68PE 5% 3000000	
1609	0160-0179	I COFAU MILA SIRE SA ROAMAN	
6010	0150-0093	C:FXD MICA 33PF 5% 300VDCW C:FXD CER 0.010F +60-20 100VDCW	
6C11	0160-0179		
6012	0150-0050	C:FXU MICA 33PF 5% 3UOVDCW	
6013	0150-0093	CIFXU CER ICOUPF 600 VUCW	
6C14	0140-0176	C:FXD CER U.01UF +80-20 100VDCW C:FXD MICA 100 PF 2% 300 VDCW	
6CR1	1901-0162	1	
6CR2	1901-0162	SEMICON DEVICE: DIOUE SULD IN SET OF 6	
6CR3	1901-0162	TO THE TAX OF THE TAX OF THE TAX OF T	
6CR4	1901-0162	SEMICON DEVICE: DIODE SOLD IN SET OF 6 SEMICON DEVICE: DIODE SOLD IN SET OF 6	
5CR5	1901-0162	SEMICON DEVICE: DIODE SOLD IN SET OF 6	
SCR6	1901-0162	SEMICON DEVICE : DIODE SOLD IN SET OF 6	
SL1	00851-8009	i I	
L2	9140-0149	COLLERF COLLERS DE L. COLLERS	
L3	9140-0146	COIL:FXD RF 1.86 UH CCIL:FXD RF 10.0 UH	
L4	9140-0152	COILIFXD RF 41.06 UH	
L5	9140-0149	COIL: FXD RF 1.86 UH	
L6	9140-0146	COLLEGE NE 10	
L7	9140-0149	COIL FXD RF 10.0 UH COIL FXD RF 1.86 UH For Service Manuals Contact	
LB	9140-0149	COIL FXD RF 1.66 DD MAUNITHON TECHNICAL SERVICES	
L9	9140-0149	COILIFAN RE 1-46 III	
L10	00851-8010	Oxon OX9 4QY Tel:- 01844-351694 Fax:- 01844-352554	
L11	9140-0235	COIL FRE TAPPED 0.95-1.5UH	
L12	9140-0159	COIL#FXD 0.478H 20%	
51	1850-0153		
		TRANSISTOR: FNP SM1642	
₹1	0683-2725	R#FXD COMP 2700 OHM 5% 1/4%	
	1		
I	[	· · · · · · · · · · · · · · · · · · ·	ļ

 $<sup>\</sup>pm$  See list of abbreviations in introduction to this section

Table 6-2. Reference Designation Index, Assembly A2 (cont'd)

A6R2 0683-6925 RIFAD COMP 6800 OHM 5% 1/4% A6R3 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A6R3 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A6R3 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A6R3 0683-1525 RIFAD COMP 6800 OHM 5% 1/4% A7C1 0150-6050 CIFAD CER 10:00F 600 VUC% A7C2 0150-6050 CIFAD CER 10:00F 600 VUC% A7C3 0140-6076 CIFAD CER 10:00F 600 VUC% A7C3 0140-6076 CIFAD CER 10:00F 600 VUC% A7C3 0150-6050 CIFAD CER 10:00F 600 VUC% A7C4 0150-6050 CIFAD CER 10:00F 600 VUC% A7C5 0150-6050 CIFAD CER 10:00F 600 VUC% A7C6 0150-6050 CIFAD CER 10:00F 600 VUC% A7C7 0140-6150 CIFAD CER 10:00F 600 VUC% A7C8 0150-6050 CIFAD CER 10:00F 600 VUC% A7C8 0150-6050 CIFAD CER 10:00F 600 VUC% A7C8 0150-6093 CIFAD CER 10:00F 600 VUC% A7C9 0150-6093 CIFAD CER 10:00F 600 VUC% A7C1 0150-6093 CIFAD CER 10:00F 600 VUC% A7C1 0150-6093 CIFAD CER 10:00F 600 VUC% A7C1 0150-6093 CIFAD CER 10:00F 600 VUC% A7C1 0150-6093 CIFAD CER 10:00F 600 VUC% A7C1 0150-6093 CIFAD CER 10:00F 600-20 100VUC% A7C1 0150-6093 CIFAD CER 10:00F 600-20 100VUC% A7C1 0150-6093 CIFAD CER 10:00F 600-20 100VUC% A7C1 0150-6093 CIFAD CER 10:00F 600-20 100VUC% A7C1 0150-6093 CIFAD CER 10:00F 600-20 100VUC% A7C1 0150-6093 CIFAD CER 10:00F 600-20 100VUC% A7C1 0150-6093 CIFAD CER 10:00F 600-20 100VUC% A7C1 1010-0013 SENICON UEVICE 10100E GERMANIUM A7CR1 1910-0011 SENICON UEVICE 10100E GERMANIUM A7CR3 1910-0011 SENICON UEVICE 10100E GERMANIUM A7CR3 1910-0011 SENICON UEVICE 10100E GERMANIUM A7CR3 1910-0011 SENICON UEVICE 10100E GERMANIUM A7CR3 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-0153 TRANSISTORIFM SMIGUZ A7A01 1850-	Reference Designation	֍ Stock No.	Description #	Note
ABR3 0683-1525 RIFXU COMP 1500 OHM 5% 174%  ABRT1 0852-0021 RITEMPERATURE COMPENSATING 100 OHM 5% 25C  ABTB1 00851-2010 BOARDICURRENT-CONTROLLED ATTEN  A7 00851-6020 ASSY1 20 MC I.F. AMPLIFIER  A7C1 0150-0093 CIFXU CER 1600FF 600 VUCW A7C2 0150-0093 CIFXU CER 1600FF 600 VUCW A7C3 0140-0170 CIFXU MICA 100 PF 2% 300 VUCW A7C3 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C4 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C5 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C6 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C8 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C8 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C8 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C10 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C11 0140-015 CIFXU CER 0.01UF +80-20 100VUCW A7C12 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C13 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.01UF +80-20 100VUCW A7C14 0150-0093 CIFXU CER 0.00F 80-20 CON EXAMILIAR A7C14 0150-0093 CIFXU CER 0.00F 80-20 CON EXAMILIAR A7C14 0150-0093 CIFXU CER 0.00F 80-20 CON EXAMILIAR A7C14 0150-0093 CIFXU CER 0.00F 80-20 CON EXAMILIAR A7C14 0803-0093 TRANSISTORIFP SHI642 A7C2 0803-0805 RIFXU COMP 680 COM EX 1/4% A7C2 0803-0805 RIFXU COMP 680 COM EX 1/4% A7C2 0803-08				
A6871 0852-0021 RIFAU COMP 1500 OHM 5% 174%  A6871 0852-0021 RIFEMPERATURE COMPENSATING 100 OHM 5% 25C  A6781 00851-2010 BOARDICURRENT-CONTROLLED ATTEN  A7 00851-6020 AS\$Y1 20 MC 1.F. AMPLIFIER  A7C1 0150-0050 CIFAD CR 1600FF 600 VDC%  A7C2 0150-0093 CIFAD CR 1600FF 600 VDC%  A7C3 0140-0176 CIFAD MICA 100 PF 2% 300 VDC%  A7C4 0150-0090 CIFAD CR 0.01UF +80-20 100VDC%  A7C5 0150-0090 CIFAD CR 0.01UF +80-20 100VDC%  A7C6 0150-0090 CIFAD CR 0.01UF +80-20 100VDC%  A7C7 0140-0190 CIFAD CR 0.01UF +80-20 100VDC%  A7C8 0150-0093 CIFAD CR 0.01UF +80-20 100VDC%  A7C8 0150-0093 CIFAD CR 0.01UF +80-20 100VDC%  A7C9 0150-0093 CIFAD CR 0.01UF +80-20 100VDC%  A7C10 0150-0093 CIFAD CR 0.01UF +80-20 100VDC%  A7C11 0150-0093 CIFAD CR 0.01UF +80-20 100VDC%  A7C12 0150-0093 CIFAD CR 0.01UF +80-20 100VDC%  A7C13 0150-0093 CIFAD CR 0.01UF +80-20 100VDC%  A7C14 0150-0093 CIFAD CR 0.01UF +80-20 100VDC%  A7C17 0150-0093 CIFAD CR 0.01UF +80-20 100VDC%  A7C18 1910-0011 SEMICON UEVICE 1010DE GERNANIUM  A7CRA 1910-0011 SEMICON UEVICE 1010DE GERNANIUM  A7CRA 1910-0011 SEMICON UEVICE 1010DE GERNANIUM  A7CRA 1910-0011 SEMICON UEVICE 1010DE GERNANIUM  A7CRA 1910-0011 SEMICON UEVICE 1010DE GERNANIUM  A7CRA 1910-0015 COLLIFAD 0.47UH 20%  COLLIFAD 1.40UH 10%  A7Q1 1850-0153 TRANSISTORIFPD SM1642  A7Q2 1850-0153 TRANSISTORIFPD SM1642  A7Q3 1850-0153 TRANSISTORIFPD SM1642  A7Q3 1850-0153 TRANSISTORIFPD SM1642  A7Q3 1850-0153 TRANSISTORIFPD SM1642  A7Q4 1853-0000 TRANSISTORIFPD SM1642  A7Q3 1850-0153 TRANSISTORIFPD SM1642  A7Q4 1853-0000 TRANSISTORIFPD SM1642  A7Q3 1850-0153 TRANSISTORIFPD SM1642  A7Q4 1853-0153 TRANSISTORIFPD SM1642  A7Q4 1853-0153 TRANSISTORIFPD SM1642  A7Q5 1850-0153 TRANSISTORIFPD SM1642  A7Q6 1850-0153 TRANSISTORIFPD SM1642  A7Q8 1850-0153 TRANSISTORIFPD SM1642  A7Q8 1850-0153 TRANSISTORIFPD SM1642  A7Q9 1850-0153 TRANSISTORIFPD SM1642  A7Q8 1850-0153 TRANSISTORIFPD SM1642  A7Q9 1850-0153 TRANSISTORIFPD SM1642  A7Q9 1850-0153 TRANSISTORIFPD SM1642  A7Q9 1850-0153 TRANSISTORIFPD SM1642  A7Q9 1850-0153 TRANSISTORIFPD SM1642  A7Q9	16R2	0683-6825	RIFXD COMP 6800 OHM 5% 1/4W	
### A6TB1	(6R3	0683-1525	RIFXU COMP 1500 OHM 5% 1/4W	
A7C1 0150-0050	16RT1	0852-0021	RITEMPERATURE COMPENSATING 100 OHM 5% 250	1
A7C1 0150-0050	6TB1	00851-2010	BOARDICURRENT-CONTROLLED ATTEN	
A7C2 0150-0093 C1FXD CER 0.01UF +80-20 100VUCC C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 600 VUCW C1FXD CER 100PF 28 300 VDCW C1FXD CER 100PF 28 300 VDCW C1FXD CER 100PF 28 300 VDCW C1FXD CER 100PF 28 300 VDCW C1FXD CER 100PF 28 300 VDCW C1FXD CER 100PF 28 300 VDCW C1FXD CER 100PF 28 300 VDCW C1FXD CER 100PF 28 300 VDCW C1FXD CER 100PF 28 300 VDCW C1FXD CER 100PF 480-20 100 VUCW C1FXD CER 100PF 480-20 100 VUCW C1FXD CER 100PF 480-20 100 VUCW C1FXD CER 100PF 480-20 100 VUCW C1FXD CER 100 VUCW C1FXD CER 100 VUCW C1FXD CER 100 VUCW C1FXD CER 100 VUCW C1FXD CER 100 VUCW C1FXD CER 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD CIPT 480-20 100 VUCW C1FXD C1PT 480-20 100 VUCW	.7	00851-6020	ASSY! 20 MC I.F. AMPLIFIER	
A7C2 0150-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C4 0150-0050 C1FXD CER LOOPF 600 VDCW   C1FXD CER LOOPF 600 VDCW   C1FXD CER LOOPF 600 VDCW   C1FXD CER LOOPF 600 VDCW   A7C7 0150-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C7 0140-0190 C1FXD MICA 39 PF 58 300 VDCW   A7C8 0150-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C9 0150-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C9 0150-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C10 0150-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C11 0140-0215 C1FXD MICA 80FF 28 300VDCW   A7C13 0150-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C13 0150-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C14 0150-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C15 0150-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C16 1910-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C17 1910-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C18 1910-0093 C1FXD CER U_01UF +80-20 100VUCC   A7C19 1910-0011 SEMICON UEVICE:010UE GERMANIUM   A7CR2 1910-0011 SEMICON UEVICE:010UE GERMANIUM   A7CR4 1910-0011 SEMICON UEVICE:010UE GERMANIUM   A7CR4 1910-0011 SEMICON UEVICE:010UE GERMANIUM   A7CR4 1910-0013 SEMICON UEVICE:010UE GERMANIUM   A7CR4 1910-0013 SEMICON UEVICE:010UE GERMANIUM   A7CR4 1910-0013 SEMICON UEVICE:010UE GERMANIUM   A7CR4 1910-0015 C01L1FXD U_04TUD 10X   A7C12 1850-0153 TRANSISTORIPMP SM1642 TRANSISTORIPMP SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1642 TRANSISTORIPM SM1644 TRANSISTO	701	0150-0050	CIEXD CER 1000PE 600 VOCH	
A7C3 0140-0170 CIFXD MICA 100 PF 2% 300 VDCW A7C4 0150-0050 CIFAD CER 1000PF 600 VDCW A7C5 0150-0093 CIFAD CER 1000PF 600 VDCW A7C6 0150-0093 CIFAD CER 0.01UF +80-20 100VDCW A7C7 0140-0190 CIFXD MICA 39 PF 5% 300 VDCW A7C8 0150-0050 CIFAD CER 0.01UF +80-20 100VDCW A7C8 0150-0053 CIFAD CER 0.01UF +80-20 100VDCW A7C10 0150-0093 CIFAD CER 0.01UF +80-20 100VDCW A7C11 0140-0215 CIFAD MICA 80PF 2% 300VDCW A7C12 0150-0093 CIFAD CER 0.01UF +80-20 100VDCW A7C13 0150-0093 CIFAD CER 0.01UF +80-20 100VDCW A7C14 0150-0093 CIFAD CER 0.01UF +80-20 100VDCW A7C15 0150-0093 CIFAD CER 0.01UF +80-20 100VDCW A7C16 1910-0011 SEMICON DEVICE:010DE GERMANIUM A7C17 1910-0011 SEMICON DEVICE:010DE GERMANIUM A7C83 1910-0011 SEMICON DEVICE:010DE GERMANIUM A7C84 1910-0011 SEMICON DEVICE:010DE GERMANIUM A7C87 1910-0011 SEMICON DEVICE:010DE GERMANIUM A7C87 1910-0011 SEMICON DEVICE:010DE GERMANIUM A7C89 1910-0013 TRANSISTORIFMP SMI642 A7C12 1850-0153 TRANSISTORIFMP SMI642 A7C2 1850-0153 TRANSISTORIFMP SMI642 A7C2 1850-0153 TRANSISTORIFMP SMI642 A7C2 1850-0153 TRANSISTORIFMP SMI642 A7C2 1850-0153 TRANSISTORIFMP SMI642 A7C2 1850-0153 TRANSISTORIFMP SMI642 A7C3 1854-0005 TRANSISTORIFMP SMI642 A7C3 1854-0005 TRANSISTORIFMP SMI642 A7C84 1854-0005 TRANSISTORIFMP SMI642 A7C85 0683-6825 RIFAD COMP 3900 OMM 5% 1/4W A7C87 0683-1825 RIFAD COMP 6800 OMM 5% 1/4W A7C87 0683-1825 RIFAD COMP 6800 OMM 5% 1/4W A7C87 0683-1825 RIFAD COMP 6800 OMM 5% 1/4W A7C87 0683-1825 RIFAD COMP 6800 OMM 5% 1/4W A7C87 0683-1825 RIFAD COMP 6800 OMM 5% 1/4W A7C87 0683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C87 0683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C8 0683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C8 0683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C8 0683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C8 0683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C8 0683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C8 OMB 683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C8 OMB 683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C8 OMB 683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C8 OMB 683-6825 RIFAD COMP 6800 OMM 5% 1/4W A7C8 OMB 683-6825 RIFAD COMP 6800 OMM 5% 1/4	702		CIEXO CER O CIVE HODOCA	i
A7C4 0150-0050 C1FAD CER 1000PF 600 VDCW A7C5 0150-0093 C1FAD CER 0.01UF +80-20 100VUCW A7C6 0150-0093 C1FAD CER 0.01UF +80-20 100VUCW A7C7 0140-0190 C1FAD CER 0.01UF +80-20 100VUCW A7C8 0150-0093 C1FAD CER 0.01UF +80-20 100VUCW A7C8 0150-0093 C1FAD CER 0.01UF +80-20 100VUCW A7C9 0150-0093 C1FAD CER 0.01UF +80-20 100VUCW A7C10 0150-0093 C1FAD CER 0.01UF +80-20 100VUCW A7C11 0150-0093 C1FAD CER 0.01UF +80-20 100VUCW A7C12 0150-0093 C1FAD CER 0.01UF +80-20 100VUCW A7C13 0150-0093 C1FAD CER 0.01UF +80-20 100VUCW A7C14 0150-0093 C1FAD CER 0.01UF +80-20 100VUCW A7C71 1910-0011 SEMICON DEVICE:010UF GERMANIUM A7C71 1910-0011 SEMICON DEVICE:010UF GERMANIUM A7C73 1910-0011 SEMICON DEVICE:010UF GERMANIUM A7C73 1910-0011 SEMICON DEVICE:010UF GERMANIUM A7C74 1910-0011 SEMICON DEVICE:010UF GERMANIUM A7C75 1910-0011 SEMICON DEVICE:010UF GERMANIUM A7C77 1910-0011 SEMICON DEVICE:010UF GERMANIUM A7C78 1910-0011 SEMICON DEVICE:010UF GERMANIUM A7C78 1910-0011 SEMICON DEVICE:010UF GERMANIUM A7C79 1850-0153 TRANSISTORIPP SMI642 A7C9 1850-0153 TRANSISTORIPP SMI642 A701 1850-0153 TRANSISTORIPP SMI642 A702 1850-0153 TRANSISTORIPP SMI642 A703 1850-0153 TRANSISTORIPP SMI642 A704 1853-0003 TRANSISTORIPP SMI642 A704 1853-0003 TRANSISTORIPP SMI642 A705 A776 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7771 0683-1825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-1825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-6825 RIFAD COMP 5800 OMM 5% 1/4W A7781 0683-		1	CIEVE MICA 100 PE ON ZO LOUGH	
A7C5 0150-0093 C1FXD CER 0.01UF +80-20 100VUCW A7C6 0150-0093 C1FXD CER 0.01UF +80-20 100VUCW A7C8 0150-0050 C1FXD CER 1600PF 600 VDCW A7C8 0150-0093 C1FXD CER 1600PF 600 VDCW A7C9 0150-0093 C1FXD CER 1600PF 600 VDCW A7C10 0150-0093 C1FXD CER 0.01UF +80-20 100VUCW A7C11 0140-0215 C1FXD CER 0.01UF +80-20 100VUCW A7C12 0150-0093 C1FXD CER 0.01UF +80-20 100VUCW A7C13 0150-0093 C1FXD CER 0.01UF +80-20 100VUCW A7C14 0150-0093 C1FXD CER 0.01UF +80-20 100VUCW A7C15 0150-0093 C1FXD CER 0.01UF +80-20 100VUCW A7C14 0150-0093 C1FXD CER 0.01UF +80-20 100VUCW A7C14 0150-0093 C1FXD CER 0.01UF +80-20 100VUCW A7C14 1910-0011 SEMICON UEVICE:1010UE GERRANIUM A7CR3 1910-0011 SEMICON UEVICE:1010UE GERRANIUM A7CR4	-	_	C.EXD CER 100 PF 28 300 VDCW	
A7C6 0150-0093			CIPAD CER 1000PF 800 VDCW	l
A7C7 A7C8 A7C9 A7C9 A7C9 A7C9 A7C9 A7C9 A7C9 A7C9	.705	0150-0093	C 1FXD CER 0.01UF +80-20 100VDCW	
A7CB A7CB A7C9 O150-0093 C1FAÜ CER 100PF 600 VDCW A7C10 O150-0093 C1FAÜ CER 0.01UF +88-20 100VUCW A7C11 O150-0093 C1FAÜ CER 0.01UF +88-20 100VUCW A7C12 O150-0093 C1FAÜ CER 0.01UF +88-20 100VUCW A7C13 A7C13 O150-0093 C1FAÜ CER 0.01UF +80-20 100VUCW A7C13 A7C14 O150-0093 C1FAÜ CER 0.01UF +80-20 100VUCW A7C13 A7C14 O150-0093 C1FAÜ CER 0.01UF +80-20 100VUCW A7C14 A7C14 O150-0093 C1FAÜ CER 0.01UF +80-20 100VUCW A7C14 A7C15 A7C14 A7C16 A7C17 A7C17 A7C17 A7C18 A7C19 A7C19 A7C19 A7C11 A7C19 A7C19 A7C11 A7C19 A7C11 A7C19 A7C11 A7C19 A7C11 A7C19 A7C11 A7C19 A7C11 A7C19 A7C11 A7C19 A7C10 A7C11 A7C19 A7C10 A7C11 A7C1		•	CIFXU CER 0.01UF +80-20 100VUC#	
A768		0140-0190	CIFXD MICA 39 PF 5% 300 VDCW	
A7C9 A7C10  A7C10  A7C10  A7C11  A7C11  A7C12  A7C12  A7C12  A7C13  A7C13  A7C13  A7C14  A7C14  A7C14  A7C15  A7C15  A7C15  A7C15  A7C15  A7C16  A7C16  A7C17  A7C17  A7C17  A7C17  A7C18  A7C18  A7C18  A7C19  A7C19  A7C19  A7C19  A7C19  A7C11  A7C19  A7C11  A7C10  A7C11  A7C10  A7C11  A7C10  A7C11  A7C1	-		CIFXD CER 1000PF 600 VDCw	
A7C11 0150-0093 C:FXD CER 0.01UF +80-20 100VDCW A7C11 0140-0215 C:FXD CER 0.01UF +80-20 100VDCW A7C12 0150-0093 C:FXD CER 0.01UF +80-20 100VDCW A7C13 0150-0093 C:FXD CER 0.01UF +80-20 100VDCW A7C14 0150-0093 C:FXD CER 0.01UF +80-20 100VDCW A7C14 0150-0093 C:FXD CER 0.01UF +80-20 100VDCW A7C14 0150-0093 C:FXD CER 0.01UF +80-20 100VDCW A7C14 1910-0011 SEMICON DEVICE:D10DE GERMANIUM A7CR2 1910-0011 SEMICON DEVICE:D10DE GERMANIUM A7CR4 1910-0011 SEMICON DEVICE:D10DE GERMANIUM A7CR4 1910-0011 SEMICON DEVICE:D10DE GERMANIUM A7CR4 1910-0011 SEMICON DEVICE:D10DE GERMANIUM A7CR4 1910-0015 COIL:FXD 0.47UH 20% COIL:FXD 1.0UH 10% A7C1 1650-0153 TRANSISTOR:FNP SM1642 A702 1850-0153 TRANSISTOR:FNP SM1642 A703 1850-0153 TRANSISTOR:FNP SM1642 A704 1853-0003 TRANSISTOR:FNP SM1642 A709 1854-0005 TRANSISTOR:FNP SM1642 A709 1854-0005 TRANSISTOR:FNP SM1642 A709 1854-0005 TRANSISTOR:FNP SM1642 A709 1854-0005 TRANSISTOR:FNP SM1642 A7R1 0683-3925 R:FXD COMP 3900 UMM 5% 1/4% A7R2 0683-6825 R:FXD COMP 6800 UMM 5% 1/4% A7R3 0683-1825 R:FXD COMP 6800 UMM 5% 1/4% A7R6 0683-6825 R:FXD COMP 680 UMM 5% 1/4% A7R6 0683-6825 R:FXD COMP 680 UMM 5% 1/4% A7R6 0683-6825 R:FXD COMP 680 UMM 5% 1/4% A7R6 0683-6825 R:FXD COMP 680 UMM 5% 1/4% A7R8 0683-1225 R:FXD COMP 680 UMM 5% 1/4% A7R8 0683-1225 R:FXD COMP 680 UMM 5% 1/4% A7R8 0683-6825 R:FXD COMP 680 UMM 5% 1/4% A7R1 0683-6825 R:FXD COMP 680 UMM 5% 1/4% A7R1 0683-6825 R:FXD COMP 680 UMM 5% 1/4% A7R1 0683-6825 R:FXD COMP 680 UMM 5% 1/4% A7R11 0683-6825 R:FXD COMP 390 UMM 5% 1/4% A7R11 0683-6825 R:FXD COMP 680 UMM 5% 1/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R11 0683-6825 R:FXD COMP 5% 0/4% A7R12 0683-6825 R:FXD COMP 5% 0/4% A7R14 U683-6825 R:FXD COMP 5% 0/4%		0150-0093	CIFXD CER 0.01UF +80-20 100VDCW	
A7C12	7010	0150-0093	C:FXD CER 0.01UF +80-20 100VDCW	
A7C12 A7C13 A7C14 D150-0093 A7C14 D150-0093 A7C14 A7C14 D150-0093 C1FXD CER U.01UF +80-20 100VUCW C1FXD CER U.01UF +80-20 100VUCW C1FXD CER U.01UF +80-20 100VUCW A7CR1 A7CR2 A7CR1 A7CR2 D1910-0011 SEMICON UEVICE:D10DE GERMANIUM A7CR3 D1910-0011 SEMICON UEVICE:D10DE GERMANIUM A7CR4 D1910-0011 SEMICON UEVICE:D10DE GERMANIUM A7CR4 D1910-0011 SEMICON UEVICE:D10DE GERMANIUM A7CR4 A7CR2 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR4 A7CR4 A7CR4 A7CR4 A7CR4 A7CR4 A7CR4 A7CR4 A7CR4 A7CR5 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR8 A7CR7 A7CR8 A7CR8 A7CR8 A7CR8 A7CR9 A7CR9 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR2 A7CR2 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR4 A7CR3 A7CR4 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR7 A7CR8 A7CR8 A7CR8 A7CR8 A7CR9 A7CR9 A7CR1	7011	0140-0215	CIFXO MICA EOPF 2% 300VDCW	
A7C14 0150-0093 C1FXD CER U.01UF +80-20 100VDCW A7C14 0150-0093 C1FXD CER U.01UF +80-20 100VDCW A7C14 0190-0011 SEMICON DEVICE:DIODE GERMANIUM A7CR3 1910-0011 SEMICON DEVICE:DIODE GERMANIUM A7CR4 1910-0011 SEMICON DEVICE:DIODE GERMANIUM A7CR4 1910-0011 SEMICON DEVICE:DIODE GERMANIUM A7CR4 1910-0011 SEMICON DEVICE:DIODE GERMANIUM A7CR4 1910-0011 SEMICON DEVICE:DIODE GERMANIUM A7CR4 1910-0015 COLLIFXD U.4TUH 20% A7CR2 9140-0158 COLLIFXD 1.0UH 10% A7CR2 1850-0153 TRANSISTORIFNP SMI642 A7CR3 1850-0153 TRANSISTORIFNP SMI642 A7CR3 1850-0153 TRANSISTORIFNP SMI642 A7CR3 1850-0153 TRANSISTORIFNP SMI642 A7CR3 1850-005 TRANSISTORIFNP SMICH2 A7CR3 1850-005 TRANSISTORIFNP SMICH2 A7CR4 0683-3925 RIFXD COMP 3900 UHM 5% 1/4W A7CR5 0683-1825 RIFXD COMP 6800 UHM 5% 1/4W A7CR6 0683-1825 RIFXD COMP 1800 UHM 5% 1/4W A7CR6 0683-6805 RIFXD COMP 6800 UHM 5% 1/4W A7CR6 0683-6805 RIFXD COMP 6800 UHM 5% 1/4W A7CR7 0683-3925 RIFXD COMP 6800 UHM 5% 1/4W A7CR8 0683-1225 RIFXD COMP 6800 UHM 5% 1/4W A7CR8 0683-1225 RIFXD COMP 6800 UHM 5% 1/4W A7CR9 0683-1225 RIFXD COMP 6800 UHM 5% 1/4W A7CR1 0683-3925 RIFXD COMP 680 UHM 5% 1/4W A7CR1 0683-3925 RIFXD COMP 680 UHM 5% 1/4W A7CR1 0683-3925 RIFXD COMP 680 UHM 5% 1/4W A7CR1 0683-3925 RIFXD COMP 680 UHM 5% 1/4W A7CR1 0683-3925 RIFXD COMP 680 UHM 5% 1/4W A7CR1 0683-3925 RIFXD COMP 680 UHM 5% 1/4W A7CR1 0683-3925 RIFXD COMP 680 UHM 5% 1/4W A7CR1 0683-3925 RIFXD COMP 680 UHM 5% 1/4W A7CR1 0683-3925 RIFXD COMP 680 UHM 5% 1/4W A7CR1 0683-3925 RIFXD COMP 680 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0683-1525 RIFXD COMP 50 UHM 5% 1/4W A7CR1 0HM 5 UND TO UHM 5% 1/4W A7CR1 0HM 5 UND TO UND TO UND TO UND TO UND TO UND T			C1FXD CER U-01UF +80-20 100VDCW	
A7C14 0150-0093 C:FXD CER C:01UF +80-20 100VUCW  A7CR1 1910-0011 SEMICON UEVICE:DIOUE GERMANIUM SEMICON UEVICE:DIOUE GEMANIUM SEMICON UEVICE:DIOUE GEMANIUM SEMICON UEVICE:DIOUE GEMANIUM SEMICON UEVICE:DIOUE GEMANIUM SEMICON UEVICE:DIOUE GERMANIUM S			C1FXD CFR 0.01UF +80-20 100VDCW	
A7CR2 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR4 A7CR3 A7CR4 A7CR5 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR7 A7CR6 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR8 A7CR8 A7CR8 A7CR8 A7CR8 A7CR9 A7CR8 A7CR9 A7CR1			C:FXD CER 0.01UF +80-20 100VUCW	
A7CR2 A7CR3 A7CR3 A7CR3 A7CR3 A7CR3 A7CR4 A7CR3 A7CR4 A7CR5 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR7 A7CR6 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR6 A7CR7 A7CR7 A7CR7 A7CR7 A7CR7 A7CR8 A7CR8 A7CR8 A7CR8 A7CR8 A7CR9 A7CR8 A7CR9 A7CR1	7CR1	1910+0011		
A7CR3 A7CR4 A7CR5 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR7 A7CR6 A7CR7 A7CR7 A7CR8 A7CR7 A7CR8 A7CR7 A7CR8 A7CR8 A7CR9 A7CR9 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR1 A7CR2 A7CR3 A7CR4 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR6 A7CR7 A7CR8 A7CR7 A7CR8 A7CR8 A7CR8 A7CR8 A7CR8 A7CR8 A7CR8 A7CR8 A7CR9 A7CR9 A7CR1			SENTION DEVICE DIOUE GERMANIUM	
A7CR4 1910-0011 SEMICON DEVICE:DIODE GERMANIUM  A7L1 9140-0159 COLL:FXD U-47UH 20%  A7Q1 1850-0153 TRANSISTOR:FNP SM1642  A7Q2 1850-0153 TRANSISTOR:FNP SM1642  A7Q3 1850-0153 TRANSISTOR:FNP SM1642  A7Q3 1850-0153 TRANSISTOR:FNP SM1642  A7Q4 1853-0003 TRANSISTOR:FNP SM1642  A7Q5 1854-0005 TRANSISTOR:PNP SILICON F SOMC MIN  A7Q5 1854-0005 TRANSISTOR:PNP SILICON F SOMC MIN  A7Q5 1854-0005 TRANSISTOR:2N7UB NPN SILICON  A7Q8 0683-6425 R:FXU COMP 6800 OHM 5% 1/4%  A7Q9 0683-6425 R:FXU COMP 1800 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 1800 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q9 A7Q9 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q1 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q1 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q1 O683-6805 R:FXU COMP 68 OHM 5% 1/4%  A7Q1 O683-6805 R:FXU COMP 5% 1/4%  A7			SEMICON DEVICE DIODE GERMANIUM	
A7L1 9140-0159 COILIFXD U-47UH 20% A7C1 1850-0153 TRANSISTORIFNP SM1642 A7G2 1850-0153 TRANSISTORIFNP SM1642 A7G3 1850-0153 TRANSISTORIFNP SM1642 A7G3 1850-0153 TRANSISTORIFNP SM1642 A7G4 1853-0003 TRANSISTORIFNP SM1642 A7G5 1854-0005 TRANSISTORIFNP SILICON F 50MC MIN A7G7 1854-0005 TRANSISTORIPNP SILICON A7G8 0683-3925 RIFXD COMP 5800 OHM 5% 1/4% A7G8 0683-1825 RIFXD COMP 1800 OHM 5% 1/4% A7G9 A7G7 0683-6825 RIFXD COMP 1800 OHM 5% 1/4% A7G7 0683-6805 RIFXD COMP 51 OHM 5% 1/4% A7G7 0683-6805 RIFXD COMP 58 0HM 5% 1/4% A7G7 0683-3925 RIFXD COMP 68 OHM 5% 1/4% A7G8 0683-3925 RIFXD COMP 5800 OHM 5% 1/4% A7G8 0683-1225 RIFXD COMP 5900 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 5900 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 5900 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 5800 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 680 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 680 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 680 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 680 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 680 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 68 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 68 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 680 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 680 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 680 OHM 5% 1/4% A7G8 0683-6825 RIFXD COMP 680 OHM 5% 1/4% A7G8 O683-6825 RIFXD COMP	•		SEMICON DEVICE DIODE GERMANIUM	
A7L2 9140-0158 COIL:FXD 1.0UH 1D%  A7Q1 1850-0153 TRANSISTOR:PNP SM1642 A7Q2 1850-0153 TRANSISTOR:PNP SM1642 A7Q3 1850-0153 TRANSISTOR:PNP SM1642 A7Q4 1853-0003 TRANSISTOR:PNP SILICON F 50MC MIN A7Q5 1854-0005 TRANSISTOR:PNP SILICON  A7R1 0683-3925 R:FXD COMP 3900 OHM 5% 1/4% A7R2 0683-6825 R:FXD COMP 6800 OHM 5% 1/4% A7R3 0683-1825 R:FXD COMP 6800 OHM 5% 1/4% A7R4 0683-5105 R:FXD COMP 51 OHM 5% 1/4% FACTURY SELECTED PARTITYPICAL VALUE GIVEN A7R6 0683-6805 R:FXD COMP 68 OHM 5% 1/4% A7R7 0683-3925 R:FXD COMP 68 OHM 5% 1/4% A7R8 0683-1225 R:FXD COMP 5900 OHM 5% 1/4% A7R8 0683-1225 R:FXD COMP 68 OHM 5% 1/4% A7R10 0683-6805 R:FXD COMP 68 OHM 5% 1/4% A7R11 0683-6805 R:FXD COMP 68 OHM 5% 1/4% A7R12 0683-6825 R:FXD COMP 68 OHM 5% 1/4% A7R11 0683-2705 R:FXD COMP 68 OHM 5% 1/4% A7R12 0683-6825 R:FXD COMP 68 OHM 5% 1/4% A7R13 0683-2705 R:FXD COMP 68 OHM 5% 1/4% A7R14 0683-1525 R:FXD COMP 68 OHM 5% 1/4% A7R15 O683-2705 R:FXD COMP 68 OHM 5% 1/4% A7R14 O683-1525 R:FXD COMP 1500 OHM 5% 1/4% A7R14 OF A7R14 R:FXD COMP 1500 OHM 5% 1/4% A7R15 OF A7R14 R:FXD COMP 1500 OHM 5% 1/4% A7R14 D:FXD COMP 1500 OHM 5% 1/4% A7R14 D:FX	76.44	1910-0011	SEMICON DEVICE FUIODE GERMANIUM	
A701 1650-0153 TRANSISTORIPNP SM1642 A702 1850-0153 TRANSISTORIPNP SM1642 A703 1850-0153 TRANSISTORIPNP SM1642 A704 1853-0003 TRANSISTORIPNP SILICON F SOME MIN A705 1854-0005 TRANSISTORIPNP SILICON  A771 0683-3925 RIFAD COMP 3900 OHM 5% 1/4% A772 0683-6825 RIFAD COMP 6800 OHM 5% 1/4% A773 0683-1825 RIFAD COMP 1800 OHM 5% 1/4% FACTORY SELECTED PARTITYPICAL VALUE GIVEN A776 0683-6805 RIFAD COMP 68 OHM 5% 1/4% FACTORY SELECTED PARTITYPICAL VALUE GIVEN A777 0683-3925 RIFAD COMP 6800 OHM 5% 1/4% A778 0683-1225 RIFAD COMP 6800 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 6800 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 5900 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 68 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 68 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 68 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 68 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 68 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 68 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 68 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 68 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 3900 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 68 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 3900 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 3900 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 3900 OHM 5% 1/4% A778 0683-3925 RIFAD COMP 3900 OHM 5% 1/4% A778 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A778 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A778 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 14 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 14 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 14 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 14 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 14 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 14 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 14 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 14 0683-1525 RIFAD COMP 1500 OHM 5% 1/4% A78 14 0683-				
A7Q2 A7Q3 A7Q3 A7Q4 A7Q4 A7Q4 A7Q5 A7Q4 A7Q5 A7Q5 A7Q6 A7Q5 A7Q6 A7Q6 A7Q6 A7Q7 A7Q6 A7Q7 A7Q6 A7Q7 A7Q7	/L2	9140-0158	COIL:FXD 1.0UH 10%	
A7Q2 A7Q3 A7Q3 A7Q4 A7Q4 A7Q5 A7Q4 A7Q5 A7Q5 A7Q6 A7Q5 A7Q6 A7Q6 A7Q6 A7Q6 A7Q7 A7Q6 A7Q7 A7Q6 A7Q7 A7Q6 A7Q7 A7Q6 A7Q7 A7Q6 A7Q7 A7Q6 A7Q7 A7Q7	701	1850-0153	TRANSISTORIENP SMIGUS	
A7Q3 A7Q4 A7Q4 A7Q5 A7Q4 A7Q5 A7Q6 A7Q5 A7Q6 A7Q6 A7Q6 A7Q6 A7Q7 A7Q6 A7Q6 A7Q7 A7Q7	702		TRANSISTORIENP SMIAU2	
A795  A795  A795  A795  A796  A796  A796  A781  A781  A782  A783  A783  A784  A784  A784  A785  A784  A786  A786  A786  A787  A788  A788  A789  A789  A7814  A7814  A7814  A7814  A7814  A7814  A7814  A7814  A7814  A7814  A7814  A7815  A7814  A7816  A7816  A7816  A7817  A7818  A7818  A7814  A7814  A7814  A7814  A7814  A7814  A7815  A7814  A7814  A7816  A7817  A7816  A7817  A7817  A7817  A7818  A7818  A7814  A7814  A7814  A7814  A7814  A7814  A7814  A7814  A7815  A7814  A7816  A7816  A7817  A7817  A7817  A7818  A7818  A7818  A7819  A7810  A7810  A7810  A7810  A7810  A7811  A7811  A7811  A7812  A7813  A7814  A7814  A7814  A7814  A7814  A7814  A7815  A7814  A7816  A7816  A7817  A7817  A7817  A7818  A7818  A7818  A7819  A7819  A7810	703		TRANSISTORIENP SMIAUS	
A7R1			TRANSISTORIPHP STLICON E SOME MIN	
A7R2 A7R3 A7R4 O683-6825 A7R4 O683-1825 A7R4 O683-5105  A7R5 A7R6 A7R7 A7R6 A7R7 O683-6805  A7R7 A7R8 A7R8 A7R9 A7R10 A7R11 O683-6805  A7R11 A7R12 A7R12 A7R13 A7R14 A7R17 A7R18 A7R14 A7R18 A7R19 A7R19 A7R19 A7R11 A7R11 A7R11 A7R11 A7R11 A7R11 A7R12 A7R11 A7R11 A7R12 A7R11 A7R12 A7R13 A7R14 A7R14 A7R14 A7R14 A7R14 A7R14 A7R14 A7R15 A7R14 A7R16 A7R17 A7R17 A7R17 A7R18 A7R18 A7R19 A7R19 A7R19 A7R19 A7R19 A7R10 A7R10 A7R10 A7R11 A7R11 A7R11 A7R11 A7R11 A7R11 A7R12 A7R11 A7R	f .		TRANSISTOR 2N7UB NPN SILICON	
A7R2 A7R3 A7R4 O683-6825 A7R4 O683-1825 A7R4 O683-5105  A7R5 A7R6 A7R7 A7R6 A7R7 O683-6805  A7R7 A7R8 A7R8 A7R9 A7R10 A7R11 O683-6805  A7R11 A7R12 A7R12 A7R13 A7R14 A7R17 A7R18 A7R14 A7R18 A7R19 A7R19 A7R19 A7R11 A7R11 A7R11 A7R11 A7R11 A7R11 A7R12 A7R11 A7R11 A7R12 A7R11 A7R12 A7R13 A7R14 A7R14 A7R14 A7R14 A7R14 A7R14 A7R14 A7R15 A7R14 A7R16 A7R17 A7R17 A7R17 A7R18 A7R18 A7R19 A7R19 A7R19 A7R19 A7R19 A7R10 A7R10 A7R10 A7R11 A7R11 A7R11 A7R11 A7R11 A7R11 A7R12 A7R11 A7R	721	0493-3026		
A7R3 A7R4 O683-1825 A7R4 O683-5105 RIFXD COMP 180U OHM 5% 1/4% FACTORY SELECTED PARTITYPICAL VALUE GIVEN A7R5 A7R6 O683-6805 RIFXD COMP 68 OHM 5% 1/4% A7R7 O683-3925 RIFXD COMP 680U OHM 5% 1/4% A7R8 O683-1225 RIFXD COMP 590U OHM 5% 1/4% A7R9 A7R10 O683-6805 RIFXD COMP 120U OHM 5% 1/4% NOT ASSIGNED RIFXD COMP 68 OHM 5% 1/4% A7R11 O683-3925 RIFXD COMP 68 OHM 5% 1/4% A7R12 O683-6425 A7R13 O683-2705 RIFXD COMP 390U OHM 5% 1/4% RIFXD COMP 680U OHM 5% 1/4% For Service Manuals Contact MAURITRON TECHNICAL SERVICES MAURITRON TECHNICAL SERVICES			RIEVI COMP 4000 OHM 5% 1/4W	
A7R4  0683-5105  RIFXD COMP 51 OHM 5% 1/4%  FACTORY SELECTED PARTITYPICAL VALUE GIVEN  RIFXD COMP 68 OHM 5% 1/4W  A7R6  0683-6805  RIFXD COMP 6800 OHM 5% 1/4W  A7R7  0683-3925  RIFXD COMP 5900 OHM 5% 1/4W  A7R8  0683-1225  RIFXD COMP 1200 OHM 5% 1/4W  NOT ASSIGNED  RIFXD COMP 68 OHM 5% 1/4W  A7R11  0683-6805  RIFXD COMP 68 OHM 5% 1/4W  A7R12  0683-6825  RIFXD COMP 6800 OHM 5% 1/4W  A7R13  0683-2705  RIFXD COMP 6800 OHM 5% 1/4W  FOR Service Manuals Contact  MAURITRON TECHNICAL SERVICES  MAURITRON TECHNICAL SERVICES			RIEVO COMP 1000 ONM 5% 1/4W	
A7R5  0683-6805  RIFXD COMP 68 OHM 5% 1/4W  A7R6  0683-6825  RIFXD COMP 6800 OHM 5% 1/4W  A7R7  0683-3925  RIFXD COMP 5900 OHM 5% 1/4W  A7R8  A7R8  A7R10  0683-6805  RIFXD COMP 1200 OHM 5% 1/4W  NOT ASSIGNED  RIFXD COMP 68 OHM 5% 1/4W  A7R11  0683-3925  RIFXD COMP 68 OHM 5% 1/4W  A7R12  0683-6825  RIFXD COMP 6800 OHM 5% 1/4W  A7R13  0683-2705  RIFXD COMP 6800 OHM 5% 1/4W  FOR Service Manuals Contact  MAURITRON TECHNICAL SERVICES  MAURITRON TECHNICAL SERVICES	1		BIEND COMP IN COM SW 1/4%	
A7R5  0683-6805  R:FXD COMP 68 OHM 5% 1/4W  A7R6  0683-6825  R:FXD COMP 6800 OHM 5% 1/4W  A7R7  0683-3925  R:FXD COMP 5900 OHM 5% 1/4W  A7R8  A7R9  A7R10  0683-6805  R:FXD COMP 1200 OHM 5% 1/4W  NOT ASSIGNED  R:FXD COMP 68 OHM 5% 1/4W  A7R11  0683-3925  R:FXD COMP 68 OHM 5% 1/4W  A7R12  0683-6825  R:FXD COMP 6800 OHM 5% 1/4W  A7R13  0683-2705  R:FXD COMP 6800 OHM 5% 1/4W  A7R14  0683-1525  R:FXD COMP 1500 OHM 5% 1/4W  A7R14  For Service Manuals Contact  MAURITRON TECHNICAL SERVICES		0883-9103	EACTORY SELECTED BARTATYPICAL MALUE CAME	
A7R6 A7R7 O683-6825 R:FXU COMP 6800 OHM 5% 1/4W A7R8 O683-1225 R:FXU COMP 5900 OHM 5% 1/4W A7R8 O683-1225 R:FXU COMP 1200 OHM 5% 1/4W NOT ASSIGNED R:FXD COMP 68 OHM 5% 1/4W A7R11 O683-6805 R:FXD COMP 68 OHM 5% 1/4W A7R12 O683-6825 R:FXD COMP 6800 OHM 5% 1/4W A7R13 O683-2705 R:FXD COMP 6800 OHM 5% 1/4W R:FXD COMP 6800 OHM 5% 1/4W A7R14 O683-1525 R:FXD COMP 1500 OHM 5% 1/4W A7R14 O683-1525 R:FXD COMP 1500 OHM 5% 1/4W A7R14 O683-1525 R:FXD COMP 1500 OHM 5% 1/4W A7R14 O683-1525	7R5	0683-6805	RIFXD COMP 68 OHM 5% 1/4W	
A7R7 A7R8 A7R9 A7R10 O683-6805 R:FXD COMP 3900 OHM 5% 1/4W NOT ASSIGNED R:FXD COMP 68 OHM 5% 1/4W A7R11 O683-6805 R:FXD COMP 68 OHM 5% 1/4W A7R12 O683-6825 A7R13 O683-2705 R:FXD COMP 6800 OHM 5% 1/4W R:FXD COMP 6800 OHM 5% 1/4W R:FXD COMP 6800 OHM 5% 1/4W R:FXD COMP 27 OHM 5% 1/4W R:FXD COMP 1500 OHM 5% 1/4W A7R14 O683-1525 R:FXD COMP 1500 OHM 5% 1/4W A7R14 O683-1525	794			
A7R8			NIFAU CUMP 6800 OHM 5% 1/4W	
A7R10 0683-6805 R:FXD COMP 68 OHM 5% 1/4w  A7R11 0683-3925 R:FXD COMP 3900 OHM 5% 1/4w  A7R12 0683-6825 R:FXD COMP 6800 OHM 5% 1/4w  A7R13 0683-2705 R:FXD COMP 27 OHM 5% 1/4w  A7R14 0683-1525 R:FXD COMP 1500 OHM 5% 1/4w  MAURITRON TECHNICAL SERVICES  A7R14 OCENTER DE CENTRE D			KIFAU COMP 3900 OHM 5% 1/4W	
A7R10 0683-6805 R:FXD COMP 68 OHM 5% 1/4w  A7R11 0683-3925 R:FXD COMP 3900 OHM 5% 1/4w  A7R12 0683-6825 R:FXD COMP 6800 OHM 5% 1/4w  A7R13 0683-2705 R:FXD COMP 27 OHM 5% 1/4w  A7R14 0683-1525 R:FXD COMP 1500 OHM 5% 1/4w  MAURITRON TECHNICAL SERVICES		0685-1225		
A7R11 0683-3925 R:FXD COMP 3900 OHM 5% 1/4W A7R12 0683-6625 R:FXD COMP 6800 OHM 5% 1/4W A7R13 0683-2705 R:FXD COMP 27 OHM 5% 1/4W A7R14 0683-1525 R:FXD COMP 1500 OHM 5% 1/4W A7R14 O683-1525 R:FXD COMP 1500 OHM 5% 1/4W A7R14 O683-1525 R:FXD COMP 1500 OHM 5% 1/4W A7R14 O683-1525 R:FXD COMP 1500 OHM 5% 1/4W ANDITION TECHNICAL SERVICES		0683-6805		
A7R12 0683-6825 R:FXU COMP 6800 0HM 5% 1/4W For Service Manuals Contact MAURITRON TECHNICAL SERVICES R:FXU COMP 1500 0HM 5% 1/4W MAURITRON TECHNICAL SERVICES		_		
A7R13 0683-2705 R:FXU COMF 27 OHM 5% 1/4W For Service Manuals Contact A7R14 0683-1525 R:FXU COMF 1500 OHM 5% 1/4W MAURITRON TECHNICAL SERVICES			R:FXD COMP 3900 OHM 5% 1/4W	
A7R13 0683-2705 R:FXU COMP 27 OHM 5% 1/4w MAURITON TECHNICAL SERVICES MAURITON TECHNICAL SERVICES	1		RIFXD COMP 6800 OHM 5% 1/4W	
A7R14 0683-1525 R:FXU COMP 1500 OHM 5% 1/4W MAURITHON TECHNICAL SERVICES		0683-2705	R:FXU COMP 27 OHM 5% 1/4# FOr Service Manuals Contact	
NTD LC NIDOTY TOO I TO TO TO TO THE STORY OF		0683-1525	RIFXU COMP 1500 OHM 5% 1/4W MAURITHON TECHNICAL SERVICES	
0683-6805 RIFXU COMP 68 OHM 5% 1/4w Oxon OX9 4QY	7R15	0683-6805	RIFXI) COMP AR OHM 5% 1/1 OCHORY FIELD TO, CHIRING	
A7R16 0683-2235 PIEXD COMP 22K OHM EX 1774	'R16	0683-2235	RIEXD COMP 22K OHM 58 1/4/4 Tel: 01844-351694 Fax: 01844-352554	
Email:- enquiries@maunitron.co.uk			ALCOMORMULE MEMBERS - HEMBERS	

 $<sup>\ \ \</sup>pm$  See list of abbreviations in introduction to this section

	Stock No.	Description #		Note
A7R17				Note
A7R18	0683-1245	RIFXU COMP 120K OHM 5% 1/4W	1	
A7R18 A7R19	0684-1021	R:FXU COMP 1000 OHM 10% 1/4W		
41113	0683-1535	R:FXD COMP 15K OHM 5% 1/4W	1	
A7RT1	0852-0020			
A7T1		RITEMPERATURE COMPENSATING 150 OHM 5% 250	1	
	9120-0090	TRANSFORMER: INTERMEDIATE FREQUENCY		
A7TB1	00851-2011	BLANK PC BUARD: IF 20MC		
CI	0150-0119	1		
C2	0150-0022	CIFAU CER 2X(0.01 UF) 20% 250VDCW	1	
C3	0150-0019		l	
C 4	0150-0019	C:FXD CER 1000FF 20% 500 VDCW	}	
C 5	0180-0076	C:FXD CER ICOUPF 20% 500 VDCW C:FXD ELECT 20UF 25VDCW	İ	
06	(1150, 001)			
	0150-0019	CIFAU CER 1000 PF 20% 500VDCW FEED-THRU		
)1	1250-0083	CONNECTOR: ENC INPUT TO A2		
.1	9140-0051			
:1		COIL:FXD 400 UH		
	0684-1021	R*FXD COMP 1000 OHM 10% 1/4W		
1	00851-6029	CADLE ASSY.		
2	00851-6030	7-INCH COAX W/MALE RF CONN CABLE ASSY.		
3	00851-6049	9-INCH COAX W/MALE RF CONN		
<b>.</b>	00851-6029	7-INCH COAX W/MALE RE COMM		
5		CABLE ASSY. 7-INCH COAX W/MALE RF CONN		
'	00851-6031	TANCE WOOL		
		22-INCH COAXEAFTEN-VERT DISPLAY SWITCH		
!	00851-6034	CABLE:5-INCH COAR WARE		
1	00851-6033	CABLE:5-INCH COAX VIDEO OUT TO VERT AM CABLE:16-1/2 INCH COAX	ł	
Pi	1250-0229			
P2	1250-0229	CONNECTORING CABLE PLUG SUB-MINIATURE	1	j
P3	1250-0229	CONNECTORING CABLE PLUG SUB-MINIATURE CONNECTORING CABLE PLUG SUB-MINIATURE	1	ł
P4	1250-0229	CONNECTORINE CABLE PLUG SUB-MINIATURE CONNECTORINE CABLE PLUG SUB-MINIATURE	}	- 1
	00851-8003	TABLE PLUG SUB-MINIATURE		
	00851-8003	FILTER:LUM-FASS	1	
	00851-8003	FILTER: LOW-PASS	ł	1
	00851-8003	FILTER: LOW-PASS	1	
<b>1</b>	00851-8003	FILTER:LOW-FASS FILTER:LOW-FASS		
		MISCELLANEOUS		
(	340-0095		1	
	-	TERMINAL : FEEDTHRU TEFLON INSULATED OUTPUT: 20 MC I.F. TO VERT AMPL		
		10 VERT AMPL		
1			1	
1				
	ĺ			
		·		
	1		1	

### Cabinet Parts, Model 851B




Table 6-3. Replaceable Parts (cont'd)

₱ Stock No.	Description #	Mfr.	Mfr. Part No.	ТС
				-
0683-4735	RIFXD COMP 47K OHM 5% 1/4W	1		ŀ
0683-5105	KIPAU COMP 51 OHM 50 1/44	01121	CB 4735	
0683-5125	RIFXD COMP 5100 OHM 59 174W	01121	CB 5105	1
0683-5625	RIFXD COMP 5600 OHM 5% 1/4W	01121	CB 5125	1
0683-6805	RIFXD COMP 68 OHM 5% 1/4W	01121	CB 5625	2
0683-6825		01121	CB 6805	3
0683-9115	RIFXD COMP 6800 OHM 5% 174W	01121	CB 6825	
0684-1001	R:FXD COMP 910 OHM 5% 1/4W	01121	CB 9115	9
0684-1011	RIFXD COMP 10 OHM 10% 1/4W	01121	CB 1001	1
0684-1021	R:FXD COMP 100 OHM 10% 1/4W R:FXD COMP 1000 OHM 10% 1/4W	01121	CB 1011	1
0684-1031	1	01121	CB 1021	2
0684-4711	RIFXD COMP 10K OHM 10% 1/4W	01121	C6 1031	
0686-1225	R:FXD COMP 470 OHM 10% 1/4#	01121	CB 4711	1
0686-1825	R:FXD COMP 1200 OHM 5% 1/2#	01121	EB 1225	1
0686-2025	RIFXD COMP 1800 OHM 5% 1/2W	01121	EB 1825	1
	R:FXD COMP 2000 OHM 5% 1/2W	01121	Lb 2025	1 1
0686 <b>-</b> 2725 0686 <b>-</b> 4725	R:FXD COMP 2700 OHM 5% 1/2W	1 1		1
0686-5115	I RIPAU COMP 4700 OHM 5% 1/2+	01121	B 2725	1
0686-5625	I KIPAU COMP 510 OHM 58 123W	01121   01121   0	8 4725	1
0686-6225	RIFXD COMP 5600 OHM 5% 1/2W	01121	5115	1
	R:FXD COMP 6200 OHM 5% 1/2W	01121 E	10 5625 18 6225	1
686-7525	RIFXD COMP 7500 OHM 5% 172W	1 1		1
0687-1011	R:FXD COMP 100 OHM 10% 172W	01121 E	.b 7525	
687-1021	RIFXD COMP 1000 OHM 10% 1/2W	01121   E	.B 1011	1
687-1031	KIFKU COMP TOK OHM TOK 1/2#	01121 E	B 1021	1 1
687-1041	RIFXD COMP 100K OHM 10x 1/2W	C1121 E	8 1031	i
687-1051	RIFXD COMP 1 MEGOHM 10% 1/2W	01121	D 1041	1
687-1231	RIFXD COMP 12K OHM 10% 1/2W	01121 E	8 1051	1.1
687-1521	R:FXD COMP 1500 OHM 10% 1/2W	01121 E	B 1231	1
687-2211	RIFXD COMP 220 OHM 10% 1/2W	01121 E	B 1521	2
687-2741	R:FXD COMP 270K OHM 10% 1/2W	01121 E	B 2211	2
687-3311	R:FXD COMP 330 OHM 10% 1/2W	01121	3 2741	1
687-3321	RIFXD COMP 330C OHM 10% 1/2W	01121 E	3 3311	1 . 1
587-3921	RIFXD COMP 3900 OHM 10% 1/2W	01121 E	3321	1
587-3931	R:FXD COMP 39K OHM 10% 172W	01121 EE	3921	2
887-4721	R:FXD COMP 4700 OHM 10% 1/2W	01121   EB	3931	1 1
87-4731		01121 28	4721	1
87-5611	R:FXD COMP 47K OHM 10% 1/2W R:FXD COMP 560 OHM 10% 1/2W	01121 EB	4731	1.1
87-5621	RIFXD COMP 5600 OHM 10% 1/2W	01121 EB	5611	1
87-6801	R:FXD COMP 68 CHM 10% 1/2W	01121 EB	5621	2
87-8221	R:FXD COMP 8200 OHM 10% 1/2W	01121 EB	6801	1
		01121 EB	8221	3
90-1231	RIFXD COMP 12K OHM 10% 1W	} }		-
90-1851	KIFXU COMP 1.8 MEGOHM LOW IW	01121 GB	1231	11
90-2721	RIFIND COMP 2700 OHM 10% 1W	01121 GB	1851	4
90-3951	KIFXU COMP 3.9 MEGOHM 10x 1w	01121 GB	2721	i
90-8241	RIFXD COMP 820K OHM 10% 1W	01121 GB 01121 GB	3951 8241	5
99-0005	R:FXD COMP 2.7 OHM 10% 1W	i i		1
27-0004	KIFXU DEPC 5 OHM 18 1/2W	01121 GB	27G1	1
27-0005	RIFXU DEPC 5.77 OHM 1/2% 1/2%	28480 072	27-0004	2
27-0008	R*FAU DEPC 11.61 OHM 1/24 1/24	28480 072	27-0005	1
27-0010	R:FXD DEPC 17.61 OHM 1/2% 1/2W	28480 072	7-0008	i
	· <del>-</del>	28480 072	.7-0010	1
	_			

Table 6-3. Replaceable Parts (cont'd)

0727-0025	480 C72' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072'	7-0025 7-0033 7-0036 7-0042 7-0048 7-0062 7-0063 7-0074 7-0094	1 2 1 2 1 2 1 2 2 2 1 2 2
0727-0025         RIFXD DEPC 51.5 OHM 1% 1/2W         268           0727-0036         RIFXD DEPC 61.11 OHM 1% 1/2W         288           0727-0042         RIFXD DEPC 71.16 OHM 1% 1/2W         288           0727-0048         RIFXD DEPC 96.25 OHM 1/2% 1/2W         289           0727-0062         RIFXD DEPC 247.5 OHM 1/2% 1/2W         289           0727-0063         RIFXD DEPC 249.5 OHM 1/2 1/2W         281           0727-0074         RIFXD DEPC 436 OHM 1/2% 1/2W         281           0727-0074         RIFXD DEPC 870 OHM 1/2% 1/2W         281           0727-0101         RIFXD DEPC 870 OHM 1/2 1/2W         281           0727-0122         RIFXD DEPC 2-51K OHM 1% 1/2W         261           0727-0123         RIFXD DEPC 3000 OHM 1% 1/2W         197           0727-0124         RIFXD DEPC 30-26 OHM 1% 1/2W         197           0727-0125         RIFXD DEPC 30-26 OHM 1% 1/2W         197           0727-0126         RIFXD DEPC 30-3K OHM 1% 1/2W         264           0727-0136         RIFXD DEPC 30-26 OHM 1% 1/2W         264           0727-0163         RIFXD DEPC 10-16 OHM 1% 1/2W         264           0727-0170         RIFXD DEPC 30-3K OHM 1% 1/2W         264           0727-0173         RIFXD DEPC 30-3K OHM 1% 1/2W         197	480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072'	7-0025 7-0033 7-0036 7-0042 7-0048 7-0062 7-0063 7-0074 7-0094	2 2 2 2 2 2
0727-0036         RIFXD DEPC 61.11 OHM 1x 1/2w         288           0727-0042         RIFXD DEPC 76.25 OHM 1/2x 1/2w         281           0727-0042         RIFXD DEPC 96.25 OHM 1/2x 1/2w         281           0727-0048         RIFXD DEPC 150 OHM 1x 1/2w         281           0727-0063         RIFXD DEPC 247.5 OHM 1x 1/2w         281           0727-0074         RIFXD DEPC 436 OHM 1/2x 1/2w         281           0727-0074         RIFXD DEPC 436 OHM 1/2x 1/2w         281           0727-0074         RIFXD DEPC 370 OHM 1x 1/2w         281           0727-0101         RIFXD DEPC 370 OHM 1x 1/2w         281           0727-0122         RIFXD DEPC 290C OHM 1x 1/2w         281           0727-0123         RIFXD DEPC 300C OHM 1x 1/2w         197           0727-0124         RIFXD DEPC 300C OHM 1x 1/2w         197           0727-0125         RIFXD DEPC 300C OHM 1x 1/2w         197           0727-0126         RIFXD DEPC 300C OHM 1x 1/2w         197           0727-0127         RIFXD DEPC 300C OHM 1x 1/2w         197           0727-0136         RIFXD DEPC 300C OHM 1x 1/2w         197           0727-0136         RIFXD DEPC 300C OHM 1x 1/2w         284           0727-0158         RIFXD DEPC 300C OHM 1x 1/2w         284	480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072:	7-0033 7-0036 7-0042 7-0048 7-0062 7-0063 7-0074 7-0094	2 2 2 2 2 2
0727-0036         RIFKO DEPC 61.11 OHM 1x 1/2w         28           0727-0042         RIFKO DEPC 76.16 OHM 1x 1/2w         28           0727-0042         RIFKO DEPC 96.25 OHM 1/2w 1/2w         28           0727-0048         RIFKO DEPC 150 OHM 1x 1/2w         28           0727-0063         RIFKO DEPC 292.5 OHM 1x 1/2w         28           0727-0074         RIFKO DEPC 436 OHM 1/2x 1/2w         28           0727-0074         RIFKO DEPC 436 OHM 1x 1/2w         28           0727-0074         RIFKO DEPC 10.3K OHM 1x 1/2w         28           0727-0107         RIFKO DEPC 25.0K OHM 1x 1/2w         28           0727-0122         RIFKO DEPC 290C OHM 1x 1/2w         28           0727-0123         RIFKO DEPC 30.0K OHM 1x 1/2w         28           0727-0124         RIFKO DEPC 30.0K OHM 1x 1/2w         19           0727-0125         RIFKO DEPC 30.0K OHM 1x 1/2w         19           0727-0126         RIFKO DEPC 30.0K OHM 1x 1/2w         19           0727-0136         RIFKO DEPC 30.0K OHM 1x 1/2w         26           0727-0136         RIFKO DEPC 30.0K OHM 1x 1/2w         28           0727-0158         RIFKO DEPC 30.0K OHM 1x 1/2w         28           0727-0170         RIFKO DEPC 30.0K OHM 1x 1/2w         28           0727-0189<	480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072: 480 072:	7-0033 7-0036 7-0042 7-0048 7-0062 7-0063 7-0074 7-0094	2 1 2 1 2 2 2
0727-0036         RifxD DEPC 71.16 OHM 1x 1/2w         28.           0727-0042         RifxD DEPC 96.25 OHM 1/2a 1/2w         28.           0727-0048         RifxD DEPC 150 OHM 1/2w 1/2w         28.           0727-0062         RifxD DEPC 247.5 OHM 1/2x 1/2w         28.           0727-0063         RifxD DEPC 243.6 OHM 1/2x 1/2w         28.           0727-0094         RifxD DEPC 436 OHM 1/2x 1/2w         28.           0727-01094         RifxD DEPC 436 OHM 1/2x 1/2w         28.           0727-0121         RifxD DEPC 436 OHM 1x 1/2w         28.           0727-0122         RifxD DEPC 290. OHM 1x 1/2w         28.           0727-0123         RifxD DEPC 290. OHM 1x 1/2w         197.           0727-0124         RifxD DEPC 30.00 OHM 1x 1/2w         197.           0727-0126         RifxD DEPC 30.00 OHM 1x 1/2w         197.           0727-0136         RifxD DEPC 3.3x OHM 1x 1/2w         284.           0727-0153         RifxD DEPC 11.6x OHM 1x 1/2w         284.           0727-0170         RifxD DEPC 11.6x OHM 1x 1/2w         197.           0727-0173         RifxD DEPC 31.5x OHM 1x 1/2w         197.           0727-0189         RifxD DEPC 3.7x OHM 1x 1/2w         197.           0727-0189         RifxD DEPC 3.7x OHM 1x 1/2w         284.	480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 480 072'	7-0036 7-0042 7-0048 7-0062 7-0063 7-0074 7-0094	2 1 2 2 2 2
0727-0042         R:FXD DEPC 96.25 OHM 1/26 1/2W         284           0727-0048         R:FXD DEPC 150 OHM 1% 1/2W         284           0727-0063         R:FXD DEPC 247.5 OHM 1/2% 1/2W         284           0727-0074         R:FXD DEPC 436 OHM 1/2% 1/2W         284           0727-0094         R:FXD DEPC 436 OHM 1/2% 1/2W         284           0727-0101         R:FXD DEPC 1.03K OHM 1% 1/2W         284           0727-0122         R:FXD DEPC 2.51K OHM 1% 1/2W         284           0727-0123         R:FXD DEPC 3.00C OHM 1% 1/2W         284           0727-0124         R:FXD DEPC 30.00 OHM 1% 1/2W         197           0727-0125         R:FXD DEPC 30.00 OHM 1% 1/2W         197           0727-0126         R:FXD DEPC 3.03K OHM 1% 1/2W         197           0727-0136         R:FXD DEPC 5.03K OHM 1% 1/2W         284           0727-0158         R:FXD DEPC 5.03K OHM 1% 1/2W         284           0727-0170         R:FXD DEPC 11.68K OHM 1% 1/2W         197           0727-0173         R:FXD DEPC 24.7K OHM 1% 1/2W         197           0727-0173         R:FXD DEPC 24.7K OHM 1% 1/2W         284           0727-0403         R:FXD DEPC 57.46K OHM 1% 1/2W         284           0727-0405         R:FXD MET FLM 75D OHM 1% 1/2W         284      <	480 072' 480 072' 480 072' 480 072' 480 072' 480 072' 701 0C1, 701 0C1, 480 072' 480 072' 480 072' 480 072'	7-0042 7-0048 7-0062 7-0063 7-0074 7-0094	2 2 2 2 2
0727-0062 0727-0063 RIFXD DEPC 247.5 OHM 1/2% 1/2% 284 0727-0063 RIFXD DEPC 292.5 OHM 1% 1/2% 284 0727-0094 RIFXD DEPC 870 OHM 1/2 1/2% 284 0727-0094 RIFXD DEPC 870 OHM 1/2 1/2% 284 0727-0101 RIFXD DEPC 1.03K OHM 1% 1/2% 284 0727-0122 RIFXD DEPC 2.51K OHM 1% 1/2% 284 0727-0123 RIFXD DEPC 2.51K OHM 1% 1/2% 284 0727-0124 RIFXD DEPC 290C OHM 1% 1/2% 197 0727-0126 RIFXD DEPC 300C OHM 1% 1/2% 197 0727-0126 RIFXD DEPC 30.0K OHM 1% 1/2% 197 0727-0136 RIFXD DEPC 3.266K OHM 1% 1/2% 284 0727-0163 RIFXD DEPC 10.1K OHM 1% 1/2% 284 0727-0163 RIFXD DEPC 10.1K OHM 1% 1/2% 284 0727-0170 RIFXD DEPC 10.K OHM 1% 1/2% 284 0727-0173 RIFXD DEPC 20K OHM 1% 1/2% 284 0727-0178 RIFXD DEPC 24.7K OHM 1% 1/2% 284 0727-0398 RIFXD DEPC 24.7K OHM 1% 1/2% 284 0727-0403 RIFXD DEPC 3.79K OHM 1/2 1/2% 197 0727-0403 RIFXD DEPC 52.3K OHM 1/2 1/2% 197 0757-0817 RIFXD DEPC 52.3K OHM 1/2 1/2% 197 0757-0831 RIFXD DEPC 52.3K OHM 1/2 1/2% 197 0757-0831 RIFXD MET FLM 75D OHM 1 1/2% 284 0757-0885 RIFXD MET FLM 15.4K OHM 1% 1/4% 284 0757-0887 RIFXD MET FLM 15.4K OHM 1% 1/4% 284 0758-0002 RIFXD MET FLM 15.4K OHM 1% 1/4% 284 0758-0003 RIFXD MET FLM 191K OHM 1% 1/4% 284 0758-0004 RIFXD MET FLM 190C OHM 5% 1/2% 071 0758-0012 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0012 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0024 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0034 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0034 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0034 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 2000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0047 RIFXD MET FLM 2000 OHM 5% 1/2% 071 0758-0047 RIFXD MET FLM 2000 OHM 5% 1/2% 071 0758-0047 RIFXD MET FLM 2000 OHM 5% 1/2% 071	480 072' 480 072' 480 072' 480 072' 480 072' 701 DC1, 701 DC1, 701 DC1, 480 072' 480 072'	7-0062 7-0063 7-0074 7-0094	2 2 2
0727-0062 0727-0063 RIFXD DEPC 247.5 OHM 1/2% 1/2% 284 0727-0074 RIFXD DEPC 292.5 OHM 1% 1/2% 284 0727-0074 RIFXD DEPC 496 OHM 1/2% 1/2% 284 0727-0094 RIFXD DEPC 870 OHM 1/2% 1/2% 284 0727-0101 RIFXD DEPC 1.03K OHM 1% 1/2% 284 0727-0122 RIFXD DEPC 2.51K OHM 1% 1/2% 284 0727-0123 RIFXD DEPC 2.51K OHM 1% 1/2% 197 0727-0124 RIFXD DEPC 290C OHM 1% 1/2% 197 0727-0126 RIFXD DEPC 300C OHM 1% 1/2% 197 0727-0126 RIFXD DEPC 300C OHM 1% 1/2% 197 0727-0136 RIFXD DEPC 3.0266K OHM 1% 1/2% 284 0727-0138 RIFXD DEPC 3.0266K OHM 1% 1/2% 284 0727-0163 RIFXD DEPC 10.1K OHM 1% 1/2% 284 0727-0163 RIFXD DEPC 10.1K OHM 1% 1/2% 284 0727-0170 RIFXD DEPC 10.K OHM 1% 1/2% 284 0727-0173 RIFXD DEPC 10.K OHM 1% 1/2% 284 0727-0179 RIFXD DEPC 10.K OHM 1% 1/2% 284 0727-0179 RIFXD DEPC 24.7K OHM 1% 1/2% 284 0727-0398 RIFXD DEPC 24.7K OHM 1% 1/2% 284 0727-0398 RIFXD DEPC 3.79K OHM 1/2% 1/2% 284 0727-0403 RIFXD DEPC 52.3K OHM 1/2% 1/2% 197 0757-0817 RIFXD DEPC 52.3K OHM 1/2% 1/2% 197 0757-0831 RIFXD DEPC 52.3K OHM 1/2% 1/2% 197 0757-0831 RIFXD DEPC 52.3K OHM 1/2% 1/2% 197 0757-0831 RIFXD MET FLM 75D OHM 1% 1/4% 284 0757-0889 RIFXD MET FLM 15.4K OHM 1% 1/4% 284 0757-0889 RIFXD MET FLM 13.24K OHM 1% 1/4% 284 0758-0002 RIFXD MET FLM 191K OHM 1% 1/4% 284 0758-0003 RIFXD MET FLM 191K OHM 1% 1/4% 284 0758-0004 RIFXD MET FLM 1900 OHM 5% 1/2% 071 0758-0012 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0024 RIFXD MET FLM 1200 OHM 5% 1/2% 071 0758-0038 RIFXD MET FLM 1200 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 2400 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 1000 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 2400 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 2400 OHM 5% 1/2% 071 0758-0044 RIFXD MET FLM 2400 OHM 5% 1/2% 071 0758-0047 RIFXD MET FLM 2400 OHM 5% 1/2% 071 0758-0047	480 072' 480 072' 480 072' 480 072' 480 072' 701 DC1, 701 DC1, 701 DC1, 480 072' 480 072'	7-0062 7-0063 7-0074 7-0094	2 2 2
0727-0063         RIFXD DEPC 292.5 OHF 1% 1/2%         284           0727-0074         RIFXD DEPC 436 OHM 1/2% 1/2%         284           0727-0094         RIFXD DEPC 870 OHM 1/2% 1/2%         284           0727-0103         RIFXD DEPC 1.03K OHM 1% 1/2%         264           0727-0123         RIFXD DEPC 290C OHM 1% 1/2%         197           0727-0124         RIFXD DEPC 300C OHM 1% 1/2%         197           0727-0126         RIFXD DEPC 300C OHM 1% 1/2%         197           0727-0136         RIFXD DEPC 300K OHM 1% 1/2%         284           0727-0136         RIFXD DEPC 10.1K OHM 1% 1/2%         284           0727-0158         RIFXD DEPC 11.68K OHM 1% 1/2%         284           0727-0163         RIFXD DEPC 11.68K OHM 1% 1/2%         284           0727-0170         RIFXD DEPC 11.68K OHM 1% 1/2%         197           0727-0173         RIFXD DEPC 24.7K OHM 1% 1/2%         284           0727-0174         RIFXD DEPC 24.7K OHM 1% 1/2%         284           0727-0189         RIFXD DEPC 35.79K OHM 1/2% 1/2%         284           0727-0403         RIFXD DEPC 53.79K OHM 1/2% 1/2%         284           0727-0405         RIFXD MET FLM 750 OHM 1% 1/2% 1/2%         284           0757-0887         RIFXD MET FLM 750 OHM 1% 1/2% 1/2%         284	480 072' 480 072' 480 072' 480 072' 701 001, 701 001, 480 072' 480 072'	7-0063 7-0074 7-0094 7-0101	2 2 2
0727-0074         RIFXD DEPC 870 OHM 1/2x 1/2w         284           0727-0094         RIFXD DEPC 870 OHM 1/2x 1/2w         284           0727-0101         RIFXD DEPC 1.03K OHM 1x 1/2w         264           0727-0122         RIFXD DEPC 2.51K OHM 1x 1/2w         297           0727-0123         RIFXD DEPC 290C OHM 1x 1/2w         197           0727-0124         RIFXD DEPC 300C OHM 1x 1/2w         197           0727-0126         RIFXD DEPC 3.2c6K OHM 1x 1/2w         197           0727-0136         RIFXD DEPC 5.03K OHM 1x 1/2w         264           0727-0163         RIFXD DEPC 10.1K OHM 1x 1/2w         284           0727-0163         RIFXD DEPC 11c8K OHM 1x 1/2w         197           0727-0170         RIFXD DEPC 20K OHM 1x 1/2w         197           0727-0173         RIFXD DEPC 24-7K OHM 1x 1/2w         284           0727-0178         RIFXD DEPC 24-7K OHM 1x 1/2w         284           0727-0189         RIFXD DEPC 37-7K OHM 1x 1/2w         284           0727-0403         RIFXD DEPC 57-46K OHM 1x 1/2w         284           0757-0807         RIFXD MET FLM 75D OHM 1x 1/2w         284           0757-0807         RIFXD MET FLM 15-4K OHM 1x 1/2w         284           0757-0808         RIFXD MET FLM 15-4K OHM 1x 1/2w         284      <	480 072' 480 072' 480 072' 480 072' 701 DC1, 701 DC1, 480 072' 480 072'	7-0074 7-0094 7-0101	2 2
0727-0094         R:FXD DEPC 870 OHM 1/2* 1/2*         284           0727-0101         R:FXD DEPC 1.03K OHM 1% 1/2*         264           0727-0122         R:FXD DEPC 2.51K OHM 1% 1/2*         284           0727-0123         R:FXD DEPC 290C OHM 1% 1/2*         197           0727-0124         R:FXD DEPC 300C OHM 1% 1/2*         197           0727-0126         R:FXD DEPC 300C OHM 1% 1/2*         197           0727-0136         R:FXD DEPC 3.266K OHM 1% 1/2*         284           0727-0158         R:FXD DEPC 10.1K OHM 1% 1/2*         284           0727-0163         R:FXD DEPC 11.68K OHM 1% 1/2*         284           0727-0173         R:FXD DEPC 12.68K OHM 1% 1/2*         284           0727-0173         R:FXD DEPC 20K OHM 1% 1/2*         197           0727-0189         R:FXD DEPC 24.7K OHM 1% 1/2*         197           0727-0189         R:FXD DEPC 3.79K OHM 1/2* 1/2*         284           0727-0403         R:FXD DEPC 3.79K OHM 1/2* 1/2*         284           0727-0403         R:FXD MET FLM 75D OHM 1* 1/2*         197           0757-0817         R:FXD MET FLM 75D OHM 1* 1/2*         284           0757-0887         R:FXD MET FLM 15.4K OHM 1* 1/2*         284           0757-0889         R:FXD MET FLM 191K OHM 1* 1/2*         284	480 072* 480 072* 480 072* 701 001, 701 001, 701 001, 480 072* 480 072*	7-0094 7-0101	2
0727-0122         R:FXD DEFC 2.51K OHM 1x 1/2W         284           0727-0123         R:FXD DEFC 290C OHM 1x 1/2W         197           0727-0124         R:FXD DEFC 300C OHM 1x 1/2W         197           0727-0126         R:FXD DEFC 3.266K OHM 1x 1/2W         197           0727-0136         R:FXD DEFC 10.1K OHM 1x 1/2W         264           0727-0158         R:FXD DEFC 11.68K OHM 1x 1/2W         284           0727-0163         R:FXD DEFC 11.68K OHM 1x 1/2W         197           0727-0170         R:FXD DEFC 11.68K OHM 1x 1/2W         197           0727-0173         R:FXD DEFC 20K OHM 1x 1/2W         197           0727-0178         R:FXD DEFC 20K OHM 1x 1/2W         284           0727-0189         R:FXD DEFC 3.7K OHM 1x 1/2W         284           0727-0189         R:FXD DEFC 3.79K OHM 1/2 1/2W         284           0727-0403         R:FXD DEFC 52.3K OHM 1/2 1/2W         197           0727-0405         R:FXD MET FLM 4.32K OHM 1/2 1/2W         197           0757-0817         R:FXD MET FLM 4.32K OHM 1x 1/4W         284           0757-0885         R:FXD MET FLM 13.4K OHM 1x 1/4W         284           0757-0887         R:FXD MET FLM 13K OHM 1x 1/4W         284           0758-0002         R:FXD MET FLM 191K OHM 1x 1/2W         284 <td>480 072 701 001, 701 001, 701 001, 480 072 480 072</td> <td></td> <td></td>	480 072 701 001, 701 001, 701 001, 480 072 480 072		
0727-0122         R:FXD DEFC 2.51K OHM 1x 1/2W         284           0727-0123         R:FXD DEFC 290C OHM 1x 1/2W         197           0727-0124         R:FXD DEFC 300C OHM 1x 1/2W         197           0727-0126         R:FXD DEFC 3.266K OHM 1x 1/2W         197           0727-0136         R:FXD DEFC 10.1K OHM 1x 1/2W         264           0727-0158         R:FXD DEFC 11.68K OHM 1x 1/2W         284           0727-0163         R:FXD DEFC 11.68K OHM 1x 1/2W         197           0727-0170         R:FXD DEFC 11.68K OHM 1x 1/2W         197           0727-0173         R:FXD DEFC 20K OHM 1x 1/2W         197           0727-0178         R:FXD DEFC 20K OHM 1x 1/2W         284           0727-0189         R:FXD DEFC 3.7K OHM 1x 1/2W         284           0727-0189         R:FXD DEFC 3.79K OHM 1/2 1/2W         284           0727-0403         R:FXD DEFC 52.3K OHM 1/2 1/2W         197           0727-0405         R:FXD MET FLM 4.32K OHM 1/2 1/2W         197           0757-0817         R:FXD MET FLM 4.32K OHM 1x 1/4W         284           0757-0885         R:FXD MET FLM 13.4K OHM 1x 1/4W         284           0757-0887         R:FXD MET FLM 13K OHM 1x 1/4W         284           0758-0002         R:FXD MET FLM 191K OHM 1x 1/2W         284 <td>480 072 701 001, 701 001, 701 001, 480 072 480 072</td> <td></td> <td>1 1</td>	480 072 701 001, 701 001, 701 001, 480 072 480 072		1 1
0727-0123 R:FXD DEPC 2900 OHM 1% 1/2% 197 0727-0124 R:FXD DEPC 3000 OHM 1% 1/2% 197 0727-0126 R:FXD DEPC 3.06K OHM 1% 1/2% 197 0727-0136 R:FXD DEPC 5.03K OHM 1% 1/2% 264 0727-0158 R:FXD DEPC 5.03K OHM 1% 1/2% 284 0727-0163 R:FXD DEPC 11.68K OHM 1% 1/2% 284 0727-0163 R:FXD DEPC 11.68K OHM 1% 1/2% 284 0727-0170 R:FXD DEPC 18K OHM 1% 1/2% 197 0727-0173 R:FXD DEPC 20K OHM 1% 1/2% 284 0727-0173 R:FXD DEPC 24.7K OHM 1% 1/2% 284 0727-0189 R:FXD DEPC 41.7K OHM 1% 1/2% 284 0727-0189 R:FXD DEPC 41.7K OHM 1% 1/2% 284 0727-0403 R:FXD DEPC 52.3K OHM 1/2% 1/2% 284 0727-0405 R:FXD DEPC 57.46K OHM 1/2% 1/2% 197 0757-0817 R:FXD MET FLM 75D OHM 1% 1/2% 284 0757-0831 R:FXD MET FLM 4.32K OHM 1% 1/4% 284 0757-0889 R:FXD MET FLM 32.4K OHM 1% 1/4% 284 0757-0889 R:FXD MET FLM 191K OHM 1% 1/4% 284 0757-0899 R:FXD MET FLM 191K OHM 1% 1/4% 284 0758-0002 R:FXD MET FLM 191K OHM 1% 1/4% 284 0758-0003 R:FXD MET FLM 191K OHM 1% 1/2% 284 0758-0004 R:FXD MET FLM 1000 OHM 5% 1/2% 071 0758-0012 R:FXD MET FLM 1000 OHM 5% 1/2% 075 0758-0012 R:FXD MET FLM 100 OHM 5% 1/2% 075 0758-0024 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0034 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0034 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0047 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0048 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0049 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0049 R:FXD MET FLM 200 OHM 5% 1/2% 071 0758-0049 R:FXD MET FLM 200 OHM 5% 1/2% 071 0758-0049 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0049 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0049 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0049 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0047 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0047 R:FXD MET FLM 100 OHM 5% 1/2% 071 0758-0047 R:FXD MET FLM 200 OHM 5% 1/2% 071	701 UC1, 701 UC1, 701 UC1, 480 U721	7-0122	1 . 1
0727-0124 0727-0126 R:FXD DEPC 3000 OHM 1% 1/2% 197 0727-0136 R:FXD DEPC 3.266K OHM 1% 1/2% 0727-0158 R:FXD DEPC 10.1K OHM 1% 1/2% 0727-0163 R:FXD DEPC 11.68K OHM 1% 1/2% 0727-0170 R:FXD DEPC 18K OHM 1% 1/2% 0727-0173 R:FXD DEPC 18K OHM 1% 1/2% 0727-0173 R:FXD DEPC 20K OHM 1% 1/2% 0727-0178 R:FXD DEPC 24.7K OHM 1% 1/2% 0727-0189 R:FXD DEPC 24.7K OHM 1% 1/2% 0727-0403 R:FXD DEPC 3.79K OHM 1/2% 1/2% 0727-0405 R:FXD DEPC 3.79K OHM 1/2% 1/2% 197 0727-0405 R:FXD DEPC 57.46K OHM 1/2% 1/2% 197 0757-0817 R:FXD MET FLM 750 OHM 1 1/2% 0757-0887 R:FXD MET FLM 15.4K OHM 1% 1/4% 284 0757-0889 R:FXD MET FLM 13.4K OHM 1% 1/4% 284 0757-0890 R:FXD MET FLM 191K OHM 1% 1/4% 284 0758-0002 R:FXD MET FLM 191K OHM 1% 1/2% 0758-0005 R:FXD MET FLM 2700 OHM 5% 1/2% 0758-0001 R:FXD MET FLM 1500 OHM 5% 1/2% 0758-0017 R:FXD MET FLM 1500 OHM 5% 1/2% 0758-0017 R:FXD MET FLM 1500 OHM 5% 1/2% 0758-0024 R:FXD MET FLM 1500 OHM 5% 1/2% 0758-0038 R:FXD MET FLM 1600 OHM 5% 1/2% 0758-0034 R:FXD MET FLM 1600 OHM 5% 1/2% 0758-0034 R:FXD MET FLM 1200 OHM 5% 1/2% 0758-0034 R:FXD MET FLM 1600 OHM 5% 1/2% 0758-0034 R:FXD MET FLM 1600 OHM 5% 1/2% 0758-0044 R:FXD MET FLM 1600 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1600 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 1800 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 2000 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 2000 OHM 5% 1/2% 0758-0047 R:FXD MET FLM 2000 OHM 5% 1/2% 0758-0047	701 DC1, 701 DC1, 480 0721		1
0727-0126         R:FXD DEPC 3.266K OHM 1% 1/2W         195           0727-0136         R:FXD DEPC 5.03K OHM 1% 1/2W         264           0727-0158         R:FXD DEPC 10.1K OHM 1% 1/2W         284           0727-0163         R:FXD DEPC 11.68K OHM 1% 1/2W         284           0727-0170         R:FXD DEPC 20K OHM 1% 1/2W         197           0727-0173         R:FXD DEPC 20K OHM 1% 1/2W         197           0727-0189         R:FXD DEPC 41.7K OHM 1% 1/2W         197           0727-0189         R:FXD DEPC 3.79K OHM 1/2 1/2W         284           0727-0403         R:FXD DEPC 5.23K OHM 1/2 1/2W         197           0727-0405         R:FXD DEPC 57.46K OHM 1/2 1/2W         197           0757-0817         R:FXD MET FLM 750 OHM 1 1/2 1/2W         197           0757-0885         R:FXD MET FLM 15.4K OHM 1% 1/4W         284           0757-0887         R:FXD MET FLM 191K OHM 1% 1/4W         284           0757-0889         R:FXD MET FLM 191K OHM 1% 1/4W         284           0758-0002         R:FXD MET FLM 191K OHM 5% 1/2W         071           0758-0003         R:FXD MET FLM 191K OHM 5% 1/2W         071           0758-0004         R:FXD MET FLM 100 OHM 5% 1/2W         071           0758-0012         R:FXD MET FLM 100 OHM 5% 1/2W         071	701 UC1, 480 072 480 072		1
0727-0136 0727-0158 0727-0158 0727-0163 0727-0163 0727-0163 0727-0170 0727-0170 0727-0170 0727-0173 0727-0173 0727-0173 0727-0173 0727-0173 0727-0173 0727-0173 0727-0173 0727-0173 0727-0173 0727-0173 0727-0173 0727-0173 0727-0178 0727-0178 0727-0189 0727-0189 0727-0403 0727-0405 0727-0405 0727-0405 0727-0405 0727-0405 0727-0887 0727-0887 0727-0887 0727-0889 0727-0	480 0721 480 0721		1
0727-0158         R:FXU DEPC 10.1K OHM 1% 1/2%         284           0727-0163         R:FXU DEPC 11.68K OHM 1% 1/2%         284           0727-0170         R:FXD DEPC 18K OHM 1% 1/2%         197           0727-0173         R:FXD DEPC 20K OHM 1% 1/2%         197           0727-0178         R:FXD DEPC 24.7K OHM 1% 1/2%         197           0727-0189         R:FXD DEPC 41.7K OHM 1% 1/2%         284           0727-0398         R:FXD DEPC 3.79K OHM 1/2% 1/2%         284           0727-0403         R:FXD DEPC 52.3K OHM 1/2% 1/2%         197           0727-0405         R:FXD MET FLM 75D OHM 1/2% 1/2%         197           0757-0817         R:FXD MET FLM 75D OHM 1/2% 1/2%         284           0757-0831         R:FXD MET FLM 15.4K OHM 1/2 1/2%         284           0757-0886         R:FXD MET FLM 15.4K OHM 1/2 1/4%         284           0757-0887         R:FXD MET FLM 13K OHM 1/2 1/4%         284           0757-0889         R:FXD MET FLM 191K OHM 1/2 1/4%         284           0758-0002         R:FXD MET FLM 191K OHM 1/2 1/4%         284           0758-0003         R:FXD MET FLM 2700 OHM 5/2 1/2%         071           0758-0004         R:FXD MET FLM 1500 OHM 5/2 1/2%         071           0758-0017         R:FXD MET FLM 1800 OHM 5/2 1/2%         <	480 072	/2A	1
0727-0163     R:FXD DEPC 11.68K OHM 1% 1/2%     0727-0170     R:FXD DEPC 12K OHM 1% 1/2%     197     0727-0173     R:FXD DEPC 20K OHM 1% 1/2%     197     0727-0178     R:FXD DEPC 20K OHM 1% 1/2%     0727-0178     R:FXD DEPC 24.7K OHM 1% 1/2%     0727-0189     R:FXD DEPC 41.7K OHM 1% 1/2%     0727-0398     R:FXD DEPC 3.79K OHM 1/2% 1/2%     0727-0403     R:FXD DEPC 52.3K OHM 1/2% 1/2%     0727-0405     R:FXD DEPC 52.3K OHM 1/2% 1/2%     197     0757-0817     R:FXD MET FLM 75D OHM 1* 1/2%     0757-0831     R:FXD MET FLM 4.32K OHM 1% 1/2%     0757-0885     R:FXD MET FLM 15.4K OHM 1% 1/4%     0757-0887     R:FXD MET FLM 32.4K OHM 1% 1/4%     284     0757-0889     R:FXD MET FLM 143K OHM 1% 1/4%     284     0757-0890     R:FXD MET FLM 191K OHM 1% 1/4%     0758-0002     R:FXD MET FLM 191K OHM 5% 1/2%     0758-0005     R:FXD MET FLM 2700 OHM 5% 1/2%     0758-0005     R:FXD MET FLM 1500 OHM 5% 1/2%     0758-0017     R:FXD MET FLM 1500 OHM 5% 1/2%     0758-0024     R:FXD MET FLM 1500 OHM 5% 1/2%     0758-0024     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0034     R:FXD MET FLM 2400 OHM 5% 1/2%     0758-0034     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0034     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0034     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0034     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0034     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0034     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0034     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0034     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0044     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0043     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0044     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0047     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0047     R:FXD MET FLM 1000 OHM 5% 1/2%     0758-0047     R:FXD MET FLM 1000 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710 OHM 5% 1/2%     0710			2
0727-0170 0727-0173 R:FXD DEPC 18K OHM 1% 1/2W 0727-0173 R:FXD DEPC 20K OHM 1% 1/2W 0727-0189 R:FXD DEPC 24.7K OHM 1% 1/2W 0727-0398 R:FXD DEPC 41.7K OHM 1% 1/2W 0727-0398 R:FXD DEPC 3.79K OHM 1/2% 1/2W 0727-0403 R:FXD DEPC 52.3K OHM 1/2% 1/2W 197 0727-0405 R:FXD DEPC 57.46K OHM 1/2% 1/2W 197 0757-0817 R:FXD MET FLM 750 OHM 1 1/2 1/2 2W 0757-0831 R:FXD MET FLM 4.32K OHM 1 1/2 1/2 2W 0757-0885 R:FXD MET FLM 32.4K OHM 1 1/2 1/4W 0757-0887 R:FXD MET FLM 191K OHM 1 1/2 1/4 2W 0757-0889 R:FXD MET FLM 143K OHM 1 1/2 1/4 2W 0758-0002 R:FXD MET FLM 191K OHM 1/2 1/2 2W 0758-0003 R:FXD MET FLM 191K OHM 1/2 1/2 2W 0758-0004 R:FXD MET FLM 1000 OHM 5/2 1/2 2W 0758-0005 R:FXD MET FLM 1000 OHM 5/2 1/2 2W 0758-0006 R:FXD MET FLM 1000 OHM 5/2 1/2 2W 0758-0007 R:FXD MET FLM 1/2 2W 0758-0008 R:FXD MET FLM 1/2 2W 0758-0009 R:FXD MET FLM 1/2 2W 0758-0004 R:FXD MET FLM 1/2 2W 0	ルロムしゅつつつ		2
0727-0173         R:FXD DEPC 20K OHM 1% 1/2W         284           0727-0178         R:FXD DEPC 24.7K OHM 1% 1/2W         197           0727-0189         R:FXD DEPC 41.7K OHM 1% 1/2W         284           0727-0398         R:FXD DEPC 3.79K OHM 1/2% 1/2W         284           0727-0403         R:FXD DEPC 52.3K OHM 1/2% 1/2W         197           0727-0405         R:FXD DEPC 57.46K OHM 1/2% 1/2W         197           0757-0817         R:FXD MET FLM 750 OHM 1 1/2% 1/2W         284           0757-0831         R:FXD MET FLM 4.32K OHM 1% 1/2W         284           0757-0885         R:FXD MET FLM 15.4K OHM 1% 1/4W         284           0757-0887         R:FXD MET FLM 12K OHM 1% 1/4W         284           0757-0889         R:FXD MET FLM 143K OHM 1% 1/4W         284           0758-0002         R:FXD MET FLM 191K OHM 5% 1/2W         071           0758-0003         R:FXD MET FLM 1000 OHM 5% 1/2W         071           0758-0005         R:FXD MET FLM 1000 OHM 5% 1/2W         071           0758-00017         R:FXD MET FLM 1500 OHM 5% 1/2W         071           0758-0012         R:FXD MET FLM 1500 OHM 5% 1/2W         071           0758-0034         R:FXD MET FLM 200 OHM 5% 1/2W         071           0758-0043         R:FXD MET FLM 1800 OHM 5% 1/2W		7-0163	1
0727-0178       R:FXD DEPC 24.7K OHM 1% 1/2W       197         0727-0189       R:FXD DEPC 41.7K OHM 1% 1/2W       284         0727-0398       R:FXD DEPC 3.79K OHM 1/2% 1/2W       284         0727-0403       R:FXD DEPC 52.3K OHM 1/2% 1/2W       197         0727-0405       R:FXD DEPC 57.46K CHM 1/2% 1/2W       197         0757-0817       R:FXD MET FLM 750 OHM 1* 1/2% 1/2W       284         0757-0831       R:FXD MET FLM 4.32K OHM 1* 1/2W       284         0757-0885       R:FXD MET FLM 15.4K OHM 1* 1/4W       284         0757-0887       R:FXD MET FLM 32.4K OHM 1* 1/4W       284         0757-0889       R:FXD MET FLM 143K OHM 1* 1/4W       284         0757-0890       R:FXD MET FLM 191K OHM 1* 1/4W       284         0758-0002       R:FXD MET FLM 560 OHM 5* 1/2W       071         0758-0003       R:FXD MET FLM 1000 OHM 5* 1/2W       071         0758-0004       R:FXD MET FLM 2700 OHM 5* 1/2W       071         0758-0005       R:FXD MET FLM 12K OHM 5* 1/2W       071         0758-0017       R:FXD MET FLM 12K OHM 5* 1/2W       071         0758-0024       R:FXD MET FLM 2400 OHM 5* 1/2W       071         0758-0038       R:FXD MET FLM 2400 OHM 5* 1/2W       071         0758-0043       R:FXD MET FLM 1800 OHM 5* 1	701 DC 17		1
0727-0189       R:FXD DEPC 41.7K OHM 1% 1/2W       284         0727-0398       R:FXD DEPC 3.79K OHM 1/2% 1/2W       284         0727-0403       R:FXD DEPC 52.3K OHM 1/2% 1/2W       197         0727-0405       R:FXD DEPC 57.46K OHM 1/2% 1/2W       197         0757-0817       R:FXD MET FLM 750 OHM 1/2 1/2W       284         0757-0831       R:FXD MET FLM 4.32K OHM 1/2 1/2W       284         0757-0885       R:FXD MET FLM 15.4K OHM 1/2 1/4W       284         0757-0887       R:FXD MET FLM 32.4K OHM 1/2 1/4W       284         0757-0889       R:FXD MET FLM 143K OHM 1/2 1/4W       284         0758-0002       R:FXD MET FLM 191K OHM 1/2 1/4W       284         0758-0003       R:FXD MET FLM 1000 OHM 5/2 1/2W       071         0758-0004       R:FXD MET FLM 2700 OHM 5/2 1/2W       071         0758-0005       R:FXD MET FLM 12K OHM 5/2 1/2W       071         0758-0012       R:FXD MET FLM 1500 OHM 5/2 1/2W       071         0758-0024       R:FXD MET FLM 8/2K OHM 5/2 1/2W       071         0758-0034       R:FXD MET FLM 2400 OHM 5/2 1/2W       071         0758-0044       R:FXD MET FLM 1800 OHM 5/2 1/2W       071         0758-0047       R:FXD MET FLM 1800 OHM 5/2 1/2W       071         0758-0047       R:FXD MET FLM	480 0727	7-0173	1
0727-0398       R:FXD DEPC 3.79k OHM 1/2% 1/2W       284         0727-0403       R:FXD DEPC 52.3k OHM 1/2% 1/2W       197         0727-0405       R:FXD DEPC 57.46k CHM 1/2% 1/2W       197         0757-0817       R:FXD MET FLM 750 OHM 1* 1/2%       284         0757-0831       R:FXD MET FLM 4.32k OHM 1* 1/2W       284         0757-0885       R:FXD MET FLM 15.4k OHM 1* 1/4W       284         0757-0887       R:FXD MET FLM 32.4k OHM 1* 1/4W       284         0757-0889       R:FXD MET FLM 191k OHM 1* 1/4W       284         0758-0002       R:FXD MET FLM 191k OHM 1* 1/4W       284         0758-0003       R:FXD MET FLM 1000 OHM 5% 1/2W       071         0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2W       071         0758-0005       R:FXD MET FLM 12K OHM 5% 1/2W       071         0758-0012       R:FXD MET FLM 1500 OHM 5% 1/2W       071         0758-0024       R:FXD MET FLM 100 OHM 5% 1/2W       071         0758-0034       R:FXD MET FLM 2400 OHM 5% 1/2W       071         0758-0038       R:FXD MET FLM 1800 OHM 5% 1/2W       071         0758-0044       R:FXD MET FLM 1200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 1800 OHM 5% 1	701 001/		1
0727-0403       R:FXD DEPC 52.3K OHM 1/2% 1/2%       197         0727-0405       R:FXD DEPC 57.46K OHM 1/2% 1/2%       197         0757-0817       R:FXD MET FLM 750 OHM 1% 1/2%       284         0757-0831       R:FXD MET FLM 4.32K OHM 1% 1/2%       284         0757-0885       R:FXD MET FLM 15.4K OHM 1% 1/4%       284         0757-0887       R:FXD MET FLM 32.4K OHM 1% 1/4%       284         0757-0889       R:FXD MET FLM 143K OHM 1% 1/4%       284         0758-0002       R:FXD MET FLM 560 OHM 5% 1/2%       071         0758-0003       R:FXD MET FLM 1000 OHM 5% 1/2%       071         0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2%       071         0758-0005       R:FXD MET FLM 12K OHM 5% 1/2%       071         0758-0012       R:FXD MET FLM 1500 OHM 5% 1/2%       071         0758-0022       R:FXD MET FLM 1500 OHM 5% 1/2%       071         0758-0034       R:FXD MET FLM 2400 OHM 5% 1/2%       071         0758-0038       R:FXD MET FLM 2400 OHM 5% 1/2%       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2%       071         0758-0043       R:FXD MET FLM 2200 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1	480 0727	7-0189	2
0727-0405       R:FXD DEPC 57.46K CHM 1/2% 1/2%       197         0757-0817       R:FXD MET FLM 750 OHM 1* 1/2%       284         0757-0831       R:FXD MET FLM 4.32K OHM 1% 1/2%       284         0757-0885       R:FXD MET FLM 15.4K OHM 1% 1/4%       284         0757-0887       R:FXD MET FLM 32.4K OHM 1% 1/4%       284         0757-0889       R:FXD MET FLM 143K OHM 1% 1/4%       284         0758-0002       R:FXD MET FLM 191K OHM 1% 1/2%       071         0758-0003       R:FXD MET FLM 560 OHM 5% 1/2%       071         0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2%       071         0758-0005       R:FXD MET FLM 12K OHM 5% 1/2%       071         0758-0012       R:FXD MET FLM 12K OHM 5% 1/2%       071         0758-0022       R:FXD MET FLM 1500 OHM 5% 1/2%       071         0758-0034       R:FXD MET FLM 2400 OHM 5% 1/2%       071         0758-0034       R:FXD MET FLM 9100 OHM 5% 1/2%       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2%       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 1800 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 1800 OHM 5% 1/	480 0727	7-0398	1
0757-0817       R:FXD MET FLM 750 OHM 1> 1/2M       284         0757-0831       R:FXD MET FLM 4.32K OHM 1% 1/2M       284         0757-0885       R:FXD MET FLM 15.4K OHM 1% 1/4W       284         0757-0887       R:FXD MET FLM 32.4K OHM 1% 1/4W       284         0757-0889       R:FXD MET FLM 143K OHM 1% 1/4W       284         0758-0002       R:FXD MET FLM 191K OHM 1% 1/4W       284         0758-0003       R:FXD MET FLM 560 OHM 5% 1/2W       071         0758-0004       R:FXD MET FLM 1000 OHM 5% 1/2W       071         0758-0005       R:FXD MET FLM 2700 OHM 5% 1/2W       071         0758-0012       R:FXD MET FLM 1500 OHM 5% 1/2W       071         0758-0022       R:FXD MET FLM 1500 OHM 5% 1/2W       071         0758-0024       R:FXD MET FLM 2400 OHM 5% 1/2W       071         0758-0034       R:FXD MET FLM 2400 OHM 5% 1/2W       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2W       071         0758-0044       R:FXD MET FLM 1800 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1	701   UC 17	/2A	1
0757-0831       R:FXD MET FLM 4.32K OHN 1% 1/2%       284         0757-0885       R:FXD MET FLM 15.4K OHM 1% 1/4%       284         0757-0887       R:FXD MET FLM 32.4K OHM 1% 1/4%       284         0757-0889       R:FXD MET FLM 191K OHM 1% 1/4%       284         0757-0890       R:FXD MET FLM 191K OHM 1% 1/4%       284         0758-0002       R:FXD MET FLM 560 OHM 5% 1/2%       071         0758-0003       R:FXD MET FLM 1000 OHM 5% 1/2%       071         0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2%       071         0758-0005       R:FXD MET FLM 12K OHM 5% 1/2%       071         0758-0012       R:FXD MET FLM 1500 OHM 5% 1/2%       071         0758-0022       R:FXD MET FLM 82K OHM 5% 1/2%       071         0758-0024       R:FXD MET FLM 2400 OHM 5% 1/2%       071         0758-0034       R:FXD MET FLM 9100 OHM 5% 1/2%       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2%       071         0758-0044       R:FXD MET FLM 2200 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2%       071	701 DC1/	/24	1
0757-0885       R:FXD MET FLM 15.4K OHM 1% 1/4%       284         0757-0887       R:FXD MET FLM 32.4K OHM 1% 1/4%       284         0757-0889       R:FXD MET FLM 143K OHM 1% 1/4%       284         0757-0890       R:FXD MET FLM 191K OHM 1% 1/4%       284         0758-0002       R:FXD MET FLM 560 OHM 5% 1/2%       071         0758-0003       R:FXD MET FLM 1000 OHM 5% 1/2%       284         0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2%       071         0758-0005       R:FXD MET FLM 2700 OHM 5% 1/2%       071         0758-0012       R:FXD MET FLM 12K OHM 5% 1/2%       071         0758-0017       R:FXD MET FLM 1500 OHM 5% 1/2%       071         0758-0024       R:FXD MET FLM 22K OHM 5% 1/2%       071         0758-0034       R:FXD MET FLM 2400 OHM 5% 1/2%       071         0758-0038       R:FXD MET FLM 9100 OHM 5% 1/2%       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2%       071         0758-0044       R:FXD MET FLM 2200 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2%       071	480 6757	7-0817	
0757-0887       R:FXD MET FLM 32.4k OHM 1% 1/4m       284         0757-0889       H:FXD MET FLM 143K OHM 1% 1/4w       284         0757-0890       R:FXD MET FLM 191K OHM 1% 1/4w       284         0758-0002       R:FXD MET FLM 560 OHM 5% 1/2w       071         0758-0003       R:FXD MET FLM 1000 OHM 5% 1/2w       201         0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2w       071         0758-0005       R:FXD MET OX 4700 OHM 5% 1/2w       071         0758-0012       R:FXD MET FLM 12K OHM 5% 1/2w       071         0758-0017       R:FXD MET FLM 1500 OHM 5% 1/2w       071         0758-0022       R:FXD MET FLM 82K OHM 5% 1/2w       071         0758-0034       R:FXD MET FLM 2400 OHM 5% 1/2w       071         0758-0038       R:FXD MET FLM 9100 OHM 5% 1/2w       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2w       071         0758-0044       R:FXD MET FLM 2200 OHM 5% 1/2w       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2w       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2w       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2w       071	480 0757	7-0831	1
0757-0889       R:FXD MET FLM 143K OHM 1% 1/4W       284         0757-0890       R:FXD MET FLM 191K OHM 1% 1/4W       284         0758-0002       R:FXD MET FLM 560 OHM 5% 1/2W       071         0758-0003       R:FXD MET FLM 1000 OHM 5% 1/2W       264         0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2W       071         0758-0005       R:FXD MET OX 4700 OHM 5% 1/2W       071         0758-0012       R:FXD MET FLM 12K OHM 5% 1/2W       071         0758-0017       R:FXD MET FLM 1500 OHM 5% 1/2W       071         0758-0022       R:FXD MET FLM 82K OHM 5% 1/2W       071         0758-0034       R:FXD MET FLM 100 OHM 5% 1/2W       071         0758-0038       R:FXD MET FLM 9100 OHM 5% 1/2W       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2W       071         0758-0044       R:FXD MET FLM 2200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2W       071	480 0757	7-0885	1
0757-C890       R:FXD MET FLM 191K OHM 1% 1/4W       284         0758-0002       R:FXD MET FLM 560 OHM 5% 1/2W       071         0758-0003       R:FXD MET FLM 1000 OHM 5% 1/2W       284         0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2W       071         0758-0005       R:FXD MET OX 4700 OHM 5% 1/2W       071         0758-0012       R:FXD MET FLM 12K OHM 5% 1/2W       071         0758-0017       R:FXD MET FLM 1500 OHM 5% 1/2W       071         0758-0022       R:FXD MET FLM 82K OHM 5% 1/2W       071         0758-0034       R:FXD MET FLM 100 OHM 5% 1/2W       071         0758-0038       R:FXD MET FLM 9100 OHM 5% 1/2W       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2W       071         0758-0044       R:FXD MET FLM 2200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2W       071	480 0757	7-0887	1
0758-0002       R:FXD MET FLM 560 OHM 5% 1/2%       071         0758-0003       R:FXD MET FLM 1000 OHM 5% 1/2%       264         0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2%       071         0758-0005       R:FXD MET FLM 2700 OHM 5% 1/2%       071         0758-0012       R:FXD MET FLM 12K OHM 5% 1/2%       071         0758-0017       R:FXD MET FLM 1500 OHM 5% 1/2%       071         0758-0022       R:FXD MET FLM 82K OHM 5% 1/2%       071         0758-0024       R:FXD MET FLM 100 OHM 5% 1/2%       071         0758-0034       R:FXD MET FLM 2400 OHM 5% 1/2%       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2%       071         0758-0044       R:FXD MET FLM 1800 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 2200 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2%       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2%       071	480 0757	7-0889	1
0758-0003       R:FXD MET FLM 1000 OHM 5% 1/2W       264         0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2W       071         0758-0005       R:FXD MET OX 4700 OHM 5% 1/2W       284         0758-0012       R:FXD MET FLM 12K OHM 5% 1/2W       071         0758-0017       R:FXD MET FLM 1500 OHM 5% 1/2W       071         0758-0022       R:FXD MET FLM 82K OHM 5% 1/2W       284         0758-0034       R:FXD MET FLM 100 OHM 5% 1/2W       071         0758-0038       R:FXD MET FLM 2400 OHM 5% 1/2W       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2W       071         0758-0044       R:FXD MET FLM 2200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2W       071	480 0757	7-0890	1
0758-0004 R:FXD MET FLM 2700 OHM 5% 1/2% 284  0758-0005 R:FXD MET FLM 12K OHM 5% 1/2% 071  0758-0012 R:FXD MET FLM 12K OHM 5% 1/2% 071  0758-0017 R:FXD MET FLM 1500 OHM 5% 1/2% 071  0758-0022 R:FXD MET FLM 82K OHM 5% 1/2% 071  0758-0024 R:FXD MET FLM 100 OHM 5% 1/2% 071  0758-0034 R:FXD MET FLM 2400 OHM 5% 1/2% 071  0758-0038 R:FXD MET FLM 9100 OHM 5% 1/2% 071  0758-0043 R:FXD MET FLM 1800 OHM 5% 1/2% 071  0758-0044 R:FXD MET FLM 2200 OHM 5% 1/2% 071  0758-0047 R:FXD MET FLM 2200 OHM 5% 1/2% 071	115 C 20		11
0758-0004       R:FXD MET FLM 2700 OHM 5% 1/2W       071         0758-0005       R:FXD MET OX 4700 OHM 5% 1/2W       084         0758-0012       R:FXD MET FLM 12K OHM 5% 1/2W       071         0758-0017       R:FXD MET FLM 1500 OHM 5% 1/2W       071         0758-0022       R:FXD MET FLM 82K OHM 5% 1/2W       071         0758-0034       R:FXD MET FLM 100 OHM 5% 1/2W       071         0758-0038       R:FXD MET FLM 9100 OHM 5% 1/2W       071         0758-0043       R:FXD MET FLM 1800 OHM 5% 1/2W       071         0758-0044       R:FXD MET FLM 2200 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2W       071         0758-0047       R:FXD MET FLM 7500 OHM 5% 1/2W       071	480 J758		Ī
0758-0005       R:FXD MET 0X 4700 0HM 5% 1/2W       284         0758-0012       R:FXD MET FLM 12K 0HM 5% 1/2W       071         0758-0017       R:FXD MET FLM 1500 0HM 5% 1/2W       071         0758-0022       R:FXD MET FLM 82K 0HM 5% 1/2W       284         0758-0034       R:FXD MET FLM 100 0HM 5% 1/2W       071         0758-0038       R:FXD MET FLM 2400 0HM 5% 1/2W       071         0758-0043       R:FXD MET FLM 1800 0HM 5% 1/2W       071         0758-0044       R:FXD MET FLM 2200 0HM 5% 1/2W       071         0758-0047       R:FXD MET FLM 2200 0HM 5% 1/2W       071         0758-0047       R:FXD MET FLM 7500 0HM 5% 1/2W       071         0758-0047       R:FXD MET FLM 7500 0HM 5% 1/2W       071	115 0 20		1
0758-0017 0758-0022 R:FXD MET FLM 1500 OHM 5% 1/2% 0758-0024 R:FXD MET FLM 82K OHM 5% 1/2% 0758-0034 R:FXD MET FLM 100 OHM 5% 1/2% 0758-0038 R:FXD MET FLM 2400 OHM 5% 1/2% 0758-0043 R:FXD MET FLM 9100 OHM 5% 1/2% 0758-0044 R:FXD MET FLM 1800 OHM 5% 1/2% 071 0758-0047 R:FXD MET FLM 2200 OHM 5% 1/2% 071 0758-0047 R:FXD MET FLM 2200 OHM 5% 1/2% 071 0758-0047	480 0758		4
0758-0017 0758-0022 0758-0022 0758-0024 0758-0024 0758-0034 0758-0038 0758-0038 0758-0043 0758-0043 0758-0043 0758-0047	115 6 20	)	2
0758-0022 R:FXD MET FLM 82K OHM 5% 1/2W 284 0758-0024 R:FXD MET FLM 100 OHM 5% 1/2W 071 0758-0034 R:FXD MET FLM 2400 OHM 5% 1/2W 071 0758-0038 R:FXD MET FLM 9100 OHM 5% 1/2W 071 0758-0043 R:FXD MET FLM 1800 OHM 5% 1/2W 071 0758-0044 R:FXD MET FLM 2200 OHM 5% 1/2W 071 0758-0047 R:FXD MET FLM 7500 OHM 5% 1/2W 071	115 C 20		1
0758-0024 R:FXD MET FLM 100 0HM 5% 1/2% 071 0758-0034 R:FXD MET FLM 2400 0HM 5% 1/2% 071 0758-0038 R:FXD MET FLM 9100 0HM 5% 1/2% 071 0758-0043 R:FXD MET FLM 1800 0HM 5% 1/2% 071 0758-0044 R:FXD MET FLM 2200 0HM 5% 1/2% 071 0758-0047 R:FXD MET FLM 7500 0HM 5% 1/2% 071	180 0758		
0758-0034 R:FXD MET FLM 2400 OHM 5% 1/2W 071 0758-0038 R:FXD MET FLM 9100 OHM 5% 1/2W 071 0758-0043 R:FXD MET FLM 1800 OHM 5% 1/2W 071 0758-0044 R:FXD MET FLM 2200 OHM 5% 1/2W 071 0758-0047 R:FXD MET FLM 7500 OHM 5% 1/2W 071	115 0 20		2
0758-0043 R:FXD MET FLM 1800 OHM 5% 1/2W 071 0758-0044 R:FXD MET FLM 2200 OHM 5% 1/2W 071 0758-0047 R:FXD MET FLM 7500 OHM 5% 1/2W 071	115 6 20		2
0758-0043 R:FXD MET FLM 1800 OHM 5% 1/2W 071 0758-0044 R:FXD MET FLM 2200 OHM 5% 1/2W 071 0758-0047 R:FXD MET FLM 7500 OHM 5% 1/2W 071	L15 C 20	)	2
0758-0044 R:FXD MET FLM 2200 OHM 5% 1/2W 071 0758-0047 R:FXD MET FLM 7500 OHM 5% 1/2W 071	15 C 20		1
0758-0047 R:FXD MET FLM 7500 OHM 5% 1/2W U71	15 C 20		2
	15 6 20		2
	15 0 20		1
0758-0057 R:FXD MET FLM 5600 OHM 5% 1/2W 071	15 C 20	1	
	15 C 32		
	180 0761		2
	37 RS5	. +917	1
	,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,,	3	2
For Service Manuals Contact MAURITRON TECHNICAL SERVICES 8 Cherry Tree Rd, Chinnor Oxio A/Q AQY Tet-01844-351694 Farr 01844-352554 Email: engulació mauritron.co.uk	134 VTA		

₩ Stock No.	Description #	Mfr.	Men D	1
		WIII.	Mfr. Part No.	TQ
0852-0020	ALTENDED	1		
0852-0021	RITEMPERATURE COMPENSATING 150 CHM 5% 250	0120	5 TM1 40 150 000	
	I THIS TENT ENATURE COMPENSATING TOO OUM FOR THE		5 TM1/8 150 OHM-5%	1
0905-0050	I SASKETTELT DEALK SZID INCH THICK		5 TM 1/8 100 OHM-5%	2
1200-0037	SUCKETICRY TURE	8547	1   OBD#	1
1200-0043	INSULATOR: TRANSISTOR ANODIZED ALUMINUM	7282	5 97094	1
	ALIAN AND TEED MEDITION	76530	294457	2
1200-0050	PINICRT SOCKET	ĺ		-
1200-0076	INSULATOR: TRANISTOR	28480	1200-0050	14
1200-0081	BUSHING: INSULATOR NYLON	02735	DF14AC	1 1
1200-0085	COVER BLATE COT CONST.	26365	974SPECIAL	4
1200-0087	COVER PLATE:CRT SOCKET 14-CONTACT CLAMP:TRANSISTOR	72825	9109-1	3
	CLAM • IRANSISTOR	02735	DF-13-A	1
1200-0092	Distriction	1	, 13-R	4
	BUSHINGITRANSISTOR	02735	495334-1	1 1
1250-0053	CAPICONNECTOR BNC WITH CHAIN	02/33	495534-1	8
1250-0083	CONNECTOR:BNC	91737	CW123A/U	1
1250-0149	CONNECTORIRE JACK CHASSIS RIGHT ANGLE	91737	UG-1094/U	2
1250-0171	CONNECTOR BNC JACK	91737	UG-1174/U	1
	The second secon	91737	11246	4
1250-0212	CONNECTOR JACK CHASSIS BNC	1	}	-
1250-0228	CONNECTOR RE JACK CHASSIS	91737	11656	1.1
1250-0229	CONNECTOR RE CARLE DIVISION	94735	RF 6650	1
1251-0148	CONNECTOR RE CABLE PLUG SUB-MINIATURE	94735	RF 6621-27	4
1400-0084	CONNECTORIPOWER 3 PIN MALE	60427	H-10611G-3L	4
	FUSEHOLDER EXTRACTOR POST TYPE	75015	342014	1
1520-000+	1	12472	272014	1
1520-0001	PLATE : MOUNTING ELECTROLYTIC CAPACITOR	20000	1520 00	
1850-0040	IMANGISTORIGERMANIUM 2N3RT PNP	28480	1520-0001	3
1850-0062	TRANSISTOR:GERMANIUM ALLAY JUNICTION	94154	2N383	2
1850-0064	I TANDIDIORIGERMANIUM DNITER DND	28480	1850-0062	5
1850-0065	TRANSISTOR: GERMANIUM 2N1370	02735	2N1183	2
	TOTAL SECURISTOR ENTOYO	01295	2N1370	4
1850-0090	TRANSISTOR:GERMANIUM 2N11838 PNP	1		7
1850-0098	TRANSISTORICE PRANTILLA DUE ATTENTA	86684	2N1183B	
1850-0153	TRANSISTOR GERMANIUM PNP SELECTED	28480	1850-0098	2
851-0017	TRANSISTORIPHP SM1642	28480	1850-0153	2
1853-0003	TRANSISTOR: 2N1304	01205	2N1304	7
.000	TRANSISTORIPHP SILICON F 50MC MIN	28480	1853-0003	6
1854-0003				3
1854-0005	TRANSISTOR:NPN SILICON	28/190	1854-0003	
	TRANSISTOR: 2N708 NPN STI TOON	07263	1034-0003	7
854-0022	IRANSISTORINEN SILITON	0/203	∠N/U8	6
854-0033	TRANSISTORISILICON NPN 2N3301	28480	1854-0022	7
901-0025	SEMICON DEVICE: DIODE JUNCTION	03508	2N3391	i
		28480	1901-0025	13
901-0029	SEMICON DEVICE:DIODE SI 600V	1		
901-0033	SEMICON DEVICE DIONE ST 1000	28480	1901-0029	
901-0045	SEMICON DEVICE: DIODE SILICON 1N485 B	07910	1N485B	1
901-0047	SEMICON DEVICE DIODE SILICON	28480	1901-0045	2
901-0049	SEMICON DEVICE DIODE JUNCTION	28480	1901-0047	2
,-1 0079	SEMICON DEVICE DIODE SILICON	28480	1901-0047	4
201-0050	· ·	-5460	·	4
901-0059	SEMICON DEVICE:DIODE 1N629	03077	Nego	
901-0096	SEMICON DEVICE DIODE SILICON	03877	N029	1
01-0142	RECTIFIER:SILICON	28480 1	901-0096	4
901-0162	SEMICON DEVICE: DIODE SOLD IN SET OF C	28480 1	901-0142	4
902-0017	SEMICON DEVICE DIODE SI	28480 1	901-0162	6
ļ		28480 1	902-0017	4
002-0025	SEMICON DEVICE DIODE SILICON			7
02-0050	SEMICON DEVICE DIOUR SILICON	28480 1	902-0025	
02-0241	SEMICON DEVICE: DIODE SI JUNG 8.66V 5%	28480 1	902-0050	2
10-0011	SEMICON DEVICE DIODE SILICON 100V 5%	28480 1	902-0241	1
00-0019	SEMICON DEVICE:DIODE GERMANTUM	28480 1	910-0011	1
0017	RIVAR COMP 500 OHM 10% LIN 1/2W	28480 2	100-0019	4
			-	1

Table 6-3. Replaceable Parts (cont'd)

The state of the s

⊕ Stock No.	Description#	Mfr.	Mfr. Part No.	TQ
	DAVAR COMP. 10K. OUM. 10K. LTM. 3H.	29/190	2100-0027	.
2100-0027	R:VAR COMP 10K OHM 10% LÎN 2W R:VAR COMP 1000 OHM LÎN		2100-0027	1
2100-0036	RIVAR COMP 1000 OHM LIN RIVAR COMP 2500 OHM 10% LIN 172W		2100-0056	1
2100-0067	ŘIVAR COMP 2500 OHM 10% LIN 1/2W		2100-0092	1
2100-0092 2100-0095	RIVAR COMP TOK OHM 20% CIN 1/5W		2100-0092	1 1
2100-0075				•
2100-0107	REVAR COMP 50K OHM 30% 173W VERNIER		2100-0107	1
2100-0144	REVAR COMP 250K OHM 30% LIN 2/5W		2100-0144	1
2100-0150	R:VAR 2-SECT 10K OHM 20% LIN 1/4W		2100-0150	1
2100-0154	R:VAR COMP 1000 OHM 30% LIN 0+15#		2100-0154	3
2100-0189	RIVAR COMP 1 MEGOHM 30% LIN 1/4W	28480	2100-0189	1
2100-0218	RIVAR COMP 1.2 MEGOHM 20% LIN 2W	28480	2100-0218	1
2100-0893	RIVAR 2K(FRONT) 750K(REAR) OHM LIN 1/2W		2100-0893	2
2100-0910	RIVAR COMP 2X35K OHM LIN 20% 1/4W		SERIES 5 TYPE 71-2	6
2100-0956	RIVAR COMP 500 OHM 20% LTN 1/20W	28480	2100-0956	1
2100-0957	REVAR COMP 5K OHM 20% LIN 1/20W	28480	2100-0957	2
	CANAR COMP TOK SUM COM C FW	201100	2100-0958	
2100-0958	RIVAR COMP 10K OHM 20% 0.5W		#313.600	1 1
2110-0016	FUSE:CARTRIDGE 0.6 AMP SLO-BLO		TYPE MDL	1
2110-0044	FUSE CARTRIDGE 0.3 AMP SLO-BLO		NE 2E1	1
2140-0018 3100-0812	LAMP:GLOW 1/10W ROTARY SWITCH		3100-0812	2
			7100 0015	١.
3100-0815	ROTARY SWITCH:2-SECTION 3-POSITION		3100-0815 3100-1500	
3100-1500	SWITCH:ROTARY	42190		1 1
3101-0033	SWITCH:SLIDE DPDT		3101 <b>-005</b> 2	li
3101-0052	SWITCH: PUSHBUTTON, NORMALLY OPEN		3130-0041	i
3130-0041	SWITCH SHIELD		4320-0007	∣ i
<b>43200007</b>	EXTRUSION:RUBBER COIL:BRACKET	28480	5000-0408	lî
5000-0408 5040-0218	COUPLING:MECHANICAL	28480	5040-0218	î
5040-0216 5040-0234	LAMPHOLDER		5040-0234	li
5040-0235	BASE:LAMPHOLDER	28480	5040-0235	1 1
5040-0400 5040-0400	SUPPORT: CAPACITOR	28480	5040-0400	1 2 2 1
5040-0401	SUPPORT:CAPACITOR	28480	5040-0401	∣ <u>ā</u>
5040-0418	INSULATOR:POTENTIOMETER	28480	5040-0401 5040-0418	1
5040-0421	INSULATOR : POTENTIOMETER		5040-0421	1
5060-0409	COIL ALIGNMENT		5060-0409	1
5083-0624	ELECTRON TUBE: CATHODE-RAY P-2 PHOSPHOR	. –	5083-0624	1
5083-0634	ELECTRON TUBE : CATHODE-RAY P-7 PHOSPHOR	28480	5083-0634	1
E007-045#	ELECTRON TUBE: CATHODE-RAY P-31 PHOSPHOR	28480	5083-0654	1
5083-0654 9100-0274	TRANSFORMER: POWER		6-2463	1
9110-0274	CHOKE: FILTER 70 MH 1.0 AMP 1.5 OHM		9110-0042	1
9120-0090	TRANSFORMER: INTERMEDIATE FREQUENCY		9120-0090	1
9120-0092	TRANSFORMER : AUDIO		9120-0092	1
	0071.4570.400.441	29#90	9140:0051	2
9140-0051	COIL:FXD 400 UH		9140-0082	2
9140-0082	COIL:FXD RF 15 UH COIL:FXD 500 UH 5%		2500-14	1
9140-0118	COIL: FXD RF 1 MH		9140-0137	3
9140-0137 9140-0146	COIL: FXD RF 10.0 UH		9140:0146	12
. 1 . 5 . 6 5			0.110.01.110	5
9140-0149	COIL:FXD RF 1.86 UH	28480	9140:0149 9140:0150	1
9140-0150 9140-0152	COIL:FXD RF 2.7 UH COIL:FXD RF 41.06 UH	28480	9140:0152	1
9140-0158	CATLATUD & OUL SOM	99800	1025-20	1
9140-0159	COIL:FXD 0.47UH 20% For Service Manuals Contact MAURITRON TECHNICAL SERVICES 8 Cherry Tree Rd. Chinnor Oxon OX9 40Y	99800	1025-SERIES	2
	Oxon CX11 4-07  Tet- 01844-351624 F 1x- 01844-352554  Email: enquies@mauritron.co.uk			
			}	

Table 6-3. Replaceable Parts (cont'd)

⊕ Stock No.	Description #	Mfr.	Mfr. Part No.	7
			1 1 1 10.	
9140-0232	COILIRF TAPPED 0-254UH-0-50UH			
9140-0235	COIL RF TAPPED 0.95-1.8UH	28480	9140-0232	1
08551-2083	BUSHING	28480	9140-0235	
00140-61606	CABLE HIGH VOLTAGE, INCLUDES 270K RES	28480	08551-2083	1
00851-0006	BRACKET: POWER SUPPLY	28480	00140-61606	- 1
	SHARK SUPPLY	28480	00851-0006	
00851-0007	SHIELD:HIGH VOLTAGE		00031-0006	- 1
00851-0008	COVER SOCKET	28480	00851-0007	- 1
00851-0009	BOACKETACHER	20400	00851-0008	
00851-0013	BRACKET: SWEEP AND HORIZ AMPL PCBD A6	28480	00851-0008	
00851-0014	I COAFLIAGE IN THE SEA TO MUNICIPALE	20400	00851-0009	1
	COVER SWITCH IF GAIN	20400	00851-0013	1
00851-0015	D 475	20400	00851-0014	1 2
00851-0016	PLATE COVER IF GAIN	20,00		- }
00851-0017	BRACKETTIF GAIN (DB) SWITCH	28480	00851-0015	2
00851-2004	COVER INPUL BP FILTER ALS	28480	00851-0016	1
00851-2005	PLANK PC BOARD LOW VOLT AGED SHOPE	28480	00851-0017	i
00931-5003	BLANK PC BOARD . SWEEP & HORIZ AMPL	28480	00851-2004	Ī
00051 2004		28480	00851-2005	li
00851-2006	BLANK PC BOARD HV POWER SUPPLY			
00851-2007 00851-2008	PLANN PU BUARDIAMBITETED 13444	28480	00851-2006	1
	I DENIN TO FURRUIFIRST, 1-10KC, SP EV	28480	00851-2007	li
00851-2009		28480	00851-2008	li
00851-2010	BOARD CURRENT-CONTROLLED ATTEN	28480	00851-2nna	li
		28480	00851-2010	li
0851-2011	BLANK PC BOARD: IF 20MC			1
0851-2013	BLANK PC BOARD: VERT. AMP.	28480	00851-2011	١.
0851-2014	BLANK PC BOARD: INPUT SWITCHING CIRCUIT	28480 (	00851-2013	1
0851-2015	BLANK PC BOARD-OUTPUT SWITCHING CIRCUIT	28480 (	00851-2014	1
0851-2016	BOARD INPUT B.P. FILTER	28480	00851-2015	1
	TOTAL BAP FILIER	28480	00851-2016	1
0851-2022	CAVITY:FILTER		2016	1
0851-2026	FILTER:CRT LT. BLUE	28480 0	0851-2022	1 1
0851-2027	KNOB: IF GAIN 0-70 DB	28480 0	0851-2026	2
0851-2028	KNOB: IF GAIN 0-10 DB	28480 0	0851-2027	1
0851-6001	HV POWER SUPPLY ASSY	28480 0	0851-2028	11
	TOWER SUPPLY ASSY	28480	0851-6001	1
0851-6002	SWITCH ACCYATE CARL	1	0051-6001	1
0851-6003	SWITCH ASSY: IF GAIN (DE) ASSY: RF CIRCUIT	28480 0	0851-6002	1 1
851-6006	SWITCH ACCY AND TO	28480 0	0851-6003	1
851-6007	SWITCH ASSY .: VERT . DISPLAY	28/180	0051-6003	1 1
851-6008	SWITCH ASSY .: I .F. BANDWIOTH	28480 0	0851-6006	1
7031 0000	CRT. SHIELD ASSEMBLY	20430 0	0851-6007	1
851-6013	CADIC	20400 0	0851-6008	1
851-6014	CABLE ASSY:ATTEN OUTPUT	20,000	N = = .	1 -
851-6015	CABLE ASSYLATTEN INPUT	28480 00	851-6013	1
021-0012	CADLE ASSY. ISWEED OUTPUT	28480 00	0851-6014	i
851-6016	CABLE ASSY.ISYNC INPLIT	28480 00	851-6015	i
851-6017	LOW VOLTAGE POWER SUPPLY ASSY.	28480 00	851-6016	i
OP4 45:-		28480 00	851-6017	1
851-6019	BOARD ASSYIVERT. AMPL. ASSY.	1 1		*
851-6020	ASSI 20 MC I.F. AMPLIFIED	28480 00	851-6019	
851-6021	ASSY CURRENT-CONTROLLED ATTEN.	28480 00	851-6020	1
351-6022	20MC AMPLIFIER ASSY	28480 00	851-6021	1
351-6023	FIRST 1=10KC BANDOASC CT.	28480 00	851-6022	1
]	FIRST 1-10KC BANDPASS FILTER ASSY.	28480 00	851-6023	1
351-6024	SECOND 1-10 KC SANDONES TO	1	-0. 0023	1
351-6025	SECOND 1-10 KC BANDPASS FILTER ASSY.	28480 00	851-602"	
351-6026	AND OF SWITCHING CIRCUIT ACCO	28480 00	951-4025	1
151-6027	DUTPUT SWITCHING CIRCUIT ASSY.	28480 00	951-4024	1
51-6028	CADLE ASSTOR IF INDICT	28480 000	351-6020	1
	FILTER ASSY. : 100KC BAND-PASS	28480 001	221-6027	1
1		28480 008	151-6028	2
1	For Service Manuals Contact MAURITRON TECHNICAL CONTACT			-1
	MAURITRON TECHNICAL SERVICES  8 Cherry Trop Did SERVICES			ŀ
ł				
	Oxon OX9 4QY		į	
	196 V 1044-35160/ E	1	j	1
	Email:- enquiries@mauritron.co.uk		j	]
1		· ·		,

<sup>#</sup> See list of abbreviations in introduction to this section

Table 6-3. Replaceable Parts (cont'd)

A CONTROL OF THE PROPERTY OF T

₩ Stock No.	Description#	Mfr.	Mfr. Part No.	TQ
00851-6029	CABLE ASSY.	28480	00851-6029	3
00851-6030	CABLE ASSY.	28480	00851-6030	1 1
00851-6031	CABLE ASSY.		00851-6031	1 1
00851-6032	CABLE ASSY.		00851-6032	1 1
00851-6033	CABLE ASSY.	28480	00851-6033	2
00851-6034	CABLE 5-INCH COAX VIDEO OUT TO VERT AMP		00851-6034	1
00851-6035	PC BOARD ASSY: INPUT BANDPASS FILTER		00851-6035	
00851-6036	CABLE ASSY:HORIZ. OUTPUT TO CRT		00851 <b>-6036</b> 00851 <b>-603</b> 7	1 1
00851-6037 00851-6038	CABLE ASSEMBLY: VERTICAL OUTPUT TO CRT BOARD ASSY: SWEEP & HORIZ AMPL		00851-6038	i
	CHITCH ACCY, CHEED TIME	28480	00851-6039	1
00851-6039	SWITCH ASSY: SWEEP TIME SWITCH ASSY: SYNC		00851-6040	i
00851-6040 00851-8001	COIL:RF FXD 0.3UH		00851-8001	i
00851-8002	COIL ORF VAR 0.3UH MAX	1 -	00851-8002	1
00851-8003	FILTER:LOW-PASS	1	00851-8003	5
00851-8004	COIL:RF	28480	00851-8004	2
00851-8005	COIL*RF		00851-8005	1
00851-8006	COIL*RF		00851-8006	1
00851-8008	COIL:RF VARIABLE	1	00851-8008	2
00851-8009	COILERF	28480	00851-8009	1
00851-8010	COIL*RF		00851-8010	1
08551-2083	BUSHING		08551-2083	11
120A-20	BEZEL*CRT.		120A-20	1
120A-83A	LIGHT FILTER: AMBER	-	120A-83A 120A-83G	1 1
120A-83G	LIGHT FILTER GREEN			
175A-83A	RETAINER : CRT. SHIELD	28480	175A-83A	1
	For Service Manuals C MAURITRON TECHNICAL SE 8 Cherry Tree Rd. Chi Oxon OX9 4QY	ERVICES innor		
	Tel:- 01844-351694 Fax:- 0184 Email:- enquiries@mauritro	14-352554	1	
	Elian:- andrues@maduso	i.co.ur		
			1	
			,	
		ŀ		
			1	
	}			
		1		L_L_